DIESE Diagnostica Senese SpA

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

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In 1932 Paul and Bunnell discovered the presence of IgM-class antibodies which agglutinate sheep and horse erythrocytes, in the serum of patients affected by mononucleosis. Later, Davidsohn differentiated the heterophile antibodies of serum sickness from those produced during mononucleosis, by introducing differential absorption with extracts of guinea pig or horse kidney (1).

The reagent contained in the kit is a suspension of latex particles coated with antigen extract obtained from bovine erythrocyte stroma. The specificity of the test is achieved by using purified antigen, thus avoiding the differential absorption of the interfering antibodies (2,3). The test is simple to perform and the results are easy to read.

3. PRINCIPLE OF THE METHOD

Latex particles are coated with the antigen which is specific for mononucleosis, obtained from bovine erythrocyte stroma. If the sample of serum or plasma contains heterophile antibodies specific for infectious mononucleosis, macroscopic agglutination will appear on mixing with the latex reagent.

4. KIT CONTENTS AND REAGENT PREPARATION:

- Reagents are sufficient for 50 determinations.
- **Bring to room temperature before use.**

LATEX LATEX 1 vial x 2.5 mL

Contents: latex particles coated with heterophile antigens extracted from bovine erythrocyte stroma. Preservative: sodium azide 0.09%.

Use: Bring to room temperature before use, and carefully mix the suspension.

CONTROL + POSITIVE CONTROL 1 vial x 0.5 mL

Contents: human serum positive for heterophile antibodies, containing sodium azide 0.09% as preservative.

CONTROL - NEGATIVE CONTROL 1 vial x 0.5 mL

Contents: human serum negative for the presence of heterophile antibodies, containing sodium azide 0.09% as preservative.

OTHER MATERIALS SUPPLIED

- Plastic-coated cardboard slides
- Mixing Sticks.

OTHER MATERIALS REQUIRED BUT NOT PROVIDED

- 40 µL pipettes.
- Normal laboratory glassware.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation**REAGENT****CONDITIONS**

LATEX

3 months at 2/8°C. **DO NOT FREEZE**

POSITIVE CONTROL SERUM	3 months at 2/8°C.
NEGATIVE CONTROL SERUM	3 months at 2/8°C.

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY. STORE AT 2-8°C

Caution: *This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.*

Waste disposal: *serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.*

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The reagents contain 0.9% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.
If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable.
4. The reaction can be performed using a rotating mixer. Regulate the speed at 100 ± 10 rpm.
5. The agglutination reaction is best interpreted by observing the slide under a lamp.

7. TYPE AND STORAGE OF SAMPLE

Fresh or defrosted serum or plasma obtained on heparin or EDTA. Samples may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers.

8. PROCEDURE

Dispense 40 µL of serum sample in a division on the card, and add one drop of Reagent. Mix carefully with a stick, spreading the mixture uniformly. Rotate and read after one minute under a lamp. A positive sample will cause agglutination of the latex suspension.

9. TEST VALIDATION

Agglutination must be obtained when the positive control is reacted with the latex suspension. No agglutination must be found when the negative control serum is tested.

10. INTERPRETATION OF THE TEST

Evident agglutination of the latex is caused by the presence of heterophile antibodies at a sufficiently high titer to indicate a probable infection by mononucleosis in the acute or subacute phase.

11. LIMITATIONS OF THE TEST

The result of the test must always be evaluated together with other data resulting from other diagnostic procedures. It is known that in patients below 3-4 years affected by infectious mononucleosis, there may be no formation of heterophile antibodies.

12. ANALYTICAL SENSITIVITY

The latex reacts by agglutinating in the presence of samples with a titer 1/56.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

The reagent was tested on 95 samples, 64 of which were negative and 31 positive. The reference result was obtained by assay in an ELISA technique for antibodies against specific EBV proteins. One false positive and 1 false negative result were found, giving a specificity and sensitivity of respectively 98.5% and 96.9%.

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