DIESSE Diagnostica Senese SpA

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

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**1. INTENDED USE**

**STAINED BACTERIAL SUSPENSIONS OF SALMONELLA AND BRUCELLA FOR THE SEMIQUANTITATIVE DETERMINATION OF SPECIFIC AGGLUTINATING ANTIBODIES. THE TEST IS PERFORMED ON A MICRPLATE.**

**2. INTRODUCTION**

Serological studies of patients with persistent fever can be useful in detecting typhoid fever and brucellosis (1,2).

The test consists in the observation of a direct agglutination reaction between the bacterial suspension and the serum sample which is suspected of containing specific antibodies. The present method, known as the Widal-Wright reaction, allows a semiquantitative titration of the serum agglutinins.

**3. PRINCIPLE OF THE METHOD**

The test involves mixing a suspension of inactivated coloured bacteria with two fold serial dilutions of test serum in U bottom wells. After incubation overnight, a positive serum determines the agglutination of bacteria. If the bacterial cells settle in compact buttons, the result of the test is negative.

**4. KIT COMPOSITION**

- The reagents are sufficient for 500 tests.
- Bring the reagents to room temperature before use.

**BACTERIAL SUSPENSION** BACTERIAL SUSPENSION 5 vials x 5 mL

**Contents:** Suspension of stained, killed bacteria at an optimized concentration for performance of the microplate test. Contains sodium azide 0,09% as preservative.

**CONTROL +** POSITIVE CONTROL 1 x 2 mL

**Contents:** Polyvalent control serum of animal origin, containing sodium azide 0,09% as preservative. Ready for use.

**OTHER MATERIALS SUPPLIED**

- 6 96-well U-bottom microtiter plates

**OTHER MATERIALS REQUIRED BUT NOT SUPPLIED**

- Pipettes with disposable needles.
- Physiological saline.

**5. STORAGE AND STABILITY OF REAGENTS**

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

**Reagents have a limited stability after opening and/or preparation**

| REAGENT              | CONDITIONS |
|----------------------|------------|
| BACTERIAL SUSPENSION | 2 months   |
| POSITIVE CONTROL     | 2 months   |

**6. PRECAUTIONS****FOR IN VITRO DIAGNOSTIC USE ONLY.**

**Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.**

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The reagents contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.
3. If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. **Allow all reagents and samples to come to room temperature (18-30°C) before use.** After use, empty the dropper and return reagents to the recommended storage temperature (2-8°C).
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure.

**7. TYPE AND STORAGE OF SAMPLE**

Fresh serum may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Use immediately after thawing.

**8. PROCEDURE**

Bring the reagent to room temperature and mix the vial of bacterial suspension carefully but not vigorously.

Mix 10 µL of test serum with 90 µL of saline in the first well of a strip. Distribute 50 µL of saline in the other wells and make serial two-fold dilutions along the row discarding the last 50 µL. Add 1 drop (50 µL) of carefully mixed homogeneous bacterial suspension to each well (see table). Mix and incubate at 37°C for at least 18 hours, covering carefully and, if possible, placing in a humid chamber to avoid evaporation. Avoid any kind of vibration. At the end of the incubation period, leave the plate to settle for at least one hour, to improve the reading.

50 µl   50 µl   50 µl   50 µl   50 µl   50 µl   50 µl   50 µl   50 µl

| 1                                | 2              | 3              | 4              | 5              | 6              | 7              | 8              | 9              | 10             | 11                | 12              |
|----------------------------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|-------------------|-----------------|
| 90 µl<br>fisiol. +<br>10µl siero | 50µl<br>fisiol | 50µl<br>fisiol | 50µl<br>fisiol | 50µl<br>fisiol | 50µl<br>fisiol | 50µl<br>fisiol | 50µl<br>fisiol | 50µl<br>fisiol | 50µl<br>fisiol | 50µl<br>Cont.Pos. | 50µl<br>fisiol. |
| 50µl<br>sospensione              |                |                |                |                |                |                |                |                |                |                   |                 |
| Titolo<br>1:20                   | 1:40           | 1:80           | 1:160          | 1:320          | 1:640          | 1:1280         | 1:2560         | 1:5120         | 1:10240        | C.P.              | NEG.            |

**9. TEST VALIDATION**

Carry out the test using 1 drop (50 µL) of positive control in a well with 1 drop of bacterial suspension. After incubation a positive reaction should be given. To control the negative reaction use 1 drop (50 µL) of saline instead of the sample. A non-agglutinated button must be clearly visible.

**10. INTERPRETATION OF THE TEST – NORMAL VALUES**Positive reaction

The agglutinated bacteria cover the bottom of the well in a homogeneous layer. When using the “H” suspension, small buttons surrounded by a layer of agglutinated bacteria has to be considered as a positive reaction.

Negative reaction

The bacteria settle in a compact brightly coloured button.

Titers of 1/80 are considered suspect; higher titers are diagnostically significant.

**11. TEST LIMITATIONS**

In the microplate agglutination test there may be a cross-reactoin between O and H suspensions, and there is therefore generally an overestimation of the H titer.

**12. DIAGNOSTIC SENSITIVITY AND SPECIFICITY**

An experimentation performed on human sera obtained from hospital laboratories, the titers of which were taken as reference, allowed an evaluation of the sensitivity of the following suspensions: S. typhi O, S. typhi H, S. paratyphi BO, S. paratyphi BH.

When S. typhi antigen H was determined on 30 positive samples, a sensitivity of 96.7% was found, while the suspension S. typhi antigen O tested on 51 positive sera gave 1 false negative and a sensitivity of 98.1%.

S. paratyphi B antigen H and antigen O were tested with respectively 12 and 24 positive samples. Sensitivity was 100% and 96.0% respectively.

S. paratyphi A antigen O and antigen H were tested in parallel with another commercial product. 12 positive and 3 doubtful results were found. The Dienes reagent gave no false negative reactions and the sensitivity was therefore 100%. Using the same criteria, 45 positive samples for Brucellosis were tested. Dienes gave 2 false negative results and the sensitivity was therefore 95.7%. The S. paratyphi B total suspension gave a sensitivity of 97.1% with 1 false negative result.

Human sera found negative in another commercial test were also studied. The following results were obtained:

| <u>Suspension</u>        | <u>No. of samples</u> | <u>Specificity %</u> |
|--------------------------|-----------------------|----------------------|
| S. typhi antigen O       | 76                    | 96.2                 |
| S. typhi antigen H       | 58                    | 95.1                 |
| S. paratyphi B antigen O | 96                    | 98.0                 |
| S. paratyphi B antigen H | 93                    | 95.9                 |
| S. paratyphi A Total     | 82                    | 98.8                 |
| S. paratyphi B Total     | 93                    | 96.9                 |
| Brucella                 | 83                    | 97,6                 |

**13. REFERENCES**

1. P. Nicoletti in F. Pasquinelli: Diagnostica e Tecnica di Laboratorio Vol. 2 p. 319. Rosini ed. Firenze 1981.
2. Edwards P.R., Ewing W.H. Identification of Enterobacteriaceae 3rd ed. Burgess Publishing Co.

|                                                                                     |                |                                                                                                                               |                |                                                                                                                                           |
|-------------------------------------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------------------|
|    | EN<br>ES<br>IT | Date of manufacture<br>Fecha de fabricación<br>Data di fabbricazione                                                          | FR<br>GR<br>PT | Date de fabrication<br>Ημερομηνία Παραγωγής<br>Data de fabrico                                                                            |
|    | EN<br>ES<br>IT | Use By<br>Fecha de caducidad<br>Utilizzare entro                                                                              | FR<br>GR<br>PT | Utiliser jusqu'e<br>Ημερομηνία λήξης<br>Prazo de validade                                                                                 |
|    | EN<br>ES<br>IT | Caution, consult accompanying documents<br>Atención, ver instrucciones de uso<br>Attenzione, vedere le istruzioni per l'uso   | FR<br>GR<br>PT | Attention voir notice d'instructions<br>Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα<br>Atenção, consulte a documentação incluída       |
|    | EN<br>ES<br>IT | Manufacturer<br>Fabricante<br>Fabbicante                                                                                      | FR<br>GR<br>PT | Fabricant<br>Κατασκευαστής<br>Fabricante                                                                                                  |
|    | EN<br>ES<br>IT | Contains sufficient for <n> tests<br>Contenido suficiente para <n> ensayos<br>Contenuto sufficiente per "n" saggi             | FR<br>GR<br>PT | Contenu suffisant pour "n" tests<br>Περιεχόμενο επαρκές για «n» εξετάσεις<br>Conteúdo suficiente para “n” ensaios                         |
|    | EN<br>ES<br>IT | Temperature limitation<br>Límite de temperatura<br>Limiti di temperatura                                                      | FR<br>GR<br>PT | Limites de température<br>Περιορισμοί θερμοκρασίας<br>Limites de temperatura                                                              |
|    | EN<br>ES<br>IT | Consult Instructions for Use<br>Consulte las instrucciones de uso<br>Consultare le istruzioni per l'uso                       | FR<br>GR<br>PT | Consulter les instructions d'utilisation<br>Συμβουλευτείτε τις οδηγίες χρήσης<br>Consulte as instruções de utilização                     |
|    | EN<br>ES<br>IT | Biological risks<br>Riesgo biológico<br>Rischio biologico                                                                     | FR<br>GR<br>PT | Risques biologiques<br>Βιολογικοί κίνδυνοι<br>Risco biológico                                                                             |
|    | EN<br>ES<br>IT | Catalogue number<br>Número de catálogo<br>Numero di catalogo                                                                  | FR<br>GR<br>PT | Référence du catalogue<br>Αριθμός καταλόγου<br>Referência de catálogo                                                                     |
|  | EN<br>ES<br>IT | In Vitro Diagnostic Medical Device<br>Producto sanitario para diagnóstico in vitro<br>Dispositivo medico-diagnostico in vitro | FR<br>GR<br>PT | Dispositif médical de diagnostic in vitro<br>In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν<br>Dispositivo médico para diagnóstico in vitro |
|  | EN<br>ES<br>IT | Batch code<br>Código de lote<br>Codice del lotto                                                                              | FR<br>GR<br>PT | Code du lot<br>Αριθμός Παρτίδας<br>Código do lote                                                                                         |

**Diesse Diagnostica Senese**  
**Via delle Rose 10**  
**53035 Monteriggioni (Siena) – Italy**  
**Tel. 0577-587111**