**CHORUS
Promonitor
VEDOLIZUMAB**

REF 86710

DIESSE

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19. BIBLIOGRAFIA

- Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
- Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
- Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S,

Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.

4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

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VDZ is able to bind to preimmobilized anti-VDZ antibodies following incubation with the sample. After washing to remove the unbound sample, HRP-labelled anti-VDZ antibody (conjugate) is added. A second incubation step allows the

conjugate to bind to the VDZ that has become attached to the microwells. After washing to remove the excess of unbound conjugate, the peroxidase substrate is added. The enzymatic reaction is stopped by adding stop solution, which makes the solution turn to yellow. The color which develops is proportional to the amount of VDZ in the patient sample.

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Use the devices within 60 min.

1. **Discard devices which show the substrate (well 2) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.

3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the CHORUS TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
The use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
5. Check that the CHORUS TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see CHORUS TRIO Operating Manual).
9. When CHORUS Autosampler (REF 81210) is used, refer to the related Operating Manual
10. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
11. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
12. Do not use the device after the expiry date.
13. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

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CALIBRATOR CALIBRATOR 1 x 0.650 ml
Contents: Human serum containing VDZ, and preservative.
Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 1.100 ml
Contents: Human serum containing VDZ, and preservative.
Liquid, ready for use.

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100.0% CI_{95%}: 91.4-99.9

Percent Negative Agreement (~Diagnostic Specificity):

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19. BIBLIOGRAPHY

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.
4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and

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19. BIBLIOGRAFIE

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. *J Crohns Colitis*. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." *Inflamm Bowel Dis*. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. *N Engl J Med*. 2013 Aug 22;369(8):711-21.
4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics,

Pharmacodynamics, and Immunogenicity of Vedolizumab. *Clin Pharmacokinet*. 2017 Nov;56(11):1287-1301.

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- Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
- Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
- Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.
- Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

INSTRUCTIONS FOR USE

CHORUS Promonitor VEDOLIZUMAB

For the quantitative detection of Vedolizumab in human serum

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

CHORUS Promonitor VEDOLIZUMAB is an immunoassay kit for automated quantitative detection of Vedolizumab (VDZ) in human serum using a disposable device applied on the CHORUS TRIO instrument.

The product, used in conjunction with other clinical and laboratory findings, is useful as an aid in the management of patients treated with Vedolizumab (VDZ), like patients with inflammatory bowel disease (IBD) and rheumatic diseases.

It must be used by professional laboratory users only.

2. INTRODUCTION

Vedolizumab (VDZ) is a humanized IgG1 monoclonal antibody expressed in Chinese hamster ovary (CHO) cells that specifically binds the $\alpha 4\beta 7$ integrin which is expressed on a discrete subset of memory T helper lymphocytes, which preferably migrate to the gastrointestinal (GI) tract, causing inflammation characteristic of ulcerative colitis and Crohn's disease. VDZ reduces inflammation by improving symptoms, characterized by more acute phases followed by clinical remission, accompanied by abdominal pain, diarrhea and bleeding that limit the social and work activities of the patients who are affected. VDZ is a biological immunosuppressant with high selectivity for intestinal tract and is indicated for the treatment of adult patients with moderate to severe active ulcerative colitis or Crohn's disease who have had an inadequate response or been intolerant to conventional therapy, or to a tumor necrosis factor alpha (TNF α) antagonist. VDZ can potentially be immunogenic and induce the formation of anti-drug antibodies (ADA), thereby potentially reducing treatment efficacy. A gradual decrease in efficacy over time following an initial response to biologics is common.

3. PRINCIPLE OF THE METHOD

The CHORUS Promonitor VEDOLIZUMAB device is ready to use for the detection of Vedolizumab (VDZ) in human serum, in the CHORUS TRIO instrument.

The assay is based on a sandwich ELISA method.

The microwell strips are pre-coated with an anti-VDZ monoclonal antibody. Calibrator, positive control and samples are added to separate device.

VDZ is able to bind to preimmobilized anti-VDZ antibodies following incubation with the sample. After washing to remove the unbound sample, HRP-labelled anti-VDZ antibody (conjugate) is added. A second incubation step allows the conjugate to bind to the VDZ that has become attached to the microwells. After washing to remove the excess of unbound conjugate, the peroxidase substrate is added. The enzymatic reaction is stopped by adding stop solution, which makes the solution turn to yellow. The color which develops is proportional to the amount of VDZ in the patient sample.

The disposable devices contain all the reagents to perform the test when applied on the CHORUS TRIO instrument.

The results are expressed in $\mu\text{g/ml}$.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the devices within 60 min.

1. **Discard devices which show the substrate (well 2) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the CHORUS TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the CHORUS TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see CHORUS TRIO Operating Manual).
9. When CHORUS Autosampler (REF 81210) is used, refer to the related Operating Manual
10. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
11. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
12. Do not use the device after the expiry date.
13. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 7: MICROPLATE WELL

Coated with anti-VDZ monoclonal antibody

Position 6: MICROPLATE WELL

Coated with anti-VDZ monoclonal antibody

Position 5: STOP SOLUTION

Contents: 0.3 M sulfuric acid solution

Position 4: CONJUGATE

HRP-labelled anti-VDZ antibody and buffer containing preservative

Position 3: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 2: TMB SUBSTRATE

Contents: Tetramethylbenzidine stabilized in citrate buffer

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.650 ml

Contents: Human serum containing VDZ, and preservative.

Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 1.100 ml

Contents: Human serum containing VDZ, and preservative.

Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 5 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be mixed carefully before use.

It is strongly recommended to avoid any harsh treatment of the samples.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

The test CHORUS Promonitor Vedolizumab (REF. 86710) can be performed simultaneously ONLY with the other CHORUS Promonitor tests.

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense the following volumes in well no. 1 of each device:

SAMPLE	130 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

At each change of batch, use a device for the calibrator.

4. Place the devices in the CHORUS TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The CHORUS TRIO instrument expresses the result in µg/ml, calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is ≥ 1.5 µg/ml

NEGATIVE: when the result is ≤ 1.1 µg/ml

EQUIVOCAL: when the result is between 1.2 and 1.4 µg/ml

Equivocal results are obtained in case of values which ranges between Limit of detection (< 1.1 µg/ml) and Low Limit of Quantitation (1.5 µg/ml).

11. LIMITATIONS

The test measures free VDZ. It does not measure VDZ bound to anti-VDZ antibodies (immunocomplexes).

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for patient management decisions and the test result should be evaluated together with the patient history, other clinical diagnostic evaluation and laboratory findings.

This product should only be used by appropriately trained laboratory professional personnel.

12. CALIBRATION RANGE

Calibration range: 1.1-60.0 µg/ml

For samples > 60.0 µg/ml retest the diluted sample in the Negative Control/Sample Diluent (REF 83607 - not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, are negative, with values < 1.1 µg/ml.

14. ANALYTICAL SPECIFICITY

3 samples (1 Negative, 1 Low Positive and 1 High Positive) were spiked with the following potentially interfering factors and then tested:

Bilirubin (4.5 mg/dl – 40 mg/dl)
Triglycerides (250 mg/dl – 1500 mg/dl)
Hemoglobin (2.5 mg/ml – 10 mg/ml)

The presence of the above interfering substances in the test serum does not alter the test result.

2 Rheumatoid Factor (RF) positive samples (200 IU/mL) were spiked with Vedolizumab:

High level – 32.4 µg/ml
Medium level – 7.2 µg/ml
Low level – < 1.1 µg/ml

The presence of RF in the test serum does not alter the test result.

15. CROSS-REACTIONS

10 samples, not containing Vedolizumab, but containing IFX (10.8 µg/ml) - 5 samples and ADL (27.0 µg/ml) - 5 samples were tested.

No cross-reactions were found.

16. METHOD COMPARISON

In an experimentation 170 samples have been tested with Diesse kit and with a competitor kit (Promonitor VDZ ELISA kit). Data are summarized in the following table:

		Reference		
		+	-	Tot.
Diesse	+	41	0	41
	-	0	129	129
	Tot.	41	129	170

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% CI_{95%}: 91.4-99.9

Percent Negative Agreement: (~Diagnostic Specificity):

100.0% CI_{95%}: 97.1-100.0

Positive Predictive Value (PPV): 100.0% CI_{95%}: 100.0-100.0

Negative Predictive Value (NPV): 100.0% CI_{95%}: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 1.

Out of 170, for 41 positive samples the correlation between Diesse kit and a competitor kit was tested.

Data are summarized in the following table:

Correlation	r	95%CI
Pearson	0.95	0.91-0.97
Spearman	0.97	0.94-0.98

The correlation between the two methods is very high.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (µg/ml)	CV%	Mean (µg/ml)	CV%
1	27.1	7.2	24.3	9.1
2	10.7	14.4	9.7	10.4
3	6.9	10.1	7.2	10.1
4	1.5	9.3	1.9	12.1
5	<1.1	-	<1.1	-
6	<1.1	-	<1.1	-

Sample	Between lots		Between instruments	
	Mean (µg/ml)	CV%	Mean (µg/ml)	CV%
1	34.5	9.2	25.6	14.3
2	13.6	10.6	9.4	10.9
3	8.3	12.5	6.8	15.0
4	2.0	10.5	1.9	14.7
5	<1.1	-	<1.1	-
6	<1.1	-	<1.1	-

19. LINEARITY

The linearity of the test was evaluated with 3 samples with known concentration of Vedolizumab. For each sample 8 serial dilutions were done; the linearity was demonstrated in the range 1.6 – 59.3 µg/ml.

19. BIBLIOGRAPHY

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as

induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.

4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

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19. BIBLIOGRAFIE

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 22 Aug 2013;369(8):711-21.
4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

	IT	Data di fabbricazione	ES	Fecha de fabricación
	EN	Date of manufacture	FR	Date de fabrication
	CS	Datum výroby	PT	Data de fabrico
	DE	Herstellungsdatum	RO	Data fabricatiei
	EL	Ημερομηνία Παραγωγής		
	IT	Utilizzare entro	ES	Fecha de caducidad
	EN	Use By	FR	Utiliser jusqu'à
	CS	Použitelné do	PT	Prazo de validade
	DE	Verwendbar bis	RO	A se folosi pana la
	EL	Ημερομηνία λήξης		
	IT	Non riutilizzare	ES	No reutilizar
	EN	Do not reuse	FR	Ne pas réutiliser
	CS	Nepoužívejte opakovaně	PT	Não reutilizar
	DE	Nicht wieder verwenden	RO	A nu se refolosi
	EL	Μην κάνετε επαναληπτική χρήση		
	IT	Attenzione, vedere le istruzioni per l'uso	ES	Atención, ver instrucciones de uso
	EN	Caution, consult accompanying documents	FR	Attention voir notice d'instructions
	CS	Pozor, čtěte přiložené dokumenty	PT	Atenção, consulte a documentação incluída
	DE	Achtung, die Gebrauchsanleitung lesen	RO	Atentie, consultați documentele însoțitoare
	EL	Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα		
	IT	Fabbricante	ES	Fabricante
	EN	Manufacturer	FR	Fabricant
	CS	Výrobce	PT	Fabricante
	DE	Hersteller	RO	Productator
	EL	Κατασκευαστής		
	IT	Contenuto sufficiente per "n" saggi	ES	Contenido suficiente para <n> ensayos
	EN	Contains sufficient for <n> tests	FR	Contenu suffisant pour "n" tests
	CS	Obsah stačí na <n> testů	PT	Conteúdo suficiente para "n" ensaios
	DE	Inhalt reicht für „n“ Tests	RO	Continut sufficient pt <n> teste
	EL	Περιεχόμενο επαρκές για «n» εξετάσεις		
	IT	Limiti di temperatura	ES	Límite de temperatura
	EN	Temperature limitation	FR	Limites de température
	CS	Teplotní omezení	PT	Limites de temperatura
	DE	Temperaturgrenzwerte	RO	Limita de temperatura
	EL	Περιορισμοί θερμοκρασίας		
	IT	Consultare le istruzioni per l'uso	ES	Consulte las instrucciones de uso
	EN	Consult Instructions for Use	FR	Consulter les instructions d'utilisation
	CS	Čtěte návod k použití	PT	Consulte as instruções de utilização
	DE	Die Gebrauchsanleitung lesen	RO	Pentru a utiliza consultați instrucțiunile
	EL	Συμβουλευτείτε τις οδηγίες χρήσης		
	IT	Rischio biologico	ES	Riesgo biológico
	EN	Biological risks	FR	Risques biologiques
	CS	Biologická rizika	PT	Risco biológico
	DE	Biologisches Risiko	RO	Risk biologic
	EL	Βιολογικοί κίνδυνοι		
REF	IT	Numero di catalogo	ES	Número de catálogo
	EN	Catalogue number	FR	Référence du catalogue
	CS	Katalogové číslo	PT	Referência de catálogo
	DE	Katalognummer	RO	Numar de catalog
	EL	Αριθμός καταλόγου		
IVD	IT	Dispositivo medico-diagnostico in vitro	ES	Producto sanitario para diagnóstico in vitro
	EN	In Vitro Diagnostic Medical Device	FR	Dispositif médical de diagnostic in vitro
	CS	Lékařské vybavení pro diagnostiku in vitro	PT	Dispositivo médico para diagnóstico in vitro
	DE	Medizinisches In-vitro-Diagnostikum	RO	Dizpozitiv medical pentru diagnosticare in vitro
	EL	Ιn Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν		
LOT	IT	Codice del lotto	ES	Código de lote
	EN	Batch code	FR	Code du lot
	CS	Kód šarže	PT	Código do lote
	DE	Chargennummer	RO	Lot
	EL	Αριθμός Παρτίδας		
CE	IT	Marcatura CE di conformità	ES	Marcado CE de conformidad
	EN	CE marking of conformity	FR	Marquage de conformité CE
	CS	Oznakowanie zgodności CE	PT	Marcação CE de conformidade
	DE	CE-Konformität Skennezeichnung	RO	Marcajul de conformitate CE
	EL	Σημανση συμμορφωσης CE		