**CHORUS
FOLATE**

REF 86402

DIESSE

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19. BIBLIOGRAFIA

1. Nazki FH, Sameer AS, Ganaie BA. Folate: Metabolism, genes, polymorphisms and the associated diseases. Gene 2014;533(1):11-20.
2. Scaglione F, Panzavolta G. Folate, folic acid and 5-methyltetrahydrofolate are not the same thing. Xenobiotica 2014;44(5):480-488.
3. Reynolds EH. The neurology of folic acid deficiency. Handb Clin Neurol 2014;120:927-43.
4. Wick M, Pinggera W, Lehmann P. Clinical Aspects and Laboratory. Iron metabolism, Anemias. Springer Verlag, Wien, New York, 6th edition 2011:41-42.
5. NCCLS CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

INSTRUCTIONS FOR USE

CHORUS FOLATE

For the quantitative determination of folate in human serum

For *in vitro* diagnostic use only

1. USE

CHORUS FOLATE is an immunological kit for the automated quantitative determination of folate.

The test is conducted on human serum using a single-use device applied to the CHORUS TRIO instrument.

The kit, intended to monitor the level of folate, is an aid in the diagnosis of folate deficiency.

It must only be used by professional laboratory personnel.

2. INTRODUCTION

Folate, also known as vitamin B9 and folacin, is one of the B vitamins. Folate is needed by the body to make DNA and RNA and to metabolise the amino acids required for cell division. Since humans cannot produce folate, it is necessary in the diet, making it an essential nutrient. Folate in the form of folic acid is used to treat anaemia caused by folate deficiency, a type of anaemia in which red blood cells become abnormally large. Folic acid is also used by pregnant women as a supplement to reduce the risk of neural tube defects (NTDs) in the baby. Low levels in early pregnancy are thought to be the cause of over half of babies born with NTDs. Long-term folic acid supplementation can help individuals prone to cardiovascular disease and some forms of cancer. There are concerns that large amounts of supplemental folic acid may mask a vitamin B12 deficiency. Folic acid and N-methyltetrahydrofolate are two common forms, the latter is more stable and has a higher concentration in serum. Folate-binding proteins are responsible for the metabolism of folate and also have the ability to bind several folate derivatives including folic acid and N-methyltetrahydrofolate. Performing a folate assay requires an extraction to free it from the proteins to which it is normally bound.

3. METHOD PRINCIPLE

The CHORUS FOLATE device is ready to use for the determination of folate in human serum in the Chorus TRIO instrument.

The test is based on a competitive ELISA method.

The microplate wells are coated with streptavidin. The calibrator, positive control and samples are added to a device. The sample and the biotinylated folate-binding protein are added to a plate and an immunocomplex between any folate present in the serum and the biotinylated folate-binding protein is formed during a first incubation, which in turn binds to the immobilised streptavidin on the solid phase.

After a short incubation, the conjugate (containing folate marked with peroxidase) is added (this delay increases the sensitivity for low-concentration samples). The conjugate competes with the folate in the test serum for sites binding to the folate-binding

protein. After washing to eliminate the excess unbound reagents, the peroxidase substrate is added. The enzymatic reaction is stopped by adding the stop solution, which causes the solution to turn yellow. The colour that develops is inversely proportional to the amount of folate in the patient's sample.

The single-use devices contain all the reagents to perform the test when applied to the Chorus TRIO instrument.

The results are expressed in ng/ml.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin that have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can completely ensure the absence of infectious agents, any material of human origin must be considered to be potentially infectious material. Follow all the normally adopted laboratory precautions when handling material of human origin.

Disposal of residues: the serum samples, calibrators and strips used must be treated as infected residues and disposed of in accordance with the provisions of current laws.

Health and safety information

1. Do not pipette with your mouth.
2. Use disposable gloves and eye protection when handling samples.
3. Wash your hands thoroughly after inserting the devices into the Chorus TRIO instrument.
4. Refer to the Safety Data Sheet (www.diesse.it) for the safety characteristics of the reagents contained in the kit.
5. Neutralised acids and other liquid waste must be disinfected by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. An exposure of 30 minutes to 1% sodium hypochlorite may be necessary to ensure effective disinfection.
6. Any spillage of potentially infectious materials must be immediately removed with absorbent paper and the contaminated area must be decontaminated, for example with 1% sodium hypochlorite, before continuing the work. If an acid is present in the spill, dry the contaminated area before using sodium hypochlorite. All materials used to decontaminate any spill, including gloves, must be disposed of as potentially infectious waste. Do not autoclave materials containing sodium hypochlorite.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) for at least 30 minutes and use within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
2. When adding the sample to the well, check that it is perfectly distributed on the bottom.
3. Check the actual presence of the reagents in the device and the integrity of the device itself; do not use devices that lack any reagents and/or have foreign bodies in the reaction well when visually inspected.

4. The devices must be used in conjunction with the Chorus TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.

The kit can only be used with CHORUS TRIO with an updated version of the software. Make sure that the software installed on the instrument is a version (Rel.) equal to or later than the one shown in the table published on the Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the Chorus TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode on the grip of the device so that it can be correctly read by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Any defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipemic or icteric samples of serum that is not completely coagulated or samples with microbial contamination may be a source of error.
11. Do not use the device after its expiry date.
12. **Check that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packs of 6 devices each

Description:

Position 9: Space available for barcode label

Position 8: Empty

Position 7: MICROPLATE WELL

Coated with streptavidin

Position 6: Empty

Position 5: BIOTINYLATED FOLATE

Content: biotinylated folate binding protein, stabilised in buffer and preservative.

Position 4: TMB SUBSTRATE

Content: Tetramethylbenzidine and hydrogen peroxide stabilised in buffer.

Position 3: STOP SOLUTION

Content: 0.3 M sulphuric acid solution

Position 2: FOLATE CONJUGATE

Content: Folate marked with peroxidase stabilised in buffer

Position 1: EMPTY WELL

into which the sample will be transferred.

Use: bring a bag to room temperature, open the bag, take out the necessary devices; store the others in the bag containing the silica gel, let the air out and **seal** by pressing on the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.450 ml

Content: Human serum containing folate and preservative.

Liquid, ready to use.

CONTROL + POSITIVE CONTROL 1 x 0.650 ml

Content: Human serum containing folate and preservative.

Liquid, ready to use.

SAMPLE DILUENT/NEGATIVE CONTROL

SAMPLE DILUENT/NEGATIVE CONTROL 1 x 1.100 ml

Content: Human serum without folate and preservative.

Liquid, ready to use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED

- CHORUS VITAMIN EXTRACTION **REF** 86500
- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITISING SOLUTION **REF** 83604
- Chorus TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Vortex
- Micropipettes capable of accurately withdrawing volumes of 50-200 µl
- Disposable gloves
- Sodium hypochlorite solution (5%)
- Containers for the collection of potentially infectious materials

6. STORAGE METHOD AND REAGENT STABILITY

The reagents must be stored at 2-8°C. If the storage temperature is incorrect, the calibration must be repeated and the correctness of the result checked using the control serum (see section 9: Test validation).

The expiry date is printed on each component and on the kit label.

The reagents have limited stability after opening:

DEVICES	8 weeks at 2-8°C
CALIBRATOR	8 weeks at 2-8°C
POSITIVE CONTROL	8 weeks at 2-8°C
SAMPLE DILUENT	8 weeks at 2-8°C

7. SAMPLE COLLECTION AND STORAGE

The sample is serum obtained from blood collected by venous blood sampling and handled with the precautions required by standard laboratory procedures.

According to the NCCLS CLSI H18-A3 guideline, serum samples to be analysed should be coagulated before centrifugation; spontaneous and complete coagulation normally occurs within 30-60 minutes at 22°C-25°C. It is recommended to physically separate the serum, by centrifugation, from contact with the cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

Fresh serum can be kept for 2 days at 2-8°C; for longer storage periods, freeze at temperatures ≤ -20°C.

The sample can undergo a maximum of 3 freeze-thaw cycles. The consequences of using other body fluids are unknown.

Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

The 1/40 extraction (Stabilising/Releasing) solution is stable for 8 hours after preparation.

The extracted sample can be used within 8 hours from extraction if kept at 4°C.

It is strongly recommended to avoid any intensive processing of the samples.

Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

The CHORUS FOLATE test (REF. 86402) can be performed simultaneously ONLY with the CHORUS VIT B12 kit (REF.86400) and with the CHORUS PROMONITOR tests.

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in section 4 Analytical Warnings.
3. Prepare the samples, calibrator, positive control and negative control (when necessary) using the product CHORUS VITAMIN EXTRACTION (REF. 86500) according to the extraction protocol given in the relative Instructions for Use.
Pre-treatment cannot be performed in the CHORUS AUTOSAMPLER instrument (REF. 81210).
4. Dispense into well no. 1 of each device:

SAMPLE	100 µl/device
CALIBRATOR	100 µl/device
POSITIVE CONTROL	100 µl/device

For every batch change use one device for the calibrator.

5. Place the devices on the CHORUS TRIO instrument. Perform the calibration (if required) and the test as indicated in the instrument's User Manual.

9. TEST VALIDATION

Use the positive control to verify the correctness of the result obtained, processing it as indicated in the instrument's User Manual. If the instrument indicates that the positive control has a value outside the acceptability limit, the calibration must be performed again. Previous results will be corrected automatically.

If the positive control result continues to be outside the acceptance range, contact Scientific Support.

Phone: 0039 0577 319554

Fax: 0039 0577 366605

email: scientificsupport@diesse.it

10. INTERPRETATION OF RESULTS

The Chorus TRIO instrument provides the result in ng/ml calculated on the basis of a batch-dependent graph stored in the instrument.

11. LIMITATIONS

All values obtained require a careful interpretation without ever disregarding other indicators relating to the same patient. In fact, the test cannot be used on its own for a clinical diagnosis and the result obtained must always be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. CALIBRATION RANGE

Calibration range: 2.0-20.0 ng/ml.

For samples >20.0 ng/ml, repeat the test by diluting the sample with Folate Sample Diluent/Negative Control provided with the kit. It is advisable to make 1:5 and 1:10 dilutions of the off-scale samples; then proceed with extraction of the samples before analysis.

13. REFERENCE RANGE

120 samples (53 men, 67 women) from an apparently healthy population between the ages of 18 and 91 years were analysed. Pregnant women were excluded.

The reference range of the Chorus FOLATE kit is 2.8-14.9 ng/ml.

These values should only be used as guidelines.

Each laboratory should investigate whether the expected values can be transferred to their patient population and, if necessary, determine its own reference range.

14. ANALYTICAL SPECIFICITY

3 samples with increasing concentration of folate were tested, to which the following interferents were added:

Bilirubin (4.5 mg/dl - 45 mg/dl)
Triglycerides (250 - 1500 mg/dl)
Haemoglobin (2.5 - 10 mg/ml)
Rheumatoid factor (44 - 220 IU/ml)

The presence of the above listed interfering substances in the serum under examination does not alter the test result, with the exception of haemoglobin and triglycerides.

2 samples with folate concentrations in the normal range were tested, to which increasing concentrations of Biotin (0.5 - 1000.0 ng/ml) were added.

The presence of Biotin in the serum under examination does not alter the test result at concentrations up to 10 ng/ml. Higher Biotin concentrations determine incorrect results.

15. COMPARISON STUDIES

In one trial, 179 samples were analysed with a Diesse kit and with another commercially available kit: the correlation data are summarised in the following table:

Correlation	r	95% CI
Pearson	0.93	0.91-0.95
Spearman	0.90	0.86-0.92

Pearson and Spearman's correlation coefficients respectively fall within the very high correlation ranges.

16. PRECISION AND REPEATABILITY

Sample	Within the session		Between sessions	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	11.3	7.3	11.3	8.8
2	5.1	5.7	7.6	13.6
3	4.6	5.9	4.7	14.5
4	3.9	4.6	4.9	12.7
5	2.2	8.3	2.9	13.3

Sample	Between batches		Between instruments	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	12.7	9.0	11.7	8.1
2	6.8	12.9	7.6	15.0

3	4.8	9.6	4.8	12.1
4	4.9	13.6	4.7	13.2
5	2.4	10.0	2.5	13.0

17. LINEARITY

The linearity of the test was evaluated using 3 samples with known concentrations of folate. Serial dilutions were made for each sample, which demonstrated the linearity of the method between 2.4-20.0 ng/ml.

18. ACCURACY

The accuracy of the test was evaluated by measuring the International Standard Vitamin B12, Serum Folate, HOLO TC (NIBSC code 03/178). Three kit batches were tested on one instrument. The mean recovery with respect to the expected standard value (5.3 ng/ml) is 92%.

19. BIBLIOGRAPHY

1. Nazki FH, Sameer AS, Ganaie BA. Folate: Metabolism, genes, polymorphisms and the associated diseases. Gene 2014; 533(1): 11-20.
2. Scaglione F, Panzavolta G. Folate, folic acid and 5-methyltetrahydrofolate are not the same thing. Xenobiotics 2014; 44(5): 480-488.
3. Reynolds EH. The neurology of folic acid deficiency. Handb Clin Neurol 2014; 120: 927-43.
4. Wick M, Pinggera W, Lehmann P. Clinical Aspects and Laboratory. Iron metabolism, Anemias. Springer Verlag, Wien, New York, 6th edition 2011: 41-42.
5. NCCLS CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline - Third Edition. Vol. 24, No. 38

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19. BIBLIOGRAFIE

1. Nazki FH, Sameer AS, Ganaie BA. Folate: Metabolism, genes, polymorphisms and the associated diseases. *Gene* 2014;533(1):11-20.
2. Scaglione F, Panzavolta G. Folate, folic acid and 5-methyltetrahydrofolate are not the same thing. *Xenobiotica* 2014;44(5):480-488.
3. Reynolds EH. The neurology of folic acid deficiency. *Handb Clin Neurol* 2014;120:927-43.
4. Wick M, Pinggera W, Lehmann P. Clinical Aspects and Laboratory. Iron metabolism, Anemias. Springer Verlag, Wien, New York, 6th edition 2011:41-42.
5. NCCLS CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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1. Nazki FH, Sameer AS, Ganaie BA. Folate: Metabolism, genes, polymorphisms and the associated diseases. Gene 2014;533(1):11-20.
2. Scaglione F, Panzavolta G. Folate, folic acid and 5-methyltetrahydrofolate are not the same thing. Xenobiotica 2014;44(5):480-488.

3. Reynolds EH. The neurology of folic acid deficiency. Handb Clin Neurol 2014;120:927-43.
4. Wick M, Pinggera W, Lehmann P. Clinical Aspects and Laboratory. Iron metabolism, Anemias. Springer Verlag, Wien, New York, 6th edition 2011:41-42.
5. NCCLS CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38.

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19. BIBLIOGRAFIA

1. Nazki FH, Sameer AS, Ganaie BA. Folate: Metabolism, genes, polymorphisms and the associated diseases. Gene 2014;533(1):11-20.
2. Scaglione F, Panzavolta G. Folate, folic acid and 5-methyltetrahydrofolate are not the same thing. Xenobiotica 2014;44(5):480-488.
3. Reynolds EH. The neurology of folic acid deficiency. Handb Clin Neurol 2014;120:927-43.
4. Wick M, Pinggera W, Lehmann P. Clinical Aspects and Laboratory. Iron metabolism, Anemias. Springer Verlag, Wien, New York, 6ª edição 2011:41-42.

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19. BIBLIOGRAFIE

1. Nazki FH, Sameer AS, Ganaie BA. Folate: Metabolism, genes, polymorphisms and the associated diseases. *Gene* 2014;533(1):11-20.
2. Scaglione F, Panzavolta G. Folate, folic acid and 5-methyltetrahydrofolate are not the same thing. *Xenobiotica* 2014;44(5):480-488.
3. Reynolds EH. The neurology of folic acid deficiency. *Handb Clin Neurol* 2014;120:927-43.
4. Wick M, Pinggera W, Lehmann P. Clinical Aspects and Laboratory. Iron metabolism, Anemias. Springer Verlag, Wien, New York, 6th edition 2011:41-42.
5. NCCLS CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 nr. 38

	IT	Data di fabbricazione	ES	Fecha de fabricación
	EN	Date of manufacture	FR	Date de fabrication
	CS	Datum výroby	PT	Data de fabrico
	DE	Herstellungsdatum	RO	Data fabricatiei
	EL	Ημερομηνία Παραγωγής		
	IT	Utilizzare entro	ES	Fecha de caducidad
	EN	Use By	FR	Utiliser jusque
	CS	Použitelné do	PT	Prazo de validade
	DE	Verwendbar bis	RO	A se folosi pana la
	EL	Ημερομηνία λήξης		
	IT	Non riutilizzare	ES	No reutilizar
	EN	Do not reuse	FR	Ne pas réutiliser
	CS	Nepoužívejte opakovaně	PT	Não reutilizar
	DE	Nicht wieder verwenden	RO	A nu se refolosi
	EL	Μην κάνετε επαναληπτική χρήση		
	IT	Attenzione, vedere le istruzioni per l'uso	ES	Atención, ver instrucciones de uso
	EN	Caution, consult accompanying documents	FR	Attention voir notice d'instructions
	CS	Pozor, čtěte příložené dokumenty	PT	Atenção, consulte a documentação incluída
	DE	Achtung, die Gebrauchsanleitung lesen	RO	Atentie, consultați documentele însoțitoare
	EL	Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα		
	IT	Fabbricante	ES	Fabricante
	EN	Manufacturer	FR	Fabricant
	CS	Výrobce	PT	Fabricante
	DE	Hersteller	RO	Productator
	EL	Κατασκευαστής		
	IT	Contenuto sufficiente per "n" saggi	ES	Contenido suficiente para <n> ensayos
	EN	Contains sufficient for <n> tests	FR	Contenu suffisant pour "n" tests
	CS	Obsah stačí na <n> testů	PT	Conteúdo suficiente para "n" ensaios
	DE	Inhalt reicht für „n“ Tests	RO	Continut sufficient pt <n> teste
	EL	Περιεχόμενο επαρκές για «n» εξετάσεις		
	IT	Limiti di temperatura	ES	Límite de temperatura
	EN	Temperature limitation	FR	Limites de température
	CS	Teplotní omezení	PT	Limites de temperatura
	DE	Temperaturgrenzwerte	RO	Limita da temperatura
	EL	Περιορισμοί θερμοκρασίας		
	IT	Consultare le istruzioni per l'uso	ES	Consulte las instrucciones de uso
	EN	Consult Instructions for Use	FR	Consulter les instructions d'utilisation
	CS	Čtěte návod k použití	PT	Consulte as instruções de utilização
	DE	Die Gebrauchsanleitung lesen	RO	Pentru a utiliza consultați instrucțiunile
	EL	Συμβουλευτείτε τις οδηγίες χρήσης		
	IT	Rischio biologico	ES	Riesgo biológico
	EN	Biological risks	FR	Risques biologiques
	CS	Biologická rizika	PT	Risco biológico
	DE	Biologisches Risiko	RO	Risk biologic
	EL	Βιολογικοί κίνδυνοι		
REF	IT	Numero di catalogo	ES	Número de catálogo
	EN	Catalogue number	FR	Référence du catalogue
	CS	Katalogové číslo	PT	Referência de catálogo
	DE	Katalognummer	RO	Numar de catalog
	EL	Αριθμός καταλόγου		
IVD	IT	Dispositivo medico-diagnostico in vitro	ES	Producto sanitario para diagnóstico in vitro
	EN	In Vitro Diagnostic Medical Device	FR	Dispositif médical de diagnostic in vitro
	CS	Lékařské vybavení pro diagnostiku in vitro	PT	Dispositivo médico para diagnóstico in vitro
	DE	Medizinisches In-vitro-Diagnostikum	RO	Dizpositiv medical pentru diagnosticare in vitro
	EL	Ιn Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν		
LOT	IT	Codice del lotto	ES	Código de lote
	EN	Batch code	FR	Code du lot
	CS	Kód šarže	PT	Código do lote
	DE	Chargennummer	RO	Lot
	EL	Αριθμός Παρτίδας		
CE	IT	Marcatura CE di conformità	ES	Marcado CE de conformidad
	EN	CE marking of conformity	FR	Marquage de conformité CE
	CS	Oznakowanie zgodności CE	PT	Marcação CE de conformidade
	DE	CE-Konformität Skennezeichnung	RO	Marcajul de conformitate CE
	EL	Σημανση συμμορφωση CE		