WW-DAT BRUCELLA POSITIVE CONTROL

REF 26007

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WW-DAT BRUCELLA POSITIVE CONTROL

Positive control serum for use in the AUTO-DAT Brucella test

For In vitro Diagnostic Use Only

1. INTENDED USE

Positive control serum for use in the AUTO-DAT Brucella test

2. INTRODUCTION

Serological studies of patients with persistent fever can be useful in detecting brucellosis.

The test consists in the observation of a direct agglutination reaction between the bacterial suspension and the serum sample which is suspected of containing specific antibodies.. The method is known as Wright's reaction. In areas where the disease is endemic and there are cases of chronic brucellosis, there may be rare false negative cases due to the presence of incomplete antibodies which inhibit direct agglutination.

3. PRINCIPLE OF THE METHOD

The serum sample is mixed with a bacterial suspension at the optimal concentration to perform the test with AUTO-DAT. In case of need to confirm a negative reaction, the addition of the Positive Control to the reaction mixture allows determining the presence of incomplete antibodies.

4. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

All material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1%

sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring to room temperature (18-30°C) before use.

- Collect the reagent taking care to avoid bacterial or chemical contamination and store at 2-8°C.

- Discard the product if it appears very turbid or maleodorant.

5. KIT COMPOSITION

The Positive Control is sufficient for 40 tests.

CONTROL + Positive Control 1 x 0,4 ml

Contents: positive control serum obtained from rabbit, containing sodium azide 0.09% as preservative. The serum contains antibodies which agglutinate the Brucella bacterial suspension.

Use: The agglutinating serum is used as a positive control serum to determine the presence of incomplete antibodies with the AUTO-DAT instrument.

6. STORAGE AND STABILITY AFTER OPENING

Reagents must be stored at 2-8°C.

Stability after opening: up to the expiration date if correctly stored.

7. PROCEDURE

Refer to WW-DAT (REF 26010) Instructions for Use. In order to exclude the presence of agglutination blocking antibodies, (incomplete antibodies) add 10uL of Positive Control to each well, making sure that the entire amount is well mixed with the content.

Repeat the Brucellosis test with the AUTO-DAT instrument.

8. TEST VALIDATION

A titer $\geq 1/40$ must be obtained.

9. INTERPRETATION OF THE TEST

Results $\geq 1/40$. They indicate the absence of incomplete antibodies in the sample under examination.

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S.p.A.