CHORUS

Chikungunya IgM

REF 81134

CE
0123

DIESSE

DIESSE Diagnostica Senese
S.p.A.

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For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

The product Chorus Chikungunya IgM is an immunoassay kit for automated qualitative detection of IgM class antibodies against Chikungunya.

The test is performed in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

The kit is intended to detect the exposure to Chikungunya infection as an aid to the relative diagnosis.

It must be used by professional laboratory users only.

2. INTRODUCTION

Chikungunya is a mosquito-borne viral disease caused by an RNA virus that belongs to the alphavirus genus of the family *Togaviridae*.

Chikungunya virus is transmitted between humans by bite of *Aedes aegypti* and *Aedes albopictus* (tiger mosquito).

After the bite of an infected mosquito, onset of illness usually occurs 4-8 days later (but can range from 2-12 days). Chikungunya is characterized by an abrupt onset of fever, frequently accompanied by joint pain, which usually lasts for a few days, but may be prolonged for weeks, months or even years. The virus can cause acute, subacute or chronic disease. Other common signs and symptoms include muscle pain, joint swelling, headache, nausea, fatigue and rash.

The symptoms in infected individuals are usually mild and the infection may go unrecognized or may be misdiagnosed. Occasional cases of ophthalmological, neurological, gastrointestinal and heart complications have been reported.

The clinical similarity of Chikungunya to other arbovirus infections (particularly Dengue and Zika) makes diagnosis difficult on the basis of clinical symptoms alone. The diagnosis could be performed directly by the detection of the presence of the virus or indirectly throughout serological tests, such as enzyme-linked immunosorbent assays (ELISA) for the determination of anti-Chikungunya IgM and IgG antibodies.

Since the product Chorus Chikungunya IgM is intended to determine the presence in patient serum of specific anti-Chikungunya IgM antibodies to establish the serological state of the patient, the device is used as an aid in the diagnosis of Chikungunya infection.

3. PRINCIPLE OF THE METHOD

The Chorus Chikungunya IgM device is ready to use for the detection of IgM antibodies against Chikungunya, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgM monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on DIESSE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the device within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.


3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested (according to the indications reported in chapter "Analytical specificity")
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with Chikungunya antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine in citrate buffer

Position 3: SAMPLE DILUENT

Contents: Saline proteic solution containing IgG goat serum and Methylorange 0.0008%.

Position 2: CONJUGATE

Contents: anti-human IgM monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the sample patient has to be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing anti-Chikungunya IgM antibodies, Tween-20 0.2%, Proclin 300 0.1% and Methylorange: 0.0008%.

Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted human serum containing anti-Chikungunya IgM antibodies Tween-20 0.2%, Proclin 300 0.1% and Methylorange: 0.0008%.

Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- Chorus Instruments
- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C (for at least 60 months), and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. At each change of batch, perform the calibration.
4. Dispense the following volumes in well no. 1 of each device:

SAMPLE 50 µl/device
 CALIBRATOR 50 µl/device
 POSITIVE CONTROL 50 µl/device

5. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum as reported in the Instrument Operating Manual to check the validity of the obtained results.

The acceptability range is reported on the enclosed Quality control certificate.

If the result of the positive control is out of the acceptability range, the calibration must be repeated.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.2

NEGATIVE: when the result is < 0.8

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

The product must be used by professional laboratory users only.

The test is not applicable to sample types different from serum.

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 0.1 and 0.7 Index.

13. ANALYTICAL SPECIFICITY

5 samples (2 Negative, 1 Cut-Off and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid Factor (44 – 220 IU/ml)

Bilirubin (4.5 mg/dl -45,0 mg/dL))

Triglycerides (10 mg/dl – 150 mg/dl)

Hemoglobin (5 mg/ml – 30 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

14. CROSS-REACTIONS

79 samples, positive to Dengue (30), Cytomegalovirus (10), Toxoplasma (10), Rubella (9) ed Epstein-Barr virus (20) were tested.

No significant cross-reactions were found for Cytomegalovirus, Toxoplasma and Rubella:

Dengue virus (with 3 positive and 2 doubtful/equivocal samples) and Epstein-Barr virus (with 2 positive and 3 doubtful/equivocal samples) could give cross-reactions.

15. METHOD COMPARISON

In an internal experimentation 138 samples have been tested with Diesse kit and with another kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	37	1	38
	-	1	99	100
	Total	38	100	138

Percent Positive Agreement (~Diagnostic Sensitivity):

97.4% CI_{95%}: 86.4-99.5

Percent Negative Agreement: (~Diagnostic Specificity):

99.0% CI_{95%}: 94.5-99.8

Positive Predicted Value (PPV): 97.4% CI_{95%}: 94.5-100.0

Negative Predicted Value (NPV): 99.0% CI_{95%}: 97.3-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 0.98.

In an external experimentation 34 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	18	0	18
	-	0	16	16
	Total	18	16	34

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% CI_{95%}: 82.4 - 99.8

Percent Negative Agreement: (~Diagnostic Specificity):

100.0% CI_{95%}: 80.6 - 99.7

Positive Predicted Value: 100.0%:

CI_{95%}: 100.0-100.0

Negative Predicted Value: 100.0%:

CI_{95%}: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 1.00.

16. PRECISION AND REPEATABILITY

Sample	Within-run Precision	
	Mean (Index)	CV%
1	3.8	12.9
2	5.5	9.1
3	3.7	8.4
4	0.6	-

Sample	Between-run precision	
	Mean (Index)	CV%
1	3.4	5.9
2	5.0	9.0
3	3.3	14.2
4	0.8	10.0

Sample	Precision between batches		Precision between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	3.1	12.9	3.4	10.0
2	4.6	13.0	5.1	12.9
3	4.1	6.3	3.6	14.7
4	0.5	18.0*	0.7	14.3

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero.

17. REFERENCES

1. <https://www.who.int/news-room/fact-sheets/detail/chikungunya>
2. Johnson B.W. et All. Laboratory Diagnosis of Chikungunya Virus Infections and Commercial Sources for Diagnostic Assays, The Journal of Infectious Diseases, 2016:2146 (Suppl 5).
3. Chelluboina S. et All., Persistence of antibody response in chikungunya, VirusDis. (July–September 2019) 30(3):469–473
4. Pierro A. et All. Persistence of anti-chikungunya virus-specific antibodies in a cohort of patients followed from the acute phase of infection after the 2007 outbreak in Italy, New Microbe and New Infect 2015; 7: 23–25
5. Grivard P et All. Molecular and serological diagnosis of Chikungunya virus infection, Pathologie Biologie, Volume 55, Issue 10, December 2007, Pages 490-494.

18. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

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1. <https://www.who.int/news-room/fact-sheets/detail/chikungunya>
2. Johnson B.W. et All. Laboratory Diagnosis of Chikungunya Virus Infections and Commercial Sources for Diagnostic Assays, The Journal of Infectious Diseases, 2016:2146 (Suppl 5).

3. Chelluboina S. et All., Persistence of antibody response in chikungunya, VirusDis. (July–September 2019) 30(3):469–473
4. Pierro A. et All. Persistence of anti-chikungunya virus-specific antibodies in a cohort of patients followed from the acute phase of infection after the 2007 outbreak in Italy, New Microbe and New Infect 2015; 7: 23–25
5. Grivard P et All. Molecular and serological diagnosis of Chikungunya virus infection, Pathologie Biologie, Volume 55, Issue 10, December 2007, Pages 490-494.

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17. BIBLIOGRAFIA

1. <https://www.who.int/news-room/fact-sheets/detail/chikungunya>
2. Johnson B.W. et All. Laboratory Diagnosis of Chikungunya Virus Infections and Commercial Sources for Diagnostic Assays, The Journal of Infectious Diseases, 2016:2146 (Suppl 5).
3. Chelluboina S. et All., Persistence of antibody response in chikungunya, VirusDis. (July–September 2019) 30(3):469–473
4. Pierro A. et All. Persistence of anti-chikungunya virus-specific antibodies in a cohort of patients followed from the acute phase of infection after the 2007 outbreak in Italy, New Microbe and New Infect 2015; 7: 23–25

5. Grivard P et All. Molecular and serological diagnosis of Chikungunya virus infection, Pathologie Biologie, Volume 55, Issue 10, December 2007, Pages 490-494.

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