RPR-DAT

REF 26030

CE
0123

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16. BIBLIOGRAFIA

- Larsen S.A et All. Laboratory Diagnosis and Interpretation of Tests for Syphilis CLINICAL MICROBIOLOGY REVIEWS (1995), p. 1-21
- Antonio Fuertes Ortiz de Urbina Syphilis (Book)
- WHO guidelines for the treatment of Treponema pallidum (Syphilis) (2016)
- Morshed M.G. et All. Recent Trends in the Serologic Diagnosis of Syphilis. Clinical and Vaccine Immunology (2015) Volume 22 Number 2

5. Meredith E. et All. RPR and the Serologic Diagnosis of Syphilis. JAMA. Author manuscript (2019)
6. Tuddenham S. et All. Syphilis Laboratory Guidelines: Performance Characteristics of Nontreponemal Antibody Tests. Clinical Infectious Diseases (2020):71 (Suppl 1)

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For In Vitro Diagnostic Use Only

1. INTENDED PURPOSE

RPR-DAT (REF 26030) is an agglutination test kit for qualitative detection of reagins in human serum, to be applied on the AUTO-DAT instrument or to be used in manual technique. The test may also be used to detect reagins in samples serially diluted to establish titer information (semi-quantitative test).

Since non-treponemal tests are based on the detection of reagins, an antibody class present in Syphilis (disease caused by *Treponema pallidum* infection), the kit is intended to determine the exposure to *Treponema pallidum* infection and to be used as an aid to the related diagnosis.

2. INTRODUCTION

Syphilis, caused by the spirochete *Treponema pallidum* subsp. *pallidum*, is a chronic infection with many diverse clinical manifestations that occur in distinct stages.

This systemic infectious disease is generally acquired by direct sexual contact and features treponemes-containing lesions. The only known hosts are human beings. There are many tests for the direct and indirect diagnosis of syphilis, still there is no single optimal test.

Direct diagnostic methods include the detection of *T. pallidum* by microscopic examination of fluid or smears from lesions, histological examination of tissues or nucleic acid amplification methods such as polymerase chain reaction (PCR).

Indirect diagnosis is based on serological tests for the detection of antibodies. Serological tests fall into two categories:

1. non-treponemal tests for screening
2. treponemal tests for confirmation

Non treponemal tests are based on the detection of reagins, an antibody class present in syphilis and occasionally in other acute and chronic conditions. Reagins can be detected in the serum about 4-6 weeks after infection, or 1-3 weeks after appearance of the chancre.

Purified beef heart extract (cardiolipin) fortified with lecithin and cholesterol, is used as antigen for reagin tests. Non treponemal tests are routinely used in syphilis serology because of their practical advantages and reproducibility, even though they are not always 100% specific.

Since the product RPR-DAT is a non-treponemal test, based on macroscopic flocculation principle, it is intended to determine the exposure to *Treponema pallidum* infection as an aid to the related diagnosis.

3. PRINCIPLE OF THE METHOD

The assay is based on the flocculation principle.

The serum sample reacts in a microplate flat bottom well with the colloidal suspension of cardiolipin, lecithin and cholesterol mixed with carbon microparticles. When antibodies (reagins) agglutinating antigen are present an aggregate is formed.

The test is applied to the AUTO-DAT instrument, which captures the image of the occurred reaction and analyzes it throughout the use of a specific algorithm, reducing the variability related to the individual interpretation.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the reagents to room temperature (18-30°C) before use.

1. After use, and return reagents to 2-8°C.
2. Do not use the reagents after the expiry date.
3. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
4. Do not modify the Test Procedure. Do not substitute reagents using reagents from other suppliers or other lots, unless it is specifically indicated that the reagent is interchangeable between lots.
5. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested

(according to the indications reported in chapter "Analytical specificity")

5. KIT CONTENTS AND REAGENT PREPARATION

Reagents are sufficient for 500 determinations.

ANTIGEN

ANTIGEN

2 x 4 mL

Contents: the antigen is composed of a stabilised mixture of cardiolipin 0.003%, cholesterol 0.09%, lecithin 0.02 % and carbon microparticles 0.2 g/L

Attention: Store in an upright position when not in use. Bring to room temperature before use. Do not shake too vigorously.

CONTROL +

POSITIVE CONTROL

1 x 0.5 mL

Contents: Human serum, reactive for syphilis, diluted up to a maximum of 15% in a protein solution with sodium azide 0.09%.

Liquid, ready for use.

CONTROL -

NEGATIVE CONTROL/SAMPLE DILUENT

1 x 4 mL

Contents: 100% non reactive human serum containing Proclin 300 at 0.029% and Gentamicin at 0.05% as preservatives.

Liquid, ready for use.

MT PLATE

6 FLAT BOTTOM MICROPLATES (8x12)

MATERIALS REQUIRED BUT NOT PROVIDED

- AUTO-DAT Instrument (REF 26000)
- Rotating mixer
- Micropipette with variable volume

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C

After the first use, the microplate must be stored at room temperature.

The expiry date is printed on each component and on the kit label.

The stability of the reagents does not change after opening the bottle, provided that the user is careful to keep the product safe from possible microbial contamination.

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum from whole blood collected by venipuncture and handled following standard laboratory techniques.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be frozen at -20°C (for at least 26 months).

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

8. ASSAY PROCEDURE

INSTRUMENTAL APPLICATION WITH AUTO-DAT

QUALITATIVE TEST

1. Bring reagents and samples to room temperature. Resuspend the reagent.
2. Draw up 25 µl of serum and drop on the bottom of the microplate well. Take care not to transfer any cellular elements. **Avoid the formation of bubbles.**
3. Repeat these steps for each sample to be tested.
4. Distribute 15 µl of antigen in each well. **Avoid the formation of bubbles.**
5. Verify that the sample-antigen mix is **evenly distributed on the bottom of the well**. If not, tap the plate on the side and / or gently tap it on the work surface.
6. Place the plates in the AUTO-DAT instrument and follow the instructions on the tablet.
7. At the end of the analysis, extract the plates from the AUTO-DAT instrument and, using the adhesive film, cover the wells in which the reaction occurred. Wells not previously used and not covered can be used for further analysis.

NOTE: Leave one strip microwell for Negative Control and one for Positive Control, to be tested exactly as the analyzed samples.

SEMI-QUANTITATIVE TEST

To be carried out on the samples reactive to the screening test.

1. Add 25 µl of Negative Control/Sample Diluent to each well of the plate.
2. Add to the first well 25 µl of the sample that gave positive result.
3. Make doubling dilutions, mixing thoroughly before transferring to each step, and proceed to dispose of the excess 25 µl present in the last well, as shown in the following scheme. **Avoid the formation of bubbles.**

Strip 1x8	A	B	C	D	E	F	G
Sample Diluent (µl)	25	25	25	25	25	25	25
Serum (µl)	25	25	25	25	25	25	25
Mix and transfer (µl)		25	25	25	25	25	25
Dilution/Titer	1/2	1/4	1/8	1/16	1/32	1/64	1/128

Proceed to dispose of the excess 25 µl

4. Follow the procedure described for the screening test from step 4 to step 7.

The titer will be given by the last dilution found reactive with the instrument.

9. TEST VALIDATION ON AUTO-DAT

Use the controls provided at each session.
Proceed as described in the paragraph "INSTRUMENTAL APPLICATION WITH AUTO-DAT".

If the result is different from the expected, contact the Scientific Support:

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

INSTRUMENTAL APPLICATION WITH AUTO-DAT

At the end of the test the instrument will print the results as non-reactive (N), doubtful (D), reactive (P).
It is advisable to visually verify the obtained results for confirmation.

ATTENTION: failure to add the antigen to the sample causes false positive results.

11. LIMITATIONS

The product must be used by professional laboratory users only and it is not intended for screening blood or tissue donors. The test is not applicable to sample types different from serum.

The presence of false negative results due to prozone reaction of high positive samples cannot be excluded. For this reason, all the specimens that give doubtful or reactive (positive) results should be confirmed with the semi-quantitative procedure

Biological false-positive reactions occur occasionally with cardiolipin antigens, e.g. in cases of drug abuse, or in certain diseases such as lupus erythematosus, mononucleosis, malaria, leprosy.

All the obtained results require a careful interpretation that must consider other indicators relative to the patient. The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALYTICAL SPECIFICITY

3 samples (2 Negative and 1 Positive) were spiked with the following potentially interfering factors and then tested:

Bilirubin (45 mg/dl)
Triglycerides (1500 mg/dl)
Hemoglobin (10 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

13. CROSS-REACTIONS

42 samples positive to *Toxoplasma gondii*, Rubella, Herpes simplex virus (Type 1 and Type 2), Cytomegalovirus, Epstein-Barr Virus, Varicella zoster, Borrelia, Rheumatoid Factor and Anti-CCP were tested.

No significant cross-reactions were found.

14. METHOD COMPARISON

In an external experimentation, 218 samples were analyzed with Diesse kit.

The experimental data are summarized below:

		Reference		
		+	-	Total
		41	3	44
Diesse	-	0	174	174
	Totale	41	177	218

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% CI95%: 91.4-99.9

Percent Negative Agreement: (~Diagnostic Specificity):

98.3% CI95%: 95.1-99.4

Positive Predictive Value (PPV): 93.2%

CI95%: 89.9-96.5

Negative Predictive Value (NPV): 100.0%

CI95%: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 0.97.

15. PRECISION AND REPEATABILITY

Precision and repeatability has been evaluated on replicates of negative and positive samples, in three different batches, with different operators performing the test, in different days with different instruments.

All the obtained results were within the acceptance criteria (Agreement % $\geq$ 95% between the expected and the obtained results) demonstrating the precision and repeatability of the method.

16. REFERENCES

1. Larsen S.A et All. Laboratory Diagnosis and Interpretation of Tests for Syphilis CLINICAL MICROBIOLOGY REVIEWS (1995), p. 1–21
2. Antonio Fuertes Ortiz de Urbina Syphilis (Book)
3. WHO guidelines for the treatment of *Treponema pallidum* (Syphilis) (2016)
4. Morshed M.G. et All. Recent Trends in the Serologic Diagnosis of Syphilis. Clinical and Vaccine Immunology (2015) Volume 22 Number 2
5. Meredith E. et All. RPR and the Serologic Diagnosis of Syphilis. JAMA. Author manuscript (2019)
6. Tuddenham S. et All. Syphilis Laboratory Guidelines: Performance Characteristics of Nontreponemal Antibody Tests. Clinical Infectious Diseases (2020):71 (Suppl 1)

17. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

18. SUMMARY OF SAFETY AND PERFORMANCE

This document, available on Eudamed database, is part of Technical Documentation.

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