CHORUS

Total IgE

REF 88000

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16. BIBLIOGRAFIA

1. Geha RS. Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.

2. Gordon RR et al. Immunoglobulin E and eczema- asthma syndrome in early childhood. *Lancet* 1982; 72-74
3. Norman PS. The clinical significance of IgE. *Hospital Practice* 1975 Aug; 10(8): 41-49
4. Johansson SGO, *Lancet* 1967; 2:951
5. Zefferstron and Johansson SGO, *Allergy* 198; 36:537
6. Wathrich B. *Clin. Allergy* 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. *J. Allergy Clin. Immunology* 1981; 68: 106-111
8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. *Clin. Allergy* 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. *J. Clin. Immunoassay* 1983 Jun; 6 (2): 131-139
10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. *J. Clin. Immunoassay* 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. *Clin. Lab. Annual* 1983; 2:351-388
12. Bukley RH. *Immunopharmacology of Allergic Disease*, 1975; 253: 474
13. Capmon A, Dessaint JP, Capron M et al. *Nature* 1975; 253:474
14. Jahonsson SGO, Bennick H, Berg T. *Clin. Immunol.* 1972; p1

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CHORUS Total IgE

For the quantitative determination of total IgE antibodies

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the quantitative determination of IgE class antibodies in human serum, using a disposable device applied on the Chorus TRIO instruments.

2. INTRODUCTION

The E (IgE) class immunoglobulins contain 2 types of polypeptide chains, their molecular weight is about 200.000 Da and they bind to the surfaces of mast cells and basophil granulocytes. The binding of allergens to the IgE bound to the cells causes the release of histamine and other vasoactive substances with a consequent series of events known as allergic reaction.

It was demonstrated that most of the patients with atopic allergic pathologies, such as atopic asthma, atopic dermatitis and hay fever, show an increase of the IgE levels in their blood.

Those individuals with no allergies present widely variable concentrations; they generally increase during childhood, reaching the typical adults levels in the second decade of life.

The measuring of total IgE concentrations is important during the first investigation of childhood allergies and as a diagnostic instrument to foresee future atopic outbreaks. A significant increase can be detected not only in allergic patients, but also in case of IgE myeloma, pulmonary aspergillosis and during the active stage of parasitic infections. Increased IgE levels are detected in case of hypergammaglobulinaemia, ulcerative colitis, hepatitis, cancer and malaria.

3. PRINCIPLE OF THE METHOD

The Chorus Total IgE devices are ready to use for the detection of total IgE antibodies, in the Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay). Anti-human IgE monoclonal antibody is bound to the solid phase.

The specific immunoglobulins are bound to the antibody through incubation with diluted human serum. After washing to eliminate the proteins which have not reacted, incubation is performed with the secondary antibody constituted by biotinylated anti-human IgE immunoglobulins.

The unbound secondary biotinylated antibody is eliminated and the Streptavidin-Peroxidase conjugate is added.

The unbound conjugate is eliminated and the peroxidase substrate is added.

The enzymatic reaction is stopped by the adding of the Stop Solution, which makes the solution turn to yellow. The colour which develops is proportional to the concentration of IgE present in the serum sample.

The disposable devices contain all the reagents to perform the test in the Chorus TRIO instruments.

The results are expressed in International Units (KIU/L) calculated in reference to "WHO 3rd International Standard Immunoglobulin E (IgE), human serum 11/234".

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Calibrator and positive control contain horse serum, obtained from healthy animals. Since no method can offer complete assurance about the absence of infectious agents, it is recommended to handle with care the components containing serum.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.


1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.
The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<https://www.diesse.it/it/downloads/downloads/strumento:39/>)
5. Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested (according to the indications reported in chapter "Analytical specificity")
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: ANTI IgE-BIOTIN

Contents: human anti-IgE conjugated with biotin in TRIS buffer pH 8.1 containing Proclin 300 and Bronidox 0.02%

Position 6: MICROPLATE WELL

Coated with anti-IgE monoclonal antibodies

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: STOP SOLUTION

Contents: Sulphuric acid solution 0.3 M

Position 2: CONJUGATE

Contents: Streptavidin-Peroxidase in a red buffer solution (pH 5.5) containing Proclin 300

Position 1: EMPTY WELL

In which undiluted serum must be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.400 ml

Contents: Horse serum containing human IgE antibodies and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.850 ml

Contents: Horse serum containing human IgE antibodies and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 5 days at 2/8°C, or frozen for longer periods at -20°C (for at least 27 months), and can be thawed a maximum of 3 times. Freeze-thawing cycles have to be avoided.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 75-100 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in International Units (KIU/L), in reference to WHO 3rd International Standard Immunoglobulin E (IgE), human serum 11/234", calculated on the basis of a lot-dependent graph stored in the instrument.

The available literature gives the following information concerning total IgE concentration in human serum:

Mean: 30 KIU/L

Age (years)	KIU/L
0 – 3	<10
3 – 4	<25
4 – 7	<50
7 – 14	<100
≥ 15	<150

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range 10-1000 KIU/L.

For samples > 1000 KIU/L retest the diluted sample in the Negative Control/Sample Diluent (PF83607-not supplied with the kit).

In order to obtain a higher precision a dilution 1:10 is suggested also for samples between 500 KIU/L and 1000 KIU/L.

13. METHOD COMPARISON

In an experimentation 78 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

Correlation	r	95%CI
Pearson	0.96	0.94-0.98
Spearman	0.97	0.96-.098

The correlation between the two methods is very high.

14. ANALITICAL SPECIFICITY

3 samples were spiked with the following potentially interfering factors and then tested:

Bilirubin (4.5 mg/dl - 45 mg/dl)

Triglycerides (250 mg/dl - 1500 mg/dl)

Hemoglobin (2.5 mg/ml - 10 mg/ml)

The presence in the serum sample of the interfering substances described above (except for Hemoglobin at concentrations > 5 mg/ml) does not affect the test result.

15. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (KIU/L)	CV%	Mean (KIU/L)	CV%
1	36.8	7.8	33.0	14.5
2	518.7	5.1	486.8	12.4
3	213.8	6.4	239.4	10.4
4	173.5	5.9	170.0	3.2
5	453.0	4.5	437.6	5.8

Sample	Precision between instruments	
	Mean (KIU/L)	CV%
1	33.6	13.1
2	466.0	10.5
3	231.8	9.5
4	167.9	4.4
5	443.1	4.4

16. REFERENCES

1. Geha RS Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.
2. Gordon RR et al. Immunoglobulin E and eczema- asthma syndrome in early childhood. Lancet 1982; 72-74
3. Norman PS. The clinical significance of IgE. Hospital Practice 1975 Aug; 10(8): 41-49
4. Johansson SGO, Lancet 1967; 2:951
5. Zefferstron and Johansson SGO, Allergy 198; 36:537
6. Wathrich B. Clin. Allergy 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J. Allergy Clin. Immunology 1981; 68: 106-111
8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. Clin. Allergy 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. J. Clin. Immunoassay 1983 Jun; 6 (2): 131-139
10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J. Clin. Immunoassay 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin. Lab. Annual 1983; 2:351-388

12. Bukley RH. Immunopharmacology of Allergic Disease, 1975;
253: 474
13. Capmon A, Dessaint JP, Capron M at al. Nature 1975;
253:474
14. Jahonsson SGO, Bennick H, Berg T. Clin. Immunol. 1972;
p1

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16. BIBLIOGRAFÍA

1. Geha RS Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.

2. Gordon RR et al. Immunoglobulin E and eczema- asthma syndrome in early childhood. *Lancet* 1982; 72-74
3. Norman PS. The clinical significance of IgE. *Hospital Practice* 1975 Aug; 10(8): 41-49
4. Johansson SGO, *Lancet* 1967; 2:951
5. Zefferstron and Johansson SGO, *Allergy* 198; 36:537
6. Wathrich B. *Clin. Allergy* 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. *J. Allergy Clin. Immunology* 1981; 68: 106-111
8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. *Clin. Allergy* 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. *J. Clin. Immunoassay* 1983 Jun; 6 (2): 131-139
10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. *J. Clin. Immunoassay* 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. *Clin. Lab. Annual* 1983; 2:351-388
12. Bukley RH. Immunopharmacology of Allergic Disease, 1975; 253: 474
13. Capmon A, Dessaint JP, Capron M et al. *Nature* 1975; 253:474
14. Jahonsson SGO, Bennick H, Berg T. *Clin. Immunol.* 1972; p1

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione	FR Date de fabrication EL Ημερομηνία Παραγωγής PT Data de fabrico
	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusqu'à EL Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser EL Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions EL Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabbrikante	FR Fabricant EL Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests EL Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température EL Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation EL Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques EL Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue EL Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro EL In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot EL Αριθμός Παρτίδας PT Código do lote
	EN CE marking of conformity ES Marcado CE de conformidad IT Marcatura CE di conformità	FR Marquage de conformité CE EL Σημανση συμμορφωσης CE PT Marcação CE de conformidade