CHORUS CFT

DIESSE

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INSTRUCTIONS FOR USE

CHORUS CFT

For the semiquantitative determination of antibodies fixing the Complement

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

For the semiquantitative determination of antibodies fixing the Complement in human serum, using a disposable device applied on the chorus and Chorus TRIO instruments.

2. INTRODUCTION

In the serology of diseases caused by various antigens, the Complement Fixation Test (CFT) has always been considered the most useful diagnostic tool, chiefly because of the particular feature of the antibodies which are detected. It is in fact well known that high titers of complement-fixing antibodies are indicative of a recent or current infection.

3. PRINCIPLE OF THE METHOD

The CFT devices are ready for use in the assay of antibodies anti-antigen using the Chorus/Chorus TRIO instruments.

The method uses antigens which are incubated for a given period with the test serum, in the presence of guinea pig complement. If specific antibodies are present, they will react with the antigen to form an immunocomplex which binds the complement, blocking its function and determining the inhibition of the haemolysis. In the absence of specific antibodies the immunocomplex will not form and the free complement will cause lysis of the haemolytic mixture composed of sheep red blood cells and monoclonal haemolysin.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

The antigen (with the exception of Coxsackie B Mix, Parainfluenza Mix and Adenovirus) is inactivated and contains 0.01% merthiolate as preservative.

In any case, material of biological origin must always be treated as potentially infectious.

Waste disposal: serum samples and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.


Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. Adding the sample into the well verify that it is perfectly distributed on the bottom.
2. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
3. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.
4. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
5. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
6. Avoid using self-defrosting freezers for the storage of the samples.
7. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
8. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
9. The use of strongly hemolyzed samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
10. Eparine and some anti virals like the AZT are strong inhibitors of the Complement.
11. Do not use the device after the expiry date.

5. KIT COMPOSITION AND REAGENT PREPARATION

Description:

Position 8: Space for application of bar code label

Position 7: SHEEP ERYTHROCYTES

Contents: Sheep erythrocytes in 10% stabilizing solution

Position 6: MICROPLATE WELL for "anticomplementary power of the serum"

Position 5: MICROPLATE WELL for the reaction between the serum and the specific antigen

Position 4: ANTIGEN

Contents: Inactivated antigen (with the exception of Coxsackie B Mix, Parainfluenza Mix and Adenovirus as inactivation destroys the antigenicity). Contains merthiolate at a concentration of 0.01% as preservative.

Position 3: LYOPHYLIZED HEMOLYSIN

Contents: Freeze-dried monoclonal antibodies anti-sheep erythrocytes

Position 2: LYOPHYLIZED COMPLEMENT

Contents: Freeze-dried guinea pig complement.

Position 1: EMPTY WELL

In which undiluted serum must be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

MATERIALS REQUIRED BUT NOT PROVIDED

- CLEANING SOLUTION 2000 REF 83609
- SANITIZING SOLUTION REF 83604 - 83608
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Water bath

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the kit label.

Once opened, if properly stored at 2-8°C, the product is stable up to the expire date.

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

All the sera have to be heat inactivated for 30 minutes at 56°C and used within 8 hours; otherwise, serum samples should be decomplexed again for 10 minutes at 56°C. Samples which present anticomplementary activity (the instrument will show P.A.C. on the monitor and no result will appear) have to be discarded.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Heat inactivate the serum for 30 minutes at 56°C.
2. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
3. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
4. Dispense 120-150 µl of undiluted test serum in well no. 1 of each device.
5. Place the devices in the Chorus/Chorus TRIO instrument. Perform the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

A positive control serum can be used to check the validity of the obtained result.

If the instrument signals that the control serum has a value outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

A semiquantitative result is furnished. Each sample gives a result relative to the specific antigen, in the range NEG - 1/128. A NEG result indicates the absence of complement fixing antibodies; in the other cases, the instrument furnishes a titer from 1/8 to 1/128 in relation to the amount of antibody which is present.

The result P.A.C. means anticomplementary power of the sample.

A diagnosis of a current infection can be obtained by checking the serum-conversion (increase of two titers) in two samples taken at a distance of two weeks.

11. LIMITATIONS

The CFT method detects IgG1, IgG2, IgG3 and IgM antibodies and does not differentiate between different classes of antibodies.

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result (false positive or false negative results cannot be excluded) should be evaluated together with the patient history and other clinical diagnostic evaluation.

Titers 1/8 and 1/16 obtained with *Listeria Monocytogenes* antigen are not significant due to the known cross-reaction with other bacterial antigens.

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