**CHORUS
Promonitor
ETANERCEPT**

REF 86706

CE
0123

DIESSE

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20. BIBLIOGRAFIA

1. Zhao S, et al. "Etanercept for the treatment of rheumatoid arthritis". Immunotherapy. 2018 Mar 1;10(6):433-445.
2. Jamnitski, A., et al... "Patients non-responding to etanercept obtain lower etanercept concentrations compared with responding patients". Ann.Rheum.Dis. 2012; 71:88-91.
3. Emi Aikawa, et al. "Immunogenicity of Anti-TNF-alpha agents in autoimmune diseases". Clin Rev.Allergy Immunol. 2010;38:82-89.
4. de Vries, M. et al. "Immunogenicity does not influence treatment with etanercept in patients with ankylosing spondylitis". Ann.Rheum.Dis. 2009; 68:531-535.
5. Dore, R. K., et.al "The immunogenicity, safety, and efficacy of etanercept liquid administered once weekly in patients with rheumatoid arthritis". Clin.Exp.Rheumatol. 2007; 25:40-46.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38.

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INSTRUCTIONS FOR USE

CHORUS Promonitor ETANERCEPT

For the quantitative detection of Etanercept in human serum

For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

CHORUS Promonitor ETANERCEPT is an immunoassay kit for automated quantitative detection of Etanercept (ETN) in human serum using a disposable device applied on the CHORUS TRIO instrument.

The product, used in conjunction with other clinical and laboratory results, is useful as an aid in the management of patients treated with Etanercept (ETN), such as patients with rheumatic diseases.

It must be used by professional laboratory users only.

2. INTRODUCTION

Etanercept (ETN) is a genetically engineered dimeric fusion protein formed by the extracellular ligand binding domain of TNFR2 (type II Tumor Necrosis Factor- α Receptor) to the Fc domain of human IgG1. It binds trimeric TNF α (Tumor Necrosis factor alpha) and LT α (lymphotoxin alpha) proinflammatory cytokines, preventing binding with their receptors and consequently the inflammatory cascade. The drug is used for the treatment of several inflammatory diseases like rheumatoid arthritis (RA), polyarticular juvenile idiopathic arthritis (JIA), psoriatic arthritis (PA), ankylosing spondylitis (AS), and plaque psoriasis (PSO). Etanercept is the least immunogenic of anti-TNF α drugs. In the case of other anti-TNF treatments like infliximab (IFX) and adalimumab (ADL), it has been extensively demonstrated that low drug levels and anti-drug antibodies (ADA) formation strongly correlate with a loss of clinical response in several pathologies. All detected anti-etanercept antibodies have been shown to be non-neutralizing. Anti-etanercept antibodies are thought to be directed to the hinge region of the molecule, therefore not having an impact on TNF α binding activity and clinical efficacy. Apart from other disease-related aspects etanercept levels in circulation seem to play a relevant role in the clinical response.

3. PRINCIPLE OF THE METHOD

The CHORUS Promonitor ETANERCEPT device is ready to use for the detection of Etanercept (ETN) in human serum, in the CHORUS TRIO instrument.

The assay is based on a sandwich ELISA method.

The microwell strips are pre-coated with an anti-ETN monoclonal antibody. Calibrator, positive control and samples are added to separate device.

ETN is able to bind to preimmobilized anti-ETN antibodies following incubation with the sample. After washing to remove the unbound sample, an HRP-labelled anti-ETN antibody (conjugate) is added. A second incubation step allows the conjugate to bind to the ETN that has become attached to the

microwells. After washing to remove the excess of unbound conjugate, the peroxidase substrate is added. The enzymatic reaction is stopped by adding Stop Solution, which makes the solution turn to yellow. The color which develops is proportional to the amount of ETN in the patient sample.

The disposable devices contain all the reagents to perform the test when applied on the CHORUS TRIO instrument.

The results are expressed in $\mu\text{g/ml}$, calculated in reference to "WHO 1st International Standard for TNF receptor II Fc fusion protein (Etanercept, Human rDNA derived), (NIBSC code: 13/204)".

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on DIESSE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Before use, bring the envelopes containing the devices to be used to room temperature (18-30°C) for at least 30 minutes. Use the device within 60 minutes.

1. **Discard devices which show the substrate (well 2) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.


3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the CHORUS TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
The use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website
(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
5. Check that the CHORUS TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see CHORUS TRIO Operating Manual).
9. When CHORUS Autosampler (REF 81210) is used, refer to the related Operating Manual
10. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
11. Do not use hemolyzed, lipemic, jaundiced samples with a higher concentration of interferents than tested (according to the guidance in the "Analytical Specificity" chapter).
12. Do not use the device after the expiry date.
13. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

DD DEVICES 2 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 7: MICROPLATE WELL

Coated with anti-ETN monoclonal antibody (concentration 0.5 µg/mL)

Position 6: MICROPLATE WELL

Coated with anti-ETN monoclonal antibody (concentration 0.5 µg/mL)

Position 5: STOP SOLUTION

Contents: 0.3 M sulfuric acid solution

Position 4: CONJUGATE

HRP-labelled anti-ETN antibody and buffer containing preservative (concentration 0.04 µg/mL).

Position 3: DILUTION BUFFER

Contents: Buffer solution containing preservative.

Position 2: TMB SUBSTRATE

Contents: Tetramethylbenzidine in citrate buffer.

Position 1: EMPTY WELL in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR **1 x 0.650 ml**
Contents: Human serum containing ETN (concentration 70 ng/mL) and preservative BND (5-Bromo-5-nitro-1,3-dioxane) 1000-1150 ppm.
Liquid, ready for use.

CONTROL + POSITIVE CONTROL **1 x 1.100 ml**
Contents: Human serum containing ETN (concentration 70 ng/mL) and preservative BND (5-Bromo-5-nitro-1,3-dioxane) 1000-1150 ppm.
Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum obtained from blood drawn by venipuncture and handled according to standard laboratory procedures.

According to NCCLS CLSI guideline H18-A3, serum samples for analysis should be coagulated before centrifugation; spontaneous and complete coagulation normally occurs within 30-60 minutes at 22°C-25°C. It is recommended that serum be physically separated by centrifugation from cell contact as soon as possible with a maximum time limit of 2 hours from the time of collection.

The fresh serum may be stored for 5 days at 2/8°C and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be mixed carefully before use.

It is strongly recommended to avoid any harsh treatment of the samples.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

The test CHORUS Promonitor Etanercept (REF 86706) can be performed simultaneously **ONLY** with the other CHORUS Promonitor tests.

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense the following volumes in well no. 1 of each device:

SAMPLE	130 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

At each change of batch, use a device for the calibrator.

4. Place the devices in the CHORUS TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The CHORUS TRIO instrument expresses the result in µg/ml, calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is ≥ 0.14 µg/ml

NEGATIVE: when the result is < 0.14 µg/ml

11. LIMITATIONS

The test measures free ETN. It does not measure ETN bound to anti-ETN antibodies (immunocomplexes).

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The clinician should evaluate the results together with other clinical diagnostic evaluation such as sampling window, methodology, ETN levels and disease activity.

The characteristics of the assay have not been determined for sample matrix other than serum.

This product should only be used by appropriately trained laboratory professional personnel.

12. CALIBRATION RANGE

Calibration range: 0.14-13.0 µg/ml

For samples > 13.0 µg/ml retest the diluted sample in the Negative Control/Sample Diluent (REF 83607 - not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, are negative, with values < 0.14 µg/ml.

14. ANALYTICAL SPECIFICITY

3 samples were tested (1 Negative, 1 Low Positive and 1 High Positive) to which the following interfering agents were added:

Bilirubin (4.5 mg/dl - 45 mg/dl)
Triglycerides (250 mg/dl - 1500 mg/dl)
Hemoglobin (2.5 mg/ml - 10 mg/ml)

The presence of the above interfering substances (with the exception of triglycerides > 750 mg/dl) in the test serum does not alter the test result.

3 Rheumatoid Factor (RF) positive samples 92.6 IU/ml (1) and 100 IU/mL (2) were spiked with Etanercept:

High level - 9.0 µg/ml
Low level - 0.50 µg/ml
Negative - < 0.14 µg/ml

The presence of RF in the test sera does not alter the test result.

15. CROSS-REACTIONS

12 samples, not containing Etanercept, but containing various concentrations, of IFX (5.3-24.2 µg/ml) - 7 samples, ADL (1.4-29.2 µg/ml) - 5 samples, were tested.

No cross-reactions were found.

16. METHOD COMPARISON

In an experimentation 171 samples have been tested with Diesse kit and with Promonitor ETN ELISA kit.

Data are summarized in the following table:

		Reference		
		+	-	Tot.
Diesse	+	42	0	42
	-	0	129	129
	Tot.	42	129	171

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% IC_{95%}: 91.6-99.9

Percent Negative Agreement (~Diagnostic Specificity):

100.0% IC_{95%}: 97.1-100.0

Positive Predictive Value (PPV): 100.0% CI_{95%}: 100.0-100.0

Negative Predictive Value (NPV): 100.0% CI_{95%}: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

Out of 171, for 42 positive samples the correlation between Diesse kit and Promonitor ETN ELISA kit was tested.

Data are summarized in the following table:

Correlation	r	95%CI
Pearson	0.97	0.95-0.99
Spearman	0.98	0.95-0.99

The correlation between the two methods is very high.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (µg/ml)	CV%	Mean (µg/ml)	CV%
1	11.7	9.1	9.2	4.3
2	5.8	3.6	4.8	10.8
3	1.0	6.0	1.0	9.4
4	0.56	3.6	0.53	3.8
5	<0.14	-	<0.14	-
6	<0.14	-	<0.14	-

Sample	Between lots		Between instruments	
	Mean (µg/ml)	CV%	Mean (µg/ml)	CV%
1	9.9	9.9	8.6	9.2
2	4.7	11.7	4.6	8.5
3	0.91	12.1	0.89	11.2
4	0.71	12.7	0.49	8.2
5	<0.14	-	<0.14	-
6	<0.14	-	<0.14	-

18. LINEARITY

The linearity of the test was evaluated with 4 samples with known concentration of Etanercept. For each sample serial dilutions were done; the linearity was demonstrated in the range 0.20-13.0 µg/ml.

19. BIOSIMILARS

A bias evaluation study was performed to test the ability of the CHORUS Promonitor ETANERCEPT device to quantify some commercially available biosimilar drugs, compared to the reference drug (Etanercept Enbrel).

The data obtained are summarized in the following table:

	CHORUS Promonitor ETANERCEPT		
	Level 1	Level 2	Level 3
Concentration Etanercept Enbrel (µg/mL)	10.4	4.4	0.59
Bias Etanercept Erelzi	2.1	1.1	0.2
Bias Etanercept Benepali	-0.1	0.1	0.0

20. BIBLIOGRAPHY

1. Zhao S, et al. "Etanercept for the treatment of rheumatoid arthritis". Immunotherapy. 2018 Mar 1;10(6):433-445.
2. Jamnitski, A., et al... "Patients non-responding to etanercept obtain lower etanercept concentrations compared with responding patients". Ann.Rheum.Dis. 2012; 71:88-91.
3. Emi Aikawa, et al. "Immunogenicity of Anti-TNF-alpha agents in autoimmune diseases". Clin Rev.Allergy Immunol. 2010;38:82-89.

4. de Vries, M. et al. "Immunogenicity does not influence treatment with etanercept in patients with ankylosing spondylitis". Ann.Rheum.Dis. 2009; 68:531-535.
5. Dore, R. K., et.al "The immunogenicity, safety, and efficacy of etanercept liquid administered once weekly in patients with rheumatoid arthritis". Clin.Exp.Rheumatol. 2007; 25:40-46.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38.

21. If a serious accident has occurred in connection with this device in the market territory of the European Union, please report it without delay to the manufacturer and the competent authority of your member state.

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20. BIBLIOGRAFIE

1. Zhao S, et al. "Etanercept for the treatment of rheumatoid arthritis". Immunotherapy. 2018 Mar 1;10(6):433-445.
2. Jamnitski, A., et al. "Patients non-responding to etanercept obtain lower etanercept concentrations compared with responding patients". Ann.Rheum.Dis. 2012; 71:88-91.
3. Emi Aikawa, et al. "Immunogenicity of Anti-TNF-alpha agents in autoimmune diseases". Clin Rev.Allergy Immunol. 2010;38:82-89.
4. de Vries, M. et al. "Immunogenicity does not influence treatment with etanercept in patients with ankylosing spondylitis". Ann.Rheum.Dis. 2009; 68:531-535.
5. Dore, R. K., et.al "The immunogenicity, safety, and efficacy of etanercept liquid administered once weekly in patients with rheumatoid arthritis". Clin.Exp.Rheumatol. 2007; 25:40-46.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38.

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4. De Vries, M. et al. "Immunogenicity does not influence treatment with etanercept in patients with ankylosing spondylitis". *Ann.Rheum.Dis.* 2009; 68:531-535.
5. Dore, R. K., et.al "The immunogenicity, safety, and efficacy of etanercept liquid administered once weekly in patients with rheumatoid arthritis". *Clin.Exp.Rheumatol.* 2007; 25:40-46.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38.

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