**CHORUS
Promonitor
USTEKINUMAB**

REF 86712

DIESSE

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For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

CHORUS Promonitor USTEKINUMAB is an immunoassay kit for automated quantitative detection of Ustekinumab (UTK) in human serum using a disposable device applied on the CHORUS TRIO instrument.

The product, used in conjunction with other clinical and laboratory findings, is useful as an aid in the management of patients treated with Ustekinumab (UTK), like patients with inflammatory bowel disease (IBD), psoriasis and rheumatic diseases.

It must be used by professional laboratory users only.

2. INTRODUCTION

Ustekinumab (UTK) is a fully human IgG1 kappa monoclonal antibody that binds with high affinity to the p40 subunit of human interleukin (IL) 12 and IL23 and has been approved for the treatment of patients with moderate to severe Crohn's disease (CD), moderate to severe ulcerative colitis (UC), psoriatic arthritis and plaque psoriasis. The drug prevents the bioactivity of IL12 and IL23 by binding and neutralizing the shared p40 subunit, preventing interaction with the cell surface receptor protein IL12Rbeta1. In the context of inflammatory bowel disease (IBD), MC and UC, the treatment regimen started with a single loading dose and a maintenance regimen with standard subcutaneous administration every 8 weeks. Clinical studies show that 2.9% to 12.4% of patients develop anti-Ustekinumab antibodies (ATU) over time. Therefore, it is important to monitor trough UTK concentrations in serum to correlate drug levels with loss of response to therapy and the occurrence of adverse events such as serum sickness and hypersensitivity reactions. Currently, Ustekinumab quantification is performed in conjunction with the immunogenicity assessment for ATU.

3. PRINCIPLE OF THE METHOD

The CHORUS Promonitor USTEKINUMAB device is ready to use for the detection of Ustekinumab (UTK) in human serum, in the CHORUS TRIO instrument.

The assay is based on a bridging ELISA method.

The microwell strips are pre-coated with an anti-UTK monoclonal antibody. Calibrator, positive control and samples are added to separate device.

UTK is able to bind to preimmobilized anti-UTK antibodies following incubation with the sample. After washing to remove the unbound sample, an HRP-labelled anti-UTK antibody (conjugate) is added. A second incubation step allows the conjugate to bind to the UTK that has become attached to the microwells. After washing to remove the excess of unbound conjugate, the peroxidase substrate is added. The enzymatic reaction is stopped

by adding Stop Solution, which makes the solution turn to yellow. The color which develops is proportional to the amount of UTK in the patient sample.

The disposable devices contain all the reagents to perform the test when applied on the CHORUS TRIO instrument.

The results are expressed in µg/ml.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on DIESSE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the devices within 60 min.

1. **Discard devices which show the substrate (well 2) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

4. The devices are for use with the CHORUS TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the CHORUS TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see CHORUS TRIO Operating Manual).
9. When CHORUS Autosampler (REF 81210) is used, refer to the related Operating Manual
10. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
11. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested (according to the indications reported in chapter "Analytical specificity")
12. Do not use the device after the expiry date.
13. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 7: MICROPLATE WELL

Coated with 0.5 µg/mL anti-UTK monoclonal antibody

Position 6: MICROPLATE WELL

Coated with 0.5 µg/mL anti-UTK monoclonal antibody

Position 5: STOP SOLUTION

Contents: 0.3 M sulfuric acid solution

Position 4: CONJUGATE

HRP-labelled anti-UTK antibody (0.0625 µg/mL max concentration) and buffer containing preservative

Position 3: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 2: TMB SUBSTRATE

Contents: Tetramethylbenzidine stabilized in citrate buffer

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in

the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.650 ml

Contents: Human serum containing 30 ng/mL UTK and preservative 1000-1150 ppm BND (5-Bromo-5-nitro-1,3-dioxane) Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 1.100 ml

Contents: Human serum containing 30 ng/mL UTK and 1000-1150 ppm preservative BND (5-Bromo-5-nitro-1,3-dioxane) Liquid, ready for use.

Confidence in measurements of Calibrator and Positive control is established with traceability to measurement standards as follows.

Calibrator and Positive control are produced diluting a known concentration of human USTEKINUMAB in its own stabilizing medium. The relative exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared and characterized, checking the consensus with a reference sera panel with different USTEKINUMAB levels.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected from the vein and handled with all precautions dictated by good laboratory practice.

According to the guideline CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 5 days at 2/8°C, or frozen for longer periods at temperature $\leq -20^{\circ}\text{C}$ (for at least 5 months) and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be mixed carefully before use.

It is strongly recommended to avoid any harsh treatment of the samples.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

The test CHORUS Promonitor Ustekinumab (REF 86712) can be performed simultaneously ONLY with the other CHORUS Promonitor tests.

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense the following volumes in well no. 1 of each device:

SAMPLE	130 $\mu\text{l}/\text{device}$
CALIBRATOR	130 $\mu\text{l}/\text{device}$
POSITIVE CONTROL	130 $\mu\text{l}/\text{device}$

At each change of batch, use a device for the calibrator.

4. Place the devices in the CHORUS TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The CHORUS TRIO instrument expresses the result in $\mu\text{g}/\text{ml}$, calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is $\geq 0.19 \mu\text{g}/\text{ml}$
NEGATIVE: when the result is $\leq 0.10 \mu\text{g}/\text{ml}$

DOUBTFUL/EQUIVOCAL: when the result is between 0.11 and 0.18 $\mu\text{g}/\text{ml}$

Doubtful/equivocal results are obtained in case of values which ranges between Limit of detection ($< 0.10 \mu\text{g}/\text{ml}$) and Low Limit of Quantitation ($0.19 \mu\text{g}/\text{ml}$).

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

The test measures free UTK. It does not measure UTK bound to anti-UTK antibodies (immunocomplexes).

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The clinician should evaluate the results together with other clinical diagnostic evaluation such as sampling window, methodology, UTK levels and disease activity.

The characteristics of the assay have not been determined for sample matrix other than serum.

This product should only be used by appropriately trained laboratory professional personnel.

12. CALIBRATION RANGE

Calibration range: 0.10-33.0 $\mu\text{g}/\text{ml}$

For samples $> 33 \mu\text{g}/\text{ml}$ retest the diluted sample in the Negative Control/Sample Diluent (REF 83607 - not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, are negative, with values $< 0.10 \mu\text{g}/\text{ml}$.

14. ANALYTICAL SPECIFICITY

3 samples were tested (1 Negative, 1 Low Positive and 1 High Positive) to which the following interfering agents were added:

Bilirubin (4.5 mg/dl - 45 mg/dl)
Triglycerides (250 mg/dl - 1500 mg/dl)
Hemoglobin (2.5 mg/ml - 10 mg/ml)

The presence of the above interfering substances in the test serum does not alter the test result.

3 Rheumatoid Factor (RF) positive samples 150 IU/ml (1) and 200 IU/mL (2) were spiked with Ustekinumab:

High level - 20.0 $\mu\text{g}/\text{ml}$
Low level - 2.0 $\mu\text{g}/\text{ml}$
Negative - $< 0.10 \mu\text{g}/\text{ml}$

The presence of RF in the test sera does not alter the test result.

15. CROSS-REACTIONS

15 samples, not containing UTK, but containing various concentrations of IFX (10 $\mu\text{g}/\text{ml}$) - 5 samples, ADL (27.0 $\mu\text{g}/\text{ml}$) - 5 samples, VDZ (65.5 $\mu\text{g}/\text{ml}$) - 5 samples were tested. No cross-reactions were found.

16. METHOD COMPARISON

In an experimentation 248 samples have been tested with Diesse kit and with Promonitor UTK ELISA kit.

Data are summarized in the following table:

		Reference		
		+	-	Tot.
	+	96	0	96
	-	1	151	152
Diesse	Tot.	97	151	248

Percent Positive Agreement (~Diagnostic Sensitivity):

99.0% IC_{95%}: 94.3-99.8

Percent Negative Agreement: (~Diagnostic Specificity):

100.0% IC_{95%}: 97.5-100.0

Positive Predictive Value (PPV): 100.0% CI_{95%}: 100.0-100.0

Negative Predictive Value (NPV): 99.3% IC_{95%}: 98.3-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

Out of 248, for 97 positive samples the correlation between Diesse kit and a competitor kit was tested.

Data are summarized in the following table:

Correlation	r	95%CI
Pearson	0.95	0.92-0.97
Spearman	0.97	0.96-0.98

The correlation between the two methods is very high.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (µg/ml)	CV%	Mean (µg/ml)	CV%
1	24.9	6.2	22.3	11.3
2	7.5	10.5	9.5	13.1
3	1.4	4.3	1.7	8.2
4	0.60	5.0	0.58	6.9
5	< 0.10	-	< 0.10	-
6	< 0.10	-	< 0.10	-

Sample	Between lots		Between instruments	
	Mean (µg/ml)	CV%	Mean (µg/ml)	CV%
1	22.4	13.7	23.7	13.6
2	7.3	13.7	8.6	15.0
3	1.6	5.0	1.9	11.1
4	0.58	13.8	0.60	11.7
5	< 0.10	-	< 0.10	-
6	< 0.10	-	< 0.10	-

18. LINEARITY

The linearity of the test was evaluated with 3 samples with known concentration of Ustekinumab. For each sample 10 serial dilutions were done; the linearity was demonstrated in the range 0.11-33.0 µg/ml.

19. BIBLIOGRAPHY

- Lucas Guillo, Ferdinando D'Amico, Silvio Danese, Laurent Peyrin-Biroulet "Ustekinumab for Extra-intestinal Manifestations of Inflammatory Bowel Disease: A Systematic Literature Review" J Crohns Colitis 2021 Jul 5;15(7):1236-1243.

- Sands BE, Sandborn WJ, Panaccione R, O'Brien CD, Zhang H, Johanns J, Adedokun OJ, Li K, Peyrin-Biroulet L, Van Assche G, Danese S, Targan S, Abreu MT, Hisamatsu T, Szapary P, Marano C; "Ustekinumab as Induction and Maintenance Therapy for Ulcerative Colitis". N Engl J Med. 2019 Sep 26;381(13):1201-1214.
- Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
- Klenske E, Osaba L, Nagore D, Rath T, Neurath MF, Atreya R. Drug Levels in the Maternal Serum, Cord Blood and Breast Milk of a Ustekinumab-Treated Patient with Crohn's Disease. J Crohns Colitis. 2019 Feb 1;13(2):267-269.
- Martín-González S, Urigoitia-Ugalde P, Careaga J, Nagore D, Navarro R, Izu-Belloso R, Blanch L, González Del Tánago J, Martínez A, Ruiz Del Agua A. Optimal concentration range of ustekinumab in patients with plaque-type psoriasis. J Am Acad Dermatol. 2019 Jun;80(6):1782-1784.
- CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

20. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

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