**CHORUS
HERPES SIMPLEX
1+2 IgG**

REF 81020

DIESSE

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3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgG antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
8. J.E. Kuhn et al.: Analysis of the IgG and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).
9. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

CHORUS HERPES SIMPLEX 1+2 IgG is an Immunoassay kit for the automated qualitative determination of IgG class antibodies against Herpes Simplex Virus 1+2. The test, performed in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

The kit is intended to detect the exposure to Herpes simplex virus (Type 1 and Type 2) infection as an aid to the related diagnosis.

It must be used by professional laboratory users only.

2. INTRODUCTION

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular cheratitis, and lesions of the central nervous system. HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion.

3. PRINCIPLE OF THE METHOD

The CHORUS HERPES SIMPLEX 1+2 IgG device is ready to use for the detection of IgG antibodies against Herpes Simplex Virus 1+2, in the Chorus/ Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off)

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on DIESSE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) for at least 30 minutes before the use. Use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the

table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)


5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, jaundiced samples with a higher concentration of interferents than tested (according to the guidance in the chapter "Analytical Specificity").
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with HSV type 1 and type 2 antigens

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: anti-human IgG 18.13 monoclonal antibodies labelled with horseradish peroxidase (max concentration 1 µg/ml), in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.425 ml

Contents: Diluted human serum containing IgG antibodies anti- ... Herpes Simplex Virus 1+2, and preservative Proclin 0.2% and Gentamycin 0.05%, Ethylene Glycol:23%) Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Diluted human serum containing IgG antibodies anti-Herpes Simplex Virus 1+2, and preservative (Proclin 0.2% and Gentamycin 0.05%, Ethylene Glycol:23%) Liquid, ready for use.

Confidence in measurements of Calibrator and Positive control is established with traceability to measurement standards as follows.

Calibrator and Positive control are produced diluting a known concentration of human IgG against Herpes Simplex Virus Type 1 and 2 in its own stabilizing medium. The relative exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected from the vein and handled with all precautions dictated by good laboratory practice. According to the guideline CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at temperatures ≤-20°C for at least 69 months, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Customer care.

Tel: 0039 0577 319554
email: scientificsupport@diesse.it
customercare@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.2

NEGATIVE: when the result is < 0.8

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

Paired specimens are required for determining seroconversion. Collect the first specimen as soon as possible after onset of symptoms, and the second specimen 2-3 weeks later. If the first sample is taken too late during the course of infection, seroconversion may not be detected. If the test results obtained for paired samples are inconclusive, reassay the specimens. Results must always be interpreted together with an evaluation of the clinical picture and results of other diagnostic investigations.

12. CROSS-REACTIONS

17 samples, positive to anti Rubella, EBV, Cytomegalovirus, Mumps IgG were tested.

No significant cross-reactions were found.

13. METHOD COMPARISON

In an experimentation 46 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	40	0	40
	-	0	6	6
	Total	40	6	46

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Positive Predictive Value (PPV): 100.% CI_{95%}: 100%

Negative Predictive Value (NPV): 100.% CI_{95%}: 100%

The agreement between the two methods is excellent with a Cohen's Kappa of 1

14. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean Index,	CV%	Mean Index,	CV%
1	2.5	6	0.3	20
2	2.5	6	1.3	12
3	2.3	10	3.5	7

15. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgG antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
8. J.E. Kuhn et al.: Analysis of the IgG and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).
9. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

16. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

17. SUMMARY OF SAFETY AND PERFORMANCE

This document, available on Eudamed database (when this is fully implemented and functioning), is part of Technical Documentation and can be requested from the manufacturer.

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15. BIBLIOGRAFIE

- G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
- B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
- C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
- M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
- D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgG antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
- R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
- J.E. Kuhn et al.: Analysis of the IgG and IgG antibody response against herpes simplex virus type 1 (HSV-1)

structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

9. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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15. BIBLIOGRAFIE

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen

- from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
 6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgG antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
 7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
 8. J.E. Kuhn et al.: Analysis of the IgG and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).
 9. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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15. BIBLIOGRAFIJA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgG antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).

7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
8. J.E. Kuhn et al.: Analysis of the IgG and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).
9. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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15. BIBLIOGRAFIJA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgG antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).

8. J.E. Kuhn et al.: Analysis of the IgG and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).
9. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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