CHORUS

**MYCOPLASMA
PNEUMONIAE IgG**

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16. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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For In Vitro Diagnostic Use Only

1. INTENDED PURPOSE

CHORUS Mycoplasma pneumoniae IgG (REF 81034) is an immunoassay kit for the automated semi-quantitative determination of IgG antibodies against Mycoplasma pneumoniae.

Mycoplasma Pneumoniae is the most common etiological agent causing pneumonia acquired in community environments. The IgG appear after the IgM and IgA and they generally reach the peak level during the 5° week, therefore the kit is used as an aid in the diagnosis of pneumonia infection.

The test, performed in human serum using a disposable device applied to the CHORUS and CHORUS TRIO instruments, must be used by professional laboratory personnel only.

2. INTRODUCTION

Mycoplasma pneumoniae is the most common etiological agent causing pneumonia acquired in community environments, particularly in the age range of 5 to 30 years; it can cause epidemics which develop slowly as the incubation period varies from 10 to 14 days, and the infection involves subjects living in close contact or in segregated groups (schools, military barracks, family groups). Pneumonia caused by Mycoplasma is also known as atypical primary pneumonia or Eaton's pneumonia.

Mycoplasma pneumoniae attacks and destroys ciliated epithelial cells of the respiratory tract mucosa. Microscopic examination evidence interstitial pneumonia, bronchitis and bronchiolitis.

Due to the recurrence of symptoms caused by different etiological agents, additional diagnostic measures such as serological tests are necessary for the diagnosis of an acute pneumonia infection.

The immune response, following a Mycoplasma pneumoniae infection, depends on the type of infection: the IgM-class antibodies are more frequently found in case of primary infections as demonstrated by their presence in younger patients. The IgM are low or not found in older patients with higher probability of reinfection. On the other hand the IgA represent the highest sensitivity marker in case of an infection in progress, a reinfection or a recent infection.

Specific IgA-class antibodies appear early during the infection and reach peak levels during the first 4 weeks, then declining quite rapidly, earlier than the IgG; diagnosis of an acute infection can often be made on the basis of a single sample test. The IgG appear after the IgM and IgA and they generally reach the peak level during the 5° week. A high titer together with a significant increase between two

samples taken approximately 2 weeks apart, confirms that an infection is in progress.

3. PRINCIPLE OF THE METHOD

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The sonicated, solubilized Mycoplasma antigen is bound to the solid phase. Specific immunoglobulins bind to the antigen through incubation with diluted human serum. After washings to eliminate the proteins which have not reacted, incubation with the conjugate (anti-human IgG monoclonal antibodies conjugated with horse radish peroxidase) is performed. The conjugate which has not been bound is eliminated and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the test serum.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in AU/ml calculated in reference to an internal standard curve.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instruments.
4. Refer to the Safety Data Sheet (available on DIESSE website: www.diesse.it) for the safety characteristics of the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry.

Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) for at least 30 minutes before the use; use within 60 min.


1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
The Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a Release (date) above the one shown in the table published on the Diesse website
(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument and do not use devices with defective labels.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument.
9. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
10. Do not use hemolyzed, lipemic or , icteric samples with a concentration of interferents higher than that tested (according to the indications given in the "Analytical specificity" chapter).
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer REF 83606.**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Use: **equilibrate a package at room temperature**, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: empty

Position 6: Coated with Mycoplasma pneumoniae antigen: maximum concentration 2.51 µg/ml

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE.

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgG (maximum concentration 1.5 µg/m) labeled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum.

CALIBRATOR CALIBRATOR 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use**.

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use**.

Confidence in measurements of Calibrator and Positive control is established with traceability to measurement standards as follows.

Calibrator and Positive control are produced diluting a known concentration of human antigens in its own stabilizing medium. The relative exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared and characterized, checking the consensus with a reference sera panel with different antigens levels.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS/CHORUS TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be

repeated and the run validated using the control serum (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation

DEVICES	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected from the vein and handled with all precautions dictated by good laboratory practice. According to the guideline CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at temperature $\leq -20^{\circ}\text{C}$, a least 47 months, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 μl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument gives a semiquantitative result in AU/mL, calculated on the basis of a lot-dependent graph stored in the instrument.

The degree of immunity of the test serum can be interpreted as follows:

IMMUNE: when the IgG concentration in the sample is >18 AU/mL

NOT IMMUNE: if the concentration is < 12 AU/mL

DOUBTFUL/EQUIVOCAL: when the concentration is in the range 12-18 AU/mL

Titer	Interpretation of the results		CHORUS Result
> 40 AU/mL	High titer	Titer compatible with a recent or current M.pneumoniae infection	P
≤ 40 AU/mL >18 AU/mL	Moderate titer	The interpretation is only possible with a new collection of sample after 8-15 days.	P
≤ 18 AU/mL ≥ 12 AU/mL	Low titer		D/E
< 12 AU/mL	Not significant titer		N

When the sign $>$ appears on the instrument the AU/mL are over 100 AU/mL; in this case dilute the sample in Sample Diluent/Negative Control REF 83607 (e.g. 1:3) and repeat the test. The result must be multiplied by the dilution used. If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS OF THE PROCEDURE

A negative result does not exclude a Mycoplasma infection: in case of clinical suspect collect a new sample after 15 days.

Serum samples collected during the acute phase of the infection can give a positive result to IgM antibodies only. It is suggested to determine the Mycoplasma IgM level using the Mycoplasma IgM kit (REF 81035).

Alternatively, a second serum sample obtained 8-14 days later, should be tested in parallel to determine an increase in the IgG antibody level.

The test result should be used in conjunction with information available from the evaluation of the case history or other diagnostic procedures.

12. CALIBRATION RANGE

Calibration range 0-100.0 AU/ml.

For samples > 100.0 AU/ml repeat the test pre-diluting the sample with Negative Control/Sample Diluent (PF83607 – not supplied with the kit).

13. ANALYTICAL SPECIFICITY

The presence in the serum sample of potentially interfering substances such as Reactive Protein C, bilirubin, rheumatoid factor, haemolysis and triglycerides does not interfere with the test.

14. METHOD COMPARISON

61 samples were tested, 42 of which were positive and 19 negatives with the Chorus kit in comparison with another commercial method. The discordant samples were tested using Western Blotting considered as reference. Doubtful samples were considered positive. The Diesse kit gave 2 false positive results.

The results are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	40	2	42
	-	0	19	19
	Total	40	21	61

Sensitivity 100.0 % CI_{95%}: 91.2 to 100.0

Specificity 90.5 % CI_{95%}: 71.1 to 97.3

Positive Predicted Value (PPV) = 95.2 CI_{95%}: 89.9 – 100.0

Negative Predicted Value (NPV) = 100.0 CI_{95%}: 100.0 – 100.0

15. PRECISION AND REPEATABILITY

Table 1. Within-run precision

	No. of replicates	Mean AU/ml	S. D.	CV%
Lot. no. 475	9	61.0	2.80	4.6
Lot. no. 459	9	45.3	4.90	10.8
Lot. no. 504	9	37.4	4.84	12.9

Table 2. Between-run and between-lots precision

Sample	Mean AU/ml	Mean	S. D.	CV%
--------	------------	------	-------	-----

	Lot 475	Lot 459	Lot 504			
MYG 1	<10	<10	<10	<10	0.00	0.0
MYG 2	15.3	14.7	14	14.7	0.65	4.4
MYG 3	41.7	39.6	40.7	40.7	1.05	2.6
MYG 4	65.7	78.3	73.7	72.6	6.38	8.8
MYG 5	>100	>100	>100	>100	0.00	0.0

16. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S. Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J. Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E. Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

17. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

18. SUMMARY OF SAFETY AND PERFORMANCE

This document, which will be made available on the EUDAMED database (when this is fully implemented and functioning), is part of the Technical Documentation and can be requested from the manufacturer.

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16. BIBLIOGRAPHIE

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of *Mycoplasma pneumoniae*. Microbiology (2000), 146, 741-747
3. S.Rotten : Interaction of *Mycoplasma* With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of *Mycoplasma pneumoniae* Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of *Mycoplasma pneumoniae* by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAFIE

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).

2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAFIJA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells. Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. BIBLIOGRAFIJA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000), 146, 741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171

5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. G. Layh-Schmitt: Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology (2000),146,741-747
3. S.Rotten : Interaction of Mycoplasma With Host Cells..Physiol.Rev.83:417-432,2003
4. J.Petitjean: Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J.Clin.Microbiol. 2002 40:165-171
5. E.Jacobs : Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion

Chromatography. Biol.Chem Vol.369, pp.1295-1299, December 1988.

6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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16. VIITED

1. G. B. Wisdom, Enzyme-Immunoassay. Clin Chem, 1976, 22, 1243.
2. G. Layh-Schmitt, Proteins complexed to the P1 adhesin of Mycoplasma pneumoniae. Microbiology, 2000, 146, lk 741–747.
3. S. Rotten, Interaction of Mycoplasma With Host Cells. Physiol Rev, 2003, 83, lk 417–432.
4. J. Petitjean, Evaluation of four Commercial Immunoglobulin G and IgM-Specific Enzyme Immunoassays for Diagnosis of Mycoplasma pneumoniae Infections. J Clin Microbiol, 2002, 40, lk 165–171.
5. E. Jacobs, Isolation of the Adherence Protein of Mycoplasma pneumoniae by Fractionated Solubilization and size exclusion Chromatography. Biol Chem, detsember 1988, 369, lk 1295–1299.
6. CLSI H18-A3. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline. Kolmas väljaanne. Kõide 24, number 38.

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