**CHORUS**  
**HELICOBACTER PYLORI**  
**IgA**

**REF** 81062

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Percent Positive Agreement (~Sensibilità Diagnostica): 90.5% (CI95%: 71.0-97.2)

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**For *In Vitro* Diagnostic Use Only**

#### **1. INTENDED PURPOSE**

CHORUS HELICOBACTER PYLORI IgA (REF 81062) is an immunoassay kit for the automated qualitative determination of IgA antibodies against *Helicobacter pylori*.

*Helicobacter pylori* infection causes a systemic (IgM, IgG) and local (IgA) immune response; therefore, the kit is used as an aid in the diagnosis of *Helicobacter pylori* infection.

The test, performed in human serum using a disposable device applied to the CHORUS and CHORUS TRIO instruments, must be used by professional laboratory personnel only.

#### **2. INTRODUCTION**

In 1983, Warren and Marshall identified *Helicobacter pylori*, a new gram-negative bacterial pathogen, in patients suffering from gastritis, and this finding led to studies on the relationship between bacterial infection and chronic gastric disease. The pathogen has been shown to be associated with peptic ulcer, chronic gastritis type B and duodenitis. It has been demonstrated that in patients with gastritis, eradication of the bacteria led to healing of the anatomical lesion.

Diagnostic procedures for the detection of the organism generally involve invasive (gastroscopic) techniques for sample collection.

However, a specific immune response is seen in infected patients. The serological test thus represents a useful alternative to the invasive bioptic technique. IgG levels rise with infection and remain constantly high until the infection is eliminated. The efficacy of antimicrobial therapy can therefore be monitored via changes in specific antibodies. The determination of the IgA antibody levels is complementary to that of the IgG. As the IgA levels decrease more rapidly than the IgG in some patients undergoing treatment, this parameter can be useful in the follow-up of patients.

For unknown reasons, about 2% of serum samples are positive only for IgA.

#### **3. PRINCIPLE OF THE METHOD**

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgA monoclonal antibodies conjugated to peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The density of the colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

#### **4. WARNINGS AND PRECAUTIONS**

##### **FOR *IN VITRO* DIAGNOSTIC USE ONLY**

**This kit contains materials of human origin which have been tested and gave a negative for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.**

**Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.**

##### **Health and Safety Information**

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instruments.
4. Refer to the Safety Data Sheet (available on DIESSE website: [www.diesse.it](http://www.diesse.it)) for the safety characteristics of the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued.

Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

### Analytical Precautions

Bring the devices to room temperature (18-30°C) for at least 30 minutes before the use; use within 60 min.


1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.  
The Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a Release (date) above the one shown in the table published on the Diesse website.  
(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument and do not use devices with defective labels.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument.
9. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
10. Do not use hemolyzed, lipemic or, icteric samples with a concentration of interferents higher than that tested (according to the indications given in the "Analytical specificity" chapter).
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer REF 83606.**

### 5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

**DD** DEVICES 6 packages each containing 6 devices.  
Use: **equilibrate a package at room temperature**, open the package and remove the devices required;

replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

### Description:

Position 8: Space for application of bar code label

Position 7: empty

Position 6: MICROPLATE WELL coated with Helicobacter pylori antigen.

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

**Contents:** Tetramethylbenzidine 0.26 mg/mL and H<sub>2</sub>O<sub>2</sub> 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

**Contents:** Proteic solution containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

**Contents:** anti-human IgA monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum

### **CALIBRATOR** CALIBRATOR 1 x 0.175 mL

**Contents:** Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

### **CONTROL +** POSITIVE CONTROL 1 x 0.425 mL

**Contents:** Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

Confidence in measurements of Calibrator and Positive control is established with traceability to measurement standards as follows.

Calibrator and Positive control are produced diluting a known concentration of human antigens in its own stabilizing medium. The relative exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared and characterized, checking the consensus with a reference sera panel with different antigens levels.

### MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607

- CHORUS/CHORUS TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

## 6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated, and the run validated using the control serum (see point 9, Test validation). The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation:

|                  |                  |
|------------------|------------------|
| DEVICES          | 8 weeks at 2/8°C |
| POSITIVE CONTROL | 8 weeks at 2/8°C |
| CALIBRATOR       | 8 weeks at 2/8°C |

## 7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

## 8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

## 9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554

email: scientificsupport@diesse.it

## 10. INTERPRETATION OF THE RESULTS

The Chorus instrument gives the result as an INDEX (ratio between the O.D. value of the examined samples and the O.D. value of the Cut-off) which can be used as the quantitative measurement since it is proportional to the quantity of specific IgA.

The sample will give the following results:

POSITIVE: when the Index is >1.2

NEGATIVE: when the Index is < 0.8

DOUBTFUL: when the values are between 0.8 and 1.2

In the case of a doubtful result, repeat the test. If the test remains doubtful, collect a new sample after 1-2 weeks.

## 11. LIMITATIONS OF THE PROCEDURE

Serum samples collected during the acute phase of the infection could give negative results with this technique. Highly lipemic, hemolyzed or heat-inactivated sera may cause false results and should not be used.

As with all serological tests, the results should be interpreted in conjunction with other clinical and diagnostic data.

## 12. ANALYTICAL SPECIFICITY

The analytical specificity of the H. Pylori IgA kit was evaluated in an inhibition test. Campylobacter coli, Campylobacter fetus, Campylobacter jejuni and E. coli were added to patients' sera. No differences were observed between the reference samples and those containing the antigens reported above. It can therefore be concluded that the test is specific for anti-Helicobacter pylori IgA antibodies-

## 13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

Table 1

|        |   | REFERENCE |    |
|--------|---|-----------|----|
|        |   | +         | -  |
| DIESSE | + | 19        | 1  |
|        | - | 2         | 41 |

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Valore Predittivo Positivo (PPV): 95% CI95%: 89.6-100.0

Valore Predittivo Negativo (NPV): 95.3% CI95%: 90.1-100.0

manufacturer and competent authority of your Member State.

## 17. SUMMARY OF SAFETY AND PERFORMANCE

This document, which will be made available on the EUDAMED database (when this is fully implemented and functioning), is part of the Technical Documentation and can be requested from the manufacturer.

## 14. PRECISION AND REPEATABILITY

**Table 2 Within-run precision**

|             | No. of replicates | Mean Index | St. Dev. | CV% |
|-------------|-------------------|------------|----------|-----|
| Lot no. 023 | 9                 | 3.2        | 0.24     | 8   |
| Lot no. 024 | 9                 | 2.9        | 0.22     | 8   |
| Lot no. 025 | 9                 | 3.1        | 0.40     | 13  |

**Table 3 Between-run precision and Precision between batches**

| Sample | Mean Index |         |         | Mean | St. Dev. | CV% |
|--------|------------|---------|---------|------|----------|-----|
|        | Lot 023    | Lot 024 | Lot 025 |      |          |     |
| PYA 1  | 0.4        | 0.3     | 0.3     | 0.3  | 0.06     | 20  |
| PYA 2  | 1.2        | 1.2     | 1.2     | 1.2  | 0.00     | -   |
| PYA 3  | 2.7        | 2.8     | 2.8     | 2.8  | 0.06     | 2   |

## 15. REFERENCES

1. Marshall B.J. and Warren J.R.: Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. *Lancet* i, 1311 (1984).
2. Jones D.M., Lessels A.M., Eldridge J.: *Campylobacter*-like organisms on the gastric mucosa: culture, histological and serological studies. *J. Clin. Pathol.* 37: 1002 (1984).
3. Blaser M.J.: *H. pylori* and the pathogenesis of gastroduodenal inflammation. *J. Inf. Dis.* 161: 626 (1990).
4. Valle J., Seppälä K., Sipponen P., Kasunen T.U.: Disappearance of gastritis after eradication of *H. pylori*: a morphometric study. *Scand. J. Gastroenterology* 26: 1057 (1991).
5. G.B. Wisdom: Enzyme-Immunoassay. *Clin. Chem.* 22: 1243 (1976).
6. CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

## 16. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the

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