Pros and Cons of serological testing in syphilis diagnosis and follow up

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SUMMARY

Since a proper diagnosis of syphilis is often difficult due to the wide variability of both clinical picture and laboratory test results, early recognition of infection caused by *Treponema pallidum* is crucial for a timely and effective treatment. In most cases, definitive diagnosis relies upon serological testing. A screening ELISA test, coupled with a quantitative RPR test and specific IgM antibodies detection, is currently regarded as the basic diagnostic procedure. In addition, a quantitative particle agglutination TP-PA test, FTA-abs IgG test and, eventually, a western-blot IgG and IgM test, allow to achieve a whole serological pattern for each patient at the time of first diagnosis.

In this study, a group of serum samples (n=107) and cerebro-spinal fluid (n=3) were retrospectively analyzed using the above mentioned tests. A population of 19 patients whose clinical picture was unremarkable for syphilis, showed border-line values at screening and negative results on confirmation tests. Thirty-three out of 91 luetic patients were diagnosed as primary or early secondary syphilis, 36 as latent syphilis, 3 as neurosyphilis, and 3 were neonates with passive specific immunization.

Quantitative RPR test and detection of specific IgM antibodies exhibited extremely high values in all 33 primary syphilis patients; a whole positive luetic pattern was also obtained by confirmation tests. Searching for IgM antibodies, a capture elisa test compared with a single device rapid elisa test showed an overall concordance of 98.1%.

In luetic patients other than primary syphilis, quantitative RPR test and detection of specific IgM antibodies provided less relevant values and a low prevalence pattern, whereas TP-PA and FTA-abs tests showed persistent positives results.

In the follow up of 19 initially treated patients, quantitative RPR values and specific IgM antibodies index showed a slow, progressive decrease until negative.

Conclusion: a comprehensive initial sieroluetic framework is a relevant diagnostic tool for primary diagnosis; along the follow up, elisa screening test, TP-PA and FTA-abs IgG have a poor clinical significance; instead, quantitative RPR titers and specific IgM antibody index change over time, and as such, they are reliable tools to evaluate the therapeutic efficacy of treatment and a rational management of luetic patients.

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BIBLIOGRAFIA

- Castro R, Prieto ES, Aguas MJ, et al. Evaluation of the *Treponema pallidum* particle agglutination technique (TP-PA) in the diagnosis of neurosyphilis. *J Clin Lab Anal* 2006; 20(6): 233-8.
- Centers for Disease Control and Prevention. Sexually transmitted diseases treatment guidelines 2002. *MMWR* 2002; 51: 1-80.
- Cole MJ, Perry KR, Parry JV. Comparative evaluation of 15 serological assays for the detection of syphilis infection. *Eur J Clin Microbiol Infect Dis* 2007 Oct; 26(10): 705-13.
- Creagan L, Bauer HM, Samuel MC, Klausner J, Liska S, Bolan G. An evaluation of the relative sensitivities of the venereal disease research laboratory test and the treponema pallidum particle agglutination test among patients diagnosed with primary syphilis. *Sex Transm Dis* 2007 Dec; 34(12): 1016-8.
- Fujimori C, Yukimasa N, Miura K, et al. Evaluation of the anti-TP antibody latex agglutination immunoassay in routine testing and a clinical viewpoint. *Rinsho Byori* 2009 mar; 57(3): 206-12.
- Gesau A, Kittler H, Hein U, Dangl-Erlach E, Stingl G, Tschachler E. Biological false-positive tests comprise a high proportion of Venereal Disease Research Laboratory reactions in an analysis of 300,000 sera. *Int J STD AIDS* 2005 Nov; 16(11): 722-6.
- Goh BT, Thornton AC. Antenatal screening for syphilis. *SEX Transm Infect* 2007; 83: 345-6.
- Kingston M, French P, Goh B, et al. UK National Guidelines on the management of Syphilis 2008. *Int J STD AIDS* 2008; 19: 729-40.
- Koek AG, Bruisten SM, Dierdorff M, van Dam AP, Templeton K. Specific and sensitive diagnosis of syphilis using a real-time PCR for *Treponema pallidum*. *Clin Microbiol Infect Dis* 2006; 12: 1221-40.
- Liu H, Rodes B, Chen CY, Steiner B. New tests for syphilis: rational design of a PCR method for detection of *Treponema pallidum* in clinical specimens using unique regions of the DNA polymerase I gene. *J Clin Microbiol* 2001; 39: 1941-6.
- ManaviK, Young H, McMillan A. The sensibility of syphilis assays in

- detecting different stages of early syphilis. *Int J STD AIDS* 2006 Nov; 17(11): 768-71.
12. Marangoni A, Moroni A, Tridapalli E, et al. Antenatal syphilis serology in pregnant women and follow-up of their infants in northern Italy. *Clin Microbiol Infect* 2008; 14: 1065-8.
 13. Marangoni A, Sambri V, Accardo S, et al. Evaluation of LIAISON® *Treponema* Screen, a novel recombinant antigen-based chemiluminescence immunoassay for the laboratory diagnosis of syphilis. *Clin Diagn Lab Immunol* 2005; 12: 1231-4.
 14. Marfin AA, Liu H, Sutton MY, Steiner B, Pillay A, Markowitz LE. Amplification of the DNA polymerase I gene of *Treponema pallidum* from whole blood of persons with syphilis. *Diagn Microbiol Infect Dis* 2001; 40: 163-6.
 15. McMillan A, Young H. Qualitative and quantitative aspects of the serological diagnosis of early syphilis. *Int J STD AIDS* 2008; 19(9): 620-4.
 16. Müller I, Brade V, Hagedorn HJ, et al. Is serological testing a reliable tool in laboratory diagnosis of syphilis? Meta-analysis of eight external quality control surveys performed by the German Infection Serology Proficiency Testing Program. *J Clin Microbiol* 2006; 44(4): 1335-41.
 17. Noordhoek GT, Wolters EC, de Jonge ME, Embden JD. Detection by polymerase chain reaction of *Treponema pallidum* DNA in cerebrospinal fluid from neurosyphilis patients before and after antibiotic treatment. *J Clin Microbiol* 1991; 29: 1976-84.
 18. Palmer HM, Higgins SP, Herring AJ, Kingston MA. Use of PCR in the diagnosis of early syphilis in the United Kingdom. *Sex Transm Infect* 2003; 79: 479-83.
 19. Sambri V, Marangoni A, Eyer C, et al. Western immunoblotting with five *Treponema pallidum* recombinant antigens for serological diagnosis of syphilis. *Clin Diagn Lab Immunol* 2001; 8: 534-9.
 20. Simms I, Ward H. Congenital syphilis in the United Kingdom. *Sex Transm Infect* 2006; 82: 1.
 21. Zetola NM, Engelman J, Jensen TP, Klausner JD. Syphilis in the United States: an update for clinicians with an emphasis on HIV coinfection. *Mayo Clin Proc* 2007; 82(9): 1091-102.