CHORUS

Helicobacter pylori Ag

DIESSE

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infection. The Maastricht III Consensus Report. Gut. 2007; 56: 772-781.

5. J.G. Kusters, A.H.M. van Vilet, E.J. Kuipers. Pathogenesis of *Helicobacter pylori* Infection. Clinical Microbiology Reviews. 2006; 19(3): 449-490.
6. S.A. Chisholm, C.L. Watson, E.L. Teare, S. Saverymuttu, R.J. Owen. Non-invasive diagnosis of *Helicobacter pylori* infection in adult dyspeptic patients by stool antigen detection: does the rapid immunochromatography test provide a reliable alternative to conventional ELISA kits? Journal of Medical Microbiology. 2004; 53: 623-627.

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INSTRUCTIONS FOR USE

CHORUS Helicobacter pylori Ag

For the qualitative determination of Helicobacter pylori antigens in stools.

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of Helicobacter pylori antigens in stools, using a disposable device applied on the Chorus TRIO instruments.

2. INTRODUCTION

In 1983, Warren and Marshall identified *Helicobacter pylori*, a new gram-negative bacterial pathogen, in patients suffering from gastritis, and this finding led to studies on the relationship between bacterial infection and chronic gastric disease. The pathogen has been shown to be associated with peptic ulcer, chronic gastritis type B and duodenitis. It has been demonstrated that in patients with gastritis, eradication of the bacteria led to healing of the anatomical lesion.

Diagnostic procedures for the detection of the organism generally involve invasive (gastroscopic) techniques for sample collection.

A valid alternative method, non-invasive, is represented by the test for the determination of the antigens contained in the stool samples. The recent Maastricht 2-2000 Consensus Report recommends the use of the stools antigen test as an aid in the diagnosis of *H. pylori* infections.

3. PRINCIPLE OF THE METHOD

The Chorus *Helicobacter pylori* Ag device is ready to use for the detection of *Helicobacter pylori* antigens present in human feces, in the Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay).

The anti-*H. pylori* monoclonal capture antibodies are bound to the solid phase. The incubation is performed with the samples and the conjugate composed of anti-*H. pylori* monoclonal antibodies conjugated to horseradish peroxidase: the specific antigen is bound to the monoclonal antibody bound to the solid phase and to the conjugate through incubation. After washing to eliminate the proteins which have not reacted, the peroxidase substrate is added.

The enzymatic reaction is stopped by the adding of the Stop Solution, which makes the solution turn to yellow. The colour which develops is proportional to the concentration of specific antigens present in the sample.

The disposable devices contain all the reagents to perform the test in the Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)
5. Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).


6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use the device after the expiry date.
11. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

DD DEVICES 2 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with a monoclonal antibody specific for H.pylori

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: STOP SOLUTION

Contents: Sulfuric acid solution 0.3 M

Position 2: CONJUGATE

Contents: ready-to-use solution containing anti-H.pylori monoclonal antibodies labelled with horseradish peroxidase and preservative.

Position 1: EMPTY WELL

In which the operator must place the undiluted sample.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 2 x 0.250 ml

Contents: Diluted H.pylori inactivated cells. Freeze-dried.

Use: reconstitute with the volume of Stool Diluent/Negative Control reported on the label.

Wait at least 5 minutes without shaking; then mix by inversion avoiding froth formation.

CONTROL + POSITIVE CONTROL 3 x 0.250 ml

Contents: Diluted H.pylori inactivated cells. Freeze-dried.

Use: reconstitute with the volume of Stool Diluent/Negative Control reported on the label.

Wait at least 5 minutes without shaking; then mix by inversion avoiding froth formation.

STOOL DILUENT/NEGATIVE CONTROL

H. PYLORI STOOL DILUENT – NEGATIVE CONTROL
1 x 8 ml

Contents: ready-to-use buffer solution and preservative.

Use: Stool Diluent (see chapter 8 Assay Procedure).

To be used to reconstitute Calibrator and Positive Control (see CALIBRATOR and CONTROL+).

The diluent can also be used as negative control solution.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-1000 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	1 day at 2/8°C after reconstitution
POSITIVE CONTROL	1 day at 2/8°C after reconstitution

7. SPECIMEN COLLECTION AND STORAGE

The type of sample is composed of feces handled with all precautions dictated by good laboratory practice.

Keep the samples at 2/8°C and test within 48 hours after collection. If testing cannot be performed within 48 hours, store the samples at -20°C or lower.

Minimize specimen freezing thawing which may cause degradation/teolysis of the antigen and result in loss of activity.

Do not keep the samples in auto-defrosting freezers.

8. ASSAY PROCEDURE

1. Set up one dilution tube for each specimen to be tested. 1.5 ml tubes are recommended for this purpose. Add 500 µl of Stool Diluent/Negative Control to each tube.
2. **Formed samples:** use a wooden applicator stick or a disposable teaspoon to transfer the fecal specimen to the tube. Transfer approximately 0.1 to 0.15 g of specimen (about the size of a small pea) into 500 µl of Stool Diluent/Negative Control, removing all the material from the collector.
Liquid samples: transfer 100 µl of specimen to the tube. Make sure the liquid specimens are evenly suspended.
3. **Thoroughly mix (vortex) the fecal specimen for 15 seconds.**
4. Let the tubes stand in vertical position for 10 minutes.
5. Centrifuge the tubes at 1000 g for 30 seconds. Ensure that the formed supernatant does not contain large particulate material.

6. Store the diluted samples at 2/8°C until the test is performed.
7. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
8. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
9. Dispense 130 µl of undiluted test sample in well no. 1 of each device; at each change of batch, use a device for the calibrator.
10. Place the devices in the Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 0.5
 NEGATIVE: when the result is ≤ 0.5

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient. The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and/or other clinical diagnostic evaluation.

ATTENTION: the Chorus Helicobacter pylori Ag test (REF 81063) can only be performed in a dedicated run or together with Chorus Clostridium difficile GDH (REF 81168) and Chorus Clostridium difficile A/B Toxins (REF 81170). Otherwise none of the programmed tests will be performed.

12. METHOD COMPARISON

In an experimentation 51 samples have been tested with Diesse kit and with a competitor kit. Data are summarized in the following table:

		Reference		
		+	-	Total
		37	0	37
Diesse	+	37	0	37
	-	0	14	14
Total		37	14	51

Percent Positive Agreement (~Diagnostic Sensitivity):

100% CI_{95%}: 90.6-99.9

Percent Negative Agreement: (~Diagnostic Specificity):

100% CI_{95%}: 78.4-99.7

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

13. PRECISION AND REPEATABILITY

Sample	Within-run Precision	
	Mean (Index)	CV%
1	2.5	5.2
2	1.6	6.9
3	0.2	5.0

Sample	Between-run precision	
	Mean (Index)	CV%
1	2.5	6.8
2	1.4	11.4
3	0.2	5.0

Sample	Precision between batches		Precision between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	2.4	6.3	2.4	7.5
2	1.3	12.3	1.5	12.0
3	0.2	10.0	0.2	5.0

14. REFERENCES

1. M. J. Blaser . Who are we? Indigenous microbes and the ecology of human diseases. EMBO Reports. 2006; 7(10): 956-60.
2. P. Bytzer, J.F. Dahlerup, J.R. Eriksen, D.E. Jarbøl, S. Rosenstock, S. Wildt. Diagnosis and treatment of Helicobacter pylori infection. Danish Medical Bulletin. 2011; 58(4): C4271.
3. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of Helicobacter pylori infection. The Maastricht 2-2000 Consensus Report. Aliment Pharmacol Ther. 2002; 16: 167-80.
4. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of Helicobacter pylori infection. The Maastricht III Consensus Report. Gut. 2007; 56: 772-781.
5. J.G. Kusters, A.H.M. van Vilet, E.J. Kuipers. Pathogenesis of Helicobacter pylori Infection. Clinical Microbiology Reviews. 2006; 19(3): 449-490.
6. S.A. Chisholm, C.L. Watson, E.L. Teare, S. Savarymuttu, R.J. Owen. Non-invasive diagnosis of Helicobacter pylori infection in adult dyspeptic patients by stool antigen detection: does the rapid immunochromatography test provide a reliable alternative to conventional ELISA kits? Journal of Medical Microbiology. 2004; 53: 623-627.

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28. LITERATUR

1. M. J. Blaser . Who are we? Indigenous microbes and the ecology of human diseases. EMBO Reports. 2006; 7(10): 956-60.
2. P. Bytzer, J.F. Dahlerup, J.R. Eriksen, D.E. Jarbøl, S. Rosenstock, S. Wildt. Diagnosis and treatment of Helicobacter pylori infection. Danish Medical Bulletin. 2011; 58(4): C4271.
3. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of Helicobacter pylori infection. The Maastricht 2-2000 Consensus Report. Aliment Pharmacol Ther. 2002; 16: 167-80.
4. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of Helicobacter pylori infection. The Maastricht III Consensus Report. Gut. 2007; 56: 772-781.
5. J.G. Kusters, A.H.M. van Vilet, E.J. Kuipers. Pathogenesis of Helicobacter pylori Infection. Clinical Microbiology Reviews. 2006; 19(3): 449-490.
6. S.A. Chisholm, C.L. Watson, E.L. Teare, S. Saverymuttu, R.J. Owen. Non-invasive diagnosis of Helicobacter pylori infection in adult dyspeptic patients by stool antigen detection: does the rapid immunochromatography test provide a reliable alternative to conventional ELISA kits? Journal of Medical Microbiology. 2004; 53: 623-627.

DIESSE Diagnostica Senese S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (SIENA)
Italien

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1. M. J. Blaser . Who are we? Indigenous microbes and the ecology of human diseases. EMBO Reports. 2006; 7(10): 956-60.
2. P. Bytzer, J.F. Dahlerup, J.R. Eriksen, D.E. Jarbøl, S. Rosenstock, S. Wildt. Diagnosis and treatment of Helicobacter pylori infection. Danish Medical Bulletin. 2011; 58(4): C4271.
3. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of Helicobacter pylori infection. The Maastricht 2-2000 Consensus Report. Aliment Pharmacol Ther. 2002; 16: 167-80.
4. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of Helicobacter pylori infection. The Maastricht III Consensus Report. Gut. 2007; 56: 772-781.
5. J.G. Kusters, A.H.M. van Vilet, E.J. Kuipers. Pathogenesis of Helicobacter pylori Infection. Clinical Microbiology Reviews. 2006; 19(3): 449-490.
6. S.A. Chisholm, C.L. Watson, E.L. Teare, S. Savarymattu, R.J. Owen. Non-invasive diagnosis of Helicobacter pylori infection in adult dyspeptic patients by stool antigen detection: does the rapid immunochromatography test provide a reliable

alternative to conventional ELISA kits? Journal of Medical Microbiology. 2004; 53: 623-627.

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14. BIBLIOGRAPHIE

1. M. J. Blaser . Who are we? Indigenous microbes and the ecology of human diseases. EMBO Reports. 2006; 7(10): 956-60.

2. P. Bytzer, J.F. Dahlerup, J.R. Eriksen, D.E. Jarbøl, S. Rosenstock, S. Wildt. Diagnosis and treatment of *Helicobacter pylori* infection. Danish Medical Bulletin. 2011; 58(4): C4271.
3. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of *Helicobacter pylori* infection. The Maastricht 2-2000 Consensus Report. Aliment Pharmacol Ther. 2002; 16: 167-80.
4. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of *Helicobacter pylori* infection. The Maastricht III Consensus Report. Gut. 2007; 56: 772-781.
5. J.G. Kusters, A.H.M. van Vilet, E.J. Kuipers. Pathogenesis of *Helicobacter pylori* Infection. Clinical Microbiology Reviews. 2006; 19(3): 449-490.
6. S.A. Chisholm, C.L. Watson, E.L. Teare, S. Saverymuttu, R.J. Owen. Non-invasive diagnosis of *Helicobacter pylori* infection in adult dyspeptic patients by stool antigen detection: does the rapid immunochromatography test provide a reliable alternative to conventional ELISA kits? Journal of Medical Microbiology. 2004; 53: 623-627.

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14. BIBLIOGRAFIA

1. M. J. Blaser. Who are we? Indigenous microbes and the ecology of human diseases. EMBO Reports. 2006; 7(10): 956-60.
2. P. Bytzer, J.F. Dahlerup, J.R. Eriksen, D.E. Jarbøl, S. Rosenstock, S. Wildt. Diagnosis and treatment of Helicobacter pylori infection. Danish Medical Bulletin. 2011; 58(4): C4271.
3. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of Helicobacter pylori infection. The Maastricht 2-2000 Consensus Report. Aliment Pharmacol Ther. 2002; 16: 167-80.
4. P. Malfertheiner, F. Megraud, O. Morain. Current concepts in the management of Helicobacter pylori infection. The Maastricht III Consensus Report. Gut. 2007; 56: 772-781.
5. J.G. Kusters, A.H.M. van Vilet, E.J. Kuipers. Pathogenesis of Helicobacter pylori Infection. Clinical Microbiology Reviews. 2006; 19(3): 449-490.
6. S.A. Chisholm, C.L. Watson, E.L. Teare, S. Saverymuttu, R.J. Owen. Non-invasive diagnosis of Helicobacter pylori infection in adult dyspeptic patients by stool antigen detection: does the rapid immunochromatography test provide a reliable alternative to conventional ELISA kits? Journal of Medical Microbiology. 2004; 53: 623-627.

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5. J.G. Kusters, A.H.M. van Vilet, E.J. Kuipers. Pathogenesis of *Helicobacter pylori* Infection. *Clinical Microbiology Reviews*. 2006; 19(3): 449-490.
6. S.A. Chisholm, C.L. Watson, E.L. Teare, S. Saverymuttu, R.J. Owen. Non-invasive diagnosis of *Helicobacter pylori* infection in adult dyspeptic patients by stool antigen detection: does the rapid immunochromatography test provide a reliable alternative to conventional ELISA kits? *Journal of Medical Microbiology*. 2004; 53: 623-627.

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