CHORUS

HERPES SIMPLEX 1+2 IgM

REF 81021 (36 tests)

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3. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
4. Warford AL et al.: Comparison of two commercial enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. AM J. Clin. Pathol. 19*896 Feb; 85(2):229-32.

INSTRUCTIONS FOR USE

CHORUS HERPES SIMPLEX 1+2 IgM

REF 81021

(English)

1. INTENDED USE

QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO HERPES SIMPLEX VIRUS 1+2 IN HUMAN SERUM USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS INSTRUMENT

2. INTRODUCTION

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular cheratitis, and lesions of the central nervous system.

HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion. The assay of specific IgM is important for the diagnosis of neonatal infection and encephalitis caused by HSV. Moreover, the presence of specific IgM indicates viral activity in progress, although it is not possible to distinguish between primary infection and reactivation.

3. PRINCIPLE OF THE METHOD

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The antigen, composed of a mix of purified Herpes Simplex 1 and extracted Herpes Simplex 2, is bound to the solid phase. Through incubation with human serum diluted in a diluent blocking the IgG, the IgM bind to the antigen.

After washings to eliminate the proteins that did not react the incubation is performed with the conjugate composed of human anti IgM monoclonal antibodies conjugated with radish peroxidase.

The conjugate that did not bind is eliminated and, after washing, the substrate for peroxidase is added.


The blue colour developed is proportional to the concentration of the specific antibodies present in the serum examined.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Use: **equilibrate a package at room temperature**, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: empty

Position 6: MICROPLATE WELL coated with a mix of Herpes Simplex 1 and 2 antigens

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE.

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: DILUENT FOR THE SAMPLES

Contents: Proteic solution containing human anti IgG, containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: monoclonal anti IgM antibodies labelled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum

CALIBRATOR CALIBRATOR 1 x 0.175 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus / Chorus Trio Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

6. PRECAUTIONS

ONLY FOR DIAGNOSTIC USE IN VITRO.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.

4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to be used to room temperature (18-30°C) before use; use within 60 min.

1. **Check that the substrate contained in well 4 is colourless. Discard devices which do not have this characteristic.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (50 µl) in well 1 of the device (see figure).
12. Do not use the device after the expiry date

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the qualitative result as an INDEX (ratio between the OD value of the test sample and that of the Cut-Off).

The test serum examined can be interpreted as follows:

POSITIVE when the index is > 1.1

NEGATIVE when the index is < 0.9

DOUBTFUL index between 0.9 and 1.1

In the case of a doubtful result, repeat the test. If the result remains doubtful, repeat the sampling after 15 days at least.

11. LIMITATIONS OF THE PROCEDURE

The results of the assay must be interpreted with caution and in conjunction with information available from the clinical evaluation and other diagnostic data.

Sera from patients in an early or late stage of the disease could give a repeatedly negative result close to the Cut-off value. In such cases, a confirmation of the result is recommended, by testing a new sample drawn from the same patient at 7-10 days interval from the first sample.

12. ANALYTICAL SPECIFICITY

The presence in the serum sample of potentially interfering substances such as PCR, Bilirubin, Rheumatoid Factor, Haemolysis and Triglycerides, does not interfere with the test.

Hyper-immune sera for other herpetic pathologies (CMV, RV, VZV, EBV) do not interfere with the test.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

285 samples were tested in an experimentation both with the Diesse kits and with reference kits (Behring Enzygnost HSV1+2 IgM).

The data of the experimentation are summarized in the following table:

Table 1.

		Reference	
		+	-
	Diesse		
	+	22	10
	-	1	252

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