CHORUS

EPSTEIN-BARR EBNA IgG

REF 81057 (36 tests)

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15. BIBLIOGRAFIA

1. A. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).

2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. *J. Virol. Methods* 21: 133 (1988).
3. C. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. *J. Inf. Dis.* 151: 984 (1985).
4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. *J. Immunol. Methods* 67: 145 (1984).

INSTRUCTIONS FOR USE

CHORUS EPSTEIN-BARR EBNA IgG

REF 81057

(English)

1. INTENDED USE

QUALITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO EPSTEIN BARR NUCLEAR ANTIGEN (EBNA) IN HUMAN SERUM, USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS INSTRUMENTS

2. INTRODUCTION

Epstein Barr Virus (EBV) is a herpesvirus which causes infectious mononucleosis (IM). It is also associated with Burkitt's lymphoma, nasopharyngeal carcinoma and lymphatic proliferative syndromes in immunodepressed patients. The virus is widespread throughout the world and 80-90% of the population is serum-positive.

The laboratory diagnosis of IM is traditionally performed by detecting heterophile antibodies which develop in the serum during the course of the infection, and which agglutinate horse erythrocytes. However, these antibodies may not always be present in patients affected by IM, particularly if below 14 years of age; furthermore, they may also persist for over a year after the infection. The determination of heterophile antibodies alone may therefore lead to an erroneous diagnosis. It is therefore important to determine the presence of antibodies towards the viral antigens. In particular, the detection of antibodies directed to the "Viral Capsid Antigen" (VCA) and the nuclear antigen (EBNA) is particularly useful.

During the course of IM, the IgM- and IgG-class antibodies to VCA appear early, while the IgG to EBNA develop later during the infection. The presence of IgM against VCA in the absence of IgG against EBNA therefore indicates that there is a current infection, while the presence of IgG against both VCA and EBNA is indicative of a prior infection.

3. PRINCIPLE OF THE METHOD

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The EBNA, obtained by the recombinant technique, is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgG monoclonal antibodies conjugated to horse radish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.


The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Use: **equilibrate a package at room temperature**, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: empty

Position 6: MICROPLATE WELL coated with a pool of EBV EBNA recombinant antigens

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum

CALIBRATOR CALIBRATOR 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus / Chorus Trio Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation) .

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

6. PRECAUTIONS

ONLY FOR DIAGNOSTIC USE IN VITRO.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to be used to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, incompletely coagulated serum or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (50 µl) in well 1 of the device (see figure).
12. Do not use the device after the expiry date.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (ratio between the OD value of the test sample and that of the Cut-Off) which can be used as a quantitative measure, as it is proportional to the amount of specific IgG present.

The test serum can be interpreted as follows:

POSITIVE: when the ratio is > 1.2

NEGATIVE: when the ratio is < 0.8 .

DOUBTFUL: for all values between 0.8 and 1.2

In the case of a doubtful result, repeat the test. If the test remains doubtful, collect a new sample after 1-2 weeks.

11. LIMITATIONS OF THE PROCEDURE

The results obtained can only be used as a guide in the diagnosis and must always be evaluated together with results of other diagnostic procedures.

There is the possibility of cross reactions with samples containing anti-E. coli antibodies.

Performance characteristics have not been determined for a visual interpretation, nor have they been evaluated in patients affected by nasopharyngeal carcinoma, Burkitt's lymphoma or other lymphadenopathies associated with EBV, apart from EBV-related mononucleosis.

As anti-EBNA IgG antibodies are present in normal convalescent sera, the diagnosis cannot be established on the basis of a single parameter. An accurate interpretation of an EBV infection must be based on the results of VCA IgG, VCA IgM, EBNA IgG, EA IgG ed EA IgM results, in addition to heterophyle antibodies.

The test must be performed on serum. The use of blood or plasma has not been established.

12. ANALYTICAL SPECIFICITY

Samples which were negative for anti-EBNA e anti-VCA IgG were tested for the determination of IgG anti-CMV, anti-Herpes Simplex and anti-Varicella. In no case did the presence of these antibodies interfere with the Enzywell EBNA IgG kit.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

101 samples, the result of which was known with another commercial method, were tested. 92 were positive and 9 negative, in agreement between the two methods. The results are summarized in the following table.

		REFERENCE	
		+	-
		92	0
CHORUS	+		
	-	0	9

The kit shows a sensitivity of 100% (CI 95%: 96–100) and a specificity of 100% (CI 95%: 70–100).

14. PRECISION

Within-run precision

	No. of replicates	Mean Index	St. Dev.	CV%
Lot no. 024	9	2.7	0.12	4
Lot no. 025	9	2.6	0.15	6
Lot no. 026	9	3.0	0.16	5

Between-run precision

	DAY1		DAY2		DAY3		DAY4		DAY5		Mean	SD	CV%
	INDEX		INDEX		INDEX		INDEX		INDEX				
	1	2	1	2	1	2	1	2	1	2			
EBN1	0.2	0.1	0.3	0.2	0.2	0.1	0.1	0.1	0.1	0.1	0.2	0.07	35.0*
EBN2	0.9	0.9	0.9	1.1	1.1	0.9	1.0	1.0	0.9	0.9	1.0	0.08	8.0
EBN3	3.7	3.7	4.3	3.8	4.3	3.8	3.9	3.1	3.7	3.5	3.8	0.35	9.2

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero

15. REFERENCES

1. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
3. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).

4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).

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16. POUŽITÁ LITERATURA

- Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
- J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
- Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
- J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).

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16. LITERATUR

- A. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
- J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).

7. C. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
8. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).

DIESSE Diagnostica Senese
Strada dei Laghi 39
53035 Monteriggioni (SIENA) Italien
Tel. 0039-0577-587111

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1. A. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
3. C. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).

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15. BIBLIOGRAFIA

1. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. *J. Med. Virology* 43: 238 (1994).
2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. *J. Virol. Methods* 21: 133 (1988).
3. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. *J. Inf. Dis.* 151: 984 (1985).
4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. *J. Immunol. Methods* 67: 145 (1984).

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16. REFERINTE

1. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
3. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).

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