CHORUS

Legionella pneumophila IgM

DIESSE
REF 81093
REF 81093/12

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15. BIBLIOGRAFIA

1. G.B. Wisdom. 1976. Enzyme-Immunoassay. Clin. Chem. 22: 1243.
2. A.K. Malan, T.B. Martins, T.D. Jaskowski, H. R. Hill, C.M. Litwin. 2003 Comparaison of two Commercial Enzyme-Linked Immunosorbent Assay with an immunofluorescence assay for detection of *Legionella pneumophila* Types 1 to 6. J. Clin. Microb. 41: 3060-3063.
3. B. Neumeister, M. Susa, B. Nowak, E. Straube, G. Ruckdeschel, J. Hacher, R. Marre. 1995 Enzyme immunoassay for detection of antibodies against isolated flagella of *Legionella pneumophila*. Eur. J. Clin. Microbiol. Infect. Dis. 14 (9) : 764-7.

4. N. Barka, J.P. Tomasi, S. Stadtsbaeder 1986 ELISA using whole *Legionella pneumophila* cell as antigen J. Immunol. Methods 93: 77-81
5. M. S. Lever. 1993. Production and characterisation of a *Legionella pneumophila* specific monoclonal antibody. FEMS Microbiol. Lett. 1 : 5-9.
6. L.B. Reller, M.P. Weinstein 2003 Diagnosis of *Legionella* Infection Clin. Infect. Dis. 36:64-69

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CHORUS LEGIONELLA PNEUMOPHILA IgM

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgM-class antibodies to *Legionella pneumophila* in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

CHORUS Legionella pneumophila IgM For the Qualitative Determination of IgM Antibodies to *Legionella pneumophila*

For In Vitro Diagnostic Use Only

2. INTRODUCTION

Legionnaire's disease is a serious infectious disease with a high mortality rate. Legionellosis is the definition given to all forms of disease caused by gram-negative bacteria of the *Legionella* genus. The total mortality rate is 5-15% which can reach 30-50% in hospital populations and 70-80% in patients in a serious clinical state or if treated late. The species most frequently involved in human pathologies is *Legionella pneumophila*, although other species have been isolated. The most serious form of *Legionella* infection is the so-called Legionnaire's disease, characterized by acute pneumonia which becomes evident 2-10 days after contact with the bacteria. Another characteristic clinical picture is Pontiac fever, an acute extra-pulmonary infection. There are also forms of subclinical infection without symptoms, which can only be detected by demonstrating serum-conversion.

Legionellae are widespread in nature, where they are generally found in the presence of water (lakes, rivers, thermal springs, water-bearing layers, and damp areas in general). Starting from these sources, *Legionellae* can colonize artificial water supplies (civic distribution systems for drinking water, water systems of individual buildings, air-conditioning plants, swimming pools, fountains, etc.), which act as amplifiers and disseminators of the microorganism. The infection can affect the whole population. The subjects considered to be most at risk are elderly males who are smokers, consumers of alcohol, affected by chronic diseases (obstructive bronchopulmonary infections, cardiovascular and renal disorders, diabetes, etc.), and those with acquired immunodeficiency following therapy (organ transplant, treatment with steroids and anti-tumoral drugs, etc.), or affected by HIV.

80 % of all the cases of legionellosis reported can be attributed to *Legionella pneumophila* serotype 1, and serum conversion is

generally considered highly predictive of the infection. Serological diagnosis of *Legionella pneumophila* after the first sample tested, at the beginning of the disease, must be confirmed by an increase in the positive response in subsequent samples tested three weeks later and after 6-8 weeks, in as far as clinical investigations have shown that significant levels of specific antibodies sometimes appear late, and there is the possibility of anti-*Legionella pneumophila* antibodies in the absence of legionellosis.

The antibody response to infection is generally a mixture of IgA, IgM and IgG. Detection of the presence of IgM can support the diagnosis of an acute infection even if the IgM may remain in circulation for a long period after infection.

3. PRINCIPLE OF THE ASSAY

The test is based on the ELISA principle (Enzyme Linked ImmunoSorbent Assay).

The antigen, which is obtained by culturing *Legionella pneumophila* serotype 1 and inactivated with formaldehyde, is bound to the solid phase. The specific IgM and IgG-class antibodies bind to the antigen through incubation with diluted human serum. After washings to eliminate the proteins which have not reacted, incubation with the conjugate (anti-human IgM monoclonal antibodies conjugated with horseradish peroxidase) is performed. The conjugate which has not been bound is eliminated and the peroxidase substrate added.

The blue color which develops is proportional to the concentration of specific antibodies present in the test serum.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: the used devices containing materials of human origin must be handled as infectious material and disposed of according to the full requirements of this issue.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instrument.
4. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol

b) The substrate is acid

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minutes exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument and do not use devices with defective labels.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument.
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed samples or samples presenting microbial contamination may all constitute sources of error.
11. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
12. Resuspend the freeze-dried calibrator and control serum homogeneously according to the volumes reported on the labels.
13. Pipette the test serum (50 µl) in well 1 of the device (see figure).
14. Do not use the device after the expiry date.

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 81093).


The kit is sufficient for 12 tests (REF 81093/12).

DD DEVICES

6 packages each containing 6 devices (REF 81093).

2 packages each containing 6 devices (REF 81093/12).

Description of device:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with *Legionella pneumophila* serotype 1 antigen

Position 5: MICROPLATE WELL

Uncoated

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing anti human IgG antibody, phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgM labelled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 mL

Content: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin. Liquid ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Content: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin. Liquid ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER (REF 83606)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 – 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution

- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples cannot be used. The test cannot be applied to plasma.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 6, Analytical Warning, points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (ratio between the OD value of the sample and that of the Cut-off).

The test serum can be interpreted as follows:

POSITIVE: when the ratio is > 1.1

NEGATIVE: when the ratio is < 0.9.

DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

In the case of a doubtful/equivocal result, repeat the test. If the test remains doubtful/equivocal, collect a new sample after 1-2 weeks.

11. LIMITATIONS

The test results should be interpreted with caution and used in conjunction with information available from the clinical evaluation and other diagnostic procedures. Serum samples from patients that are in an early or late stage of this disease may give repeatedly negative results close to the base level of the cut off. In such a case it is recommendable to check the result by doing the test with a second sample obtained 7-10 days later than the first one.

12. ANALYTICAL SPECIFICITY

The following 47 samples containing potentially interfering substances were tested:

Rheumatoid factor (n=6)
Bilirubin high titer (n=3)
Triglycerides (n=12)
C-Reactive Protein (n=3)
Hemolyzed (n=23)

The following 27 sera, positive to other diseases, were tested:

Mycoplasma pneumoniae (n=11)
Chlamydia pneumoniae (n=7)
Helicobacter pylori (n=1)
Clostridium Tetani (n=2)
TBC (n=2)
Yersinia enterocolitica (n=2)
Epstein-Barr (n=2)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

The diagnostics sensitivity and specificity were tested analyzing 130 samples with DIESSE kit as well as with another ELISA commercial kit (ELISA Viracell). The discordant samples were tested using Focus IFA test.

Below are the schematized data of the trial:

		Reference		
		+	-	Total
Diesse	+	18	5	23
	-	0	107	107
	Total	18	112	130

Diagnostic Sensitivity: 100.00 % CI_{95%}: 82.41 – 100.00

Diagnostic Specificity: 95.54 % CI_{95%}: 92.42 – 99.08

14. PRECISION

WITHIN-RUN PRECISION

	Results (Index)			Mean	SD	CV%
LGM 1	2.8	3.0	3.1	2.9	0.14	4.8
LGM 2	1.1	1.0	0.8	1.0	0.15	15.0
LGM 3	0.5	0.6	0.6	0.6	0.06	10.0

BETWEEN-RUNS PRECISION

	Results (Index)			Mean	SD	CV%
	Day 1	Day 2	Day 3			
LGM 1	2.8	2.7	2.2	2.6	0.36	13.8
	3.0	2.7	2.2			
	3.1	2.7	2.0			
LGM 2	1.1	1.1	0.9	1.0	0.14	14.0
	1.0	1.3	1.0			
	0.8	1.1	1.1			
LGM 3	0.5	0.5	0.4	0.5	0.07	14.0
	0.6	0.6	0.5			
	0.6	0.5	0.6			

15. REFERENCES

1. G.B. Wisdom. 1976. Enzyme-Immunoassay. Clin. Chem. 22: 1243.
2. A.K. Malan, T.B. Martins, T.D. Jaskowski, H. R. Hill, C.M. Litwin. 2003 Comparaison of two Commercial Enzyme-Linked Immunosorbent Assay with an immunofluorescence assay for detection of *Legionella pneumophila* Types 1 to 6. J. Clin. Microb. 41: 3060-3063.
3. B. Neumeister, M. Susa, B. Nowak, E. Straube, G. Ruckdeschel, J. Hacher, R. Marre. 1995 Enzyme immunoassay for detection of antibodies against isolated flagella of *Legionella pneumophila*. Eur. J. Clin. Microbiol. Infect. Dis. 14 (9): 764-7.
4. N. Barka, J.P. Tomasi, S. Stadtsbaeder 1986 ELISA using whole *Legionella pneumophila* cell as antigen J.Immunol. Methods 93: 77-81
5. M. S. Lever. 1993. Production and characterisation of a *Legionella pneumophila* specific monoclonal antibody. FEMS Microbiol. Lett. 1: 5-9.
6. L.B. Reller, M.P. Weinstein 2003 Diagnosis of Legionella Infection Clin.Infect. Dis. 36:64-69

16. GLOSSARY OF LABELING SYMBOLS

	Date of manufacture
	Use By
	Do not reuse.
	Caution, consult accompanying documents.
	Manufacturer
	Contains sufficient for <n> tests
	Temperature limitation
	Consult instructions for use
	Biological risks
	Catalog number
	<i>In vitro</i> diagnostic medical device
	Batch code

Manufactured by
DIESSE Diagnostica Senese S.p.A.
 Strada dei Laghi 39
 53035 Monteriggioni (SIENA)
 Italy

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15. POUŽITÁ LITERATURA

1. G.B. Wisdom. 1976. Enzyme-Immunoassay. Clin. Chem. 22: 1243.
2. A.K. Malan, T.B. Martins, T.D. Jaskowski, H. R. Hill, C.M. Litwin. 2003 Comparaison of two Commercial Enzyme-Linked Immunosorbent Assay with an immunofluorescence assay for detection of *Legionella pneumophila* Types 1 to 6. J. Clin. Microb. 41: 3060-3063.
3. B. Neumeister, M. Susa, B. Nowak, E. Straube, G. Ruckdeschel, J. Hacher, R. Marre. 1995 Enzyme immunoassay for detection of antibodies against isolated flagella of *Legionella pneumophila*. Eur. J. Clin. Microbiol. Infect. Dis. 14 (9): 764-7.
4. N. Barka, J.P. Tomasi, S. Stadtsbaeder 1986 ELISA using whole *Legionella pneumophila* cell as antigen J. Immunol. Methods 93: 77-81
5. M. S. Lever. 1993. Production and characterisation of a *Legionella pneumophila* specific monoclonal antibody. FEMS Microbiol. Lett. 1: 5-9.
6. L.B. Reller, M.P. Weinstein 2003 Diagnosis of Legionella Infection Clin. Infect. Dis. 36:64-69

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15. LITERATUR

- G.B. Wisdom. 1976. Enzyme-Immunoassay. Clin. Chem. 22: 1243.
- A.K. Malan, T.B. Martins, T.D. Jaskowski, H. R. Hill, C.M. Litwin. 2003 Comparaison of two Commercial Enzyme-Linked Immunosorbent Assay with an immunofluorescence assay for detection of *Legionella pneumophila* Types 1 to 6. J. Clin. Microb. 41: 3060-3063.
- B. Neumeister, M. Susa, B. Nowak, E. Straube, G. Ruckdeschel, J. Hacher, R. Marre. 1995 Enzyme immunoassay for detection of antibodies against isolated flagella of *Legionella pneumophila*. Eur. J. Clin. Microbiol. Infect. Dis. 14 (9) : 764-7.
- N. Barka, J.P. Tomasi, S. Stadtsbaeder 1986 ELISA using whole *Legionella pneumophila* cell as antigen J. Immunol. Methods 93: 77-81
- M. S. Lever. 1993. Production and characterisation of a *Legionella pneumophila* specific monoclonal antibody. FEMS Microbiol. Lett. 1 : 5-9.
- L.B. Reller, M.P. Weinstein 2003 Diagnosis of Legionella Infection Clin. Infect. Dis. 36:64-69

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- G.B. Wisdom. 1976. Enzyme-Immunoassay. Clin. Chem. 22: 1243.
- A.K. Malan, T.B. Martins, T.D. Jaskowski, H. R. Hill, C.M. Litwin. 2003 Comparaison of two Commercial Enzyme-Linked Immunosorbent Assay with an immunofluorescence assay for detection of Legionella pneumophila Types 1 to 6. J. Clin. Microb. 41: 3060-3063.
- B. Neumeister, M. Susa, B. Nowak, E. Straube, G. Ruckdeschel, J. Hacher, R. Marre. 1995 Enzyme immunoassay for detection of antibodies against isolated flagella of Legionella pneumophila. Eur. J. Clin. Microbiol. Infect. Dis. 14 (9) : 764-7.
- N. Barka, J.P. Tomasi, S. Stadtsbaeder 1986 ELISA using whole Legionella pneumophila cell as antigen J. Immunol. Methods 93: 77-81
- M. S. Lever. 1993. Production and characterisation of a Legionella pneumophila specific monoclonal antibody. FEMS Microbiol. Lett. 1 : 5-9.
- L.B. Reller, M.P. Weinstein 2003 Diagnosis of Legionella Infection Clin. Infect. Dis. 36:64-69

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15. BIBLIOGRAFÍA

1. G.B. Wisdom. 1976. Enzyme-Immunoassay. Clin. Chem. 22: 1243.
2. A.K. Malan, T.B. Martins, T.D. Jaskowski, H. R. Hill, C.M. Litwin. 2003 Comparaison of two Commercial Enzyme-Linked Immunosorbent Assay with an immunofluorescence assay for detection of *Legionella pneumophila* Types 1 to 6. J. Clin. Microb. 41: 3060-3063.
3. B. Neumeister, M. Susa, B. Nowak, E. Straube, G. Ruckdeschel, J. Hacher, R. Marre. 1995 Enzyme immunoassay for detection of antibodies against isolated flagella of *Legionella pneumophila*. Eur. J. Clin. Microbiol. Infect. Dis. 14 (9) : 764-7.
4. N. Barka, J.P. Tomasi, S. Stadtsbaeder 1986 ELISA using whole *Legionella pneumophila* cell as antigen J.Immunol. Methods 93: 77-81
5. M. S. Lever. 1993. Production and characterisation of a *Legionella pneumophila* specific monoclonal antibody. FEMS Microbiol. Lett. 1 : 5-9.
6. L.B. Reller, M.P. Weinstein 2003 Diagnosis of Legionella Infection Clin.Infect. Dis. 36:64-69

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15. BIBLIOGRAPHIE

- G.B. Wisdom. 1976. Enzyme-Immunoassay. Clin. Chem. 22: 1243.
- A.K. Malan, T.B. Martins, T.D. Jaskowski, H. R. Hill, C.M. Litwin. 2003 Comparaison of two Commercial Enzyme-Linked Immunosorbent Assay with an immunofluorescence assay for detection of *Legionella pneumophila* Types 1 to 6. J. Clin. Microb. 41: 3060-3063.
- B. Neumeister, M. Susa, B. Nowak, E. Straube, G. Ruckdeschel, J. Hacher, R. Marre. 1995 Enzyme immunoassay for detection of antibodies against isolated flagella of *Legionella pneumophila*. Eur. J. Clin. Microbiol. Infect. Dis. 14 (9) : 764-7.
- N. Barka, J.P. Tomasi, S. Stadtsbaeder 1986 ELISA using whole *Legionella pneumophila* cell as antigen J. Immunol. Methods 93: 77-81
- M. S. Lever. 1993. Production and characterisation of a *Legionella pneumophila* specific monoclonal antibody. FEMS Microbiol. Lett. 1 : 5-9.
- L.B. Reller, M.P. Weinstein 2003 Diagnosis of *Legionella* Infection Clin. Infect. Dis. 36:64-69

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