CHORUS

TOSCANA VIRUS IgG

REF 81067

REF 81067/12

Prodotto da/Manufactured by/Fabricado por:
DIESSE Diagnostica Senese
 Strada dei Laghi 39
 53035 Monteriggioni (SIENA) – Italy

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15. BIBLIOGRAFIA

1. Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). Am. J. Trop. Med. Hyg. 45, 429-434.
2. Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. J. Medical Virology 49, 83-86.
3. Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in Escherichia coli. J. Clin. Microbiol. 36, 3170-3172.
4. Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. J. Clin. Microbiol. 37, 649-652.

INSTRUCTIONS FOR USE

CHORUS TOSCANA VIRUS IgG

(English)

1. INTENDED USE

QUALITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO TOSCANA VIRUS IN HUMAN SERUM, USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS INSTRUMENT

2. INTRODUCTION

The Toscana Virus, belonging to the Phlebovirus genus of the Bunyaviridae family, was isolated for the first time in Italy in 1971. Since then, its presence in the mediterranean area has been established and its role as an etiological agent in human disease has been demonstrated. The disease caused by Toscana Virus is characterized by fever, strong headache and meningitis or meningoencephalitis. The infection is transmitted by hematophagous insects of the Phlebotomus genus and is seasonal with a peak in the summer months, due to the life cycle of the vector. Epidemiological data indicate that in areas where the Toscana Virus is endemic, many, perhaps the majority of cases of aseptic meningitis during the summer period are caused by this virus. Until now, the diagnosis of Toscana virus infection has been mainly performed serologically. Specific antibodies of the IgM class are present during the acute phase, while the presence of specific IgG antibodies in the absence of IgM indicates a previous infection. All the serological methods adopted up to the present time require propagation of the virus in cell cultures; these methods can only be performed in specialized laboratories, while the infection cannot be diagnosed in the majority of serological laboratories. As it has been shown that the viral nucleoprotein is the chief target of the anti-Toscana virus antibodies, this has been expressed by recombinant DNA technology in E. coli and used to set up an immunoenzymatic kit for the detection of specific IgG and IgM, as an aid in the serological study of this virus.

3. PRINCIPLE OF THE METHOD

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay). The IgG are determined after incubation of the test serum with a precoated solid phase sensitized with recombinant nucleoprotein of highly purified Toscana virus. The specific immunoglobulins are bound to the solid phase and are detected through further incubation with a mouse Mab anti-human IgG conjugated to peroxidase.

The complex is determined by addition of the substrate which takes on a blue colour in the presence of peroxidase. The colour which develops is proportional to the concentration of specific antibodies present in the serum sample. The disposable devices contain all the reagents to perform the test when applied on the Chorus instrument.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 81067).


The kit is sufficient for 12 tests (REF 81067/12).

DD DEVICES

6 packages each containing 6 devices (REF 81067).

2 packages each containing 6 devices (REF 81067/12).

Use : **equilibrate a package at room temperature**, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: empty

Position 6: MICROPLATE WELL coated with recombinant nucleoprotein of Toscana Virus

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum

CALIBRATOR CALIBRATOR 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C

6. PRECAUTIONS

ONLY FOR DIAGNOSTIC USE IN VITRO.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol

b) The substrate is acid

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (50 µl) in well 1 of the device (see figure).
12. Do not use the device after the expiry date

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (ratio between the OD value of the test sample and that of the Cut-Off) which can be used as a quantitative measure, as it is proportional to the amount of specific IgG present.

The test serum can be interpreted as follows:

POSITIVE: when the ratio is > 1.2

NEGATIVE: when the ratio is < 0.8 .

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2

In the case of a doubtful/equivocal result, repeat the test. If the test remains doubtful/equivocal, collect a new sample after 1-2 weeks.

11. LIMITATIONS OF THE PROCEDURE

The test result should be used in conjunction with information available from the evaluation of other diagnostic procedures.

12. ANALYTICAL SPECIFICITY

Cross-reactions cannot be excluded with other correlated serotypes of the Phlebovirus genus such as Sicily and Naples.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

52 samples were tested in an experimentation in comparison with another commercial test: 33 were positive, 13 negative. The results are summarized in the following table, which shows 5 false negative and 1 false positive result.

		REFERENCE	
		+	-
DIESSE	+	33	1
	-	5	13

The kit shows a sensitivity of 87% (CI 95%: 72–94) and a specificity of 93% (CI 95%: 69–99).

14. PRECISION

Within-run precision

	No. of replicates	Mean Index	St. Dev.	CV%
Lot no. 017	8	1.8	0.11	6
Lot no. 018	9	2.9	0.26	9
Lot no. 019	9	3.3	0.27	8

Between-run precision and precision between batches

Sample	Mean Index			Mean	St. Dev.	CV%
	Lot 017	Lot 018	Lot 019			
TVG 1	0.7	0.7	0.5	0.6	0.12	20
TVG 2	1.2	1.2	1.0	1.1	0.12	11
TVG 3	2.2	2.0	2.3	2.2	0.15	7

15. REFERENCES

1. Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). Am. J. Trop. Med. Hyg. 45, 429-434.
2. Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. J. Medical Virology 49, 83-86.
3. Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in Escherichia coli. J. Clin. Microbiol. 36, 3170-3172.
4. Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. J. Clin. Microbiol. 37, 649-652.

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15. LITERATUR

1. Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). Am. J. Trop. Med. Hyg. 45, 429-434.
5. Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. J. Medical Virology 49, 83-86.
6. Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in Escherichia coli. J. Clin. Microbiol. 36, 3170-3172.
7. Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. J. Clin. Microbiol. 37, 649-652.

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15. BIBLIOGRAFIA

- Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). Am. J. Trop. Med. Hyg. 45, 429-434.
- Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. J. Medical Virology 49, 83-86.
- Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in Escherichia coli. J. Clin. Microbiol. 36, 3170-3172.
- Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. J. Clin. Microbiol. 37, 649-652.

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15. BIBLIOGRAPHIE

1. Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). *Am. J. Trop. Med. Hyg.* 45, 429-434.
2. Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. *J. Medical Virology* 49, 83-86.
3. Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in Escherichia coli. *J. Clin. Microbiol.* 36, 3170-3172.
4. Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. *J. Clin. Microbiol.* 37, 649-652.

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15. BIBLIOGRAFIA

1. Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). *Am. J. Trop. Med. Hyg.* 45, 429-434.
2. Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. *J. Medical Virology* 49, 83-86.
3. Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in *Escherichia coli*. *J. Clin. Microbiol.* 36, 3170-3172.
4. Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. *J. Clin. Microbiol.* 37, 649-652.

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusqu'e Ημερομηνία λήξης Prazo de validade
	EN ES IT	Do not reuse No reutilizar Non riutilizzare	FR GR PT	Ne pas réutiliser Μην κάνετε επαναληπτική χρήση Não reutilizar
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbricante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

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