CHORUS

25OH Vitamin D Total

DIESSE

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INSTRUCTIONS FOR USE

CHORUS 25OH Vitamin D Total

**For the quantitative determination of 25(OH)
Vitamin D**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the quantitative determination of 25(OH) Vitamin D in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Vitamin D is a fat-soluble vitamin that is involved in calcium-phosphorus metabolism. It exists in two forms: cholecalciferol (Vitamin D3), and ergocalciferol (Vitamin D2). In the bloodstream D2 and D3 are conjugated to carrier proteins, mainly VDBP (Vitamin D-Binding Protein). Both the forms are biologically inert and suffer, in the body, two important hydroxylations: the first in the liver, with formation of calcidiol (25(OH) Vitamin D); the second in the kidney, with the formation of calcitriol (1.25(OH)₂ Vitamin D), the biologically active metabolite.

25(OH) Vitamin D is generally considered the main form of storage and the determination of its serological concentration is considered a reliable indicator of overall Vitamin D status (more than 95% of Vitamin D in the serum is represented by 25(OH)). Vitamin D is essential for good health of bones. A deficiency of Vitamin D results in a defective mineralization of the bone structure which can lead to childhood rickets, osteoporosis and osteomalacia in adults. An excess of Vitamin D, "Hypervitaminosis D" determines pathological effects on the organism.

Vitamin D has also extra-skeletal functions, including the immune modulating and anti-proliferative, has a possible role in chronic inflammatory and neoplastic diseases, diabetes, cardiovascular and autoimmune diseases and in dermatology, where a deficiency of Vitamin D is associated to psoriasis, vitiligo, atopic dermatitis, melanoma and epithelial tumors of the skin.

Several studies have shown that, in Mediterranean countries such as Italy, despite the latitude, the population is among those with the lowest serum levels of Vitamin D in Europe.

3. PRINCIPLE OF THE METHOD

The Chorus 25OH Vitamin D Total device is ready to use for the detection of 25(OH) Vitamin D, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA (Enzyme linked ImmunoSorbent Assay) competition principle.

The monoclonal antibody is bound to the solid phase and 25(OH) Vitamin D, present in the sample, competes with biotinylated Vitamin D.

The higher the concentration of 25(OH) Vitamin D in the sample, the lower the amount of biotinylated Vitamin D that binds.

After washing to eliminate the components which have not reacted, incubation is performed with the conjugate, composed of Streptavidin conjugated to horseradish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The colour which develops is inversely proportional to the concentration of 25 (OH) Vitamin D present in the serum sample.

The competition reaction is performed using an extraction buffer (Incubation buffer), required to release 25 (OH) Vitamin D from VDBP.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in ng/ml.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.


1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.
The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: MICROPLATE WELL with BIOTYNILATED VITAMIN D

Position 6: MICROPLATE WELL

Coated with anti-Vitamin D monoclonal antibodies

Position 5: MICROPLATE WELL

Uncoated

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine and H₂O₂ stabilized in buffer.

Position 3: INCUBATION BUFFER

Contents: Proteic solution containing preservative

Position 2: CONJUGATE

Contents: Streptavidin labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

Contents: Diluted human serum containing Vitamin D and preservative. Freeze-dried.

Use: reconstitute with the volume of deionised water reported on the label, wait at least 5 minutes without shaking; then mix by inversion avoiding froth formation.

CONTROL CONTROL SERUM

Contents: Diluted human serum containing Vitamin D and preservative. Freeze-dried.

Use: reconstitute with the volume of deionised water reported on the label, wait at least 5 minutes without shaking; then mix by inversion avoiding froth formation.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at -20°C
CONTROL SERUM	8 weeks at -20°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for short periods at -20°C. Avoid the freezing of the samples for long periods.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 90 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in ng/ml, calculated on the basis of a lot-dependent graph stored in the instrument.

Attention: do not consider the qualitative report N; D/E; P on the instrument.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range: 4.0-115.0 ng/ml.

13. REFERENCE RANGE

The table below reports the expected values of 25(OH) Vitamin D, as reported in literature.

State	Range (ng/ml)
Lack	< 20
Insufficient	20 – 30
Normal	30 – 100
Toxic	> 100

The values observed in the population, determined by examining 399 sera, were between 4.0 and 46.9 ng/ml, with a median of 14.2 ng/ml.

14. ANALITICAL SPECIFICITY

3 samples were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (44 IU/ml - 220 IU/ml)
Bilirubin (4.5 mg/dl - 45 mg/dl)
Triglycerides (10 mg/dl - 250 mg/dl)
Hemoglobin (5 mg/ml - 30 mg/ml)

The presence in the serum sample of the interfering substances described above (except for hemoglobin) does not affect the test result.

15. CROSS-REACTIONS

7 samples were tested adding Vitamin D2 and Vitamin D3. No significant cross-reactions were found for Vitamin D2 up to 1000 ng/ml and Vitamin D3 up to 125 ng/ml.

16. METHOD COMPARISON

In an experimentation 217 samples have been tested with Diesse kit and with a competitor kit.

Results:

Correlation	r	95%CI
Pearson	0.92	0.89-0.94
Spearman	0.92	0.90-0.94

The correlation between the two methods is very high.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	17.5	5.8	17.3	8.3
2	23.4	7.6	27.2	7.0
3	27.6	8.6	27.9	6.5
4	31.4	3.5	33.4	4.7
5	33.2	3.1	35.4	4.2
6	37.2	4.0	37.5	2.7
7	39.4	3.8	44.9	12.6
8	38.3	1.4	41.7	10.5
9	78.1	3.1	77.9	5.3

Sample	Between lots		Between Instruments	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	12.9	14.4	12.9	3.9
2	22.5	9.4	22.5	3.6
3	21.9	8.4	21.9	2.3
4	26.4	10.3	26.4	4.6
5	30.9	4.0	30.9	2.3
6	32.6	4.5	32.6	2.2
7	36.1	3.4	36.1	2.5
8	35.8	2.7	35.8	4.8
9	73.4	7.1	73.4	9.3

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