CHORUS

SYPHILIS SCREEN RECOMBINANT

REF 81100 (36 tests)

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15. BIBLIOGRAFIA

1. S. J. Norris: Polypeptides of *Treponema pallidum*: progress toward understanding their structural, functional and immunologic roles. Microbiological Reviews, 1993.
2. P.A. Hanff et al. Humoral immune response in human syphilis to polypeptides of *Treponema pallidum*. J. Immunology 129: 1287 (1982).
3. L. Lewis et al. Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
4. H. Young et al., Enzywell recombinant enzyme immunoassay for the serological diagnosis of syphilis. Int. J. STD and AIDS. 11: 288-291 (2000).

INSTRUCTIONS FOR USE

CHORUS SYPHILIS SCREEN RECOMBINANT

REF 81100

(English)

1. INTENDED USE

QUALITATIVE DETERMINATION OF ANTI-TREPONEMA PALLIDUM IMMUNOGLOBULINS IN HUMAN SERUM, USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS INSTRUMENTS

2. INTRODUCTION

The serological diagnosis of syphilis is performed by demonstrating the presence of significant levels of specific *Treponema pallidum* (TP) antibodies in the serum sample.

The reference method used is the FTA-ABS technique but its execution is laborious and the interpretation of the results is not simple; alternative methods have therefore been introduced to simplify the procedure.

The TPHA test is the preferred technique for screening purposes, in as far as it detects specific Ig even at low titers. Unfortunately, the result of the test is determined by subjective interpretation as the test cannot be completely automated. The TPHA test is not sensitive in the primary phase.

The present kit allows screening with the ELISA method. It reveals the presence of specific antibodies of any class. The antigen used is obtained by a recombinant technique.

3. PRINCIPLE OF THE METHOD

The test is based on the “sandwich” principle, a solid phase enzyme-linked immunoassay technique. During incubation, the antibodies present in the sample are bound to the solid phase antigen as well as to the antigen labelled with horse-radish peroxidase (HRP), creating an antigen-antibody-antigen-HRP “sandwich”. The unbound conjugate is eliminated by washing and the bound enzymatic activity is determined by adding a chromogen substrate. The intensity of the colour is proportional to the anti-*Treponema pallidum* antibody concentration in the sample.


The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Use: **equilibrate a package at room temperature**, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: empty

Position 6: MICROPLATE WELL coated with recombinant anti-T. pallidum antigens

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: empty

Position 2: CONJUGATE

Contents: Mixture of T. pallidum recombinant proteins labelled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum

CALIBRATOR CALIBRATOR 1 x 0.175 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin, **liquid ready for use.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER [REF] 83606
- CLEANING SOLUTION 2000 [REF] 83609
- SANITIZING SOLUTION [REF] 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT [REF] 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C

6. PRECAUTIONS

ONLY FOR DIAGNOSTIC USE IN VITRO.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minutes exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Check that the substrate contained in well 4 is colourless. Discard devices which do not have this characteristic.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (50 µl) in well 1 of the device (see figure).
12. Do not use the device after the expiry date

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (ratio between the OD value of the sample and that of the Cut-off).

POSITIVE: when the ratio is > 1.2

NEGATIVE: when the ratio is < 0.8.

DOUBTFUL: for all values between 0.8 and 1.2

In the case of a doubtful result, repeat the test. If the test remains doubtful, collect a new sample after 1-2 weeks.

11. LIMITATIONS OF THE PROCEDURE

The test is not able to discriminate between the presence of IgG and IgM.

Highly positive samples could show a Prozone effect and give false negative results.

In case of negative result and presence of persistent clinical symptoms, it is advised to repeat the test with different methods.

However, a negative result does not exclude the possibility of a *Treponema pallidum* infection.

All the values obtained require a careful interpretation that must consider other indicators relative to the patient. The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALYTICAL SPECIFICITY

Analytical specificity was studied by testing samples containing potentially interfering factors such as Anti-Nuclear EBV antibodies, heterophile antibodies, Borrelia antibodies, Rheumatoid Factor (from 40 to 1080 IU/mL), Triglycerides (≤ 870 mg/dL), bilirubin (≤ 11 mg/dL), hemoglobin (≤ 10 mg/dL). None of these factors caused interference in the test.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

53 serum samples were tested, in comparison with another commercial method. 27 were negative and 26 positive samples were found. The samples were compared with another commercial method. The results are summarized in the following table.

Table 1

	REFERENCE	
	+	-
	26	0
CHORUS	+	0
	-	27

The kit shows a specificity of 100% (CI 95%: 88–100) and a sensitivity of 100% (CI 95%: 87–100).

14. PRECISION

Table 2 Within-run precision

	No. of replicates	Mean Index	St. Dev.	CV%
Lot no. 030	9	1.8	0.11	6
Lot no. 031	9	1.5	0.07	5
Lot no. 034	9	3.3	0.18	5

Table 3 Between-run precision and Precision between batches

Sample	Mean Index			Mean	St. Dev.	CV%
	Lot 030	Lot 031	Lot 034			
TPN 1	0.3	0.3	0.3	0.3	0.00	0
TPN 2	0.8	0.7	0.6	0.7	0.10	14
TPN 3	1.2	1.4	1.1	1.2	0.15	12
TPN 4	4.3	4.3	4.1	4.2	0.12	3

15. REFERENCES

1. S. J. Norris: Polypeptides of Treponema pallidum: progress toward understanding their structural, functional and immunologic roles. Microbiological Reviews, 1993.
2. P.A. Hanff et al. Humoral immune response in human syphilis to polypeptides of Treponema pallidum. J. Immunology 129: 1287 (1982).
3. L. Lewis et al. Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
4. H. Young et al., Enzywell recombinant enzyme immunoassay for the serological diagnosis of syphilis. Int. J. STD and AIDS. 11: 288-291 (2000).

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5. S. J. Norris: Polypeptides of *Treponema pallidum*: progress toward understanding their structural, functional and immunologic roles. Microbiological Reviews, 1993.
6. P.A. Hanff et al. Humoral immune response in human syphilis to polypeptides of *Treponema pallidum*. J. Immunology 129: 1287 (1982).
7. L. Lewis et al. Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
8. H. Young et al., Enzywell recombinant enzyme immunoassay for the serological diagnosis of syphilis. Int. J. STD and AIDS. 11: 288-291 (2000).

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15. BIBLIOGRAPHIE

1. S. J. Norris: Polypeptides of *Treponema pallidum*: progress toward understanding their structural, functional and immunologic roles. Microbiological Reviews, 1993.
2. P.A. Hanff et al. Humoral immune response in human syphilis to polypeptides of *Treponema pallidum*. J. Immunology 129: 1287 (1982).

3. L. Lewis et al. Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. *J. Clin. Microbiol.* 28: 296 (1990).
4. H. Young et al., Enzywell recombinant enzyme immunoassay for the serological diagnosis of syphilis. *Int. J. STD and AIDS.* 11: 288-291 (2000).

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