ENZY-WELL

HELICOBACTER PYLORI IgA

REF 91062 (96 tests)

Prodotto da/Manufactured by/ Hersteller: Κατασκευάζεται
 από/ Fabricado por/ Fabriqué par:
 DIESSE Diagnostica Senese
 Strada dei Laghi, 39
 53035 Monteriggioni (Siena) - Italy

	Capitolo Section Capítulo
Modifiche introdotte nella revisione corrente Changes introduced in the current revision Eingeführte Änderungen in der vorliegenden Fassung Τροποποιήσεις που εισήχθησαν στην τρέχουσα διόρθωση Cambios introducidos en la revisión actual Modifications apportées à la révision courante Alterações introduzidas na revisão atual	4 – 5

INDICE / INDEX / INHALT/ ΠΕΡΙΕΧΟΜΕΝΑ / INDICE / SOMMAIRE / CONTEÚDO

1. UTILIZZAZIONE / INTENDED USE / VERWENDUNGSZWECK/ ΥΠΟΔΕΙΚΝΥΟΜΕΝΗ ΧΡΗΣΗ / INDICACIONES DE USO/ UTILISATION/UTILIZAÇÃO PRETENDIDA
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST / PHYSIOLOGIE/ ΕΙΣΑΓΩΓΗ / RESUMEN Y EXPLICACIÓN DEL TEST / INTERET CLINIQUE / SUMÁRIO E EXPLICAÇÃO DO TESTE
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / TESTPRINZIP/ ΑΡΧΗ ΤΗΣ ΜΕΘΟΔΟΥ / PRINCIPIO DEL MÉTODO / PRINCIPE/ PRINCÍPIOS DO TESTE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / INHALT DES KITS UND VORBEREITUNG DER REAGENZIEN / ΣΥΝΘΕΣΗ ΤΟΥ ΚΙΤ ΚΑΙ ΠΡΟΕΤΟΙΜΑΣΙΑ ΤΩΝ ΑΝΤΙΑΠΑΣΤΗΡΙΩΝ / COMPONENTES DEL KIT Y PREPARACIÓN DEL REACTIVO / CONTENU DU KIT ET PREPARATION DES REACTIFS / CONTEÚDA DO KIT E PREPARAÇÃO DO REAGENTE /
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / LAGERUNG UND STABILITÄT DER REAGENZIEN / ΤΡΟΠΟΣ ΑΠΟΘΗΚΕΥΣΗΣ ΚΑΙ ΣΤΑΘΕΡΟΤΗΤΑ ΤΩΝ ΑΝΤΙΑΠΑΣΤΗΡΙΩΝ/ CONSERVACIÓN Y ESTABILIDAD DE LOS REACTIVOS / MODALITES DE CONSERVATION ET STABILITE DES REACTIFS / ARMAZENAMENTO E ESTABILIDADE
6. PRECAUZIONI / PRECAUTIONS / HINWEISE UND VORSICHTSMASSNAHMEN/ ΠΡΟΦΥΛΑΞΕΙΣ / PRECAUCIONES DE USO / PRECAUTIONS D'UTILISATION / PRECAUÇÕES
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / ART UND LAGERUNG DER PROBEN / ΕΙΔΟΣ ΔΕΙΓΜΑΤΟΣ ΚΑΙ ΑΠΟΘΗΚΕΥΣΗ / TIPO DE MUESTRA Y CONSERVACION / TYPE D'ECHANTILLON ET CONSERVATION / TIPO E ARMAZENAMENTO DAS AMOSTRAS
8. PROCEDIMENTO / TEST PROCEDURE / TESTDURCHFÜHRUNG / ΔΙΑΔΙΚΑΣΙΑ / PROCEDIMIENTO DEL TEST / PROCEDURE / PROCEDIMENTO DE TESTE /
9. SCHEMA DEL SAGGIO / SCHEME OF TEST PROCEDURE / TESTSCHEMA /ΔΙΑΓΡΑΜΜΑ ΤΗΣ ΑΝΑΛΥΣΗΣ/ ESQUEMA DEL PROCEDIMIENTO DEL TEST / SCHEMA D'ANALYSE/ ESQUEMA DE PROCEDIMENTO DE TESTE /
10. VALIDAZIONE DEL TEST / TEST VALIDATION / TESTVALIDIERUNG / ΕΓΚΥΠΟΤΗΤΑ ΤΟΥ ΤΕΣΤ VALIDACIÓN DEL TEST / VALIDATION DU TEST / VALIDAÇÃO DO TESTE /
11. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / INTERPRETATION DER ERGEBNISSE / ΕΡΜΗΝΕΙΑ ΤΩΝ ΑΠΟΤΕΛΕΣΜΑΤΩΝ / INTERPRETACIÓN DE RESULTADOS / INTERPRETATION DU TEST / INTERPRETAÇÃO DOS RESULTADOS
12. LIMITAZIONI DELLA PROCEDURA / LIMITATIONS OF THE PROCEDURE / GRENZEN DES VERFAHRENS / ΠΕΡΙΟΡΙΣΜΟΙ ΣΤΗΝ ΜΕΘΟΔΟ / LIMITACIONES DEL TEST / LIMITES DU TEST/ LIMITAÇÕES DO PROCEDIMENTO /
13. SPECIFICITA' ANALITICA / ANALYTICAL SPECIFICITY / ANALYTISCHE SPEZIFITÄT /ΑΝΑΛΥΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ/ ESPECIFICIDAD ANALÍTICA / SPECIFICITE ANALYTIQUE / ESPECIFICIDADE ANALÍTICA /
14. SENSIBILITA' E SPECIFICITA' DIAGNOSTICA / DIAGNOSTIC SENSITIVITY AND SPECIFICITY / DIAGNOSTISCHE SENSITIVITÄT UND SPEZIFITÄT /ΕΥΑΙΣΘΗΣΙΑ ΚΑΙ ΔΙΑΓΝΩΣΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ / SENSIBILIDAD Y ESPECIFICIDAD DIAGNOSTICAS / SENSIBLITE ET SPECIFICITE DIAGNOSTIQUES / ESPECIFICIDADE E SENSIBILIDADE DO DIAGNÓSTICO
15. PRECISIONE / PRECISION / PRÄZISION / ΑΚΡΙΒΕΙΑ / PRECISIÓN / PRECISION / PRECISÃO
16. GUIDA DEI PROBLEMI DI UTILIZZO / “TROUBLE SHOOTING” / “TROUBLE SHOOTING” (PROBLEMLÖSUNGEN) /ΟΔΗΓΟΣ ΕΠΙΛΥΣΗΣ ΠΡΟΒΛΗΜΑΤΩΝ /GUIA DE RESOLUCION DE PROBLEMAS / RESOLUTION DES PROBLEMES / RESOLUÇÃO DE PROBLEMAS
17. BIBLIOGRAFIA / REFERENCES / LITERATUR / ΒΙΒΛΙΟΓΡΑΦΙΑ / BIBLIOGRAFÍA / BIBLIOGRAPHIE / REFERÊNCIAS

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ENZY-WELL HELICOBACTER PYLORI IgA

REF 91062

(Italiano)

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(English)

1. INTENDED USE

IMMUNOENZYMATIC KIT FOR THE QUALITATIVE DETERMINATION OF IgA-CLASS ANTIBODIES TO HELICOBACTER PYLORI IN HUMAN SERUM

2. SUMMARY AND EXPLANATION OF THE TEST

In 1983, Warren and Marshall (1) identified *Helicobacter pylori*, a new gram-negative bacterial pathogen, in patients suffering from gastritis, and this finding led to studies on the relationship between bacterial infection and chronic gastric disease (2). The pathogen has been shown to be It has been demonstrated that in patients with gastritis, eradication of the bacteria led to healing of the anatomical lesion (4).

Diagnostic procedures for the detection of the organism generally involve invasive (gastroscopic) techniques for sample collection.

However, a specific immune response is seen in infected patients. The serological test thus represents a useful alternative to the invasive bioptic technique. IgG levels rise with infection and remain constantly high until the infection is eliminated. The efficacy of antimicrobial therapy can therefore be monitored via changes in specific IgG antibody. The determination of the IgA antibody levels is complementary to that of the IgG. As the IgA levels decrease more rapidly than the IgG in some patients undergoing treatment, this parameter can be useful in the follow-up of patients.

For the same reason, a low IgA level in IgG-positive patients may indicate a past infection.

For unknown reasons, about 2% of serum samples are positive only for IgA (5,6).

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay) (5).

The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgA monoclonal antibodies labelled with peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The density of the colour which develops is proportional to the concentration of specific antibodies present in the serum sample. When the enzymatic reaction is interrupted by the addition of a sulfuric acid solution, the resulting yellow colour can be easily read on an ELISA microplate reader.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.

- **Bring reagents to room temperature before use.**

MT PLATE MICROPLATE 12 x 8 wells coated with *Helicobacter pylori*.

Use: open the package at the opposite end from the code (P, followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONTROL + POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum, at known concentration of anti-H. Pylori IgA, in Phosphate buffer 0.01 mol/L containing BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONTROL CUT-OFF CUT-OFF CONTROL 1 x 2 mL

Contents: Diluted human serum, at known concentration of anti-H. Pylori IgA, in Phosphate buffer 0.01 mol/L containing BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONJ CONJUGATE 1x16 mL

Contents: anti-human IgA monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

IgA CONTROL - IgA NEGATIVE CONTROL (PF93910) 1 x 1.6 mL

Contents: Human serum in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

IgA SAMP DIL 50x IgA DILUENT 50X (PF30267) 1x4.5 mL For dilution of serum samples.

Contents: Proteic solution with Phenol 0.05% and Bronidox 0.02%.

Preparation: Dilute the required volume 1:50 in the washing buffer to obtain the diluent ready for use.

SUBS TMB SUBSTRATE (PF93619) 1x12 mL Ready for use. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader, wavelength 450 nm or at 450/620 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes ranging between 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	6 weeks at 2/8°C, polythene bag
Controls	6 weeks at 2/8°C
Conjugate	6 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15/30°C; store in the dark
IgA Sample Diluent	1 week at 15/30°, 8 weeks at 2/8°C
Wash Buffer	ready for use, 2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS

For in vitro diagnostic use only.

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be carefully

mixed before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided.

The test is not applicable to human plasma.

8. TEST PROCEDURE

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).
- Prepare the required diluent for samples by adding 1 part of IgA Diluent 50x to 49 parts of the diluted washing buffer (example: 2 mL + 98 mL of diluted washing buffer).

Leave one well empty for the blank (only the substrate will be used). Place UNDILUTED controls in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control. Dilute the samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended).

The wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µL), add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed (100 µL/well).

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or at 450/620 nm within 30 min.

9. SCHEME OF TEST PROCEDURE

STEP 1	Place 100 µL of diluted sample//control sera in the wells of the strips.
--------	--

Incubate for 45 min. at 37°C

Wash 4 times (300 µL)

STEP 2	Add 100 µL of conjugate to each well
--------	--------------------------------------

Incubate for 45 min. at 37°C

Wash 4 times (300 µL)

STEP 3	Add 100 µL of Substrate to each well
--------	--------------------------------------

Incubate for 15 min. at R.T.

STEP 4	Add 100 µL of Stop Solution
--------	-----------------------------

Read absorbance at 450 nm within 30 min

10. TEST VALIDATION

Subtract the value of the blank (≤ 0.150) from all the other readings. The O.D. values of the control Cut-off serum when tested in triplicate must be within 25% of the mean value. Disregard any abnormal value and recalculate the mean. The Positive control must have an O.D. at least 1.5 times that of the Cut-Off serum. The ratio between Negative Control and Cut-off must be ≤ 0.6 . The O.D. of the cut-off must be ≥ 0.2 at 450 nm and ≥ 0.16 at 450/620 nm.

11. INTERPRETATION OF THE RESULTS

Calculate the ratio between the average O.D. value of the sample and that of the Cut-off (INDEX). The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the Cut-Off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE PROCEDURE

Highly lipemic, icteric, hemolysed or heat-inactivated sera may cause false results and should not be used.

As with all serological tests, the results should be interpreted in conjunction with other clinical and diagnostic data.

13. ANALYTICAL SPECIFICITY

The analytical specificity of the Enzywell H. Pylori IgA kit was evaluated in an inhibition test. *Campylobacter coli*, *Campylobacter fetus*, *Campylobacter jejuni* and *E. coli* were added to patients' serum samples. No differences were observed

between the reference samples and those containing the antigens reported above. It can therefore be concluded that the test is specific for anti-Helicobacter pylori IgA.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In a clinical evaluation, 171 samples were tested after characterization as positive or negative by endoscopic examination: 110 proved positive for the presence of H. Pylori in the gastric biopsy.

		RIFERIMENTO	
		+	-
DIESSE	+	95	2
	-	15	59

The sensitivity and specificity were 86.4% and 96.7% respectively.

15. PRECISION

Within run Precision on three different batches

Cut Off $n=16$	Batch. N. 31	Batch. N. 32	Batch. N. 33
O.D.	0.267	0.472	0.494
CV%	10	14	12

Between run Precision

Campione	AU/ml				
	I run	II run	III run	Media	CV%
PYA1	0.2	0.3	0.3	0.3	17
PYA2	1.2	1.0	1.4	1.2	18
PYA3	2.5	2.0	2.8	2.4	18

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert, paragraph 4, for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow). Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

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