CHORUS

**Toxoplasma IgA
Capture**

DIESSE

REF 81044

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14. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. P. Herbrink, A. van Loon et al.: Interlaboratory evaluation of indirect enzyme-linked immunosorbent assay, antibody capture enzyme-linked immunosorbent assay, and immunoblotting for detection of immunoglobulin M antibodies to *Toxoplasma gondii*. J. Clin. Microbiol. 25: 100 (1987).
3. F. Wielaard, H. van Gruijthuisen et al.: Diagnosis of acute toxoplasmosis by an enzyme immunoassay for specific immunoglobulin M antibodies. J. Clin. Microbiol. 17: 981 (1983).
4. J. Fung, A. Clogston et al.: Serologic diagnosis of toxoplasmosis with emphasis on the detection of toxoplasma-specific immunoglobulin M antibodies. Am. J. Clin. Pathol. 83: 196 (1985).
5. Remington J.S., Miller M.J., Brownlee I.: IgM antibodies in acute toxoplasmosis. I. Diagnostic significance in congenital cases and a method for their rapid demonstration. Pediatrics 41: 1082 (1968).
6. Pomares C., Montoya J.G. Laboratory Diagnosis of Congenital Toxoplasmosis. Journal of Clinical Microbiology Volume 54, Number 50 (2016): 2448-2454
7. Remington J.S., Thulliez P. and Montoya J.G. Recent Developments for Diagnosis of Toxoplasmosis JOURNAL OF CLINICAL MICROBIOLOGY (2004): 941-945
8. Bessières M.H., Roques C., Berrebi A., Barre V., Cazaux M., Séguéla J.P. IgA antibody response during acquired and congenital toxoplasmosis. J Clin Pathol (1992); 45: 605-608

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INSTRUCTIONS FOR USE

CHORUS Toxoplasma IgA Capture

**For the qualitative determination of Toxoplasma
IgA gondii antibodies**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic capture method for the qualitative determination of IgA class antibodies against Toxoplasma gondii in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Toxoplasma gondii is an ubiquitous protozoa which causes infection in all mammalian species. In man, serological studies have shown that a high percentage of the adult population has been infected by this parasite. The disease is particularly dangerous in immunodepressed patients and during pregnancy, as it is transmitted from mother to foetus.

For this reason, it is important to determine the antibody titer before the beginning of pregnancy; those subjects who do not present immunity should be followed up with frequent analyses to check for serum conversion.

IgA antibodies can persist up to 6-9 months from the beginning of the infection.

3. PRINCIPLE OF THE METHOD

The Chorus Toxoplasma IgA Capture device is ready to use for the detection of IgA antibodies against Toxoplasma gondii, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay).

Monoclonal anti-human IgA antibodies are bound to the solid phase. The IgA immunoglobulins bind to anti-IgA antibodies through incubation with diluted sample.

After washing to eliminate the proteins which have not reacted, incubation is performed with the antigen bound to specific anti-Toxoplasma gondii monoclonal antibodies, conjugated to horseradish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The enzymatic reaction is subsequently blocked by adding the Stop Solution which makes the solution turn into yellow. The colour which develops is proportional to the concentration of specific antibodies present in the examined serum.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

BEFORE STARTING THE TEST, MAKE SURE THAT THE TEMPERATURE OF THE STRIP, VISIBLE ON THE INSTRUMENT DISPLAY (PAGE CHECK), HAS REACHED 30°C ± 2°C.

Before use bring the package containing the devices to room temperature (18-30°C) for 30 minutes at least.

Use within 60 minutes.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)

5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use strongly hemolyzed, lipemic, icteric samples, serum not completely coagulated or samples presenting microbial contamination.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**


5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 81044).

DD DEVICES

6 packages each containing 6 devices (REF 81044).

Description:

Position 9: Space for application of bar code label

Position 8: FREEZE-DRIED ANTIGEN

Contents: purified and inactivated Toxoplasma gondii

Position 7: MICROPLATE WELL

Coated with anti-human IgA monoclonal antibodies

Position 6: Uncoated empty MICROPLATE WELL

Position 5: STOP SOLUTION

0.3 mol/L H₂SO₄ solution, ready for use.

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution in phosphate buffer containing sodium azide 0.09% and dye.

Position 2: CONJUGATE

Contents: anti-Toxoplasma gondii monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.1%.

Position 1: EMPTY WELL

In which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in

the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR

1 x 0.450 mL

Contents: Human serum containing anti-Toxoplasma gondii IgA antibodies, diluted in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%. Liquid, ready for use

CONTROL +

POSITIVE CONTROL

1 x 0.900 mL

Contents: Human serum containing anti-Toxoplasma gondii IgA antibodies, diluted in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%. Liquid, ready for use

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER [REF] 83606
- CLEANING SOLUTION 2000 [REF] 83609
- SANITIZING SOLUTION [REF] 83604 - 83608
- Chorus/Chorus TRIO Instrument [REF] 81000-81200
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 7 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. **BEFORE STARTING THE TEST, MAKE SURE THAT THE TEMPERATURE OF THE STRIP, VISIBLE ON THE INSTRUMENT DISPLAY (PAGE CHECK), HAS REACHED 30°C ± 2°C.**
2. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.

3. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
4. Dispense the following volumes in well no. 1 of each device:

SAMPLE	50 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

At each change of batch, use a device for the calibrator.

5. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.1
 NEGATIVE: when the result is < 0.9
 DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient. The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. METHOD COMPARISON

In an experimentation 103 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	23	1	24
	-	2	77	79
	Total	25	78	103

Percent Positive Agreement (~Diagnostic Sensitivity):

92.0% CI_{95%}: 74.9-97.7

Percent Negative Agreement: (~Diagnostic Specificity):

98.7% CI_{95%}: 93.0-99.8

The agreement between the two methods is excellent with a Cohen's Kappa of 0.91.

13. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (Index)	CV%	Mean (Index)	CV%
1	2.4	9.6	2.5	4.0
2	1.8	11.1	1.8	11.7
3	0.4	10.0	0.4	-

Sample	Precision between batches		Precision between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	2.5	8.8	2.4	11.3
2	1.9	10.0	1.9	14.2
3	0.5	10.0	0.4	12.5

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- G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- P. Herbrink, A. van Loon et al.: Interlaboratory evaluation of indirect enzyme-linked immunosorbent assay, antibody capture enzyme-linked immunosorbent assay, and immunoblotting for detection of immunoglobulin M antibodies to Toxoplasma gondii. J. Clin. Microbiol. 25: 100 (1987).
- F. Wielaard, H. van Gruijthuisen et al.: Diagnosis of acute toxoplasmosis by an enzyme immunoassay for specific immunoglobulin M antibodies. J. Clin. Microbiol. 17: 981 (1983).
- J. Fung, A. Clogston et al.: Serologic diagnosis of toxoplasmosis with emphasis on the detection of toxoplasma-specific immunoglobulin M antibodies. Am. J. Clin. Pathol. 83: 196 (1985).
- Remington J.S., Miller M.J., Brownlee I.: IgM antibodies in acute toxoplasmosis. I. Diagnostic significance in congenital cases and a method for their rapid demonstration. Pediatrics 41: 1082 (1968).
- Pomares C., Montoya J.G. Laboratory Diagnosis of Congenital Toxoplasmosis. Journal of Clinical Microbiology Volume 54, Number 50 (2016): 2448-2454
- Remington J.S., Thulliez P. and Montoya J.G. Recent Developments for Diagnosis of Toxoplasmosis JOURNAL OF CLINICAL MICROBIOLOGY (2004): 941-945
- Bessières M.H., Roques C., Berberi A., Barre V., Cazaux M., Séguéla J.P. IgA antibody response during acquired and congenital toxoplasmosis. J Clin Pathol (1992); 45: 605-608

DIESE Diagnostica Senese
 S.p.A.
 Strada dei Laghi 39
 53035 Monteriggioni (Siena)
 Italy

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22. LITERATUR

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976)
2. P. Herbrink, A. van Loon et al.: Interlaboratory evaluation of indirect enzyme-linked immunosorbent assay, antibody capture enzyme-linked immunosorbent assay, and immunoblotting for detection of immunoglobulin M antibodies to *Toxoplasma gondii*. J. Clin. Microbiol. 25: 100 (1987).
3. F. Wielaard, H. van Gruijthuisen et al.: Diagnosis of acute toxoplasmosis by an enzyme immunoassay for specific immunoglobulin M antibodies. J. Clin. Microbiol. 17: 981 (1983).
4. J. Fung, A. Clogston et al.: Serologic diagnosis of toxoplasmosis with emphasis on the detection of toxoplasma-specific immunoglobulin M antibodies. Am. J. Clin. Pathol. 83: 196 (1985).
5. Remington J.S., Miller M.J., Brownlee I.: IgM antibodies in acute toxoplasmosis. I. Diagnostic significance in congenital cases and a method for their rapid demonstration. Pediatrics 41: 1082 (1968).
6. Pomares C., Montoya J.G. Laboratory Diagnosis of Congenital Toxoplasmosis. Journal of Clinical Microbiology Volume 54, Number 50 (2016): 2448-2454
7. Remington J.S., Thulliez P. and Montoya J.G. Recent Developments for Diagnosis of Toxoplasmosis JOURNAL OF CLINICAL MICROBIOLOGY (2004): 941-945
8. Bessi eres M.H., Roques C., Berrebi A., Barre V., Cazaux M., S egu ela J.P. IgA antibody response during acquired and congenital toxoplasmosis. J Clin Pathol (1992); 45: 605-608

DIESE Diagnostica Senese
S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (Siena)
Italy

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1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. P. Herbrink, A. van Loon et al.: Interlaboratory evaluation of indirect enzyme-linked immunosorbent assay, antibody capture enzyme-linked immunosorbent assay, and immunoblotting for detection of immunoglobulin M antibodies to *Toxoplasma gondii*. J. Clin. Microbiol. 25: 100 (1987).
3. F. Wielaard, H. van Gruijthuijsen et al.: Diagnosis of acute toxoplasmosis by an enzyme immunoassay for specific immunoglobulin M antibodies. J. Clin. Microbiol. 17: 981 (1983).
4. J. Fung, A. Clogston et al.: Serologic diagnosis of toxoplasmosis with emphasis on the detection of toxoplasma-specific immunoglobulin M antibodies. Am. J. Clin. Pathol. 83: 196 (1985).
5. Remington J.S., Miller M.J., Brownlee I.: IgM antibodies in acute toxoplasmosis. I. Diagnostic significance in congenital cases and a method for their rapid demonstration. Pediatrics 41: 1082 (1968).
6. Pomares C., Montoya J.G. Laboratory Diagnosis of Congenital Toxoplasmosis. Journal of Clinical Microbiology Volume 54, Number 50 (2016): 2448-2454
7. Remington J.S., Thulliez P. and Montoya J.G. Recent Developments for Diagnosis of Toxoplasmosis JOURNAL OF CLINICAL MICROBIOLOGY (2004): 941-945
8. Bessi eres M.H., Roques C., Berrebi A., Barre V., Cazaux M., S egu ela J.P. IgA antibody response during acquired and congenital toxoplasmosis. J Clin Pathol (1992); 45: 605-608

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- G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- P. Herbrink, A. van Loon et al.: Interlaboratory evaluation of indirect enzyme-linked immunosorbent assay, antibody capture enzyme-linked immunosorbent assay, and immunoblotting for detection of immunoglobulin M antibodies to *Toxoplasma gondii*. J. Clin. Microbiol. 25: 100 (1987).
- F. Welaard, H. van Gruijthuisen et al.: Diagnosis of acute toxoplasmosis by an enzyme immunoassay for specific immunoglobulin M antibodies. J. Clin. Microbiol. 17: 981 (1983).
- J. Fung, A. Clogston et al.: Serologic diagnosis of toxoplasmosis with emphasis on the detection of toxoplasma-specific immunoglobulin M antibodies. Am. J. Clin. Pathol. 83: 196 (1985).
- Remington J.S., Miller M.J., Brownlee I.: IgM antibodies in acute toxoplasmosis. I. Diagnostic significance in congenital cases and a method for their rapid demonstration. Pediatrics 41: 1082 (1968).
- Pomares C., Montoya J.G. Laboratory Diagnosis of Congenital Toxoplasmosis. Journal of Clinical Microbiology Volume 54, Number 50 (2016): 2448-2454
- Remington J.S., Thulliez P. and Montoya J.G. Recent Developments for Diagnosis of Toxoplasmosis JOURNAL OF CLINICAL MICROBIOLOGY (2004): 941-945

8. Bessi res M.H., Roques C., Berrebi A., Barre V., Cazaux M., S gu la J.P. IgA antibody response during acquired and congenital toxoplasmosis. J Clin Pathol (1992); 45: 605-608

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- toxoplasma-specific immunoglobulin M antibodies. Am. J. Clin. Pathol. 83: 196 (1985).
21. Remington J.S., Miller M.J., Brownlee I.: IgM antibodies in acute toxoplasmosis. I. Diagnostic significance in congenital cases and a method for their rapid demonstration. Pediatrics 41: 1082 (1968).
 22. Pomares C., Montoya J.G. Laboratory Diagnosis of Congenital Toxoplasmosis. Journal of Clinical Microbiology Volume 54, Number 50 (2016): 2448-2454
 23. Remington J.S., Thulliez P. and Montoya J.G. Recent Developments for Diagnosis of Toxoplasmosis JOURNAL OF CLINICAL MICROBIOLOGY (2004): 941–945
 24. Bessi eres M.H., Roques C., Berrebi A., Barre V., Cazaux M., S egu ela J.P. IgA antibody response during acquired and congenital toxoplasmosis. J Clin Pathol (1992); 45: 605-608

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14. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. P. Herbrink, A. van Loon et al.: Interlaboratory evaluation of indirect enzyme-linked immunosorbent assay, antibody capture enzyme-linked immunosorbent assay, and immunoblotting for detection of immunoglobulin M antibodies to *Toxoplasma gondii*. J. Clin. Microbiol. 25: 100 (1987).
3. F. Wielaard, H. van Gruijthuisen et al.: Diagnosis of acute toxoplasmosis by an enzyme immunoassay for specific immunoglobulin M antibodies. J. Clin. Microbiol. 17: 981 (1983).
4. J. Fung, A. Clogston et al.: Serologic diagnosis of toxoplasmosis with emphasis on the detection of toxoplasma-specific immunoglobulin M antibodies. Am. J. Clin. Pathol. 83: 196 (1985).
5. Remington J.S., Miller M.J., Brownlee I.: IgM antibodies in acute toxoplasmosis. I. Diagnostic significance in congenital cases and a method for their rapid demonstration. Pediatrics 41: 1082 (1968).
6. Pomares C., Montoya J.G. Laboratory Diagnosis of Congenital Toxoplasmosis. Journal of Clinical Microbiology Volume 54, Number 50 (2016): 2448-2454

7. Remington J.S., Thulliez P. and Montoya J.G. Recent Developments for Diagnosis of Toxoplasmosis JOURNAL OF CLINICAL MICROBIOLOGY (2004): 941-945
8. Bessières M.H., Roques C., Berrebi A., Barre V., Cazaux M., Séguéla J.P. IgA antibody response during acquired and congenital toxoplasmosis. J Clin Pathol (1992); 45: 605-608

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