CHORUS

DIESSE

**Cytomegalovirus
IgM Capture**

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15. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to Cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to Cytomegalovirus. J. Clin. Microbiol. 13: 416
5. M. Musiani et al.: Rapid detection of antibodies against Cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of Cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
7. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
8. Munro S.C et al.: Diagnosis of and Screening for Cytomegalovirus Infection in Pregnant Women. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2005, Vol 43, No. 9 : 4713-4718
9. Dollard S.C. et al.: National Prevalence Estimates for Cytomegalovirus IgM and IgG Avidity and Association between High IgM Antibody Titer and Low IgG Avidity. CLINICAL AND VACCINE IMMUNOLOGY, Nov. 2011, Vol. 18, No. 11 : 1895-1899

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0123

INSTRUCTIONS FOR USE

CHORUS Cytomegalovirus IgM Capture

**For the qualitative determination of IgM
antibodies anti-Cytomegalovirus**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic capture method for the qualitative determination of IgM class antibodies against Cytomegalovirus in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Cytomegalovirus is a herpes virus transmitted by close human contact. No symptoms of infection are apparent in the majority of cases.

However, the virus is very dangerous and may be fatal in immunodepressed patients.

Serum-negative female patients who become infected during pregnancy may transmit the disease to the fetus. In 95% of cases this occurs without symptoms, but some neonates may present jaundice, hepato-splenomegaly and retarded psycho-motorial development. For this reason it is of great importance to determine the immunitary state of the patient before the start of pregnancy and check for serum conversion.

The assay of specific IgM is of great importance in the diagnosis of primary infection.

3. PRINCIPLE OF THE METHOD

The Chorus Cytomegalovirus IgM Capture device is ready to use for the detection of IgM antibodies against Cytomegalovirus, in the Chorus/ Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay).

Monoclonal anti-human IgM antibodies are bound to the solid phase. The IgM immunoglobulins bind to anti-IgM antibodies through incubation with diluted sample.

After washing to eliminate the proteins which have not reacted, incubation is performed with the antigen bound to specific anti-Cytomegalovirus monoclonal antibodies, conjugated to horseradish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The enzymatic reaction is subsequently blocked by adding the Stop Solution which makes the solution turn into yellow. The colour which develops is proportional to the concentration of specific antibodies present in the examined sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

BEFORE STARTING THE TEST, MAKE SURE THAT THE TEMPERATURE OF THE STRIP, VISIBLE ON THE INSTRUMENT DISPLAY (PAGE CHECK), HAS REACHED 30°C ± 2°C.

Before use bring the package containing the devices to room temperature (18-30°C) for 30 minutes at least.

Use within 60 minutes.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)

5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use strongly hemolyzed, lipemic, icteric samples, serum not completely coagulated or samples presenting microbial contamination.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES

6 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: FREEZE-DRIED ANTIGEN

Contents: purified and inactivated Cytomegalovirus

Position 7: MICROPLATE WELL

Coated with anti-human IgM monoclonal antibodies

Position 6: Uncoated empty MICROPLATE WELL

Position 5: STOP SOLUTION

0.3 mol/L H₂SO₄ solution, ready for use.

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution in phosphate buffer containing sodium azide 0.09% and dye.

Position 2: CONJUGATE

Contents: anti-Cytomegalovirus monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and preservative.

Position 1: EMPTY WELL

In which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in

the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR

1 x 0.450 ml

Contents: Human serum containing anti-Cytomegalovirus IgM antibodies, diluted in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%.

Liquid, ready for use

CONTROL +

POSITIVE CONTROL

1 x 0.900 ml

Contents: Human serum containing anti-Cytomegalovirus IgM antibodies, diluted in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%.

Liquid, ready for use

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus/Chorus TRIO Instrument **REF** 81000-81200
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. **BEFORE STARTING THE TEST, MAKE SURE THAT THE TEMPERATURE OF THE STRIP, VISIBLE ON THE INSTRUMENT DISPLAY (PAGE CHECK), HAS REACHED 30°C ± 2°C.**

- Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
- Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
- Dispense the following volumes in well no. 1 of each device:

SAMPLE	50 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

At each change of batch, use a device for the calibrator.

- Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.1
NEGATIVE: when the result is < 0.9
DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. METHOD COMPARISON

In an experimentation 140 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	47	1	48
	-	0	92	92
	Total	47	93	140

Percent Positive Agreement (~Diagnostic Sensitivity):
100% CI_{95%}: 92.4-99.9

Percent Negative Agreement: (~Diagnostic Specificity):
98.9% CI_{95%}: 94.1-99.8

The agreement between the two methods is excellent with a Cohen's Kappa of 0.98.

13. CROSS-REACTIONS

44 samples, positive to Rheumatoid Factor, Varicella zoster, Epstein-Barr virus and Herpes simplex virus were tested.

No significant cross-reactions were found (except for Epstein-Barr virus and Herpes simplex virus).

14. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (Index)	CV%	Mean (Index)	CV%
1	2.3	12.6	2.6	5.8
2	1.1	12.7	1.2	3.3
3	0.5	10.0	0.5	10.0

Sample	Precision between batches		Precision between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	3.1	10.3	2.4	9.2
2	1.1	10.9	1.1	11.8
3	0.4	12.5	0.4	12.5

15. REFERENCES

- G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
- Grint P.C.A. et al.: Screening tests for antibodies to Cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
- Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to Cytomegalovirus. J. Clin. Microbiol. 13: 416
- M. Musiani et al.: Rapid detection of antibodies against Cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
- F. de Ory et al.: Serological diagnosis of Cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
- P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
- Munro S.C et al.: Diagnosis of and Screening for Cytomegalovirus Infection in Pregnant Women. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2005, Vol 43, No. 9 : 4713-4718
- Dollard S.C. et al.: National Prevalence Estimates for Cytomegalovirus IgM and IgG Avidity and Association between High IgM Antibody Titer and Low IgG Avidity. CLINICAL AND VACCINE IMMUNOLOGY, Nov. 2011, Vol. 18, No. 11 : 1895-1899

DIESE Diagnostica Senese
S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (Siena)
Italy

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- G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
- Grint P.C.A. et al.: Screening tests for antibodies to Cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
- Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to Cytomegalovirus. J. Clin. Microbiol. 13: 416
- M. Musiani et al.: Rapid detection of antibodies against Cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
- F. de Ory et al.: Serological diagnosis of Cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
- P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
- Munro S.C. et al.: Diagnosis of and Screening for Cytomegalovirus Infection in Pregnant Women. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2005, Vol 43, No. 9 : 4713-4718
- Dollard S.C. et al.: National Prevalence Estimates for Cytomegalovirus IgM and IgG Avidity and Association between High IgM Antibody Titer and Low IgG Avidity.

CLINICAL AND VACCINE IMMUNOLOGY, Nov. 2011, Vol.
18, No. 11 : 1895–1899

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23. LITERATUR

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to Cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to Cytomegalovirus. J. Clin. Microbiol. 13: 416
5. M. Musiani et al.: Rapid detection of antibodies against Cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of Cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
7. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
8. Munro S.C et al.: Diagnosis of and Screening for Cytomegalovirus Infection in Pregnant Women. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2005, Vol 43, No. 9 : 4713–4718
9. Dollard S.C. et al.: National Prevalence Estimates for Cytomegalovirus IgM and IgG Avidity and Association between High IgM Antibody Titer and Low IgG Avidity. CLINICAL AND VACCINE IMMUNOLOGY, Nov. 2011, Vol. 18, No. 11 : 1895–1899

DIESSE Diagnostica Senese
S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (Siena)
Italien

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15. BIBLIOGRAFÍA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to Cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to Cytomegalovirus. J. Clin. Microbiol. 13: 416
5. M. Musiani et al.: Rapid detection of antibodies against Cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of Cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).

7. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
8. Munro S.C et al.: Diagnosis of and Screening for Cytomegalovirus Infection in Pregnant Women. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2005, Vol 43, No. 9 : 4713–4718
9. Dollard S.C. et al.: National Prevalence Estimates for Cytomegalovirus IgM and IgG Avidity and Association between High IgM Antibody Titer and Low IgG Avidity. CLINICAL AND VACCINE IMMUNOLOGY, Nov. 2011, Vol. 18, No. 11 : 1895–1899

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15. BIBLIOGRAPHIE

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).

3. Grint P.C.A. et al.: Screening tests for antibodies to Cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to Cytomegalovirus. J. Clin. Microbiol. 13: 416
5. M. Musiani et al.: Rapid detection of antibodies against Cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of Cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
7. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
8. Munro S.C et al.: Diagnosis of and Screening for Cytomegalovirus Infection in Pregnant Women. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2005, Vol 43, No. 9 : 4713–4718
9. Dollard S.C. et al.: National Prevalence Estimates for Cytomegalovirus IgM and IgG Avidity and Association between High IgM Antibody Titer and Low IgG Avidity. CLINICAL AND VACCINE IMMUNOLOGY, Nov. 2011, Vol. 18, No. 11 : 1895–1899

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15. BIBLIOGRAFIE

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to Cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to Cytomegalovirus. J. Clin. Microbiol. 13: 416
5. M. Musiani et al.: Rapid detection of antibodies against Cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of Cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
7. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).

8. Munro S.C et al.: Diagnosis of and Screening for Cytomegalovirus Infection in Pregnant Women. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2005, Vol 43, No. 9 : 4713–4718
9. Dollard S.C. et al.: National Prevalence Estimates for Cytomegalovirus IgM and IgG Avidity and Association between High IgM Antibody Titer and Low IgG Avidity. CLINICAL AND VACCINE IMMUNOLOGY, Nov. 2011, Vol. 18, No. 11 : 1895–1899

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