CHORUS

Borrelia IgM

DIESSE

REF 81047/C03

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

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17. BIBLIOGRAFIA

1. L.A. Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickstein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi Antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444
7. Lindgren E., Jaenson T.G.T Lyme borreliosis in Europe: influences of climate and climate change, epidemiology, ecology and adaptation measures. World Health Organization 2006.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

INSTRUCTIONS FOR USE

CHORUS Borrelia IgM

For the qualitative determination of anti-Borrelia burgdorferi sensu latu IgM antibodies

For In Vitro Diagnostic Use Only

1. INTENDED USE

Third generation immunoenzymatic method for the qualitative determination of IgM class antibodies against *Borrelia burgdorferi sensu latu* in human serum and citrate, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Lyme disease (Borreliosis) is a multi-systemic infectious disease caused by *Borrelia burgdorferi* bacteria belonging to the Spirochaete family.

The infection is transmitted to humans by the bite of ticks of the Ixodes genus and it can be divided in three clinical phases: early localized infection, early disseminated infection, late disseminated infection. During the first phase about two-thirds of patients develop a skin erythema together with flu-like symptoms. During the second phase, which occurs after weeks or months, the infection spreads through the bloodstream and the lymphatic system causing muscle-skeletal and neurological disorders. The late phase occurs after months or years and can provoke acrodermatitis chronica atrophicans (ACA), chronic neuroborreliosis, polyarthrits.

The diagnosis of Borreliosis is based on the patient history, clinical picture and results of laboratory tests.

The antibody production can be extremely slow during the early phase of the disease; on the other hand the antibodies may remain for more than ten years.

The acute phase is generally characterized by high titers IgM antibodies. High IgG titers with a low or absent IgM concentration appear during the phase of recovery or during the chronic stage.

3. PRINCIPLE OF THE METHOD

The Chorus Borrelia IgM device is a third generation immunoenzymatic method ready to use for the detection of IgM antibodies against *Borrelia burgdorferi sensu latu* with high diagnostic sensitivity and specificity, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted sample.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgM antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase

substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Before use bring the package containing the devices to room temperature (18-30°C) for 30 minutes at least.

Use within 60 minutes.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

4. The devices are for use with the Chorus/Chorus TRIO instrument, the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)

5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed and lipemic samples, of plasma with anticoagulants different from citrate or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**


5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 30 tests.

DD DEVICES

- 1 package each containing 6 devices for calibration
- 4 packages each containing 6 devices for sample

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Combination of recombinant antigens: OspC (B. afzelii, B. garinii, B. burgdorferi sensu stricto, B. spielmanii), VlsE, internal flagelin - p41i, p39, p17 and OspE of Borrelia burgdorferi sensu latu species.

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic buffer solution at pH 7.1

Position 2: CONJUGATE

Contents: anti-human IgM antibodies labelled with peroxidase.

Position 1: EMPTY WELL

In which the operator must place the undiluted sample

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in

the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.150 mL

Contents: Buffer containing anti-Borrelia burgdorferi IgM antibodies. Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.250 mL

Contents: Buffer containing anti-Borrelia burgdorferi IgM antibodies. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER [REF] 83606
- CLEANING SOLUTION 2000 [REF] 83609
- SANITIZING SOLUTION [REF] 83604 – 83608
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum and citrate collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh sample may be stored for a week at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted calibrator in well no. 1 of a device for calibration (yellow label) and 50 µl of undiluted positive control in well no. 1 of another device for calibration (yellow label); repeat at each change of batch.

4. Dispense 50 µl of undiluted test sample in well no. 1 of each device for sample (white label).
5. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 1.2

NEGATIVE: when the result is < 0.8

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CROSS-REACTIONS

The serological diagnosis is difficult because of the wide genetic range of the *Borrelia burgdorferi* species, the possible cross-reactions with extraneous antigens and the abundance of *Borrelia* heat shock proteins.

46 samples, positive to Cytomegalovirus, Treponema, Herpes simplex virus, Epstein-Barr virus, Rheumatoid Factor and ANA have been tested.

It is not possible to exclude cross-reactions of antibodies directed against Cytomegalovirus, Treponema, Epstein-Barr virus in acute phase and Rheumatoid Factor with titer superior to 140 IU/ml.

13. METHOD COMPARISON

In an experimentation 91 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	49	0	49
	-	1	41	42
	Total	50	41	91

Percent Positive Agreement (~Diagnostic Sensitivity):

98.0% CI_{95%}: 89.4-99.6

Percent Negative Agreement (~Diagnostic Specificity):

100.0% CI_{95%}: 91.4-99.9

The agreement between the two methods is excellent with a Cohen's Kappa of 0.98.

14. DIAGNOSTIC SENSITIVITY

97 samples from patients at different stages of the disease have been analysed.

The samples were from subjects clinically classified in the different stages of the disease, subsequently characterized with the kit.

The obtained data are summarized in the table below:

Stage of the disease	Samples	Diagnostic Sensitivity (%)	CI _{95%}
Erythema migrans (Stage I)	38	42.1	27.8-57.7
Neuroborreliosis (Stage II)	29	100	88.3-99.8
Lyme Arthritis (Stage III)	30	26.7	14.2-44.3

15. DIAGNOSTIC SPECIFICITY

100 samples from healthy donors have been analysed.

The samples were from European subjects, living in a not endemic geographic area (central Italy) with low seroprevalence (< 20 cases per year).

The obtained data are summarized in the table below:

Diagnostic Specificity (%)	CI _{95%}
96.0	90.1-98.4

16. PRECISION AND REPEATABILITY

Sample	Within run	
	Mean (Index)	CV%
1	0.3	-
2	0.1	40.0*
3	1.5	-
4	1.4	2.9

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero.

Sample	Between run	
	Mean (Index)	CV%
1	0.3	20.0*
2	0.2	25.0*
3	1.4	10.7
4	1.4	13.6

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero.

Sample	Between lots		Between Instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	0.4	22.5*	0.4	25.0*
2	0.2	25.0*	0.2	25.0*
3	1.8	8.3	1.6	13.8
4	1.6	9.4	1.4	12.9

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero.

17. REFERENCES

1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to *Borrelia burgdorferi* 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235–60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to *Borrelia burgdorferi*. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to *Borrelia burgdorferi* Antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437–7444
7. Lindgren E., Jaenson T.G.T Lyme borreliosis in Europe: influences of climate and climate change, epidemiology, ecology and adaptation measures. World Health Organization 2006.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235–60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi Antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437–7444
7. Lindgren E., Jaenson T.G.T Lyme borreliosis in Europe: influences of climate and climate change, epidemiology, ecology and adaptation measures. World Health Organization 2006.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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17. LITERATUR

1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickstein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi Antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444
7. Austin Lindgren E., Jaenson T.G.T Lyme borreliosis in Europe: influences of climate and climate change, epidemiology, ecology and adaptation measures. World Health Organization 2006.

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1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious Disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi Antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444
7. Lindgren E., Jaenson T.G.T Lyme borreliosis in Europe: influences of climate and climate change, epidemiology, ecology and adaptation measures. World Health Organization 2006.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34

6. Austin Vaz, Lisa Glickstein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi Antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444
7. Lindgren E., Jaenson T.G.T Lyme borreliosis in Europe: influences of climate and climate change, epidemiology, ecology and adaptation measures. World Health Organization 2006.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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17. BIBLIOGRAFIA

1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, nº 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi Antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444
7. Lindgren E., Jaenson T.G.T Lyme borreliosis in Europe: influences of climate and climate change, epidemiology, ecology and adaptation measures. World Health Organization 2006.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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17. BIBLIOGRAFIE

1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi Antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444
7. Lindgren E., Jaenson T.G.T Lyme borreliosis in Europe: influences of climate and climate change, epidemiology, ecology and adaptation measures. World Health Organization 2006.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione	FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico
	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabbrikante	FR Fabricant GR Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote