TPHA-DAT

REF 26035

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INSTRUCTIONS FOR USE

TPHA-DAT

Test for the semiquantitative determination of antibodies against *Treponema pallidum* in human serum, by passive haemoagglutination, applied to AUTO-DAT instrument.

For In Vitro Diagnostic Use Only

1. INTENDED USE

Test, by passive haemoagglutination, for the semiquantitative determination of antibodies against *Treponema pallidum* in human serum.

2. INTRODUCTION

Syphilis is a chronic infection that evolves into distinct stages of primary, secondary, latent and tertiary infection. These stages determine different clinical symptoms, generally producing an initial lesion known as syphiloma, which then evolves into a syphilitic eruption followed by long latency periods. If left untreated, the infection can cause cardiovascular problems and neurosyphilis.

The infection is caused by a bacterium of the spirochetes order, *Treponema pallidum*; it is, generally, sexually transmitted, although the disease can also be transmitted by infected blood transfusions. Intrauterine infections are also possible.

The diagnosis is based on the use of four test categories: direct microscopic examination, direct antigen test, non-treponemic antibody test and treponemic antibody test. TPHA test, which determines the presence of specific anti-treponemic antibodies, is widely used for the screening of syphilis.

3. PRINCIPLE OF THE METHOD

The specific antibodies, present in serum of the patient, bind to the treponemic antigens (Nichols strain), coated to avian erythrocytes.

In case of cell agglutination, the reaction is positive. If the cells form a compact button, the reaction is negative.

The results are printed by the Auto-DAT instrument. In any case the titer can be estimated visually.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as

potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples

Waste disposal: samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the device to room temperature (18-30°C) before use.

1. After use, and return reagents to 2-8°C.
2. Do not use the reagents after the expiry date.
3. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
4. Do not modify the Test Procedure. Do not substitute reagents using reagents from other suppliers or other lots, unless it is specifically indicated that the reagent is interchangeable between lots.
5. Do not extend the incubation further than 24 hours.
6. Do not use hemolyzed, lipemic, icteric samples, serum not completely coagulated or samples presenting microbial contamination.

5. KIT CONTENTS AND REAGENT PREPARATION

Reagents are sufficient for 50 determinations.

TEST CELLS TES CELLS 2 x 4 mL

Contents: avian erythrocytes, coated with *Treponema pallidum* antigens.

CONTROL CELLS CONTROL CELLS 1 x 4 mL

Contents: uncoated avian erythrocytes.

SAMPLE DILUENT SAMPLE DILUENT 1 x 10 mL

Contents: saline solution containing substances which inhibit non-specific reactions.

CONTROL + POSITIVE CONTROL 1 x 0.3 ml

Contents: diluted human serum containing anti-Treponema pallidum antibodies. Liquid, ready to use.

CONTROL - NEGATIVE CONTROL 1 x 0.3 mL

Contents: diluted human serum not containing anti-Treponema pallidum antibodies. Liquid, ready to use.

96 U-WELLS MT MICROPLATES: 3

BRING THE REAGENTS TO ROOM TEMPERATURE BEFORE USE.

MATERIALS REQUIRED BUT NOT PROVIDED:

- Pipettes with disposable tips
- AUTO-DAT instrument (REF 26000)

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

Do not freeze. Keep the vials in vertical position.

The expiry date is printed on each component and on the kit label.

If stored and used according to the Instructions, the reagents are stable up to the expiration date.

7. SPECIMEN COLLECTION AND STORAGE

Fresh serum stored at 2-8°C for a maximum of 4 days. For longer storage, freeze at -20°C.

Two subsequent freezing and thawing can be performed.

Thawed samples must be well mixed before the test.

Do not keep the samples in auto-defrosting freezers.

8. ASSAY PROCEDURE

Follow the indications reported on the User Manual to perform the test with the AUTO-DAT instrument.

8.1 Assay with AUTO-DAT

Assay range: the instrument expresses the results from negative (NEG) to $\geq 1/2560$. The evaluation is performed by observing the sample agglutination at the dilutions of 1/80 and 1/640.

In order to dilute the samples 1/80, use the wells of rows A and E.

- Add 95 µl of Sample diluent to well A1 and well E1 and then 5 µl of the first and second samples respectively.
- Mix thoroughly.
- Collect 25 µl of solution present in A1 and transfer to well B1.
- Collect 25 µl of solution present in A1 and transfer to well C1.
- Collect 25 µl of solution present in E1 and transfer to well F1.
- Collect 25 µl of solution present in E1 and transfer to well G1.
- Resuspend thoroughly by inversion both the Test Cells and the Control Cells.
- Add 75 µl of Control Cells to wells B1 and F1.

- Add 75 µl of Test Cells to wells C1 and G1.

In order to dilute the samples 1:640, use the wells of rows D and H.

- Collect 3.5 µl of the solutions in A1 and E1, adding them to 25 µl of diluent.
- Mix thoroughly and discard 3.5 µl of solution
- Add 75 µl of Test Cells to each well

	1	2
A	Sample 1 5 µl + 95 µl sample diluent	Positive Control
B	25 µl solution present in A1 + 75 µl Control Cells	25 µl Positive Control + 75 µl Control Cells
C	25 µl solution present in A1 + 75 µl Test Cells	25 µl Positive Control + 75 µl Test Cells
D	3.5 µl solution present in A1 + 25 µl sample diluent + 75 µl Test Cells	25 µl Positive Control + 75 µl Test cells
E	Sample 2 5 µl + 95 µl sample diluent	Negative Control
F	25 µl solution present in E1 + 75 µl Control Cells	25 µl Negative Control + 75 µl Control Cells
G	25 µl solution present in E1 + 75 µl Test Cells	25 µl Negative Control + 75 µl Test Cells
H	3.5 µl solution present in E1 + 25 µl sample diluent + 75 µl Test Cells	25 µl Negative Control + 75 µl Test Cells

Repeat the operation for all the subsequent pairs of samples using the columns of plate from 2 to 12.

Once the dilutions have been prepared and all the reagents added:

- Mix thoroughly avoiding the formation of bubbles.
- Place the plate inside the AUTO-DAT to perform the incubation for 60 minutes and the relative reading.

Note: the positive and negative controls are supplied ready for use and they do not require dilution.

Positive and negative controls should be used following the samples scheme; 4 wells are used for the positive control and 4 wells are used for the negative control.

The first well is left empty (row A and E); in the 3rd and 4th well the test with test cells (in double) is performed.

If the instrument is already in use or if more than 24 samples have to be tested, the incubation can be performed on the lab bench.

In this case, insert the plate into the AUTO-DAT, with high caution, avoiding accidental knocks.

In case the AUTO-DAT instrument cannot be used, the test can be performed in manual mode with a doubling dilution and a visual reading of the result.

In this case follow the procedure described at paragraph 8.2.

8.2 Assay with manual procedure

In order to dilute the samples, use the wells of rows A and E.

The test is carried out using the wells of the plate column.

The wells in row A are used for the dilution of the sample, those in row B are used to perform the test with the Control Cells.

- Dilute the sample in well A1 as described at paragraph 8.1.
- Leave both the wells B1 and C1 empty.
- Add 25 µl of sample diluent to all the subsequent wells in column 1.
- Add 25 µl of sample to the second, third and fourth well. From the latter make a doubling dilution, up to the eighth well. Discard 25 µl from this last one.
- Add 75 µl of Control Cells to well 2 and 75 µl of Test Cells to all the others.

The titer of the sample can vary from 1/80 to 1/2560.

	1
A	Sample 1 5 µl + 95 µl sample diluent
B	25 µl sample present in A1 + 75 µl Control Cells
C	25 µl sample present in A1 + 75 µl Test Cells
D	25 µl sample diluent + 25 µl sample present in A1 + 75 µl Test Cells
E	25 µl sample diluent + 25 µl sample present in D1 + 75 µl Test Cells
F	25 µl sample diluent + 25 µl sample present in E1 + 75 µl Test Cells
G	25 µl sample diluent + 25 µl sample present in F1 + 75 µl Test Cells
H	25 µl sample diluent + 25 µl sample present in G1; discard 25 µl of solution; add 75 µl Test Cells

- Wait for 60 minutes on a vibration-free bench and observe the formation of agglutinate or erythrocytes compact button.

9. TEST VALIDATION

Perform the test with positive and negative control.

If the result is different from the expected, contact the Scientific Support:

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

In absence of agglutination the erythrocytes sediment, forming a compact button in the center of the well.

In case of agglutination the erythrocytes sediment, forming a homogeneously diffused layer on the bottom of the well and/or forming a characteristic ring structure.

The agglutination of the Test Cells and the simultaneous absence of agglutination of the Control Cells is an index of the presence of anti-Treponema pallidum antibodies.

The absence of agglutination indicates that the antibodies, if present, are at a concentration lower than the test sensitivity limit.

The Control Cells sediment forming a particularly compact button; in the event of their agglutination, the test is not to be considered valid.

The semi-quantitative titer is given by the last dilution able to make agglutination.

The presence of an erythrocytes compact button on the bottom is interpreted as a negative result.

The titer 1/80 is considered the limit between non-reactive and reactive samples.

11. LIMITATIONS

The antibodies detected by TPHA develop 4-5 weeks after the onset of the disease.

No tests in haemoagglutination on serum can distinguish between antibodies from Treponema pallidum infections and antibodies from infections caused by other pathogen treponemas.

Lepromatous leprosy and infectious mononucleosis can show false positive reactions with this type of test.

The test result should be always evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALYTICAL SENSITIVITY

Analytical sensitivity: 1/80.

13. METHOD COMPARISON

In an experimentation 419 samples have been tested.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	207	0	207
	-	3	209	212
	Total	210	209	419

Percent Positive Agreement (~Diagnostic Sensitivity):

98.6% CI_{95%}: 95.9.-99.5

Percent Negative Agreement: (~Diagnostic Specificity):

100% CI_{95%}: 98.2-100

The agreement between the two methods is excellent with a Cohen's Kappa of 0.98.

14. REFERENCES

- Lesiński J., Krach J. and Kadziewicz E. Specificity, sensitivity and diagnostic value of the TPHA test. British Journal of Venereal Diseases (1974) 50,334.
- Rathlev T. Haemagglutination test utilizing pathogenic Treponema pallidum for the sero-diagnosis of Syphilis. British Journal of Venereal Diseases (1967), 43, 181.
- Uete T., Fukazawa S., Ogi K. and Takeuchi Y. Clinical evaluation of T. pallidum haemagglutination test. British Journal of Venereal Diseases (1971) 47, 73.
- Ratnam S. The laboratory diagnosis of syphilis. Can J Infect Dis Med Microbiol (2005); 16(1):45-51.
- Laboratory diagnosis of sexually transmitted infections, including human immunodeficiency virus. WHO (2013).

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusqu'à Ημερομηνία λήξης Prazo de validade
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbrikante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote
	EN ES IT	CE marking of conformity Marcado CE de conformidad Marcatura CE di conformità	FR GR PT	Marquage de conformité CE Σημανση συμμορφωσης CE Marcação CE de conformidade