CHORUS

PSA

DIESSE

REF 86550

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

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CHORUS PSA

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18. BIBLIOGRAFIA

- Christensson A., Laurell C.B., Lilja H., Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine proteinase inhibitors. Eur. J. Biochem., 194, 755-763 (1990).
- Watt KW et al., Human prostate-specific antigen: Structural and functional similarity with serine proteases. Proc. Natl. Acad. Sci. USA, 83, 3166-3170 (1986).
- Chen Z., Prestigiacomo A., Stamey T., Purification and characterization of Prostate-Specific Antigen (PSA) Complexed to α_1 -Anticymotrypsin: Potential Reference Material for International Standardization of PSA Immunoassays. Clin. Chem., 41/9, 1273-1282 (1995).
- Junker R., Brandt B., Zechel C., Assmann G., Comparison of prostate-specific antigen (PSA) measured by four combinations of free PSA and total PSA assays. Clin. Chem., 43, 1588-1594 (1997).
- Stenman UH., Leinonen J., Alfthan H., Rannikko S., Tuhkanen K., Alfthan O., A complex between prostate-specific antigen and α_1 -anticymotrypsin is the major form of prostate-specific antigen in serum of patients with prostatic cancer: assay of the complex improves clinical sensitivity for cancer. Cancer Research, 51, 222-226 (1991).

6. Park S1, Wians FH Jr, Cadeddu JA, Spurious prostate-specific antigen (PSA) recurrence after radical prostatectomy: interference by human antimouse heterophile antibodies. *Int J Urol.* (2007) Mar; 14(3):251-3.
7. El Ezzi A.A. and El-Saidi M.A. Stability of total and free Prostate Specific Antigen in serum submitted to intermittent cold storage conditions. *Indian Journal of Clinical Biochemistry* (2009), 24 (2) 166-174.

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

0123

INSTRUCTIONS FOR USE

CHORUS PSA

For the quantitative determination of Total Prostate Specific Antigen (tPSA)

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the quantitative determination of Total Prostate Specific Antigen (tPSA) in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Prostate Specific Antigen (PSA) is a serine protease with chymotrypsin-like activity. This single chain glycoprotein with a molecular weight of 28.4 kDa has the function to avoid the formation of aggregates in the seminal fluid.

Small quantities of PSA are present in the serum, the amount increases in case of manipulation of the prostate gland or in case of pathologies such as prostatitis, hypertrophy or prostate cancer.

Since prostate cancer is the second most prevalent form of male neoplasia, the detection of elevated PSA levels plays an important role in the early diagnosis.

3. PRINCIPLE OF THE METHOD

The Chorus PSA device is ready to use for the detection of Total Prostate Specific Antigen (tPSA), in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The streptavidin is bound to the solid phase. The addition of the sample and of the enzymatic reagent, constituted by an anti-PSA biotinylated monoclonal antibody and anti-PSA conjugated monoclonal antibody, directed towards distinct epitopes, allows the formation of a sandwich complex that is bound to the solid phase.

After washing to eliminate the proteins which have not reacted, the peroxidase substrate is added. The blue colour which develops is proportional to the concentration of tPSA present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in ng/ml.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can

offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the devices within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>).
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.


8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested (see the "Analytical specificity" paragraph), serum not completely coagulated or samples presenting microbial contamination to avoid unattended results.
11. Do not use the product after the expiry date reported in the box label.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with Streptavidin

Position 5: Empty

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: Empty

Position 2: ENZYMATIC REAGENT

Contents: anti-PSA biotinylated monoclonal antibodies and anti-PSA monoclonal antibodies conjugated with horseradish peroxidase, in buffer solution with preservatives.

Position 1: EMPTY WELL

In which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.250 mL

Contents: Diluted human serum containing tPSA and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.500 mL

Contents: Diluted human serum containing tPSA and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus/Chorus TRIO Instrument **REF** 81000 – 81200
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.

- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

CALIBRATOR 8 weeks at 2/8°C

POSITIVE CONTROL 8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 5 days at 2/8°C, or frozen at -20°C for longer periods (for at least 20 months), and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 75 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554

Fax: 0039 0577 366605

email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in ng/ml, calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

Healthy males should have values < 4 ng/ml

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range 0.10-100 ng/ml.

Functional sensitivity: 0.30 ng/ml.

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 0.11 and 3.8 ng/ml.

Normal range for a population is dependent upon a lot of factors like the population tested and the specificity and sensitivity of the used method.

For this reason, each laboratory should determinate an internal range of expected values using the test with the local population.

14. ANALYTICAL SPECIFICITY

A study has been carried out using the reagent of the kit in presence of the following substances.

SUBSTANCE	CONCENTRATION
ACETYSALICYLIC ACID	100 µg/ml
ASCORBIC ACID	100 µg/ml
CAFFEINE	100 µg/ml
CEA	10 µg/ml
AFP	10 µg/ml
CA-125	10.000 U/ml
hCG	1000 IU/ml
hLH	10 IU/ml
hTSH	100 ng/ml
hPRL	100 µg/ml

No interferences have been detected.

Moreover, 5 samples were spiked with the following potentially interfering factors and then tested:

Rheumatoid Factor (44 – 220 IU/ml)
 Bilirubin (4.5 mg/dl – 18 mg/dl)
 Triglycerides (10 mg/dl – 150 mg/dl)
 Hemoglobin (5 mg/ml – 15 mg/ml)

The presence in the serum sample of the interfering substances described above (up to the tested concentration) does not affect the test result.

15. CROSS-REACTIONS

The interference of human anti-mouse heterophile antibodies (HAMA) is known.

16. METHOD COMPARISON

In an experimentation 167 samples have been tested with Diesse kit and with a competitor kit.

Results of the analysis:

Correlation	r	95%CI
Pearson	0.99	0.98-0.99
Spearman	0.98	0.98-0.99

The correlation between the two methods is very high.

17. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	20.6	2.6	21.4	7.1
2	5.2	4.2	1.7	14.7
3	2.0	3.0	4.4	14.5

Sample	Precision between batches		Precision between instruments	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	23.9	13.0	21.7	6.8
2	1.7	12.9	1.6	13.1
3	4.6	3.5	4.2	11.4

18. REFERENCES

- Christensson A., Laurell C.B., Lilja H., Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine proteinase inhibitors. Eur. J. Biochem., 194, 755-763 (1990).
- Watt KW et al., Human prostate-specific antigen: Structural and functional similarity with serine proteases. Proc. Natl. Acad. Sci. USA, 83, 3166-3170 (1986).
- Chen Z., Prestigiacomo A., Stamey T., Purification and characterization of Prostate-Specific Antigen (PSA) Complexed to α_1 -Anticymotrypsin: Potential Reference Material for International Standardization of PSA Immunoassays. Clin. Chem., 41/9, 1273-1282 (1995).
- Junker R., Brandt B., Zechel C., Assmann G., Comparison of prostate-specific antigen (PSA) measured by four combinations of free PSA and total PSA assays. Clin. Chem., 43, 1588-1594 (1997).
- Stenman UH., Leinonen J., Alfthan H., Rannikko S., Tuhkanen K., Alfthan O., A complex between prostate-specific antigen and α_1 -anticymotrypsin is the major form of prostate-specific antigen in serum of patients with prostatic cancer: assay of the complex improves clinical sensitivity for cancer. Cancer Research, 51, 222-226 (1991).
- Park S1, Wians FH Jr, Cadeddu JA, Spurious prostate-specific antigen (PSA) recurrence after radical prostatectomy: interference by human antimouse heterophile antibodies. Int J Urol. (2007) Mar; 14(3):251-3.
- El Ezzi A.A. and El-Saidi M.A. Stability of total and free Prostate Specific Antigen in serum submitted to intermittent cold storage conditions. Indian Journal of Clinical Biochemistry (2009), 24 (2) 166-174.

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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- Christensson A., Laurell C.B., Lilja H., Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine proteinase inhibitors. Eur. J. Biochem., 194, 755-763 (1990).
- Watt KW et al., Human prostate-specific antigen: Structural and functional similarity with serine proteases. Proc. Natl. Acad. Sci. USA, 83, 3166-3170 (1986).
- Chen Z., Prestigiacomo A., Stamey T., Purification and characterization of Prostate-Specific Antigen (PSA) Complexed to α_1 -Anticymotrypsin: Potential Reference Material for International Standardization of PSA Immunoassays. Clin. Chem., 41/9, 1273-1282 (1995).

4. Junker R., Brandt B., Zechel C., Assmann G., Comparison of prostate-specific antigen (PSA) measured by four combinations of free PSA and total PSA assays. Clin. Chem., 43, 1588-1594 (1997).
5. Stenman UH., Leinonen J., Alfthan H., Rannikko S., Tuhkanen K., Alfthan O., A complex between prostate-specific antigen and α_1 -antitrypsin is the major form of prostate-specific antigen in serum of patients with prostatic cancer: assay of the complex improves clinical sensitivity for cancer. Cancer Research, 51, 222-226 (1991).
6. Park S1, Wians FH Jr, Cadeddu JA, Spurious prostate-specific antigen (PSA) recurrence after radical prostatectomy: interference by human antimouse heterophile antibodies. Int J Urol. (2007) Mar; 14(3):251-3.
7. El Ezzi A.A. and El-Saidi M.A. Stability of total and free Prostate Specific Antigen in serum submitted to intermittent cold storage conditions. Indian Journal of Clinical Biochemistry (2009), 24 (2) 166-174.

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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18. BIBLIOGRAPHIE

1. Christensson A., Laurell C.B., Lilja H., Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine proteinase inhibitors. Eur. J. Biochem., 194, 755-763 (1990).

2. Watt KW et al., Human prostate-specific antigen: Structural and functional similarity with serine proteases. Proc. Natl. Acad. Sci. USA, 83, 3166-3170 (1986).
3. Chen Z., Prestigiacomo A., Stamey T., Purification and characterization of Prostate-Specific Antigen (PSA) Complexed to α_1 -Anticymotrypsin: Potential Reference Material for International Standardization of PSA Immunoassays. Clin. Chem., 41/9, 1273-1282 (1995).
4. Junker R., Brandt B., Zechel C., Assmann G., Comparison of prostate-specific antigen (PSA) measured by four combinations of free PSA and total PSA assays. Clin. Chem., 43, 1588-1594 (1997).
5. Stenman UH., Leinonen J., Alfthan H., Rannikko S., Tuhkanen K., Alfthan O., A complex between prostate-specific antigen and α_1 -anticymotrypsin is the major form of prostate-specific antigen in serum of patients with prostatic cancer: assay of the complex improves clinical sensitivity for cancer. Cancer Research, 51, 222-226 (1991).
6. Park S1, Wians FH Jr, Cadeddu JA, Spurious prostate-specific antigen (PSA) recurrence after radical prostatectomy: interference by human antimouse heterophile antibodies. Int J Urol. (2007) Mar; 14(3):251-3.
7. El Ezzi A.A. and El-Saidi M.A. Stability of total and free Prostate Specific Antigen in serum submitted to intermittent cold storage conditions. Indian Journal of Clinical Biochemistry (2009), 24 (2) 166-174.

DIESE Diagnostica Senese
S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

0123

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione	FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico
	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabbicante	FR Fabricant GR Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote