**iRapid
SARS-CoV-2 Quant
"Neutralizing" Ab**

REF 70100

CE

DIESSE

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17. BIBLIOGRAFIA

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

INSTRUCTIONS FOR USE

iRapid SARS-CoV-2 Quant "Neutralizing" Ab

For the semi-quantitative determination of IgG antibodies anti-S1(RBD) of SARS-CoV-2 virus as an aid in the determination of the immune reaction to SARS-CoV-2 or to the vaccine

**For *In Vitro* Diagnostic Use Only
For Professional Use Only**

1. INTENDED PURPOSE

Immunochromatographic test for the semi-quantitative determination of anti-S1(RBD) SARS-CoV-2 virus IgG antibodies in human serum, plasma and fingerstick whole blood.

The test is for professional use only and is intended to be used as an aid in the determination of the immune status after infection or vaccination.

The results can be interpreted as Binding Antibody Units (BAU/ml) referred to the first international standard WHO 20/136 for anti-SARS-CoV-2.

The product may be used in any laboratory and can not be used as the sole basis for disease management.

2. INTRODUCTION

Coronaviruses (CoV) are positive single-stranded RNA viruses belonging to the Coronaviridae family, which are divided into four genera: Alpha, Beta, Delta and Gammacoronavirus. CoVs are commonly found in many animal species. Occasionally, animal CoVs may acquire genetic mutations by errors during genome replication, which can further expand their tropism to humans.

A total of six human CoV types have been identified as responsible for human respiratory disorders, which include two CoV alpha and four CoV beta (the two most recent being SARS-CoV and MERS-CoV). Typically, these CoVs cause asymptomatic infection or severe acute respiratory illness, including fever, cough, and shortness of breath. However, other symptoms such as gastroenteritis and neurological diseases of varying severity have also been reported. Since December 2019, the World Health Organization has reported the seventh human betacoronavirus (2019-nCoV) that caused a pneumonia epidemic. This new virus is linked to an outbreak of febrile respiratory disease in Wuhan city, Hubei province, China. Some patients with 2019-nCoV developed severe pneumonia, pulmonary edema, ARDS, or multiple organ failure and died.

A serological assay helps to identify individuals exposed to the virus and will be useful to decide the extent of containment measures in COVID-19 disease.

The determination of IgG anti-S1(RBD) antibody titer can be used as an indication of immune reaction in post-COVID-19


patients or in vaccinated individuals also indicating the presence of neutralizing antibodies.

3. PRINCIPLE OF THE METHOD

The iRapid SARS-CoV-2 Quant "Neutralizing" Ab is a semi-quantitative membrane-based rapid immunoassay for the detection of IgG antibodies directed against the receptor binding domain (RBD) of SARS-CoV-2 virus in whole blood, serum and plasma samples.

Antibodies detecting human IgG are coated on the test line region (T) of the membrane while an anti-SARS-CoV-2 RBD is coated on the control line region (C) of the membrane.

During testing, the specimen is added to the sample pad of the strip (S) and two drops of buffer are added. The antibodies against SARS-CoV-2 RBD in the sample will bind to RBD protein coated to colloidal gold present in the conjugate paper.

The mixture then migrates upward on the membrane chromatographically by capillary action and reacts with the monoclonal anti-human IgG in the test line region. Later during migration, the RBD conjugated with colloidal gold will bind to the anti-SARS-CoV-2 RBD present on the control line.

If the specimen contains IgG antibodies against SARS-CoV-2 RBD at a concentration above limit of detection, a colored line will appear in the test line region (T).

If there are no IgG antibodies against SARS-CoV-2 reception binding domain or to a lower concentration than limit of detection, only the control line will appear. The test will be validated only if the control line appears, otherwise the test is invalid and cannot be interpreted.

The results are shown as coloured lines: the colour intensity is proportional to antibodies concentration, measured in Binding Antibody Units (BAU/ml), calculated with reference to the first international standard WHO 20/136 for anti-SARS-CoV-2.

4. WARNINGS AND PRECAUTIONS

**For *In Vitro* Diagnostic Use Only
For Professional Use Only**

Waste disposal: samples, cassettes and any other materials once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. The lancet has not to be reused.
2. Do not pipette by mouth.
3. Handle all samples as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
4. Wear disposable gloves and eye protection while handling specimens.
5. Wash hands thoroughly after using the cassette.
6. Consult the relative Material Safety Data Sheet (available on DIESE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
7. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the

contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

8. The used test should be discarded according to local regulations.
9. Do not eat, drink or smoke in the area where the specimens or kits are handled.

Analytical precautions

1. Do not use test if pouch is torn or damaged.
2. Do not freeze the product
3. Do not use expired test.
4. In case the pouch has been stored at 4-8°C, allow at least 10 minutes for the device to come to room temperature.
5. The test device should remain in the sealed pouch until use.
6. Perform the test quickly after opening the aluminum pouch.
7. The test must be placed on a flat surface while waiting for the results. The test should never be oriented upwards.
8. Read carefully the Instruction for use to guarantee the correct execution and interpretation of the test.
9. The test device should not be reused.
10. Please ensure that an appropriate amount of samples are used for testing. Too much or too little sample size may lead to deviation of results.
11. Humidity and temperature can adversely affect results.

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 20 tests.

DD DEVICES 20 packages each containing 1 test cassette
Each device is sealed in foil pouch with dessicant.

BUFFER BUFFER SOLUTION 1 x 3.0 ml
Contents: Saline buffer containing detergent and Sodium Azide (< 0.1%) in plastic dropper bottle
Liquid, ready for use (sufficient for 20 tests).

The kit contains also:

- 20 calibrated micropipettes (capillary tubes) capable delivering 20 µL accurately.
- 20 lancets (for finger prick whole blood)
- 20 alcohol prep swabs

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer
- Disposable gloves
- Specimen collection containers
- Calibrated micropipette for 10µL

6. STORAGE AND STABILITY OF REAGENTS

Store as packaged in the sealed pouch between 4 and 30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use.

The expiry date is printed on the components and on the kit label.

7. SPECIMEN COLLECTION

The sample is composed of:

- fresh serum or plasma collected with the normal venipuncture and handled with all precautions dictated by good laboratory practice.
- whole blood collected from fingerprick using accessories (lancet and micropipette for sample collection) provided in the kit.

Possible consequences, in case of use of other biological liquids, are not known.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

iRapid SARS-CoV-2 Quant "Neutralizing" Ab can be performed using whole blood (from fingerprick), serum or plasma.

Whole Blood – Capillary Collection:

1. Wear protective gloves.
2. Clean the patient's finger with the alcohol swab. Allow to dry.
3. Bring the kit components to room temperature before testing.
4. Open the pouch and remove the device. Once opened, the device must be used immediately.
5. Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
6. Open the lancet by completely unscrewing the blue cap (do not pull on it).
7. Prick the lateral side of the patient's fingertip to obtain a drop of blood (20 µL). Using gentle pressure, massage the finger towards the fingertip to encourage a drop of blood to form.
8. Hold the capillary micropipette (included in the kit) horizontally, and touch with the tip of the micropipette the blood sample.
9. To expel the sample, place the tip of the micropipette on the S well and squeeze the micropipette bulb. Transfer to the sample well must be immediate in order to avoid sample clotting.
10. Add 2 drops of buffer solution (2x40 µL) into the S well while holding the bottle vertically. Start the timer.
11. Interpret the results at 10 minutes.

Plasma & Serum (from a prepared sample in a laboratory):

1. Collect a 10 µL serum or plasma sample using a calibrated micropipette (not included in the kit). Only samples that are clean, clear and with good fluidity can be used for the assay (do not use samples that are viscous or have high lipid levels).
2. Dispense the sample on the S well. Add 2 drops of buffer solution (2x40 µL) into the same S well while holding the bottle vertically.
3. Interpret the results at 10 minutes.

9. TEST VALIDATION

Internal procedural control is included in the test (coloured line appearing in the control region (C)). It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is

recommended the use of known positive and negative controls as a good laboratory practice to confirm the test procedure and to verify proper test performance.

In case of not valid test, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

IgG POSITIVE

Two coloured lines appear: one coloured line appears in the control line region (C) and another line is in the test line region (T).

The colour intensity is proportional to the antibodies concentration (BAU/ml) as follows:

>1000 BAU/mL	
800 BAU/mL	
600 BAU/mL	
450 BAU/mL	
<300 BAU/mL	

IgG NEGATIVE

Only one coloured line appears in the control line region (C). No line appears in the test line region (T).

INVALID TEST

No coloured line appears at control line (C).

11. LIMITATIONS

1. The iRapid SARS-CoV-2 Quant "Neutralizing" Ab is for professional use only. This test should be used for detection of IgG antibodies to SARS-CoV-2 receptor binding domain in whole blood from fingerprick, plasma and serum.
2. The iRapid SARS-CoV-2 Quant "Neutralizing" Ab will only indicate the presence of IgG antibodies to SARS-CoV-2 RBD in the specimen and should not be used as the sole criteria for the diagnosis of SARS-CoV-2 infections.
3. As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
4. In case of suspected false negative result, retest a new collected sample or perform additional follow-up using other methods. A negative result at any time does not rule out the possibility of presence of IgG antibodies anti-SARS-CoV-2.

5. The test will show negative results under the following conditions: the concentration of the specific antibodies in the sample is lower than the minimum detection limit of the test or the specific antibodies have not appeared at the time of sample collection.

12. CROSS-REACTIONS

212 samples, positive to Adenovirus IgG (5), Adenovirus IgM (5), Antinuclear Antibodies (ANA) (10), Bordetella pertussis IgG (5), Borrelia IgG (10), Chlamydomphila pneumoniae IgG/IgA (10), Common Coronavirus NL63 IgG (10), Coronavirus 229E IgG (5), Coronavirus HKU1 IgG (5), Coronavirus OC43 IgG (5), Cytomegalovirus (CMV) IgG (10), Enterovirus IgG (1), Epstein-Barr Virus, Epstein-Barr Nuclear Antigen (EBV: EBNA) IgG (10), Haemophilus influenzae IgG (5), Hepatitis C Virus (HCV) Ig (10), Herpes Simplex Virus (HSV) IgG (10), HIV Ig (11), Influenza A IgG (5), Influenza B IgG (5), Legionella IgG (5), MERS IgG (5), Mycoplasma pneumoniae IgG/IgA (10), Mycoplasma pneumoniae IgM (5), Parainfluenza IgG (5), Plasmodium falciparum (5), Rheumatoid Factor (RF) IgG (10), RSV IgG (5), RSV IgM (5), Rubella IgG (5), Rubella IgM (5), SARS IgG (5), Toxoplasma IgG (5) were tested.

Possible cross-reactions can be observed with Rheumatoid Factor (RF) (2 samples on 10), Adenovirus IgM (1 sample on 5) Adenovirus IgG (1 sample on 5), Hepatitis C Virus (HCV) Ig (1 sample on 10), Antinuclear Antibodies (ANA) (1 sample on 10) and Plasmodium falciparum (1 sample on 5).

13. ANALYTICAL SPECIFICITY

2 samples (1 Negative and 1 Positive) were spiked with the following potentially interfering factors and then tested:

Bilirubin – 18 mg/dL
Triglycerides – 150 mg/dL
Hemoglobin – 15 mg/mL
Oximetatazoline – 1.25 mg/dL
Acetaminophen – 1.56 mg/dL
Ibuprofen – 2.19 mg/dL

No interferences were found.

14. METHOD COMPARISON

745 samples have been tested with iRapid SARS-CoV-2 Quant "Neutralizing" Ab:

- 515 pre-pandemic samples
- 230 positive (100 pandemic samples and 130 vaccinated samples). Pandemic samples were collected from 15 days to 6 months after symptoms onset. Vaccinated samples were collected from 3 days to 3 months after the end of vaccination.

Positive samples were tested in comparison to a reference test based on a different test, which detected neutralizing antibodies. 70% of positive samples, included the discordant ones, were also analyzed with the in-house seroneutralization test using UK variant.

Results are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	221	10	231
	-	9	505	514
	Total	230	515	745

Percent Positive Agreement (~Diagnostic Sensitivity):

96.1% CI_{95%}: 92.7-97.9

Percent Negative Agreement: (~Diagnostic Specificity):

98.1% CI_{95%}: 96.5-98.9

Positive Predictive Value: 95.7% CI_{95%}: 94.2-97.2

Negative Predictive Value: 98.2% CI_{95%}: 97.2 - 99.2

The agreement between the two methods is excellent with a Cohen's Kappa of 0.93.

15. LIMIT OF DETECTION

The limit of detection based on the WHO international standard 20/136 was found at 125 BAU/ml.

16. PRECISION AND REPEATABILITY

Accuracy and repeatability were evaluated on replicates of positive and negative samples, on three lots.

All samples always give comparable results, without interpretation difficulties, regardless of the lot used and the test replicates.

17. BIBLIOGRAPHY

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

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1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
 2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
 3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
 4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
- Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology -vol 42 n5 (2306-2309); May 2004

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