CHORUS

dsDNA-G

REF 86032

REF 86032/12

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INSTRUCTIONS FOR USE

CHORUS dsDNA-G

For the quantitative determination of IgG-class antibodies against dsDNA

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the quantitative determination of IgG-class antibodies against dsDNA in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Antibodies binding to DNA belong to the group of anti-nuclear antibodies (ANA) that have been observed in several autoimmune diseases. Antibodies reacting with native double-stranded (ds) DNA are regarded as being specific for systemic lupus erythematosus (SLE) and have been observed in approximately 50-80% of the patients.

Antibodies against dsDNA are found during active phases of SLE. The amount of the serum concentration is positively correlated with the severity of the disease. Thus, detection of these autoantibodies is important for the diagnosis and the clinical monitoring of SLE. Consequently it has been established as one of the 11 ACR-criteria for the diagnosis of SLE.

Most patients with SLE display IgG class antibodies against dsDNA. These autoantibodies are associated with lupus nephritis. Approximately 30% of the SLE patients develop IgA class anti-dsDNA antibodies, additionally. There have been suggestions that the presence of these IgA class anti-dsDNA antibodies may define a certain subset of SLE patients. Indeed studies demonstrated the association of this subclass with certain parameters of the disease activity, such as elevated erythrocyte sedimentation rate, or the consumption of complement component C3, as well as the clinical parameters of cutaneous vasculitis, acral necrosis and erythema, while no association was found for nephritis and arthritis.

IgM class anti-dsDNA antibodies were found in 52% of the sera from patients with SLE. In contrast to IgG and IgA class autoantibodies, the subclass IgM antibodies do not correlate with disease activity. However, a highly significant negative correlation between IgM anti-dsDNA antibodies and lupus nephritis, including its laboratory parameters was demonstrated. Therefore IgM class anti-dsDNA antibodies may indicate a subset of lupus patients being protected against the risk of developing nephritis.

3. PRINCIPLE OF THE METHOD

The dsDNA-G devices are ready to use for the assay of IgG-class antibodies against dsDNA in the CHORUS instruments. The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The antigen, highly purified, is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of monoclonal antibodies, highly specific for anti-human immunoglobulins conjugated to horse radish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

The Chorus dsDNA-G is calibrated against reference sera from the CDC (Atlanta). The results are expressed in IU/ml.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and found negative for both HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. All reagents and samples must be handled in accordance with the safety standards normally adopted in the laboratory.

Disposal of residues: the serum samples, calibrators and the strips used must be treated as infected residues, therefore disposed of in accordance with the provisions of current laws.

Personal safety warnings

1. Do not pipette with your mouth.
2. Use single-use gloves and eye protection when handling samples.
3. Wash your hands thoroughly after inserting the devices into the CHORUS TRIO instrument.
4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
5. Neutralised acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. Exposure to 1% hypochlorite sodium for 30 minutes should be sufficient to ensure effective disinfection.
6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% hypochlorite sodium, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried.

All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste.

Do not autoclave materials containing hypochlorite sodium.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) and use within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
2. When adding the sample to the well, check that it is perfectly distributed at the bottom.
3. Check the actual presence of reagents in the device and the integrity of the device itself. Do not use devices that lack any reagents and/or that have any foreign bodies in the reaction well when visually inspected.
4. The devices must be used in conjunction with the CHORUS TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.

Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a Release (date) above the one shown in the table published on the Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)

5. Check that the CHORUS TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode placed on the handle of the device in order to allow it to be read correctly by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipemic, icteric, non-fully coagulated serum samples or of samples with microbial contamination may be a source of errors.
11. Do not use the device after its expiry date
12. **Check that the instrument is connected to the Washing Buffer Autoimmunity REF 86004.**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86032).


The kit is sufficient for 12 tests. (REF 86032/12).

DD DEVICES

6 packages each containing 6 devices (REF 86032).

2 packages each containing 6 devices (REF 86032/12).

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with highly purified dsDNA

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing Tween-20 0.2% and Proclin 0.1%.

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labeled with horse radish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.175 ml

Contents: Protein solution containing specific antibodies capable of binding the microplate and preservative. Liquid, ready to use.

CONTROL + POSITIVE CONTROL

1 x 0.425 ml

Contents: Protein solution containing specific antibodies capable of binding the microplate and preservative. Liquid, ready to use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C

POSITIVE CONTROL 8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples cannot be used.

The test cannot be applied to plasma.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result in IU/ml (WHO/80) calculated on the basis of a batch-dependent curve stored in the instrument.

The test serum can be interpreted as follows:

POSITIVE: when the result is > 30.0 IU/ml

NEGATIVE: when the result is < 20.0 IU/ml

DOUBTFUL/EQUIVOCAL: for all values between 20.0 and 30.0 IU/ml

In the case of a doubtful result, repeat the test. If the test remains doubtful, collect a new sample.

11. LIMITATIONS

All the positive results require a careful interpretation.

The test cannot be used as only method for a clinical diagnosis.

Negative results may not exclude an eventual infection.

The test results should be interpreted in conjunction with information available from the clinical evaluation and other diagnostic procedures.

12. CALIBRATION RANGE

Calibration range: 10.0-150.0 IU/ml.

For samples > 150.0 IU/ml repeat the test pre-diluting the sample with Negative Control/Sample Diluent (PF83607 – **not supplied with the kit**).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 10.0 and 13.2 IU/ml.

14. ANALYTICAL SPECIFICITY

4 samples (2 negative and 2 Positive) containing the following interfering substances were tested:

Rheumatoid Factor (110 UI/ml)
Bilirubin (45 mg/dl)
Triglycerides (1500 mg/dl)
Hemoglobin (10 mg/ml)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test.

15. CROSS-REACTIONS

28 Samples positive to tTG-A, tTG-G, Gliadin-G, MPO, PR3, CCP, RF-M and ASCA were tested.

No significant cross-reactions were found.

16. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an experimentation 157 samples were tested with the Diesse kit and with another commercial kit.

The results are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	39	3	42
	-	1	114	115
	Total	40	117	157

Percent Positive Agreement (~Diagnostic Sensitivity):

97.5% CI_{95%}: 87.0-99.5

Percent Negative Agreement: (~Diagnostic

Specificity): 97.4% CI_{95%}: 92.7-99.1

Positive Predictive Value (PPV): 92.9% IC_{95%}: 88.9-96.9

Negative Predictive Value (NPV): 99.1% IC_{95%}: 97.6-100.0

17. PRECISION

Sample	Within-run Precision		Between-run precision	
	Mean (IU/ml)	CV%	Mean (IU/ml)	CV%
1	12.8	14.0	11.1	10.5
2	66.9	12.8	77.7	10.4
3	107.4	12.9	115.2	7.3
4	136.4	12.3	135.8	10.2
5	139.0	6.8	107.6	13.8
6	11.8	10.6	11.2	8.5

Sample	Precision between batches		Precision between instruments	
	Mean (IU/ml)	CV%	Mean (IU/ml)	CV%
1	11.4	9.8	11.4	12.8
2	66.8	12.3	72.0	13.2
3	120.4	12.6	119.9	11.9
4	140.8	5.3	141.5	6.8
5	121.5	13.2	114.0	13.4
6	12.2	8.4	11.2	7.2

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