CHORUS COPROCOLLECT

REF 86602

DIESSE

REF 86602/40

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INSTRUCTIONS FOR USE

CHORUS Coprocollect

**Device for stools collection, for use with the
CHORUS Calprotectin test**

For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

Single-use device for manual collection of the amount of faeces required to perform the CHORUS Calprotectin test (REF. 86600-86600/C03-86601).

It must only be used by professional laboratory personnel.

2. INTRODUCTION

The device for stools collection allows to collect the amount necessary to perform the CHORUS Calprotectin test directly from the primary container without weighing the sample.

3. PRINCIPLE OF THE METHOD

The device consists of a tube containing the extraction solution and a stick shaped with seven radial grooves for the sample collection. The upper end of the device consists of two removable parts having opposite rotations. The screw cap (white), connected to the stick, is removed by turning it in anti-clockwise direction; the funnel (green), used to retain the coarsest stool impurities, can be removed with the white screw cap by turning it clockwise.

After performing this procedure, the sample can be transferred into the collection tube.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

Wait until stools samples reach room temperature (18-30°C) before collecting the samples.

The product, once used or expired or in case of deterioration, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. All material of human origin must be handled as potentially infectious. All the samples must be handled according to the safety precautions normally adopted in laboratory practice.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after use.

4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

1. Place the unused devices in the primary package.
2. Make sure that the product packaging is intact.
3. Verify that the device is correctly closed.
4. Do not use broken or damaged devices.
5. Do not use the device in case of liquid leakage.
6. Do not use the devices after the expiry date.
7. Adopt all the precautions to avoid contamination of the product during its use.

5. COMPOSITION

The package is sufficient for 20 determinations (REF 86602).

The package is sufficient for 40 determinations (REF 86602/40).

MATERIALS REQUIRED BUT NOT PROVIDED

- Vortex
- Shaker for tubes
- Centrifugal machine
- Micropipettes for the accurate collection of 50-100 µl volumes
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on the external label of the package.

7. SPECIMEN COLLECTION AND STORAGE

1. The collection device containing the stools sample must be stored at 2-8°C and processed within 48 hours after collection (see chapter 8 Procedure).
2. In case of transfer of the device containing the stools sample from the collection center to the laboratory, the sample needs to be stored in a refrigerated thermal container only if the delivery requires more than 60 minutes.
3. The faecal extract may be stored for a maximum of 4 days at 2/8°C or frozen at -20°C in accordance with the Instructions for Use of the CHORUS Calprotectin kit.
4. In case of temporary unavailability of the device, maintain the stools sample in a sterile primary container at 2-8°C for a maximum of 7 days. For longer periods freeze at -20°C.

8. PROCEDURE

Procedure for stools collection

Sampling procedure for faeces of normal consistency:

1. To collect a stools sample pull out the shaped stick by turning the white screw-cap in anti-clockwise direction.
2. Dip the shaped stick in the stools sample.
3. Repeatedly rotate the stick in the stools sample until all the grooves are filled with faecal material.
4. Before inserting the stick into the device, remove the faecal material in excess by rotating the stick on the internal wall of the primary container.
5. Fully insert the shaped stick, with faecal material, into the tube containing the extraction solution and turn the cap clockwise until it is completely closed, keeping blocked the funnel (green) to avoid its opening.

Sampling procedure for faeces of hard consistency (the shaped stick can hardly be dipped into the sample):

1. Pipette 50-100 µl of saline solution inside the stools container.
2. Keep the container with the stools and the saline solution at room temperature for 60 minutes at least.
3. Proceed following the instructions for "Sampling procedure for faeces of normal consistency" from point 2.

Sampling procedure for liquid faeces or anyhow having a consistency that does not allow the use of the shaped stick:

1. A laboratory pipette is necessary.
2. By using the pipette, aspirate 56 µl of liquid faeces.
3. Pull out the shaped stick and the funnel (green) by turning the funnel clockwise.
4. Pipette the liquid faeces sample directly into the collection device.
5. Fully insert the shaped stick back into the tube and turn the funnel (green) in anti-clockwise direction, until it is completely closed.

Sampling procedure for frozen faeces:

1. In case the sample is frozen, let the sample reach room temperature before collecting it with the device.
2. Proceed following the instructions for "Sampling procedure for faeces of normal consistency" from point 2.

Laboratory procedure

1. Follow the procedure described below for all the samples, whichever collection procedure has been used.
2. Shake the device by vortex for 60 seconds in order to properly homogenate the content.
3. Place the device on a roller shaker for tubes and shake for 20 minutes.
4. Repeat steps 2 and 3 for all the samples to be processed.
5. After the stirring, remove the screw-cap (white) and the funnel (green) by turning the funnel clockwise.

6. Close the tube with the extra screw-cap (present in the package) and centrifugate it for at least 10 minutes at 5000 rpm.
7. After the centrifuge, transfer the supernatant into a clean tube making sure to not touch the pellet on the bottom of the tube.
8. The faecal extract is now ready to be tested following the Instructions for Use of the CHORUS Calprotectin kits (REF 86600-86600/C03-86601).

In case of any problem contact the Scientific Support:

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

9. INCIDENT REPORTING

If a serious incident has occurred in relation to this device in the European Union market territory, please report it to the manufacturer and the competent authority of its Member State without delay.

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