SEDICHECK CUBE 200/80

REF 10438

DIESSE
DIESSE

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SEDICHECK CUBE 200/80

REF 10438 4 x 3 ml (2 x Level I + 2 x Level II)

For the precision monitoring of the VES-Matic Cube 200, VES-Matic Cube 80 and VES-Matic 5 instruments

For in vitro diagnostic use only

1. INTENDED PURPOSE

Double level control material to be used for the precision monitoring of the VES-Matic Cube 200, VES-Matic Cube 80 and VES-Matic 5 instruments, used for the determination of Erythrocyte Sedimentation Rate (ESR).

It must be used by professional laboratory users only.

2. INTRODUCTION AND PRINCIPLE OF THE METHOD

It is good laboratory practice to use materials of control to check the precision of the methods used.

The product SediCheck Cube 200/80 is a material of control to be used to monitor the precision of the Ves-Matic Cube 200, Ves-Matic Cube 80 and VES-Matic 5 instruments.

3. WARNINGS

1. Check that the batch number on the tube corresponds to the lot number reported on the Certificate of Analysis (CoA) enclosed with the present Instructions for Use.
2. DO NOT FREEZE AND DO NOT EXPOSE TO TEMPERATURES OVER 30°C.
3. Good function of the product is guaranteed only if it is stored and used as described in the present Instructions for Use.
4. To resuspend the product turn the tube upside down repeatedly.

4. REAGENTS

Contents:

Level I - stabilized suspension of photo-absorbing particles.

Preservative: glycerol > 40%.

Level II - stabilized suspension of photo-absorbing particles.

Preservative: ethanol 20%.

Preparation: ready for use

Stability: Store at 2 - 30°C.

The unopened product is stable up to the expiry date.

5. PROCEDURE

1. During first use turn the tube upside down repeatedly until complete resuspension. Check that the bottom

and the walls of the tube do not contain residuals of the product that are not resuspended. Keep shaking for 2 minutes at least. DO NOT USE THE VORTEX.

2. Insert the tube in the instrument like any other sample (see the operating manual of the instrument).
3. Compare the result found with the value reported on the CoA. If the result obtained is within the range reported on the CoA and on the barcode the instrument will print the data without any alarm. Out of range values will be reported with an asterisk.
4. Do not use the same tube for more than 14 consecutive days.

6. EXPECTED RESULTS

If the results of the control are not within the range reported on the CoA but are reproducible, the relative working conditions may be the cause of the values not falling within the given range. In this case it is advisable to perform several runs following the laboratory's usual method and thus determine the relative acceptable range.

If the results of the control are out of the range reported in the CoA and are not reproducible, it is advisable to repeat the control using a new batch.

It must be remembered that SediCheck Cube 200/80 is used to check the precision of the system.

7. LIMITATIONS

SediCheck Cube 200/80 is not intended for use as a reference standard.

SediCheck Cube 200/80 cannot be used to control any hematological parameters other than the ESR.

8. REFERENCES

1. D. Chesher Evaluating assay precision. Clin Biochem Rev Vol 29 Suppl (i) August 2008; S23-S26
2. J. Ezzelle, I. R. Rodriguez-Chavez, J. M. Dardenb, M. Stirewalta, N. Kunwarb, R. Hitchcocka, T. Waltera, and M. P. D'Souzac Guidelines on Good Clinical Laboratory Practice: Bridging Operations between Research and Clinical Research Laboratories. J Pharm Biomed Anal. 2008 January 7; 46(1): 18-29.

9. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

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