CHORUS
WEST NILE VIRUS
IgG

DIESSE

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INSTRUCTIONS FOR USE

CHORUS WEST NILE VIRUS IgG

For the qualitative determination of anti-West Nile Virus IgG antibodies

For in vitro diagnostic use only

1. USE

Enzyme immunoassay method for the qualitative determination of anti-West Nile Virus (WNV) IgG class antibodies in human serum with a single-use device applied to Chorus and Chorus TRIO instruments.

2. INTRODUCTION

The West Nile virus belongs to the Flaviviridae family. The virus reservoirs are wild birds and mosquitoes (most frequently type Culex), whose bites are the main means of transmission to humans.

The incubation period from the time of the bite of the infected mosquito varies between 2 and 14 days, but may also be 21 days in those with immune deficiency.

West Nile infections are asymptomatic in 80% of people but may cause fever, headache and nausea; in some cases meningitis and/or encephalitis may occur, which can cause serious brain damage or death. There is no specific therapy for West Nile fever. In most cases, symptoms go away on their own after a few days or may last for a few weeks. In the most severe cases, hospitalization is required, where treatments given include intravenous fluids and assisted breathing.

Diagnosis is mainly performed by serological testing for the determination of WNV-specific IgG and IgM antibodies. Alternatively, it may also be performed by PCR or viral culture on serum and cerebrospinal fluid specimens.

3. PRINCIPLE OF THE METHOD

The Chorus West Nile Virus IgG device is ready to use for the determination of anti-WNV IgG antibodies in Chorus/Chorus TRIO instruments.

The test is based on the ELISA (Enzyme Linked ImmunoSorbent Assay) principle. The antigen is bound to the solid phase. Specific immunoglobulins bind to the antigen following incubation with diluted human serum. After washing to remove proteins that have not reacted, incubation is carried out with the conjugate consisting of human anti-IgG monoclonal antibodies conjugated with horseradish peroxidase. The unbound conjugate is removed and the substrate for the peroxidase is added. The blue colour that develops is proportional to the concentration of specific antibodies in the serum under examination.

The single-use devices contain all the reagents to perform the test when applied to the Chorus/Chorus TRIO instruments. The results are expressed in Index (sample OD/cut-off OD).

4. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

This kit contains materials of human origin which have been tested with FDA-approved tests and found negative for both HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. All reagents and samples must be handled in accordance with the safety standards normally adopted in the laboratory.

Disposal of residues: the serum samples, calibrators and strips used must be treated as infected residues, therefore disposed of in accordance with the provisions of current laws.

Personal safety warnings

1. Do not pipette with your mouth.
 2. Use single-use gloves and eye protection when handling samples.
 3. Wash your hands thoroughly after inserting the devices into the Chorus/Chorus TRIO instrument.
 4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
 5. Neutralised acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. Exposure to 1% sodium hypochlorite for 30 minutes should be sufficient to ensure effective disinfection.
 6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% sodium hypochlorite, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried.
- All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste.
- Do not autoclave materials containing sodium hypochlorite.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) and use within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
 2. When adding the sample to the well, check that it is perfectly distributed at the bottom.
 3. Check the actual presence of the reagents in the device and the integrity of the device itself. Do not use devices that lack any reagents and/or that have any foreign bodies in the reaction well when visually inspected.
 4. The devices must be used in conjunction with the Chorus/Chorus TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.
- Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a release (Rel.) date later than the one shown in the table published on the Diesse website**

(<http://www.diesse.it/it/Support/Download/strumento:39/>)


5. Check that the Chorus/Chorus TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode placed on the handle of the device in order to allow it to be read correctly by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipemic, icteric, non-fully coagulated serum samples or of samples with microbial contamination may be a source of errors.
11. Do not use the device after its expiry date.
12. **Check that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests

DD DEVICES packs of 6 devices each

Description:

Position 8: Space available for barcode label

Position 7: Empty

Position 6: MICROPLATE WELL

Sensitized with WNV NS1 recombinant protein antigen **Position**

5: MICROPLATE WELL

Not sensitized.

Position 4: TMB SUBSTRATE

Content: 0.26 mg/mL tetramethylbenzidine and 0.01% H₂O₂ stabilized in 0.05 mol/L citrate buffer (pH 3.8).

Position 3: DILUENT FOR SAMPLES

Content: Protein solution containing 0.05% phenol, 0.02% Bronidox and an indicator to detect the presence of serum.

Position 2: CONJUGATE

Content: Peroxidase-labelled human anti-IgG monoclonal antibodies in phosphate buffer solution containing 0.05% phenol and 0.02% Bronidox.

Position 1: EMPTY WELL

Where the sample is transferred.

Usage: bring a bag to room temperature, open the bag, take out the necessary devices; store the others in the bag containing the silica gel, let the air out and **seal** by pressing on the closure. Store at 2-8°C.

CALIBRATOR 1 x 0.175 ml

Content: Buffer containing a preservative and specific antibodies capable of binding the antigen present on the microplate. Liquid, ready to use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Content: Buffer containing a preservative and specific antibodies capable of binding the antigen present on the microplate. Liquid, ready to use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED:

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus/Chorus TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Micropipettes capable of accurately withdrawing volumes of 50-200µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials

6. STORAGE METHOD AND REAGENT STABILITY

The reagents must be stored at 2-8°C. In the event of an incorrect storage temperature, the calibration must be repeated and the correctness of the result checked using the control serum (see section 9: Test validation).

The expiry date is printed on each component and on the outer label of the package.

The reagents have limited stability after opening and/or preparation:

DEVICES	8 weeks at 2-8°C
CALIBRATOR	8 weeks at 2-8°C
POSITIVE CONTROL	8 weeks at 2-8°C

7. SAMPLE TYPE AND STORAGE

The sample type is serum obtained from blood collected by normal venous blood sampling and handled as required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

Fresh serum can be kept for 5 days at 2-8°C; for longer storage periods, freeze at -20°C.

The sample can undergo up to a maximum of 3 thaws.

Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

Heat inactivation may give erroneous results. Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in section 4 Analytical Warnings.
3. Dispense into well no. 1 of each device:
 - SAMPLE: 50 µl/device
 - CALIBRRATOR: 50 µl/device
 - POSITIVE CONTROL: 50 µl/device

At each batch change use a device for the calibrator.

4. Introduce the devices on the Chorus/Chorus TRIO instrument. Perform the calibration (if required) and the test as indicated in the User Manual of the instrument.

9. TEST VALIDATION

Use the positive control serum to verify the correctness of the result obtained, processing it as indicated in the User Manual of the instrument. If the instrument indicates that the control serum has a value outside the acceptability limit, the calibration must be performed again. Previous results will be corrected automatically.

If the control serum result continues to be outside the acceptability range, contact Scientific Support.

Phone: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

The Chorus/Chorus TRIO instrument provides the result in Index (sample OD/cut-off OD).

The test on the serum under examination can be interpreted as follows:

POSITIVE: when the result is >1.1

NEGATIVE: when the result is <0.9

DOUBTFUL/UNCERTAIN: when the result is between 0.9 and 1.1

In case of a doubtful/uncertain result, repeat the test. If the result remains doubtful/uncertain, repeat the sampling.

11. TEST LIMITATIONS

All values obtained require a careful interpretation, without ever disregarding other indicators relating to the same patient.

The test, in fact, cannot be used on its own for a clinical diagnosis and the result obtained must always be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. ANALYTICAL SPECIFICITY

5 samples were tested (2 Negative, 1 at Cut-Off and 2 Positive) to which the following interfering agents were added:

Rheumatoid factor (44 IU/ml - 220 IU/ml)

Bilirubin (4.5 mg/dl - 45 mg/dl)

Lipoprotein (10 mg/dl - 250 mg/dl)

Haemoglobin (5 mg/ml - 30 mg/ml)

The presence of the above interfering substances in the test serum does not alter the test result.

13. CROSS-REACTIVITY

21 samples positive for Dengue (9) and Zika (12) were tested.

Cross-reactions with Dengue and Zika were detected.

14. COMPARISON STUDIES

In one trial, 172 samples were analysed with a Diesse kit and with another commercially available kit.

The trial data are outlined below:

		Reference		Total
		+	-	
Diesse	+	69	2	71
	-	0	101	101
	Total	69	103	172

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% CI_{95%}: 94.7-99.9

Percent Negative Agreement (~Diagnostic Specificity):

98.1% CI_{95%}: 93.2-99.4

The degree of agreement between the two methods is excellent with a K value (Cohen's coefficient) of 0.98.

15. ACCURACY AND REPEATABILITY

Sample	Within the session		Between sessions	
	Mean (Index)	CV%	Mean (Index)	CV%
1	1.0	6.0	1.0	5.0
2	1.4	3.6	1.4	9.3
3	2.0	1.0	1.9	8.9
4	0.6	8.3	0.5	12.0
5	0.4	5.0	0.4	7.5

Sample	Between batches		Between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	1.0	7.0	1.2	14.2
2	1.4	5.0	1.6	13.1
3	1.9	4.7	2.2	15.0
4	0.6	5.0	0.7	14.3
5	0.5	6.0	0.6	15.0

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- Epidemiology and Transmission Dynamics of West Nile Virus Disease
Emerg Infect Dis. 2005 Aug; 11(8): 1167-1173
2. Centers for Disease Control and Prevention
Outbreak of West Nile-Like Viral Encephalitis – New York, 1999
 3. E. McDonald, S. Mathis, S.W. Martin, J.E. Staples, M. Fischer, N.P. Lindsey
Surveillance for West Nile Virus Disease – United States, 2009-2018
MMWR Surveill Summ. 2021 Mar 5;70(1):1-15
 4. Y. Lustig, D. Sofer, E.D. Bucris, E. Mendelson
Surveillance and Diagnosis of West Nile Virus in the Face of Flavivirus Cross-Reactivity
Front Microbiol. 2018 Oct; 9: 2421

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Percent Positive Agreement (~Sensibilitatea diagnosticului):

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Epidemiology and Transmission Dynamics of West Nile Virus Disease

Emerg Infect Dis. 2005 Aug; 11(8): 1167-1173

18. Centers for Disease Control and Prevention

Outbreak of West Nile-Like Viral Encephalitis – New York, 1999

19. E. McDonald, S. Mathis, S.W. Martin, J.E. Staples, M. Fischer, N.P. Lindsey

Surveillance for West Nile Virus Disease – United States, 2009-2018

MMWR Surveill Summ. 2021 Mar 5;70(1):1-15

20. Y. Lustig, D. Sofer, E.D. Bucris, E. Mendelson
Surveillance and Diagnosis of West Nile Virus in the Face of
Flavivirus Cross-Reactivity
Front Microbiol. 2018 Oct; 9: 2421

	IT	Data di fabbricazione	ES	Fecha de fabricación
	EN	Date of manufacture	FR	Date de fabrication
	CS	Datum výroby	PT	Data de fabrico
	DE	Herstellungsdatum	RO	Data fabricatiei
	EL	Ημερομηνία Παραγωγής		
	IT	Utilizzare entro	ES	Fecha de caducidad
	EN	Use By	FR	Utiliser jusqu'à
	CZ	Použitelné do	PT	Prazo de validade
	DE	Verwendbar bis	RO	A se folosi pana la
	EL	Ημερομηνία λήξης		
	IT	Non riutilizzare	ES	No reutilizar
	EN	Do not reuse	FR	Ne pas réutiliser
	CS	Nepoužívejte opakovaně	PT	Não reutilizar
	DE	Nicht wieder verwenden	RO	A nu se refolosi
	EL	Μην κάνετε επαναληπτική χρήση		
	IT	Attenzione, vedere le istruzioni per l'uso	ES	Atención, ver instrucciones de uso
	EN	Caution, consult accompanying documents	FR	Attention voir notice d'instructions
	CZ	Pozor, čtěte příložené dokumenty	PT	Atenção, consulte a documentação incluída
	DE	Achtung, die Gebrauchsanleitung lesen	RO	Atentie, consultați documentele însoțitoare
	EL	Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα		
	IT	Fabbricante	ES	Fabricante
	EN	Manufacturer	FR	Fabricant
	CS	Výrobce	PT	Fabricante
	DE	Hersteller	RO	Productator
	EL	Κατασκευαστής		
	IT	Contenuto sufficiente per "n" saggi	ES	Contenido suficiente para <n> ensayos
	EN	Contains sufficient for <n> tests	FR	Contenu suffisant pour "n" tests
	CZ	Obsah stačí na <n> testů	PT	Conteúdo suficiente para "n" ensaios
	DE	Inhalt reicht für „n“ Tests	RO	Continut suficient pt <n> teste
	EL	Περιεχόμενο επαρκές για «n» εξετάσεις		
	IT	Limiti di temperatura	ES	Límite de temperatura
	EN	Temperature limitation	FR	Limites de température
	CS	Teplotní omezení	PT	Limites de temperatura
	DE	Temperaturgrenzwerte	RO	Limita da temperatura
	EL	Περιορισμοί θερμοκρασίας		
	IT	Consultare le istruzioni per l'uso	ES	Consulte las instrucciones de uso
	EN	Consult Instructions for Use	FR	Consulter les instructions d'utilisation
	CS	Čtěte návod k použití	PT	Consulte as instruções de utilização
	DE	Die Gebrauchsanleitung lesen	RO	Pentru utilizarea consultați instrucțiunile
	EL	Συμβουλευτείτε τις οδηγίες χρήσης		
	IT	Rischio biologico	ES	Riesgo biológico
	EN	Biological risks	FR	Risques biologiques
	CS	Biologická rizika	PT	Risco biológico
	DE	Biologisches Risiko	RO	Risk biologic
	EL	Βιολογικοί κίνδυνοι		
REF	IT	Numero di catalogo	ES	Número de catálogo
	EN	Catalogue number	FR	Référence du catalogue
	CS	Katalogové číslo	PT	Referência de catálogo
	DE	Katalognummer	RO	Numar de catalog
	EL	Αριθμός καταλόγου		
IVD	IT	Dispositivo medico-diagnostico in vitro	ES	Producto sanitario para diagnóstico in vitro
	EN	In Vitro Diagnostic Medical Device	FR	Dispositif médical de diagnostic in vitro
	CS	Lékařské vybavení pro diagnostiku in vitro	PT	Dispositivo médico para diagnóstico in vitro
	DE	Medizinisches In-vitro-Diagnostikum	RO	Dizpozitiv medical pentru diagnosticare in vitro
	EL	Ιn Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν		
LOT	IT	Codice del lotto	ES	Código de lote
	EN	Batch code	FR	Code du lot
	CS	Kód šarže	PT	Código do lote
	DE	Chargennummer	RO	Lot
	EL	Αριθμός Παρτίδας		
CE	IT	Marcatura CE di conformità	ES	Marcado CE de conformidad
	EN	CE marking of conformity	FR	Marquage de conformité CE
	CS	Označení shody CE	PT	Marcação CE de conformidade
	DE	CE-Konformität Skennezeichnung	RO	Marcajul de conformitate CE
	EL	Σήμανση συμμορφωσης CE		