**CHORUS
25OH Vitamin D
Total**

REF 86900

DIESSE

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18. BIBLIOGRAFIA

1. Lappe J.M. et All., The Role of Vitamin D in Human Health: A Paradigm Shift. Journal of Evidence-Based Complementary & Alternative Medicine (2010); 16(1)

2. Su Z. et All., 25-Hydroxyvitamin D: Analysis and clinical application. Clinica Chimica Acta 433 (2014) 200–205
3. Matyjaszek-Matuszek B. et All., Clinical implications of vitamin D deficiency. Prz Menopauzalny (2015); 14(2): 75-81
4. Gisondi P., Romanelli M., Dini V., Prignano F., Vinceti M., Pellacani G., Girolimoni G. Deficit di Vitamin D in dermatologia – Linee guida e raccomandazioni SIdMaST (2015).

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For In Vitro Diagnostic Use Only

1. INTENDED USE

The product Chorus 25OH Vitamin D Total (REF 86900) is an immunoassay kit for automated quantitative detection of 25(OH) Vitamin D, performed in human serum using a disposable device applied on the Chorus and Chorus TRIO instruments.

Since 25(OH) Vitamin D is considered the main form of storage of Vitamin D, the kit is intended to monitor the overall Vitamin D status and the relative therapeutic measures.

It must be used by professional laboratory users only.

2. INTRODUCTION

Vitamin D is a fat-soluble vitamin that is involved in calcium-phosphorus metabolism. There are three ways to maintain an adequate Vitamin D status: sunshine exposure, dietary intake and pharmaceutical supplementation.

There are two types of physiologically important Vitamins D: Cholecalciferol (Vitamin D3), mainly synthesized in the skin upon exposure to ultraviolet B (UV-B) radiation and Ergocalciferol (Vitamin D2), obtained from diet.

Both the types of Vitamin D are transported in the bloodstream conjugated to Vitamin D-binding proteins (VDBP).

In the liver, Vitamin D is metabolized to 25-hydroxyvitamin D-25(OH)D; then circulates to the kidneys where it undergoes a second hydroxylation to the biologically active form of Vitamin D: 1.25(OH)2D (calcitriol).

Vitamin D plays an important role in bone growth and bone remodeling and rickets, osteomalacia and osteoporosis are common manifestations of Vitamin D deficiency. Moreover Vitamin D deficiency may be associated with increased risk of a wide range of diseases, such as cardiovascular diseases, specific infections (active tuberculosis, upper respiratory tract infection and influenza), cancer and diabetes.

25(OH)D is the most plentiful and stable metabolite of vitamin D in the human bloodstream, so it represents the functional indicator of vitamin D status.

Since the product Chorus 25OH Vitamin D Total is intended to detect the concentration of 25(OH) Vitamin D in patient serum, it can be used to monitor the overall Vitamin D status and the relative therapeutic measures.

3. PRINCIPLE OF THE METHOD

The Chorus 25OH Vitamin D Total device is ready to use for the detection of 25(OH) Vitamin D, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA (Enzyme linked ImmunoSorbent Assay) competition principle.

The monoclonal antibody is bound to the solid phase and 25(OH) Vitamin D, present in the sample, competes with biotinylated Vitamin D.

The higher the concentration of 25(OH) Vitamin D in the sample, the lower the amount of biotinylated Vitamin D that binds.

After washing to eliminate the components which have not reacted, incubation is performed with the conjugate, composed of Streptavidin conjugated to horseradish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The colour which develops is inversely proportional to the concentration of 25 (OH) Vitamin D present in the serum sample. The competition reaction is performed using an extraction buffer (Incubation buffer), required to release 25 (OH) Vitamin D from VDBP.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in ng/ml.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on DIESSE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the device within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (Errore. Riferimento a collegamento ipertestuale non valido.)


5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested (according to the indications reported in chapter "Analytical specificity")
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: MICROPLATE WELL with BIOTINYLATED VITAMIN D (final concentration: 0.55 µg/ml)

Position 6: MICROPLATE WELL

Coated with anti-Vitamin D monoclonal antibodies (final concentration: 3 µg/well)

Position 5: MICROPLATE WELL

Uncoated

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine and H₂O₂ stabilized in buffer.

Position 3: INCUBATION BUFFER

Contents: Proteic solution (Bovine Casein: 5 g/l) containing preservative (Proclin 300: 0.1%)

Position 2: CONJUGATE

Contents: Streptavidin labelled with horseradish peroxidase (maximum concentration: 0.3 µg/ml), in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the patient sample has to be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

Contents: 100% serum containing Vitamin D, Proclin 300 0.2% and Gentamicin 0.05%. Freeze-dried.

Use: reconstitute with the volume of deionised water reported on the label, wait at least 5 minutes without shaking; then mix by inversion avoiding froth formation.

CONTROL CONTROL SERUM

Contents: 100% serum containing Vitamin D, Proclin 300 0.2% and Gentamicin 0.05%. Freeze-dried.

Use: reconstitute with the volume of deionised water reported on the label, wait at least 5 minutes without shaking; then mix by inversion avoiding froth formation.

Calibrator and Positive control exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared using reference materials previously calibrated to the ID-LC-MS/MS reference method, as encompassed in the Vitamin D Standardization Program initiative.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus Instruments
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at -20°C
CONTROL SERUM	8 weeks at -20°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum from whole blood collected by venipuncture and handled following standard laboratory techniques.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored up to the 7th day at 2/8°C, or frozen for longer periods at -20°C (at least 23 months) and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers.

Heat-inactivation can rise to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. At each change of batch, perform the calibration.
4. Dispense the following volumes in well no. 1 of each device:

SAMPLE	90 µl/device
CALIBRATOR	90 µl/device
POSITIVE CONTROL	90 µl/device

5. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum as reported in the Instrument Operating Manual to check the validity of the obtained results.

The acceptability range is reported on the enclosed Quality control certificate.

If the result of the positive control is out of the acceptability range, the calibration must be repeated.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in ng/ml, calculated on the basis of a lot-dependent graph stored in the instrument.

Attention: do not consider the qualitative report N; D/E; P on the instrument.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range: 4.0-115.0 ng/ml.

13. REFERENCE RANGE

The table below reports the expected values of 25(OH) Vitamin D, as reported in literature.

State	Range (ng/ml)
Lack	< 20
Insufficient	20 – 30
Normal	30 – 100
Toxic	> 100

The values observed in the population, determined by examining 399 sera, were between 4.0 and 46.9 ng/ml, with a median of 14.2 ng/ml.

14. ANALYTICAL SPECIFICITY

3 samples were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (44 IU/ml - 220 IU/ml)

Bilirubin (4.5 mg/dl - 45 mg/dl)

Triglycerides (10 mg/dl - 250 mg/dl)

Hemoglobin (5 mg/ml - 30 mg/ml)

The presence in the serum sample of the interfering substances described above (except for hemoglobin) does not affect the test result.

15. CROSS-REACTIONS

7 samples were tested adding Vitamin D2 and Vitamin D3.

No significant cross-reactions were found for Vitamin D2 up to 1000 ng/ml and Vitamin D3 up to 125 ng/ml.

16. METHOD COMPARISON

In an experimentation 217 samples have been tested with Diesse kit and with a competitor kit.

Results:

Correlation	r	95%CI
Pearson	0.92	0.89-0.94
Spearman	0.92	0.90-0.94

The correlation between the two methods is very high.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	17.5	5.8	17.3	8.3
2	23.4	7.6	27.2	7.0
3	27.6	8.6	27.9	6.5
4	31.4	3.5	33.4	4.7
5	33.2	3.1	35.4	4.2
6	37.2	4.0	37.5	2.7
7	39.4	3.8	44.9	12.6
8	38.3	1.4	41.7	10.5
9	78.1	3.1	77.9	5.3

Sample	Between lots		Between Instruments	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	12.9	14.4	12.9	3.9
2	22.5	9.4	22.5	3.6
3	21.9	8.4	21.9	2.3
4	26.4	10.3	26.4	4.6
5	30.9	4.0	30.9	2.3
6	32.6	4.5	32.6	2.2
7	36.1	3.4	36.1	2.5
8	35.8	2.7	35.8	4.8
9	73.4	7.1	73.4	9.3

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19. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

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