CHORUS

**CHLAMYDOPHILA
PNEUMONIAE IgA**

REF 81252

DIESSE

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

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INSTRUCTIONS FOR USE

CHORUS CHLAMYDOPHILA PNEUMONIAE IgA

For the semiquantitative determination of IgA antibodies anti- Chlamydomphila pneumoniae

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the semiquantitative determination of IgA class antibodies against Chlamydomphila pneumoniae in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Chlamydiae are gram-negative bacteria, with an obligatory intracellular energetic parasitism.

Chlamydomphila pneumoniae is an infective agent which causes different respiratory diseases.

The Chlamydomphila pneumoniae infections usually evolve in an asymptomatic way or with mild symptoms; however it is possible to encounter more severe diseases, like community acquired pneumonia, bronchitis, sinusitis and pharyngitis.

The disease can have also chronic character. Possible associations with asthma, lung cancer, atherosclerosis, myocardial infarction and Alzheimer's disease have been described. Infections can be primary, chronic or recurrent.

The primary infection leads to IgM class antibodies response (within about 3 weeks) and to a delayed response of IgA and IgG class antibodies (about 6 weeks). IgM and IgA antibodies titres rapidly decrease (within 4-6 months), while the IgG levels decrease slower. During chronic or recurrent infections mostly IgA and IgG antibodies rapidly increase; a high titre of IgA-antibodies is a good indication of a chronic infection; presence of IgG indicate a past infection and is detectable in 30-50 % of blood donors sera.

3. PRINCIPLE OF THE METHOD

The Chorus Chlamydomphila pneumoniae IgA device is ready to use for the detection of IgA antibodies against Chlamydomphila pneumoniae, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgA monoclonal antibodies conjugated to horseradish

peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Arbitrary Units (AU/ml).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.


3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with Chlamydomydia pneumoniae antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgA labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing IgA antibodies anti-Chlamydomydia pneumoniae and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted human serum containing IgA antibodies anti-Chlamydomydia pneumoniae and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Arbitrary Units (AU/ml) calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 18.0

NEGATIVE: when the result is < 12.0

DOUBTFUL/EQUIVOCAL: for all values between 12.0 and 18.0

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

Interpretation of the results, based on the combined assay of IgG, IgM and IgA, is reported in the table below:

Chlamydomphila pneumoniae antibody level			
IgG	IgM	IgA	Interpretation
Negative	negative	negative	No indication of infection
negative or positive	positive	negative or positive	Indication of infection in progress
Positive	negative	negative	Indication of past infection
negative or positive	negative	positive	Indication of infection or chronic or recurrent in progress

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range 10.0-100.0 AU/ml.

For samples > 100.0 AU/ml repeat the test pre-diluting the sample with Negative Control/Sample Diluent (PF83607 – not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 10.0 and 11.8 AU/ml.

14. ANALYTICAL SPECIFICITY

5 samples (2 Negative, 1 Cut-Off and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid Factor (44 – 220 UI/ml)
Bilirubin (4.5 mg/dl – 45 mg/dl)
Triglycerides (10 mg/dl – 250 mg/dl)
Hemoglobin (5 mg/ml – 30 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

15. METHOD COMPARISON

In an experimentation 157 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table :

		Reference		
		+	-	Total
Diesse	+	64	2	66
	-	6	85	91
	Total	70	87	157

Percent Positive Agreement (~Diagnostic Sensitivity):

91.4% CI_{95%}: 82.5-96.0

Percent Negative Agreement (~Diagnostic Specificity):

97.7% CI_{95%}: 92.0-99.3

16. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean AU/ml	CV%	Mean AU/ml	CV%
1	10.6	7.4	10.7	3.7
2	12.5	5.5	12.9	7.9
3	13.2	7.7	13.6	9.6
4	19.7	9.6	18.4	6.7
5	20.5	6.0	18.3	6.0
6	40.8	4.0	42.7	9.0
7	46.4	6.1	45.5	9.6
8	76.2	6.1	81.3	8.8
9	102.3	4.0	97.6	7.1

Sample	Precision between batches		Precision between instruments	
	Mean AU/ml	CV%	Mean AU/ml	CV%
1	10.2	3.1	10.2	3.1
2	14.5	3.8	14.5	6.6
3	16.7	4.2	16.6	4.0
4	18.5	3.5	18.5	5.6
5	22.1	4.7	22.4	3.4
6	42.3	5.3	42.3	5.2
7	45.5	4.0	45.5	7.1
8	83.9	3.3	83.9	2.4
9	92.6	4.9	92.6	1.2

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DIESSE Diagnostica Senese S.p.A.