CHORUS

PR3

REF 86066

REF 86066/12

Prodotto da/
Manufactured by/
Fabricado por:

DIESSE

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INSTRUCTIONS FOR USE

CHORUS PR3

For the semiquantitative determination of anti-PR3 antibodies

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the semiquantitative determination of anti-proteinase 3 (PR3) antibodies in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Antibodies against proteinase 3 (PR3) belong to the group of anti-neutrophil cytoplasmic antibodies (ANCA) which are directed against cytoplasmic components of neutrophilic granulocytes and monocytes. Indirect immunofluorescence tests on ethanol-fixed neutrophils has been the established method for the detection of ANCAs. It became apparent that some ANCAs create a cytoplasmic fluorescence pattern (thus called cANCA) while others create a perinuclear pattern (pANCA). As both patterns may cover multiple antigens, immunofluorescence is not suitable for a satisfying differential diagnosis of vasculitis; Thus each IFT should be verified with specific ELISA tests.

PR3 is a serine protease from the azurophilic granules (lysosomes) of neutrophils with a molecular weight of 29 kDa. The cationic protein has a proteolytic activity to elastin, hemoglobin and collagen VII. In addition, PR3 promotes the activation of platelets by cathepsin G and inactivates the C1-inhibitor.

ANCA are important markers for the differential diagnosis of autoimmune vasculitis. Autoantibodies to PR3 are a specific serological marker for Wegener's granulomatosis (WG). Although the etiology and pathogenesis is still obscure, autoantibodies to PR3 are thought to play an active role in the pathogenesis of WG. Anti-PR3 titers are strongly associated with disease activity and inhibit the proteolytic activity of PR3.

3. PRINCIPLE OF THE METHOD

The PR3 devices are ready to use for the assay of antibodies against PR3 in the CHORUS instruments. The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The antigen, highly purified, is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed

of monoclonal antibodies, highly specific for anti-human immunoglobulins conjugated to horse radish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

The Chorus PR3 is calibrated against reference sera from the CDC (Atlanta). The results are expressed in AU/ml.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instrument.
4. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument.
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed samples or samples presenting microbial contamination may all constitute sources of error.
11. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
12. Pipette the test serum (50 µl) in well 1 of the device (see figure).
13. Do not use the device after the expiry date.
14. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity REF 86004.**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86066)


The kit is sufficient for 12 tests (REF 86066/12).

DD DEVICES

6 packages each containing 6 devices (REF 86066).

2 packages each containing 6 devices (REF 86066/12).

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with highly purified proteinase 3 (PR3)

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing Tween-20 0.2% and Proclin 0.1%.

Positon 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labeled with horse radish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing IgG antibodies anti-PR3 and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted human serum containing IgG antibodies anti-PR3 and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can

be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results. Strongly lipemic, icteric, hemolyzed or contaminated samples cannot be used.

The test cannot be applied to plasma.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result in AU/ml calculated on the basis of a batch-dependent curve stored in the instrument.

The test serum can be interpreted as follows:

POSITIVE: when the result is > 18.0 AU/ml

NEGATIVE: when the result is < 12.0 AU/ml

DOUBTFUL/EQUIVOCAL: for all values between 12.0 and 18.0 AU/ml

In the case of a doubtful result, repeat the test. If the test remains doubtful, collect a new sample.

11. LIMITATIONS

All the positive results require a careful interpretation. The test cannot be used as only method for a clinical diagnosis. Negative results may not exclude an eventual infection. The test results should be interpreted in conjunction with information available from the clinical evaluation and other diagnostic procedures.

12. CALIBRATION RANGE

Calibration range 3.0-100.0 AU/ml.

For samples > 100.0 AU/ml repeat the test pre-diluting the sample with Negative Control/Sample Diluent (PF83607 – not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 3.0 and 11.2 AU/ml.

14. ANALYTICAL SPECIFICITY

5 samples (2 negative, 1 Cut-Off and 2 Positive) containing the following interfering substances were tested:

Rheumatoid Factor (44 – 220 IU/ml)
Bilirubin (4.5 mg/dl – 45 mg/dl)
Triglycerides (10 mg/dl – 250 mg/dl)
Hemoglobin (5 mg/ml – 30 mg/ml)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test.

15. CROSS-REACTIONS

Samples positive to MPO, Cardiolipin, Gliadin, AMA-M2, Transglutaminase, SS-A, SS-B, Sm, Cern-B, Scl-70, Jo-1, U1-70, snRNP, dsDNA IgG and ASCA were tested.

No significant cross-reactions were found.

16. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an experimentation 118 samples were tested with the Diesse kit and with another commercial kit.

The results are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	12	0	12
	-	1	105	106
	Total	13	105	118

Percent Positive Agreement (~Diagnostic Sensitivity):

92.3% CI_{95%}: 66.5-98.5

Percent Negative Agreement (~Diagnostic Specificity):

100.0% CI_{95%}: 96.5-100.0

17. PRECISION

Sample	Within-run Precision		Between-run precision	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	7.1	6.2	8.3	8.4
2	9.9	7.9	10.5	9.4
3	8.2	7.2	9.9	8.5
4	15.9	5.7	14.5	5.8
5	20.9	6.5	20.9	9.5
6	29.7	6.1	29.3	7.1
7	72.8	5.9	76.9	8.3
8	55.6	6.9	66.2	5.9

Sample	Precision between batches		Precision between instruments	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	8.2	4.4	8.2	1.2
2	9.2	4.6	9.2	5.4
3	8.6	7.9	8.6	6.7
4	14.9	1.9	14.9	5.6
5	21.1	3.8	21.0	2.1
6	27.1	2.8	27.1	3.1
7	74.9	6.3	74.9	1.9
8	66.6	7.4	66.6	3.3

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DIESSE Diagnostica Senese S.p.A.