CHORUS

ANA-8

REF 86010

REF 86010/12

DIESSE
DIESSE

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INSTRUCTIONS FOR USE

CHORUS ANA-8

**For the qualitative determination of IgG
antibodies anti-Sm, U1-snRNP, SS-A, SS-B,
snRNP/Sm, Scl-70, Jo-1 and CenpB**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgG class antibodies against 8 cellular and nuclear antigens (Sm, U1-snRNP, SS-A, SS-B, snRNP/Sm, Scl-70, Jo-1 and CenpB) in human serum, using a disposable device applied on the chorus and Chorus TRIO instruments.

2. INTRODUCTION

Anti-nuclear antibodies (ANA) are an important tool for the differential diagnosis of systemic rheumatic diseases. Indirect immunofluorescence test (IFT) on eukaryotic cells like HeLa has been the established method for the detection of ANAs. Single antibody specificities are distinguished by fluorescence patterns but more specific testing by ELISAs employing the target antigens are available too for a simple and reliable differentiation of ANAs.

ANAs are especially found in active and inactive systemic lupus erythematosus (SLE), mixed connective tissue diseases (MCTD), scleroderma, Sjögren's syndrome, polymyositis.

Antibodies against:

- Sm (Smith antigen) are directed against core proteins (B,B', D1-D3, E,F,G) of small nuclear ribonucleoproteins (snRNPs). Anti-Sm as well as antibodies against double stranded DNA (dsDNA) are highly specific for SLE and thus are included in diagnostic and classification criteria for SLE.
- U1 snRNP is directed to the 70 kDa protein of U1 snRNP. They are pathognomic for MCTD but do also occur in SLE. A high titer of antibodies against this antigen is typical for the Sharp-Syndrome.
- snRNP/Sm are directed against Sm and U1-snRNP proteins (70 kDa, A and C). They occur in SLE, Sjogren's syndrome, scleroderma and polymyositis.
- SS-A (Ro; soluble cytoplasmic and/or nuclear ribonucleoproteins of 52 kDa and 60 kDa) and antibodies against SS-B (La; 48 kDa protein associated with RNA polymerase III) are mainly found in high titers for primary and secondary Sjogren's

syndrome but also in SLE, congenital heartblock and neonatal lupus.

- Scl-70 are directed against DNA-topoisomerase I. They are highly specific for systemic scleroderma and give a hint for a severe course.
- CenpB (80 kDa centromere protein B) are typical for the CREST-Syndrome (69% of CREST.patients) which is a more protracted type of systemic sclerosis.
- Jo-1 are directed against histidyl-tRNA synthetase (cytoplasmic protein involved in protein biosynthesis) and are found in 20-40% of patients with polymyositis and dermatomyositis.

3. PRINCIPLE OF THE METHOD

The Chorus ANA-8 device is ready to use for the detection of IgG antibodies against 8 cellular and nuclear antigens, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay). The antigens are bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum. After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human immunoglobulins antibodies conjugated to horse radish peroxidase.

The unbound conjugate is eliminated and the peroxidase substrate is added. The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test in the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off) calculated in reference to the CDC Atlanta.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.

IO-09/268-C IFU 86010 – 86010/12 – Ed. 01.09.2015

4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86010).


The kit is sufficient for 12 tests. (REF 86010/12).

DD DEVICES

6 packages each containing 6 devices (REF 86010).

2 packages each containing 6 devices (REF 86010/12).

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with recombinant and/or highly purified antigens

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: saline proteic solution with Proclin (0.1%)

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labeled with horse radish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

in which undiluted serum must be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing IgG antibodies anti-Sm, U1-snRNP, SS-A, SS-B, snRNP/Sm, Scl-70, Jo-1 and CenpB and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted human serum containing IgG antibodies anti-Sm, U1-snRNP, SS-A, SS-B, snRNP/Sm, Scl-70, Jo-1 and CenpB and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.2

NEGATIVE: when the result is < 0.8

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALITICAL SPECIFICITY

5 samples (2 Negative, 1 Cut-Off and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (44-220 IU/ml)

Bilirubin (4.5 mg/dl - 45 mg/dl)

Triglycerides (10 mg/dl - 250 mg/dl)

Hemoglobin (5 mg/ml - 30 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

13. CROSS-REACTIONS

24 samples, positive to PR-3, MPO, Tg, a-Tg, a-TPO, Cardiolipin, Gliadin, AMA-M2, SS-A, SS-B, Sm, Jo-1, RF-G, RF-M, GBM, RA-CP, Intrinsic Factor, tTgG, tTgA and ASCA were tested.

No significant cross-reactions were found.

14. METHOD COMPARISON

In an experimentation 56 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table :

		Reference		
		+	-	Tot.
Diesse	+	17	0	17
	-	3	36	39
	Tot.	20	36	56

Percent Positive Agreement (~Diagnostic Sensitivity):

85.0% CI_{95%}: 63.9- 94.5

Percent Negative Agreement: (~Diagnostic Specificity):

100.0% CI_{95%}: 90.3- 99.9

The agreement between the two methods is excellent with a Cohen's Kappa of 0.89.

15. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean Index	CV%	Mean Index	CV%
1	0.2	25.0*	0.1	-
2	0.4	12.5	0.4	12.5
3	0.5	6.0	0.5	10.0
4	0.8	6.3	0.9	5.6
5	1.1	8.2	1.1	4.5
6	1.3	4.6	1.3	6.2
7	1.5	4.7	1.6	5.6
8	1.8	4.4	1.9	3.7

Sample	Between lots		Between Instruments	
	Mean Index	CV%	Mean Index	CV%
1	0.2	30.0*	0.2	30.0*
2	0.4	-	0.4	-
3	0.6	10.0	0.6	10.0
4	0.8	-	0.8	-
5	1.1	-	1.1	5.5
6	1.3	-	1.3	7.7
7	1.6	-	1.6	6.3
8	1.9	-	1.9	6.3

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero

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1. Peter J.B., Shoenfeld Y (1996) Sciences B.V., Amsterdamas.