ENZYWELL

TREPONEMA IgG

REF 91050 (96 tests)

Prodotto da/Manufactured by/ Hersteller:
Κατασκευάζεται από/ Fabricado por/ Fabriqué par:
DIESSE Diagnostica Senese
Strada dei Laghi, 39
53035 Monteriggioni (Siena) - Italy

INDICE / INDEX / INHALT/ ΠΕΡΙΕΧΟΜΕΝΑ / INDICE / SOMMAIRE / CONTEÚDO

1. UTILIZZAZIONE / INTENDED USE / VERWENDUNGSZWECK/ ΥΠΟΔΕΙΚΝΥΟΜΕΝΗ ΧΡΗΣΗ / INDICACIONES DE USO/ UTILISATION/UTILIZAÇÃO PRETENDIDA
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST / PHYSIOLOGIE/ ΕΙΣΑΓΩΓΗ / RESUMEN Y EXPLICACIÓN DEL TEST / INTERET CLINIQUE / SUMÁRIO E EXPLICAÇÃO DO TESTE
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / TESTPRINZIP/ ΑΡΧΗ ΤΗΣ ΜΕΘΟΔΟΥ / PRINCIPIO DEL MÉTODO / PRINCIPE/ PRINCÍPIOS DO TESTE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / INHALT DES KITS UND VORBEREITUNG DER REAGENZIEN / ΣΥΝΘΕΣΗ ΤΟΥ ΚΙΤ ΚΑΙ ΠΡΟΕΤΟΙΜΑΣΙΑ ΤΩΝ ΑΝΤΙΑΠΑΣΤΗΡΙΩΝ / COMPONENTES DEL KIT Y PREPARACIÓN DEL REACTIVO / CONTENU DU KIT ET PREPARATION DES REACTIFS / CONTEÚDA DO KIT E PREPARAÇÃO DO REAGENTE /
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / LAGERUNG UND STABILITÄT DER REAGENZIEN / ΤΡΟΠΟΣ ΑΠΟΘΗΚΕΥΣΗΣ ΚΑΙ ΣΤΑΘΕΡΟΤΗΤΑ ΤΩΝ ΑΝΤΙΑΠΑΣΤΗΡΙΩΝ/ CONSERVACIÓN Y ESTABILIDAD DE LOS REACTIVOS / MODALITES DE CONSERVATION ET STABILITE DES REACTIFS / ARMAZENAMENTO E ESTABILIDADE
6. PRECAUZIONI / PRECAUTIONS / HINWEISE UND VORSICHTSMASSNAHMEN/ ΠΡΟΦΥΛΑΞΕΙΣ / PRECAUCIONES DE USO / PRECAUTIONS D'UTILISATION / PRECAUÇÕES
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / ART UND LAGERUNG DER PROBEN / ΕΙΔΟΣ ΔΕΙΓΜΑΤΟΣ ΚΑΙ ΑΠΟΘΗΚΕΥΣΗ / TIPO DE MUESTRA Y CONSERVACION / TYPE D'ECHANTILLON ET CONSERVATION / TIPO E ARMAZENAMENTO DAS AMOSTRAS
8. PROCEDIMENTO / TEST PROCEDURE / TESTDURCHFÜHRUNG / ΔΙΑΔΙΚΑΣΙΑ / PROCEDIMIENTO DEL TEST / PROCEDURE / PROCEDIMENTO DE TESTE /
9. SCHEMA DEL SAGGIO / SCHEME OF TEST PROCEDURE / TESTSCHEMA /ΔΙΑΓΡΑΜΜΑ ΤΗΣ ΑΝΑΛΥΣΗΣ/ ESQUEMA DEL PROCEDIMIENTO DEL TEST / SCHEMA D'ANALYSE/ ESQUEMA DE PROCEDIMENTO DE TESTE /
10. VALIDAZIONE DEL TEST / TEST VALIDATION / TESTVALIDIERUNG / ΕΓΚΥΡΟΤΗΤΑ ΤΟΥ ΤΕΣΤ / VALIDACIÓN DEL TEST / VALIDATION DU TEST / VALIDAÇÃO DO TESTE /
11. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / INTERPRETATION DER ERGEBNISSE / ΕΡΜΗΝΕΙΑ ΤΩΝ ΑΠΟΤΕΛΕΣΜΑΤΩΝ / INTERPRETACIÓN DE RESULTADOS / INTERPRETATION DU TEST / INTERPRETAÇÃO DOS RESULTADOS
12. LIMITAZIONI DELLA PROCEDURA / LIMITATIONS OF THE PROCEDURE / GRENZEN DES VERFAHRENS / ΠΕΡΙΟΡΙΣΜΟΙ ΣΤΗΝ ΜΕΘΟΔΟ / LIMITACIONES DEL TEST / LIMITES DU TEST/ LIMITAÇÕES DO PROCEDIMENTO /
13. SPECIFICITA' ANALITICA / ANALYTICAL SPECIFICITY / ANALYTISCHE SPEZIFITÄT /ΑΝΑΛΥΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ/ ESPECIFICIDAD ANALÍTICA / SPECIFICITE ANALYTIQUE / ESPECIFICIDADE ANALÍTICA /
14. SENSIBILITA' E SPECIFICITA' DIAGNOSTICA / DIAGNOSTIC SENSITIVITY AND SPECIFICITY / DIAGNOSTISCHE SENSITIVITÄT UND SPEZIFITÄT /ΕΥΑΙΣΘΗΣΙΑ ΚΑΙ ΔΙΑΓΝΩΣΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ / SENSIBILIDAD Y ESPECIFICIDAD DIAGNOSTICAS / SENSIBLITE ET SPECIFICITE DIAGNOSTIQUES / ESPECIFICIDADE E SENSIBILIDADE DO DIAGNÓSTICO
15. PRECISIONE / PRECISION / PRÄZISION / ΑΚΡΙΒΕΙΑ / PRECISIÓN / PRECISION / PRECISÃO
16. GUIDA DEI PROBLEMI DI UTILIZZO / “TROUBLE SHOOTING” / “TROUBLE SHOOTING” (PROBLEMLÖSUNGEN) /ΟΔΗΓΟΣ ΕΠΙΛΥΣΗΣ ΠΡΟΒΛΗΜΑΤΩΝ /GUIA DE RESOLUCION DE PROBLEMAS / RESOLUTION DES PROBLEMES / RESOLUÇÃO DE PROBLEMAS
17. BIBLIOGRAFIA / REFERENCES / LITERATUR / ΒΙΒΛΙΟΓΡΑΦΙΑ / BIBLIOGRAFÍA / BIBLIOGRAPHIE / REFERÊNCIAS

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ENZYWELL TREPONEMA IgG

REF 91050

(Italiano)

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17. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S.A. Lukehart et al. Characterization of monoclonal antibodies to Treponema Pallidum. J. Immunology 134: 585 (1985).
3. S.A. Lukehart: Identification and characterization of Treponema Pallidum antigens by monoclonal antibodies. in: Monoclonal Antibodies, Academic Press 1986, p. 1.
4. L. Lewis et al.: Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
5. M. Norgard et al.: Sensitivity and specificity of monoclonal antibodies directed against antigenic determinants of Treponema Pallidum Nichols in the diagnosis of syphilis. J. Clin. Microbiol. 20: 711 (1984).
6. M.J. Bailey et al.: Production of murine monoclonal antibodies to the major axial filament polypeptide of Treponema pallidum. J. Gen. Microbiol. 133: 1805 (1987).

INSTRUCTIONS FOR USE

ENZYWELL TREPONEMA IgG

REF 91050

(English)

1. INTENDED USE

IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO TREPONEMA PALLIDUM IN HUMAN SERUM.

2. SUMMARY AND EXPLANATION OF THE TEST

The serological diagnosis of Syphilis is performed by demonstrating the presence of significant levels of specific Treponema pallidum antibodies in the serum sample.

The reference method used is the FTA-ABS technique which allows the detection of both IgG and IgM. However, its execution is laborious and the interpretation of the results is not simple; alternative methods have therefore been introduced to simplify the procedure. ELISA gives completely automated, reproducible and specific results, and for this reason is favoured by many operators.

The specific assay of IgM as opposed to IgG is of particular importance in the diagnosis of congenital Syphilis.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay) (1-6).

A pool of T. pallidum antigens obtained by recombinant technique in E. Coli is bound to the solid phase (8-well strips). The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgG monoclonal antibodies labelled with peroxidase.

The unbound conjugate is eliminated and the peroxidase substrate is added.

The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.

- **Bring to room temperature before use.**

MT PLATE MICROPLATE 12x8 wells coated with a pool of recombinant proteins of Treponema pallidum.

Use: open the package at the end opposite the code (SR followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONTROL + POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum, at known concentration of anti-Treponema IgG, in Phosphate buffer 0.01 mol/L containing BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONTROL CUT-OFF CUT OFF CONTROL 1 x 2 mL

Contents: Diluted human serum, at known concentration of anti-Treponema IgG, in Phosphate buffer 0.01 mol/L containing BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONJ CONJUGATE 1x16 mL

Contents: monoclonal antibodies labelled with Peroxidase, in phosphate buffer with phenol 0.05% and Bronidox 0,02%. Ready for use without further dilution.

CONTROL IgG - IgG NEGATIVE CONTROL (PF93910) 1 x 1.6 mL **INTERCHANGEABLE BETWEEN**

LOTS

Contents: Human serum in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 10 DILUENT 10 (PF93621) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

Contents: PBS containing proteins 10% w/v and sodium azide 0,09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619) 1x12 mL Ready for use. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED

- Incubator at 37°C
- Microplate reader, wavelength 450 or 450/620 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	8 weeks at 2/8°C, polythene bag
Controls	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	up to the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS**FOR IN VITRO DIAGNOSTIC USE ONLY.**

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. Use of the kit with automatic instruments must be validated by the user.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted

samples must be carefully shaken before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results. Strongly lipemic, icteric or contaminated samples should be avoided.

The test cannot be performed on human plasma.

8. TEST PROCEDURE

Manual Technique

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).

Dilute samples 1:51 distributing 20 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control.

Leave one well for the blank, performed using 100 µL of the substrate mixture.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µL), add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well. After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min.

9. SCHEME OF TEST PROCEDURE

Manual Technique

STEP 1 Place 100 µL of diluted sample//controls to the wells of the strips.

-

Incubate for 45 min. at 37°C

-

Wash 4 times (300 µL)

STEP 2 Add 100 µL of conjugate to each well

-

Incubate for 45 min. at 37°C

-

Wash 4 times (300 µL)

STEP 3 Add 100 µL of Substrate to each well

-

Incubate for 15 min. at R.T.

STEP 4 Add 100 µL of Stop Solution

-

Read absorbance at 450 nm within 30 min

10. VALIDATION OF THE TEST

Subtract the value of the blank (≤ 0.150) from all the other readings. The Positive control must have an O.D. at least 1.5 times that of the Cut-off serum. The ratio between Negative Control and Cut-off must be ≤ 0.6 .

The O.D. of the cut-off must be ≥ 0.2 at 450 nm and ≥ 0.16 at 450/620 nm.

11. INTERPRETATION OF THE RESULTS

Qualitative results

If the adsorbance of the sample is higher than that of the Cut-off, the sample is positive for the presence of specific IgG. Calculate the ratio between the average O.D. value of the sample and that of the Cut-off (INDEX). The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the Cut-off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE PROCEDURE

A serum sample obtained during the acute phase of infection, when only IgM antibodies are present, may be negative by this procedure.

The *Treponema pallidum* IgM level should be determined using Enzywell *Treponema pallidum* IgM kit. Alternatively, a second serum sample obtained 8-14 days later, should be tested in parallel to determine an increase in the IgG antibody level.

The test result should be used in conjunction with information available from the evaluation of other diagnostic procedures.

There is the possibility of cross reactions if samples contain anti-*E. coli* antibodies.

Although the control sera used in the kit are calibrated to WHO reference serum, certain discordances of results may be observed when the same serum is tested by different serological techniques.

13. ANALYTICAL SPECIFICITY

25 patient's sera containing anti-*Borrelia burgdorferi* antibodies, 12 containing anti-nuclear (ANA) antibodies, and 12 containing anti-EBV antibodies were tested. The same samples were also tested with two other commercial kits. In only one sample (EBV) was a positive reaction found with one of the other kits used.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In a clinical trial performed in a hospital laboratory, 205 samples from either healthy donors or from the serum collection were analysed. The TPHA method was used for comparison, while the reference method was the FTA-ABS technique. 51 samples were positive to the TPHA test (titre $\geq 1/160$), and 100% agreement was found on these samples between the two methods.

Of 24 samples with titre of 1/80, the ELISA method gave a positive result in 19 cases; two were doubtful and 3 were negative. The reference method confirmed the 3 negative samples found with ELISA. The results are summarized in the following table:

		REFERENCE	
		+	-
DIESSE	+	74	0
	-	2	129

The Enzywell kit offers 97.4% sensitivity and 100% specificity when compared with the FTA-ABS reference method and considering doubtful results as positive.

15. PRECISION

Within run Precision

Sample	TPG1 (Negative < Cut off)	TPG 2 (Positive >Cut off)	TPG 3 (Strongly positive)	Cut Off	Positive Control
N. Replicates	24	24	24	12	12
Media O.D.	0.18	0.50	1.14	0.28	1.69
CV %	6.4	14.0	6.2	3.8	3.6

Between run Precision

Sample	INDEX		
	Media	Standard Deviation	CV%
Pos. Control.	6.52	0.81	12.4
TPG 1	0.58	0.04	6.9
TPG 2	1.73	0.13	7.5
TPG 3	4.53	0.64	14.1

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing

		the plate (see instructions for use, section 4 , for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S.A. Lukehart et al. Characterization of monoclonal antibodies to Treponema Pallidum. J. Immunology 134: 585 (1985).
3. S.A. Lukehart: Identification and characterization of Treponema Pallidum antigens by monoclonal antibodies. in: Monoclonal Antibodies, Academic Press 1986, p. 1.
4. L. Lewis et al.: Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
5. M. Norgard et al.: Sensitivity and specificity of monoclonal antibodies directed against antigenic determinants of Treponema Pallidum Nichols in the diagnosis of syphilis. J. Clin. Microbiol. 20: 711 (1984).
6. M.J. Bailey et al.: Production of murine monoclonal antibodies to the major axial filament polypeptide of Treponema pallidum. J. Gen. Microbiol. 133: 1805 (1987).

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17. LITERATUR

- G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- S.A. Lukehart et al. Characterization of monoclonal antibodies to Treponema Pallidum. J. Immunology 134: 585 (1985).
- S.A. Lukehart: Identification and characterization of Treponema Pallidum antigens by monoclonal antibodies. in: Monoclonal Antibodies, Academic Press 1986, p. 1.
- L. Lewis et al.: Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
- M. Norgard et al.: Sensitivity and specificity of monoclonal antibodies directed against antigenic determinants of Treponema Pallidum Nichols in the diagnosis of syphilis. J. Clin. Microbiol. 20: 711 (1984).
- M.J. Bailey et al.: Production of murine monoclonal antibodies to the major axial filament polypeptide of Treponema pallidum. J. Gen. Microbiol. 133: 1805 (1987).

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17. REFERENCIAS

- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2.S.A. Lukehart et al. Characterization of monoclonal antibodies to Treponema Pallidum. J. Immunology 134: 585 (1985).
- 3.S.A. Lukehart: Identification and characterization of Treponema Pallidum antigens by monoclonal antibodies. in: Monoclonal Antibodies, Academic Press 1986, p. 1.
- 4.L. Lewis et al.: Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
- 5.M. Norgard et al.: Sensitivity and specificity of monoclonal antibodies directed against antigenic determinants of Treponema Pallidum Nichols in the diagnosis of syphilis. J. Clin. Microbiol. 20: 711 (1984).
- 6.M.J. Bailey et al.: Production of murine monoclonal antibodies to the major axial filament polypeptide of Treponema pallidum. J. Gen. Microbiol. 133: 1805 (1987).

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