ENZYWELL

CYTOMEGALOVIRUS IgG

REF 91010 (96 tests)

Prodotto da/Manufactured by/ Κατασκευάζεται από/
Fabricado por/ Fabriqué par:
DIESSE Diagnostica Senese
Strada dei Laghi 39
53035 Monteriggioni (Siena) – Italy

INDICE / INDEX / ΠΕΡΙΕΧΟΜΕΝΑ / INDICE / SOMMAIRE / CONTEÚDO

1. UTILIZZAZIONE / INTENDED USE / ΥΠΟΔΕΙΚΝΥΟΜΕΝΗ ΧΡΗΣΗ / INDICACIONES DE USO/ UTILISATION/UTILIZAÇÃO PRETENDIDA
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST / ΕΙΣΑΓΩΓΗ / RESUMEN Y EXPLICACIÓN DEL TEST / INTERET CLINIQUE / SUMÁRIO E EXPLICAÇÃO DO TESTE
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / ΑΡΧΗ ΤΗΣ ΜΕΘΟΔΟΥ / PRINCIPIO DEL MÉTODO / PRINCIPE/ PRINCÍPIOS DO TESTE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / ΣΥΝΘΕΣΗ ΤΟΥ ΚΙΤ ΚΑΙ ΠΡΟΕΤΟΙΜΑΣΙΑ ΤΩΝ ΑΝΤΙΑΔΡΑΣΤΗΡΙΩΝ / COMPONENTES DEL KIT Y PREPARACIÓN DEL REACTIVO / CONTENU DU KIT ET PREPARATION DES REACTIFS / CONTEÚDA DO KIT E PREPARAÇÃO DO REAGENTE /
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / ΤΡΟΠΟΣ ΑΠΟΘΗΚΕΥΣΗΣ ΚΑΙ ΣΤΑΘΕΡΟΤΗΤΑ ΤΩΝ ΑΝΤΙΑΔΡΑΣΤΗΡΙΩΝ/ CONSERVACIÓN Y ESTABILIDAD DE LOS REACTIVOS / MODALITES DE CONSERVATION ET STABILITE DES REACTIFS / ARMAZENAMENTO E ESTABILIDADE
6. PRECAUZIONI / PRECAUTIONS / ΠΡΟΦΥΛΑΞΕΙΣ / PRECAUCIONES DE USO / PRECAUTIONS D'UTILISATION / PRECAUÇÕES
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / ΕΙΔΟΣ ΔΕΙΓΜΑΤΟΣ ΚΑΙ ΑΠΟΘΗΚΕΥΣΗ / TIPO DE MUESTRA Y CONSERVACION / TYPE D'ECHANTILLON ET CONSERVATION / TIPO E ARMAZENAMENTO DAS AMOSTRAS
8. PROCEDIMENTO / TEST PROCEDURE / ΔΙΑΔΙΚΑΣΙΑ / PROCEDIMIENTO DEL TEST / PROCEDURE / PROCEDIMENTO DE TESTE /
9. SCHEMA DEL SAGGIO / SCHEME OF TEST PROCEDURE / ΔΙΑΓΡΑΜΜΑ ΤΗΣ ΑΝΑΛΥΣΗΣ/ ESQUEMA DEL PROCEDIMIENTO DEL TEST / SCHEMA D'ANALYSE/ ESQUEMA DE PROCEDIMENTO DE TESTE /
10. VALIDAZIONE DEL TEST / TEST VALIDATION / ΕΓΚΥΡΟΤΗΤΑ ΤΟΥ TEST VALIDACIÓN DEL TEST / VALIDATION DU TEST / VALIDAÇÃO DO TESTE /
11. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / ΕΡΜΗΝΕΙΑ ΤΩΝ ΑΠΟΤΕΛΕΣΜΑΤΩΝ / INTERPRETACIÓN DE RESULTADOS / INTERPRETATION DU TEST / INTERPRETAÇÃO DOS RESULTADOS
12. LIMITAZIONI DELLA PROCEDURA / LIMITATIONS OF THE PROCEDURE / ΠΕΡΙΟΡΙΣΜΟΙ ΣΤΗΝ ΜΕΘΟΔΟ / LIMITACIONES DEL TEST / LIMITES DU TEST/ LIMITAÇÕES DO PROCEDIMENTO /
13. SPECIFICITA' ANALITICA / ANALYTICAL SPECIFICITY / ΑΝΑΛΥΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ/ ESPECIFICIDAD ANALÍTICA / SPECIFICITE ANALYTIQUE / ESPECIFICIDADE ANALÍTICA /
14. CROSS-REATTIVI /CROSS REACTION / ΔΙΑΣΤΑΥΡΟΥΜΕΝΕΣ ΑΝΤΙΔΡΑΣΕΙΣ / REACCIONES CRUZADAS / RÉACTIONS CROISÉES / REAÇÕES CRUZADAS
15. SENSIBILITA' E SPECIFICITA' DIAGNOSTICA / DIAGNOSTIC SENSITIVITY AND SPECIFICITY / ΕΥΑΙΣΘΗΣΙΑ ΚΑΙ ΔΙΑΓΝΩΣΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ / SENSIBILIDAD Y ESPECIFICIDAD DIAGNOSTICAS / SENSIBLITE ET SPECIFICITE DIAGNOSTIQUES / ESPECIFICIDADE E SENSIBILIDADE DO DIAGNÓSTICO
16. PRECISIONE / PRECISION / ΑΚΡΙΒΕΙΑ / PRECISIÓN / PRECISION / PRECISÃO
17. GUIDA DEI PROBLEMI DI UTILIZZO / “TROUBLE SHOOTING” / ΟΔΗΓΟΣ ΕΠΙΛΥΣΗΣ ΠΡΟΒΛΗΜΑΤΩΝ / GUIA DE RESOLUCION DE PROBLEMAS / RESOLUTION DES PROBLEMES / RESOLUÇÃO DE PROBLEMAS
18. BIBLIOGRAFIA / REFERENCES / ΒΙΒΛΙΟΓΡΑΦΙΑ / BIBLIOGRAFÍA / BIBLIOGRAPHIE / REFERÊNCIAS

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ENZYWELL CYTOMEGALOVIRUS IgG

REF 91010

(Italiano)

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18. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
5. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
6. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
7. R. Ziegalmaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

INSTRUCTIONS FOR USE

ENZYWELL CYTOMEGALOVIRUS IgG

REF 91010

(English)

1. INTENDED USE

IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE AND QUANTITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO CYTOMEGALOVIRUS IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF CYTOMEGALOVIRUS INFECTION.

2. SUMMARY AND EXPLANATION OF THE TEST

Cytomegalovirus is a herpes virus transmitted by close human contact. No symptoms of infection are apparent in the majority of cases. However, the virus is very dangerous and may be fatal in immunodepressed patients. Serum-negative female patients who become infected during pregnancy may transmit the disease to the fetus. In 95% of cases this occurs without symptoms, but some neonates may present jaundice, hepato-splenomegaly and retarded psycho-motorial development. For this reason it is of great importance to determine the immunitary state of the patient and check for serum conversion. A significant increase in the anti-Cytomegalovirus IgG titer is indicative of a recent infection or reactivation of a latent infection.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay) (1-7).

The antigen, composed of partially purified and inactivated Cytomegalovirus, is bound to the solid phase (8-well strips). The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgG monoclonal antibodies conjugated to horse radish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample. When the enzymatic reaction is interrupted by the addition of a sulphuric acid solution, the yellow colour which develops can be easily read using a microplate reader.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.

Bring to room temperature before use.

MT PLATE MICROPLATE. 12x8 wells coated with Cytomegalovirus.

Use: open the package on the opposite side from the code (C plus lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expell the air and seal by pressing the closure.

CAL CALIBRATORS 5 x 1.6 mL

Contents: Diluted human serum, containing known concentrations of anti-CMV IgG antibodies, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution. The values can be expressed in IU/ml or EU/ml (see section 11 for conversion table). Calibrator 2 (1 IU/ml or 10 EU/mL) corresponds to the cut-off value and can be used in the qualitative test. The calibrators have the following values: 0.5, 1, 5, 10, 20 IU/mL, obtained by titration against the **“Proposed International Standard WHO”**.

Colour: the colour of the calibrators is proportional to the relative antibody titer.

CONJ CONJUGATE. 1x16 mL.

Contents: monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

CONTROL IgG - IgG NEGATIVE CONTROL (PF93910) 1 x 1.6 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Human serum not containing anti-CMV IgG antibodies, diluted in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution (Calibrator 0).

Stability: The product is stable up to the expiry date when stored unopened at 2/8°C.

WASH BUF 10x WASHING BUFFER 10X (PF93603). 1 x 100 mL. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 2 DILUENT 2 (PF93611). 1 x 100 mL. **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

Contents: Proteic solution in phosphate buffer with sodium azide 0,09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619). 1x12 mL. Ready for use. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602). 1x16 mL. **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader (wave length 450 or 450/620 nm and 405 nm, with linearity up to OD >= 2000)
- Microplate washer (preferable) able to dispense volumes ranging between 225-375 µl
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	6 weeks at 2/8°C, polythene bag
Calibrators	6 weeks at 2/8°C
Conjugate	6 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	up to the expiry date at 2/8°C
Rinsing Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Rinsing Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. The use of the kit with automated equipment has to be validated by the user.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted

samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided.

Human plasma cannot be used in the test.

8. TEST PROCEDURE

Manual Technique

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Rinsing Buffer 10x (100 mL + 900 mL H₂O).

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED calibrators (if possible, in duplicate) in a strip (100 µL in each well). Leave one well for the blank, performed using 100 µL of the substrate mixture.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µl), add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well. After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min. Take the reading again at 405 nm if there are OD higher than 2,000.

9. SCHEME OF TEST PROCEDURE

Manual Technique

STEP 1 Place 100 µL of diluted sample/controls to the wells of the strips.

-
- Incubate for 45 min. at 37°C
-
- Wash 4 times (300 µl)
-

STEP 2 Add 100 µL of conjugate to each well

-
- Incubate for 45 min. at 37°C
-
- Wash 4 times (300 µl)
-

STEP 3 Add 100 µL of Substrate to each well

-
- Incubate for 15 min. at R.T.
-

STEP 4 Add 100 µL of Stop Solution

-
- Read absorbance at 450 nm within 30 min

10. TEST VALIDATION

1. **CONTROL IgG** - (Calibrator 0): the O.D. of the negative control must be ≤ 0.6 times the O.D. of Calibrator 2.
2. The O.D. of Calibrator 2 must be ≥ 0.2 at 450 nm; ≥ 0.16 at 450/620 nm.
3. The O.D. of Calibrator 5 must be ≥ 1.5 .

11. INTERPRETATION OF THE RESULTS

QUANTITATIVE RESULTS

Report the OD of the calibrators on a graph after subtracting the OD of the blank. The corresponding titer of the test sample can be obtained by extrapolation.

NOTE: A standard curve must be performed for each run.

If the O.D. of any sample or calibrator is over 2.0, take the reading at 405 nm and multiply the value by 3.

Conversion of the O.D. into units/mL

The anti-Cytomegalovirus IgG can be expressed in EU/mL (arbitrary units) or in IU/mL (proposed anti CMV IgG WHO standard), by interpolating the results of the 6 calibrators and comparing the OD of the sample with the curve obtained.

Table 1: Calibrator Values in arbitrary and international units

	EU/mL	IU/mL
Calibrator 0	0	0
Calibrator 1	5	0.5
Calibrator 2	10	1
Calibrator 3	50	5
Calibrator 4	100	10
Calibrator 5	200	20

The degree of immunity can be interpreted as follows:

IMMUNE: when the anti-Cytomegalovirus IgG concentration in the sample is > 12 EU/mL or 1.2 IU/mL

NON IMMUNE: when the anti-Cytomegalovirus concentration is < 8 EU/mL or < 0.8 IU/mL

DOUBTFUL : if the result is between the two values. In this case it is advisable to repeat the test in duplicate.

QUALITATIVE RESULTS

Calculate the ratio between the O.D. value of the sample and that of the Cut-off (Calibrator 2). The sample is considered:

Positive: if the ratio is > 1.2

Doubtful: >= 0.8 - <= 1.2

Negative: if the ratio is < 0.8

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE PROCEDURE

A serum sample obtained during the acute phase of infection, when only IgM antibodies are present, may be negative by this procedure.

Heat-inactivation of the test serum does not give rise to erroneous results in the determination of anti-CMV IgG, as in the case of anti-CMV IgM.

The Cytomegalovirus IgM level should be determined using Enzywell Cytomegalovirus IgM kit. Alternatively, a second serum sample obtained 8-14 days later, should be tested in parallel to determine an increase in the IgG antibody level.

The test result should be used in conjunction with information available from the evaluation of the case history or other diagnostic procedures.

13. ANALYTICAL SPECIFICITY

34 serum samples containing potential interfering substances were tested:

- Rheumatoid factor (n=4) < 800 UI/ml
- Heterophile antibodies (n=6)
- Bilirubin (n=10) < 11mg/dl
- Triglycerides (n=10) < 1100 mg/dl
- Haemolysis (n=4) > 10 mg/dl

No interference was noted in any case.

14. CROSS REACTION

23 samples which were negative for CMV, but containing IgG antibodies against viruses such as Rubella, Epstein Barr, Herpes Simplex and Mumps, were tested. In no case did the presence of these antibodies influence the test.

15. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In a clinical trial performed in a hospital laboratory, 264 samples were analysed, 47 of which proved negative and 217 positive. The samples were analysed with another commercial immunoenzymatic method: 100% agreement was found between the two methods, both in the positive and the negative samples.

The Enzywell kit offers 100% sensitivity and specificity.

16. PRECISION**Table 1 Within-run Precision performed on 3 different lots**

Cut Off n=21	Batch. N.104	Batch. N. 105	Batch. N. 106
O.D.	0.453	0.374	0.433
CV%	4.62	4.73	4.07

Table 2 “Between-run” Precision

Sample	EU/ml			
	I run	II run	III run	CV%
n.1	12.6	17.5	17.6	18
n.2	78.5	87.5	96.3	10

17. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert paragraph 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

18. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
5. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
6. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
7. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2.H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
- 3.Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
- 4.M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
- 5.F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
- 6.P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
- 7.R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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- 1.G.B. Wisdom: Enzyme-Immunoassay.
- 2.Clin. Chem. 22: 1243 (1976).
- 3.H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
- 4.Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
- 5.M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
- 6.F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
- 7.P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
- 8.R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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17. BIBLIOGRAPHIE

8. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
9. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
10. Grint P.C.A. et él.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
11. M. Musiani et él.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
12. F. de Ory et él.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
13. P. Dal Monte et él.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
14. R. Zieglermaier et él.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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18. REFERÊNCIAS

1. G.B. Wisdom: Enzyme-Immunoassay.
2. Clin. Chem. 22: 1243 (1976).

3. H.O. Kangro: Cytomegalovirus serology: does it give the answer? *Serodiagnosis and Immunotherapy* 1: 91 (1987).
4. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. *J. Clin. Pathol.* 38: 1059 (1985).
5. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. *J. Clin. Pathol.* 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. *Serodiagnosis and Immunotherapy* 2: 423 (1988).
7. P. Dal Monte et al.: Cytomegalovirus umano. *Diagnosis* 2: 67 (1990).
8. R. Zieglermaier et al.: ELISA. *la Ricerca Clin. Lab.* 10: 83 (1980).

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusqu'e Ημερομηνία λήξης Prazo de validade
	EN ES IT	Caution,consult accompanying documents Atención,ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbicante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

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