DIESSE Diagnostica Senese SpA

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

ISTRUZIONI D'USO**W. W. TUBE TEST****(Italiano)**

REF 21120	S. typhi antigene O (gruppo D)	REF 21125	S. paratyphi A totale
REF 21121	S. typhi antigene H	REF 21126	S. paratyphi B antigene O (gruppo B)
REF 21122	S. typhi totale	REF 21127	S. paratyphi B antigene H
REF 21123	S. paratyphi A antigene O (gruppo A)	REF 21128	S. paratyphi B totale
REF 21124	S. paratyphi A antigene H	REF 21129	Brucella

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1. INTENDED USE

SALMONELLA AND BRUCELLA SUSPENSIONS FOR THE DETERMINATION OF SPECIFIC AGGLUTINATING ANTIBODIES IN SERUM USING A TUBE TEST.

2. INTRODUCTION

Serological studies of patients with persistent fever can be useful in detecting typhoid fever and brucellosis (1,2). The test consists in the observation of a direct agglutination reaction between the bacterial suspension and the serum sample which is suspected of containing specific antibodies. The present method, known as the Widal-Wright reaction, allows a semiquantitative titration of the serum agglutinins.

3. PRINCIPLE OF THE METHOD

A colourless suspension of inactivated bacteria is mixed with serially diluted serum to be tested. After incubation the agglutination pattern is read: the homogeneous suspension indicates a negative reaction while agglutinated cells on the bottom of the tube and a clear supernatant indicate a positive reaction.

4. KIT COMPOSITION

- The reagents are sufficient for **40 tests**.
- **Bring the reagents to room temperature before use.**

BACTERIAL SUSPENSION	BACTERIAL SUSPENSION 1 vial x 20 mL
-----------------------------	-------------------------------------

Contents: Colourless suspension of killed bacteria at an optimized concentration for performance of the tube test. Contains sodium azide 0,09% as preservative.

OTHER MATERIALS REQUIRED BUT NOT SUPPLIED

- Polyvalent Positive Control **REF** 21250
- Pipettes with disposable needles.
- Physiological saline.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
BACTERIAL SUSPENSIONS	2 months at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

All materials of biological origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The reagents contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.
3. If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. **Allow all reagents and samples to come to room temperature (18-30°C) before use.** After use, empty the dropper and return reagents to the recommended storage temperature (2-8°C).
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure.

7. TYPE AND STORAGE OF SAMPLE

Fresh serum may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Use immediately after thawing.

8. PROCEDURE

Bring the reagent to room temperature and mix the vial of bacterial suspension carefully but not vigorously.

Dilute the sample 1:20 with physiological saline (0,05 mL serum + 0,950 mL saline). Take 0,5 mL from the first tube and make serial dilutions with 0,5 mL of saline, as shown in the scheme. Discard 0,5 mL of diluted sample from the last tube. Add 0,5 mL of antigen suspension to each tube and incubate overnight at 37°C. See the scheme shown below. The blank is performed by mixing 0,5 mL of suspension with 0,5 mL of saline.

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9. TEST VALIDATION

Perform the tube test using one drop of **Positive Control** (available separately); evident agglutination must be found.

10. INTERPRETATION OF THE TEST – NORMAL VALUES

The titer is given by the highest dilution able to cause agglutination of at least 50% of the bacteria present.

Positive reaction:

The bacteria are agglutinated in a granular (antigen O) or flocculant (antigen H) manner; the fluid supernatant is clear, in comparison with the blank.

Negative reaction:

The suspension appears homogeneous as in the control tube.

NORMAL VALUES

Titers of 1/80 are considered suspect. Higher titers are considered diagnostically significant.

With the exception of Brucella, almost all the Widal Wright titers are below 1/640.

11. TEST LIMITATIONS

The titers against O antigens are not specific for the different Salmonelle. It is known that there are cross-reactions between *S. typhi*, *S. paratyphi B* and *S. paratyphi A* because they have antigen 12 (Kauffman-White classification) in common.

12. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

An experimentation performed on human sera obtained from hospital laboratories, the titers of which were taken as reference, allowed an evaluation of the sensitivity of the following suspensions: *S. typhi* O, *S. typhi* H, *S. paratyphi* BO, *S. paratyphi* BH.

When *S. typhi* antigen H was determined on 30 positive samples, a sensitivity of 96.7% was found, while the suspension *S. typhi* antigen O tested on 51 positive sera gave 3 false negatives and a sensitivity of 94.4%.

S. paratyphi B antigen H and antigen O were tested with respectively 12 and 24 positive samples. Sensitivity was 92.3% for both the B total suspensions had a slightly lower sensitivity of 91.9%

S. paratyphi A antigen O and antigen H were tested in parallel with another commercial product. 12 positive and 3 doubtful results were found. The Diesse reagent gave 1 false negative and the sensitivity was therefore 93.7%, considering the doubtful results as positive. The same sensitivity was found with the *S. paratyphi A* total suspension. Using the same criteria, 48 positive samples for Brucellosis were tested. Diesse gave 2 false negative results and the sensitivity was therefore 96%.

Human sera found negative in another commercial test were also studied. The following results were obtained:

<u>Suspension</u>	<u>No. of samples</u>	<u>Specificity %</u>
<i>S. typhi</i> antigen O	76	97.4
<i>S. typhi</i> antigen H	58	99.3
<i>S. paratyphi B</i> antigen O	96	100
<i>S. paratyphi B</i> antigen H	93	100
<i>S. paratyphi A</i> antigen O	80	98.8
<i>S. paratyphi A</i> antigen H	80	100
<i>S. paratyphi A</i> total	80	98.8
<i>S. paratyphi B</i> total	93	98.9
Brucella	83	100

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