**DIESSE Diagnostica Senese SpA**

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

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### **MONO SLIDE**

**REF** 02190

(Italiano)

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(English)

#### 1. INTENDED USE

**SLIDE TEST FOR THE QUALITATIVE DETERMINATION OF HETEROPHILE ANTIBODIES IN HEMOAGGLUTINATION, FOR THE DIAGNOSIS OF INFECTIOUS MONONUCLEOSIS**

#### 2. INTRODUCTION

In 1932, Paul and Bunnell detected IgM-class antibodies in the serum of patients affected by mononucleosis, which caused agglutination of sheep, horse and bovine erythrocytes. Later Davidsohn succeeded in distinguishing the heterophile antibodies found in serum sickness from those found in mononucleosis, by selective absorptions with homogenates of guinea pig or horse kidney (1).

The kit contains the reagents to perform the screening test with ease and rapidity. The stabilized horse erythrocytes reveal the presence of heterophile antibodies of mononucleosis with specificity and sensitivity.

#### 3. PRINCIPLE OF THE METHOD

The specific heterophile antibodies found in mononucleosis give an evident agglutination when reacted with horse or sheep erythrocytes.

#### 4. KIT CONTENTS AND REAGENT PREPARATION:

- The reagents are sufficient for 90 tests.
- Bring reagents to room temperature before use.

**GR REAGENT** HORSE ERYTHROCYTES 1 vial x 3.5 mL

Contents: Horse erythrocytes stabilized in PBS with 0.01% merthiolate as preservative.

**CONTROL +** POSITIVE CONTROL 1 x 0.5 mL

Contents: Human serum containing heterophile antibodies, with sodium azide 0.09%.

**CONTROL -** NEGATIVE CONTROL 1 x 0.5 mL

Contents: Human serum negative for the presence of heterophile antibodies. Contains sodium azide 0.09% as preservative.

#### OTHER MATERIALS SUPPLIED

- Disposable slides. Only use new, or well washed and dried slides.
- Mixing Sticks.

#### OTHER MATERIALS REQUIRED BUT NOT PROVIDED

- 40 µL pipettes with disposable tips
- Timer
- Saline solution.

#### ADDITIONAL MATERIAL FOR THE DIFFERENTIAL TEST

REF 02191 PAUL BUNNELL KIT DIESSE.

#### 5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

#### **Reagents have a limited stability after opening and/or preparation**

REAGENT	CONDITIONS
GR REAGENT	3 months

POSITIVE CONTROL	3 months
NEGATIVE CONTROL	3 months.

## **6. PRECAUTIONS**

**FOR IN VITRO DIAGNOSTIC USE ONLY. STORE AT 2-8°C**

**Caution:** This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

*Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.*

### Health and Safety Information

- Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- The reagents contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.  
If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
- Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

### Analytical precautions

- Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature.**
- Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable.
- The reaction can be performed using a rotating mixer. Regulate the speed at  $100 \pm 10$  rpm.
- The agglutination reaction is best interpreted by observing the slide under a lamp.

## **7. TYPE AND STORAGE OF SAMPLE**

Fresh or defrosted serum. Samples may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers.

## **8. PROCEDURE**

### Screening Test

Dispense one drop (40 µl) of sample in a circle on the slide. Add one drop of GR Reagent, mix and rotate the slide for 1 minute. Observe the agglutination. Positive samples will cause agglutination of the erythrocytes. These samples can be confirmed with the differential test and the heterophile antibodies titrated with the quantitative test. In both cases, the Diesse Paul Bunnell kit REF 02191 is required.

### Differential test

Dispense one drop of sample (50 µL) in two circles on the slide. Add one drop of Sorbent K to one, and one drop of Sorbent S to the other. Mix well with two mixing sticks and add one drop of revealer (GR Reagent) to each mixture. Rotate the slide and observe the agglutination after 1 minute. Compare the results in the two circles. A clear agglutination in the circle containing Sorbent K, and a negative result with Sorbent S in the other circle is considered positive for infectious mononucleosis. Results are negative when there is no agglutination, or if there is a more marked agglutination in the circle containing Sorbent S. If a strong agglutination is found in both circles, dilute the samples 1:10 with saline and repeat the test.

Titration of the heterophile antibodies: see instructions for use of the Paul Bunnell kit **REF** 02190.

**9. TEST VALIDATION**

Agglutination must be obtained when the positive control is reacted with the latex suspension. No agglutination must be found when the negative control serum is tested.

**10. INTERPRETATION OF THE TEST**

Evident agglutination of the erythrocytes is caused by the presence of heterophile antibodies at a sufficiently high titer to indicate a probable infection by mononucleosis in the acute or subacute phase.

**11. LIMITATIONS OF THE TEST**

The result of the test must always be evaluated together with other data resulting from other diagnostic procedures. It is known that in patients below 3-4 years affected by infectious mononucleosis, there may be no formation of heterophile antibodies.

**12. ANALYTICAL SENSITIVITY**

The test detects antibody titers of 1/56; according to the Paul Bunnell method, such titers are indicative of a current or recent mononucleosis infection.

**13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY**

The reagent was tested on 95 samples, 64 of which were negative and 31 positive. The reference result was obtained by assay in an ELISA technique for antibodies against specific EBV proteins. One false positive and 1 false negative result were found, giving a specificity and sensitivity of respectively 98.5% and 96.9%.

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