DIESSE Diagnostica Senese SpA

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

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STAINED PROTEUS SUSPENSIONS TO DETECT SPECIFIC AGGLUTINATING ANTIBODIES IN SERUM. THE SCREENING TEST IN SLIDE AND THE SEMIQUANTITATIVE TEST IN TUBE CAN BE PERFORMED.

2. INTRODUCTION

Stained suspensions of Proteus are used for the diagnosis of rickettsiae infection (Weil-Felix test). The test is based on the fact that the somatic antigens of some strains of Proteus are identical to those of some species of Rickettsie. The serum from subjects infected by some forms of Rickettsie agglutinate suspensions of Proteus strains.

3. PRINCIPLE OF THE METHOD

The serum sample is mixed with a bacterial suspension at the optimal concentration for the slide test. After 1 minute a visible agglutination will appear if the serum contains specific antibodies against the bacteria in the suspension. Positive samples must then be tested in the tube test, using one drop of the same suspension.

4. KIT COMPOSITION

- The reagents are sufficient for 60 tests.
- Bring the reagents to room temperature before use.

BACTERIAL SUSPENSION *BACTERIAL SUSPENSION* 1 vial x 3 mL

Contents: Suspension of stained, killed bacteria at an optimized concentration for performance of the slide test. Contains sodium azide 0,09% as preservative.

OTHER MATERIALS REQUIRED BUT NOT SUPPLIED

- **REF** 21350 Febrile Antigen Polyvalent Positive Control (2 mL)
- Slides. Only use new slides.
- Mixing sticks
- Pipettes with disposable needles.
- Physiological saline.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
BACTERIAL SUSPENSION	2 months at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

All reagents must be handled in accordance with the safety precautions adopted in the laboratory. All samples must be considered potentially infectious.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The reagents contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.
3. If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. **Allow all reagents and samples to come to room temperature (18-30°C) before use.** After use, empty the dropper and return reagents to the recommended storage temperature (2-8°C).
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure.

7. TYPE AND STORAGE OF SAMPLE

Fresh serum may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Use immediately after thawing.

8. PROCEDURE

Bring the reagent to room temperature and mix the vial of bacterial suspension carefully but not vigorously.

* SCREENING TEST

Distribute 40 µL of serum in a circle on the slide. Add one drop of suspension and mix with the serum. Rotate the slide. Read the results after exactly 1 minute. Positive samples will give agglutination while negative samples will remain homogeneous.

* TITRATION ON SLIDE

Dispense decreasing volumes of sample (40-20-10-5 µl) in four circles on the slide. Add one drop of suspension and mix rapidly starting from the circle containing the smallest sample volume. Turn the slide around and observe the agglutination after exactly one minute. The serum titer can be determined using the following table:

µl serum	titer
40	1:20 - 1:40
20	1:40 - 1:80
10	1:80 - 1:160
5	1:160 - 1:320

* TUBE TITRATION

Dilute the sample 1:20 with physiological saline (0.1 mL serum + 1.9 mL saline). Take 1 mL from the first tube and make serial dilutions with 1 mL of saline, as shown in the scheme. Add one drop of antigen suspension to each tube and incubate overnight at 37°C. See the scheme shown below. The blank is performed by mixing the suspension with saline

Tube	1	2	3	4	5	6	7
Saline mL	1.9	1	1	1	1	1	1
Serum ml	0,1						
Mix and transfer ml							
Dilution / Titer	1/20	1/40	1/80	1/160	1/320	1/640	1/1280

9. TEST VALIDATION

Perform the slide test using one drop of **Positive Control Serum** (available separately). Evident agglutination must be seen.

10. INTERPRETATION OF THE TEST – NORMAL VALUES

Titers of 1/80 are considered suspect; higher titers are diagnostically significant.

<i>Disease</i>	<i>Agent</i>	<i>Weil-Felix reaction</i>		
		<i>OX19</i>	<i>OX2</i>	<i>OXK</i>
<i>Endemic Typhoid fever</i>	<i>R. prowazekii</i>	++	+	-
<i>Murine typhoid fever</i>	<i>R. mooseri</i>	++	+	-
<i>Rocky Mountain spotted fever</i>	<i>R. rickettsii</i>	+ / ++	++	- / +
<i>Boutonneuse fever</i>	<i>R. conori</i>	+ / ++	++	-
<i>Tsutsuga-mushi fever</i>	<i>R. orientalis</i>	-	-	++

11. TEST LIMITATIONS

Not all rickettsiae give a cross-reaction with *Proteus*. *R. Burnetii*, the etiological agent of Q-Fever, does not react.

12. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

The three reagents were compared with another commercial product, using sera from a hospital routine practice (89), 16 of which had a request for the Weil Felix test. Only 1 sample was positive, in agreement with the other method; all the others were negative in both methods.

Diagnostic specificity was therefore 100%. The sensitivity was also 100%, but this data is only of relative value due to the low number of positive samples found during the experimentation.

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