DIESSE Diagnostica Senese SpA

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

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1. INTENDED USE

STAINED BACTERIAL SUSPENSIONS OF SALMONELLA AND BRUCELLA FOR THE SEMIQUANTITATIVE DETERMINATION OF SPECIFIC AGGLUTINATING ANTIBODIES. THE TEST IS PERFORMED ON A MICRPLATE.

2. INTRODUCTION

Serological studies of patients with persistent fever can be useful in detecting typhoid fever and brucellosis (1,2).

The test consists in the observation of a direct agglutination reaction between the bacterial suspension and the serum sample which is suspected of containing specific antibodies. The present method, known as the Widal-Wright reaction, allows a semiquantitative titration of the serum agglutinins.

3. PRINCIPLE OF THE METHOD

The test involves mixing a suspension of inactivated coloured bacteria with two fold serial dilutions of test serum in U bottom wells. After incubation overnight, a positive serum determines the agglutination of bacteria. If the bacterial cells settle in compact buttons, the result of the test is negative.

4. KIT COMPOSITION

- The reagents are sufficient for 500 tests.
- Bring the reagents to room temperature before use.

BACTERIAL SUSPENSION BACTERIAL SUSPENSION 5 vials x 5 mL

Contents: Suspension of stained, killed bacteria at an optimized concentration for performance of the microplate test. Contains sodium azide 0,09% as preservative.

CONTROL + POSITIVE CONTROL 1 x 2 mL

Contents: Polyvalent control serum of animal origin, containing sodium azide 0,09% as preservative. Ready for use.

OTHER MATERIALS SUPPLIED

- 6 96-well U-bottom microtiter plates

OTHER MATERIALS REQUIRED BUT NOT SUPPLIED

- Pipettes with disposable needles.
- Physiological saline.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
BACTERIAL SUSPENSION	2 months
POSITIVE CONTROL	2 months

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

- Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- The reagents contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.
- If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
- Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

- Allow all reagents and samples to come to room temperature (18-30°C) before use.** After use, empty the dropper and return reagents to the recommended storage temperature (2-8°C).
- Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- Do not modify the Test Procedure.


7. TYPE AND STORAGE OF SAMPLE

Fresh serum may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Use immediately after thawing.

8. PROCEDURE

Bring the reagent to room temperature and mix the vial of bacterial suspension carefully but not vigorously.

Mix 10 µL of test serum with 90 µL of saline in the first well of a strip. Distribute 50 µL of saline in the other wells and make serial two-fold dilutions along the row discarding the last 50 µL. Add 1 drop (50 µL) of carefully mixed homogeneous bacterial suspension to each well (see table). Mix and incubate at 37°C for at least 18 hours, covering carefully and, if possible, placing in a humid chamber to avoid evaporation. Avoid any kind of vibration. At the end of the incubation period, leave the plate to settle for at least one hour, to improve the reading.

	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl		
											
1	2	3	4	5	6	7	8	9	10	11	12
90 µl fisiol. + 10µl siero	50µl fisiol	50µl fisiol	50µl fisiol	50µl fisiol	50µl fisiol	50µl fisiol	50µl fisiol	50µl fisiol	50µl fisiol	50µl Cont.Pos.	50µl fisiol.
50µl sospensione											
Titolo 1:20	1:40	1:80	1:160	1:320	1:640	1:1280	1:2560	1:5120	1:10240	C.P.	NEG.

9. TEST VALIDATION

Carry out the test using 1 drop (50 µL) of positive control in a well with 1 drop of bacterial suspension. After incubation a positive reaction should be given. To control the negative reaction use 1 drop (50 µL) of saline instead of the sample. A non-agglutinated button must be clearly visible.

10. INTERPRETATION OF THE TEST – NORMAL VALUES

Positive reaction

The agglutinated bacteria cover the bottom of the well in a homogeneous layer. When using the “H” suspension, small buttons surrounded by a layer of agglutinated bacteria has to be considered as a positive reaction.

Negative reaction

The bacteria settle in a compact brightly coloured button.

Titers of 1/80 are considered suspect; higher titers are diagnostically significant.

11. TEST LIMITATIONS

In the microplate agglutination test there may be a cross-reactoin between O and H suspensions, and there is therefore generally an overestimation of the H titer.

12. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

An experimentation performed on human sera obtained from hospital laboratories, the titers of which were taken as reference, allowed an evaluation of the sensitivity of the following suspensions: S. typhi O, S. typhi H, S. paratyphi BO, S. paratyphi BH.

When S. typhi antigen H was determined on 30 positive samples, a sensitivity of 96.7% was found, while the suspension S. typhi antigen O tested on 51 positive sera gave 1 false negative and a sensitivity of 98.1%.

S. paratyphi B antigen H and antigen O were tested with respectively 12 and 24 positive samples. Sensitivity was 100% and 96.0% respectively.

S. paratyphi A antigen O and antigen H were tested in parallel with another commercial product. 12 positive and 3 doubtful results were found. The Dienes reagent gave no false negative reactions and the sensitivity was therefore 100%. Using the same criteria, 45 positive samples for Brucellosis were tested. Dienes gave 2 false negative results and the sensitivity was therefore 95.7%. The S. paratyphi B total suspension gave a sensitivity of 97.1% with 1 false negative result.

Human sera found negative in another commercial test were also studied. The following results were obtained:

<u>Suspension</u>	<u>No. of samples</u>	<u>Specificity %</u>
S. typhi antigen O	76	96.2
S. typhi antigen H	58	95.1
S. paratyphi B antigen O	96	98.0
S. paratyphi B antigen H	93	95.9
S. paratyphi A Total	82	98.8
S. paratyphi B Total	93	96.9
Brucella	83	97,6

13. REFERENCES

1. P. Nicoletti in F. Pasquinelli: Diagnostica e Tecnica di Laboratorio Vol. 2 p. 319. Rosini ed. Firenze 1981.
2. Edwards P.R., Ewing W.H. Identification of Enterobacteriaceae 3rd ed. Burgess Publishing Co.

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusqu'e Ημερομηνία λήξης Prazo de validade
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbicante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

Diesse Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) – Italy
Tel. 0577-587111