DIESSE Diagnostica Senese SpA

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DIESSE Diagnostica Senese SpA

Via delle Rose 10, 53035 Monteriggioni SI
Tel. ++39/0577/ 587111 - Fax ++39/0577/318690**INSTRUCTIONS FOR USE****STREPTOSLIDE EXTRACTING REAGENTS****REF** 0217**(English)****1. INTENDED USE****CHEMICAL EXTRACTANTS FOR USE IN THE SLIDE TEST FOR THE IDENTIFICATION OF STREPTOCOCCAL GROUPS A, B, C, D, F AND G.****2. INTRODUCTION**

Lancefield demonstrated that streptococci have group specific antigens on their cellular wall. These antigens (polysaccharides) have been extracted and recognised by specific homologous antisera. Most of the human pathogenic streptococci belong to the groups A,B,C,D,F and G (1). Group E streptococci are rarely isolated. The extraction procedures have been improved and several reagents based on chemical as well as enzymatic reactions are now available (1,2,3,4).

In the STREPTOSLIDE kit, the identification of Streptococci group A, B, C, F and G is performed without incubation by a simple chemical extraction (micromethod with nitrous acid) using latex with group specific antibodies. Group D is usually identified by direct reaction with the specific latex reagent for Group D. The isolated colonies negative to all group specific reagents and suspected to belong to the Group D, are tested again with the specific Group D latex reagent after enzymatic digestion.

3. PRINCIPLE OF THE METHOD

A few isolated colonies of streptococci are suspended in the reaction mixture for extraction of the group antigen. An aliquot of the suspension is mixed with latex reagent coated with immunoglobulins specific for the group antigens. Agglutination allows identification of the group.

4. KIT CONTENTS**The reagents are sufficient for 50 tests.****Bring to room temperature before use.****REAGENT 1** REAGENT 1 1 x 3 mL

Contents: Sodium nitrite solution.

REAGENT 2 REAGENT 2 1 x 3 mL

Contents: Acetic acid diluted in a solution of methylorange.

REAGENT 3 REAGENT 3 2 x 4 mL

Contents: Ammonium carbonate solution, contains sodium azide 0,09%.

LYTIC ENZYME LYTIC ENZYME freeze-dried 2 x 2 mL

Contents: Enzymatic extracting reagent.

Preparation: Freeze-dried, to be reconstituted with 2.0 mL of distilled water.

WARNING AND PRECAUTIONS

Reagents 1, Reagent 2, Reagent 3 and Lytic Enzyme are classified and labelled as dangerous in accordance with European legislation. All safety data are reported in the Material Safety Data Sheet.

CONTROL + POSITIVE CONTROL 1 vial x 0.5 mL

Contents: Polyvalent cell extract containing sodium azide 0,09% as preservative.

	Section
Changes introduced in the current revision	4, 6, Glossary

OTHER MATERIALS

- Slides with 6 reaction divisions, black background.
- 200 mixing sticks.

OTHER MATERIALS REQUIRED BUT NOT SUPPLIED

- Normal laboratory glassware
- Physiological Solution
- 25 µl Pipette

REF	0211	STREPTOSLIDE A
REF	0212	STREPTOSLIDE B
REF	0213	STREPTOSLIDE C
REF	0214	STREPTOSLIDE D
REF	0215	STREPTOSLIDE F
REF	0216	STREPTOSLIDE G

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
REAGENT 1	3 months
REAGENT 2	3 months
REAGENT 3	3 months
LYTIC ENZYME	5 weeks
POSITIVE CONTROL	3 months

6. PRECAUTIONS**FOR IN VITRO DIAGNOSTIC USE ONLY. STORE AT 2-8°C**

All reagents must be handled in accordance with the safety precautions adopted in the laboratory. All samples must be considered potentially infectious.

Waste disposal: samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

- ATTENTION: Some of the reagents necessitate particular attention, as reported on the relative labels.**

Reagent 1:

Signal Word: Danger

Contains: Sodium Nitrate

INDEX No. 007-010-00-4

Hazard Statements: H272 – May intensify fire; oxidiser; H301 – Toxic if swallowed; H400 – Very toxic to aquatic life

Precaution statements: P210 – Keep away from heat/sparks/open flames/hot surfaces. — No smoking; P220 – Keep/Store away from clothing/.../combustible materials; P221 – Take any precaution to avoid mixing with combustibles...; P264 – Wash ... thoroughly after handling; P270 – Do not eat, drink or smoke when using this product; P273 – Avoid release to the environment; P280 – Wear protective gloves/protective clothing/eye protection/face protection; P301+P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician; P330 – Rinse mouth; P370+P378 – In case of fire: Use water, foam, CO₂ for extinction; P391 Collect spillage; P405 Store locked up; P501 Dispose of contents/container according to law

	Section
Changes introduced in the current revision	4, 6, Glossary

Reagent 2:Signal Word: CautionContains: Acetic Acid

INDEX No. 607-002-00-6

Hazard Statements: H315 – Causes skin irritation; H319 – Causes serious eye irritationPrecaution statements: P264 – Wash ... thoroughly after handling; P280 – Wear protective gloves/protective clothing/eye protection/face protection; P302+P352 – IF ON SKIN: Wash with plenty of soap and water; P305+P351+P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.Reagent 3:Signal Word: CautionHazard Statements: H315 – Causes skin irritation; H319 – Causes serious eye irritationPrecaution statements: P264 – Wash ... thoroughly after handling; P280 – Wear protective gloves/protective clothing/eye protection/face protection; P302+P352 – IF ON SKIN: Wash with plenty of soap and water; P305+P351+P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.Lytic Enzyme:Signal Word: DangerHazard Statements: H300 – Fatal if swallowed; H315 – Causes skin irritation; H319 – Causes serious eye irritation; H335 – May cause respiratory irritation; H411 – Toxic to aquatic life with long lasting effects.Precaution statements: P273 – Avoid release to the environment; P280 – Wear protective gloves/protective clothing/eye protection/face protection; P301+P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician; P302+P352 – IF ON SKIN: Wash with plenty of soap and water; P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing; P305+P351+P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
EUH032 – Contact with acids liberates very toxic gas.

- Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- The reagents contain 0.9% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.
- If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
- Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.
- Reagent 3 is alkaline (pH>9) and can develop ammoniac fumes. **Avoid contact with the mucosa.**

Analytical precautions

- Allow all reagents and samples to come to room temperature (18-30°C) before use.** Immediately after use return reagents to the recommended storage temperature.
- Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable.
- The reaction can be performed on a rotating mixer. Regulate the speed at 100 ± 10 rpm.
- The agglutination reaction is best interpreted by examining the slide under a lamp.
- Make a Gram coloration to confirm the presence of Streptococci.

7. SAMPLES

Streptococci from primary cultures on solid media can be used, if present in well isolated colonies. Blood agar is the best medium, allowing the haemolytic reaction to be observed before performing the group identification procedure. Streptococci from broth culture may also be used.

	Section
Changes introduced in the current revision	4, 6, Glossary

8. PROCEDURE**- REACTION FOR IDENTIFICATION OF GROUPS A, B, C, F and G:**

Dispense one drop of Reagent 1 and one drop of Reagent 2 in a small pipette, add 1-5 well isolated colonies (or 25 µL of broth culture) and mix until a homogeneous suspension is obtained. Add 2 drops of Reagent 3 and mix. Shake the Grouping latexes A, B, C, F and G carefully. Dispense one drop of the first latex on the first reaction circle and subsequently one drop of all the others on the next circles. Place 25 µL extract from the tube on each circle, mix carefully with a mixing stick and rotate the slide for 1 minute. Watch each circle for agglutination (reaction showing positive identification). If the presence of Group D Streptococci is suspected, perform the direct test for suspected colonies with the Grouping latex D. If the result is negative, perform the enzymatic extraction.

- DIRECT REACTION FOR GROUP D:

Dilute 4-5 isolated colonies together with one drop of Physiological Solution on a reaction circle. Add one drop of Grouping latex D, mix with a stick, rotate and watch for the possible agglutination within one minute. A clear agglutinate indicates colonies of Group D Streptococci. This procedure identifies most of the Group D strains.

- ENZYMATIC EXTRACTION FOR GROUP D:

Using the Lytic Enzyme, the highest percentage of correct identification of strains belonging to Group D is attained. Reconstitute the freeze-dried reagent with 2 mL of distilled water and dispense 100 µL (2 drops) in a small tube. Dilute 1-5 colonies (collected with the ansa or a stick) in the enzymatic solution and incubate at 37° C for 10 minutes. Dispense one drop of Grouping latex D and 25 µL of extract on a reaction circle, mix and rotate the slide. A clear agglutinate obtained within one minute indicates a positive reaction for the identification of the Group D Streptococci.

9. TEST VALIDATION

Distribute 25 µL of Polyvalent Positive Control on each test card circle and add one drop of each latex; positive reactions should result with each Grouping latex (evident agglutination) within one minute.

10. INTERPRETATION

The group to which the streptococcal strain under test belongs is identified by the latex which gives a clear agglutination according to the method described above.

11. LIMITATIONS OF THE TEST

If the operator is certain of the identification of the family, the absence of agglutination in the conditions described in the procedure could be due to an insufficient sample. Non specific agglutination of all the latexes may be due to the presence of Fc receptors on the bacterial wall.

12. PERFORMANCE

Using strains obtained both from clinical isolates and from strain collections, it has been shown that the specificity and sensitivity of the test is 100%.

13. REFERENCES

1. Topley and Wilson's VII Ed. Vol. 2 pag. 173-217 (1983)
2. Slifkin M., Gil G.M. J. Clin. Microbiol. 15, 187-189 (1982)
3. Calandra G.B., Henson C.L, J. Clin. Microbiol. 15, 184-189 (1982)
- 4 Slifkin M., Cumbie R.J. Clin. Microbiol. 25, 1555-1556 (1987).

	Section
Changes introduced in the current revision	4, 6, Glossary

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