DIESSE Diagnostica Senese SpA

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

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DIESSE Diagnostica Senese SpA

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

INSTRUCTIONS FOR USE**WW SLIDE****REF** 021**(English)****1. INTENDED USE****STAINED SALMONELLA AND BRUCELLA SUSPENSIONS TO DETECT SPECIFIC AGGLUTINATING ANTIBODIES IN SERUM.****2. INTRODUCTION**

Serological studies of patients with persistent fever can be useful in detecting typhoid fever and brucellosis (1,2).

The test consists in the observation of a direct agglutination reaction between the bacterial suspension and the serum sample which is suspected of containing specific antibodies. The present method, known as the Widal-Wright reaction, allows a semiquantitative titration of the serum agglutinins.

3. PRINCIPLE OF THE METHOD

The serum sample is mixed with a bacterial suspension at the optimal concentration for the slide test. After 1 minute a visible agglutination will appear if the serum contains specific antibodies against the bacteria in the suspension. Positive samples must then be tested in the tube test, using one drop of the sale suspension.

4. KIT COMPOSITION

- The reagents are sufficient for 100 tests.
- Bring the reagents to room temperature before use.

S. TYPHI O**S. TYPHI H****S. PARATYPHI AO****S. PARATYPHI AH****S. PARATYPHI BO****S. PARATYPHI BH****BRUCELLA****BACTERIAL SUSPENSION**

BACTERIAL SUSPENSION 1 vial x 5 mL of each suspension

Contents: Suspension of stained, killed bacteria at an optimized concentration for performance of the slide test. Contains sodium azide 0,09% as preservative.**CONTROL +** POSITIVE CONTROL 1 vial x 1 mL**Contents:** Polyvalent positive control serum of animal origin, containing sodium azide 0,09% as preservative.**CONTROL -** NEGATIVE CONTROL 1 vial x 1 mL**Contents:** Negative control serum of animal origin, containing sodium azide 0,09% as preservative.**OTHER MATERIALS**

- 15 Slides. Only use new slides.
- Mixing sticks

OTHER MATERIALS REQUIRED BUT NOT SUPPLIED

- Pipettes with disposable needles.
- Physiological saline.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation**REAGENT****CONDITIONS**

BACTERIAL SUSPENSION

2 months at 2/8°C

CONTROLS

2 months

	Section
Changes introduced in the current revision	8

6. PRECAUTIONS**FOR IN VITRO DIAGNOSTIC USE ONLY.**

All materials of biological origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The reagents contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.
3. If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. **Allow all reagents and samples to come to room temperature (18-30°C) before use.** After use, empty the dropper and return reagents to the recommended storage temperature (2-8°C).
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure.

7. TYPE AND STORAGE OF SAMPLE

Fresh serum may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Use immediately after thawing.

8. PROCEDURE

Bring the reagent to room temperature and mix the vial of bacterial suspension carefully but not vigorously.

*** SCREENING TEST**

Distribute 20 µL of serum in a circle on the slide. Add one drop of suspension and mix with the serum. Rotate the slide. Read the results after exactly 1 minute. Positive samples will give agglutination while negative samples will remain homogeneous.

*** TUBE TITRATION**

Dilute the sample 1:20 with physiological saline (0.1 mL serum + 1.9 mL saline). Take 1 mL from the first tube and make serial dilutions with 1 mL of saline, as shown in the scheme. Add one drop of antigen suspension to each tube and incubate overnight at 37°C. See the scheme shown below. The blank is performed by mixing the suspension with saline

Tube	1	2	3	4	5	6	7
Saline mL	1.9	1	1	1	1	1	1
Serum ml	0,1						
Mix and transfer ml							
Dilution / Titer	1/20	1/40	1/80	1/160	1/320	1/640	1/1280

9. SCREENING TEST VALIDATION

Perform the slide test using 20 µL of **Positive Control Serum**. Evident agglutination must be seen.

	Section
Changes introduced in the current revision	8

10. INTERPRETATION OF THE TUBE TEST – NORMAL VALUES

Positive reaction

The bacteria are agglutinated in a granular (antigen O) or flocculant (antigen H) manner; the fluid supernatant is clear, in comparison with the blank.

Negative reaction

The suspension appears homogeneous as in the control tube.

The titer is given by the highest dilution able to cause agglutination of at least 50% of the bacteria present.

Titers of 1/80 are considered suspect; higher titers are diagnostically significant.

11. TEST LIMITATIONS

In slide agglutination using undiluted serum, the prozone phenomenon may occur. If the sample is highly suspected of being positive, repeat the test diluting 1:8 in saline. The brucella suspension does not agglutinate in the presence of incomplete antibodies.

In the slide test, agglutination of H antigens may also be caused by antibodies which react with the lipopolysaccharide. Differentiation is best obtained by examining the flaky agglutination in the tube.

12. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

An experimentation performed on human sera obtained from hospital laboratories, the titers of which were taken as reference, allowed an evaluation of the sensitivity of the following suspensions: *S. typhi* O, *S. typhi* H, *S. paratyphi* BO, *S. paratyphi* BH.

When *S. typhi* antigen H was determined on 30 positive samples, a sensitivity of 96.8% was found, while the suspension *S. typhi* antigen O tested on 51 positive sera gave 2 false negatives and a sensitivity of 96.2%.

S. paratyphi B antigen H and antigen O were tested with respectively 12 and 24 positive samples. Sensitivity was 100% and 91.7% respectively.

S. paratyphi A antigen O and antigen H were tested in parallel with another commercial product. 12 positive and 3 doubtful results were found. The Diesse reagent gave 1 false negative and the sensitivity was therefore 95.0%, considering the doubtful results as positive. Using the same criteria, 45 positive samples for Brucellosis were tested. Diesse gave 2 false negative results and the sensitivity was therefore 95.7%.

Human sera found negative in another commercial test were also studied. The following results were obtained:

<u>Suspension</u>	<u>No. of samples</u>	<u>Specificity %</u>
<i>S. typhi</i> antigen O	93	97.4
<i>S. typhi</i> antigen H	58	95.1
<i>S. paratyphi</i> B antigen O	96	97.0
<i>S. paratyphi</i> B antigen H	93	95.9
<i>S. paratyphi</i> A antigen O	82	98.8
<i>S. paratyphi</i> A antigen H	80	97.6
Brucella	83	97.6

13. REFERENCES

1. P. Nicoletti in F. Pasquinelli: Diagnostica e Tecnica di Laboratorio Vol. 2 p. 319. Rosini ed. Firenze 1981.
2. Edwards P.R., Ewing W.H. Identification of Enterobacteriaceae 3rd ed. Burgess Publishing Co.
3. Cruickshank R. Medical Microbiology 11th Ed., p. 907 (1965).
4. Huddleson I.F. and Abell E. J. Inf. Dis. 42: 242 (1928).

	Section
Changes introduced in the current revision	8

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