DIESESE Diagnostica Senese SpA

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(English)

1. INTENDED USE

NON-TREPONEMAL FLOCCULATION TEST FOR THE QUALITATIVE AND SEMIQUANTITATIVE DETERMINATION OF REAGINS IN SERUM OR PLASMA

2. INTRODUCTION

Infection with *Treponema pallidum* gives rise to a complex antibody response. Non treponemal tests are based on the detection of reagins, an antibody class present in syphilis and occasionally in other acute and chronic conditions (1). Reagins can be detected in the serum about 4-6 weeks after infection, or 1-3 weeks after appearance of the chancre. Purified beef heart extract (cardiolipin) fortified with lecithin and cholesterol, is used as antigen for reagin tests. Non treponemal tests are routinely used in syphilis serology because of their practical advantages and reproducibility, even though they are not always 100% specific.

3. PRINCIPLE OF THE METHOD

The test is performed on serum or plasma. The antigen is composed of a colloidal suspension of cardiolipin, lecithin and cholesterol mixed with carbon microparticles. When a reactive serum is mixed with the antigen on the test card, a macroscopical flocculant is formed.

4. KIT CONTENTS AND REAGENT PREPARATION

- Bring reagents to room temperature before use.

ANTIGEN ANTIGEN 2 vials x 3.2 mL

Contents: The antigen is composed of a stabilised mixture of cardiolipin, cholesterol, lecithin and carbon microparticles. Attention: Store in an upright position when not in use. Bring to room temperature before use. Do not shake too vigorously.

CONTROL + POSITIVE CONTROL SERUM 0.5 mL

Contents: Human serum, reactive for syphilis, diluted in a solution of BSA with sodium azide 0,09%; liquid, ready for use without further dilution.

CONTROL - NEGATIVE CONTROL SERUM 0.5 mL

Contents: Human serum, non reactive for syphilis, diluted in a solution of BSA with sodium azide 0,09%; liquid, ready for use without further dilution.

OTHER MATERIALS

- No.6 Flat-bottomed microtiter plates.

MATERIALS REQUIRED BUT NOT PROVIDED

- Mechanical rotor for agglutination reactions, adjustable to 100 ± 4 rpm.
- Light source
- 50 and 10 μ L micropipette (variable volume)
- 0.9% Physiological Saline solution

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT CONDITIONS

ANTIGEN	2 months
POSITIVE CONTROL	4 months
NEGATIVE CONTROL	4 months

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY. STORE AT 2-8°C

Caution: *This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.*

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The reagents contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.
If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. **Allow all reagents and samples to come to room temperature (18-30°C) before use.** Immediately after use return reagents to the recommended storage temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable.

7. TYPE AND STORAGE OF SAMPLE

Fresh or defrosted serum. Fresh samples may be stored at 2/8°C for 5 days. Plasma collected with citrate, EDTA or heparin may be used and stored at 2/8°C up to 2 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers.

8. PROCEDURE

SCREENING TEST

1. Bring reagents and samples to room temperature. Resuspend the reagent.
2. Draw up 25 µL of test serum and drop onto a microplate well. Take care not to transfer any cellular elements.
3. Repeat these steps for each sample to be tested.
4. Distribute 12,5 µL of reagent in each well.
5. Mix the reagent and sample for 1 min using a plate shaker at max. speed.
6. Place the plate on the mechanical rotor and rotate at 100 ± 4 rpm for 10 min.
7. Examine the result under a good light source. Reactions must be read within 10 minutes.

QUANTITATIVE TEST

Make at least 5 doubling dilutions of the test sample using physiological saline. If a microplate is used, 25 µL of sample can be diluted in progressive wells after placing 25 µL of saline in each well. Discard 25 µL of sample diluted from the last well. The method is analogous to the screening method starting from point 4.

9. TEST VALIDATION

Use the controls in each session, dispensing 25 µl of each control into a well of the microplate. Proceed as described in the section "Test Procedure" below.

10. INTERPRETATION OF THE TEST

1. Qualitative results of the screening test

a) REACTIVE:

The presence of clearly visible black aggregates is typical of a positive sample. The positive control must give evident agglutination.

b) NON REACTIVE:

No aggregate formation.

The Negative Control must not give aggregate formation.

NOTE: Any sample giving a reactive result should be tested in the quantitative method.

2. Quantitative Test

The result is given by the highest dilution showing any reactivity.

NOTE: In cases where the patient's history or clinical signs contrast with the diagnosis, a confirmatory treponemal test should be performed.

N.B. False-negative or false-positive results may occur if the described procedure is not followed very carefully.

11. LIMITATIONS OF THE TEST

Biological false-positive reactions occur occasionally with cardiolipin antigens, e.g. in cases of drug abuse, or in certain diseases such as lupus erythematosus, mononucleosis, malaria, leprosy.

A pro-zone reaction is sometimes found: complete or partial inhibition of reactivity occurs with undiluted serum and maximum reactivity is found only with diluted serum. This can lead to a weakly positive result in the qualitative test, whereas a strong reaction will be found when the sample is diluted. For this reason, all specimens showing any reactivity should be tested in the quantitative procedure.

12. DIAGNOSTIC SENSITIVITY AND SPECIFICITY


272 samples were tested, 200 of which were negative and 72 positive. Another commercial FDA-approved method was used in comparison. In the trial, results which were +/- between the two methods were considered in agreement. 1 because they came from healthy donors who were negative to TPHA and to an ELISA screening test. The product was compared with another commercial kit. The results can be summarized as follows:

	Pos.	Neg.	FP	FN
Diesse	71	199	1	1
Other method	72	200		

The sensitivity was 98.6% while the specificity was 99,5% in comparison with the other commercial kit.

13. REFERENCES

1. Manual of Tests for Syphilis. Public Health Service, Publication no. 411 (1969).
2. Portnoy J., Bossak H.W. et al. Public Health Reports 76, 933 (1961).
3. March R.W. et al. Abstracts of Ann. Meeting of Am. Soc. Microbiol. p. 93 (1974).
4. Bossak H.N., Harris A.D. and Olansky S. Brit. J. Ven. Dis. 231, 33 (1955).
5. Dyckman J.D., Storms S. and Huber T.W. J. Clin. Microbiol. 12, 629 (1980).
6. Fiumara N.J. JAMA 243, 2500 (1980).
7. Garner M.F., Backhouse J.L. et al. J. Clin. Path. 26, 258 (1973).
8. Hicks C.B. et al. Ann. Int. Med. 107, 492 (1987).
9. Luger A.F.H. in "Immunological Diagnosis of sexually transmitted disease (Clinical and Biochemical analysis; 23)". New York, N.Y. Dekker, 1988, 249.

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusque Ημερομηνία λήξης Prazo de validade
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbicante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

Diessa Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) – Italy
Tel. 0577-587111