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## 19. BIBLIOGRAFIA

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.
4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics,

## INSTRUCTIONS FOR USE

### CHORUS Promonitor ANTI-VEDOLIZUMAB

**For the semiquantitative detection of anti-Vedolizumab antibodies in human serum**

**For *In Vitro* Diagnostic Use Only**

#### 1. INTENDED USE

CHORUS Promonitor ANTI-VEDOLIZUMAB is an immunoassay kit for automated semiquantitative detection of anti-Vedolizumab antibodies in human serum using a disposable device applied on the CHORUS TRIO instrument.

The product, used in conjunction with other clinical and laboratory findings, is useful as an aid in the management of patients treated with Vedolizumab (VDZ), like patients with inflammatory bowel disease (IBD) and rheumatic diseases.

It must be used by professional laboratory users only.

#### 2. INTRODUCTION

Vedolizumab (VDZ) is a humanized IgG1 monoclonal antibody expressed in Chinese hamster ovary (CHO) cells that specifically binds the  $\alpha 4\beta 7$  integrin which is expressed on a discrete subset of memory T helper lymphocytes, which preferably migrate to the gastrointestinal (GI) tract, causing inflammation characteristic of ulcerative colitis and Crohn's disease. VDZ reduces inflammation by improving symptoms, characterized by more acute phases followed by clinical remission, accompanied by abdominal pain, diarrhea and bleeding that limit the social and work activities of the patients who are affected.

VDZ is a biological immunosuppressant with high selectivity for intestinal tract and is indicated for the treatment of adult patients with moderate to severe active ulcerative colitis or Crohn's disease who have had an inadequate response or been intolerant to conventional therapy, or to a tumor necrosis factor alpha (TNF $\alpha$ ) antagonist.

VDZ can potentially be immunogenic and induce the formation of anti-drug antibodies (ADA), thereby potentially reducing treatment efficacy. A gradual decrease in efficacy over time following an initial response to biologics is common.

#### 3. PRINCIPLE OF THE METHOD

The CHORUS Promonitor ANTI-VEDOLIZUMAB device is ready to use for the detection of anti-Vedolizumab antibodies in human serum, in the CHORUS TRIO instrument.

The assay is based on a bridging ELISA method.

The bridging antibody ELISA method is widely used for anti-drug antibody screening and is well accepted as a standard for immunogenicity screening.

The microwell strips are pre-coated with Vedolizumab (VDZ). Calibrator, positive control and samples are added to separate device.

Anti-VDZ antibodies are able to bind to preimmobilized VDZ following incubation with the sample. After washing to remove the unbound sample, HRP-labelled VDZ (conjugate) is added. A second incubation step allows the conjugate to bind to the anti-VDZ antibodies, previously attached to the microwells.

After washing to remove the unbounded conjugate, the peroxidase substrate is added. The enzymatic reaction is stopped by adding Stop Solution, which makes the solution turns to yellow. The color which develops is proportional to the amount of anti-VDZ antibodies in the patient sample.

The disposable devices contain all the reagents to perform the test when applied on the CHORUS TRIO instrument.

The results are expressed in Arbitrary Units (AU/ml).

#### 4. WARNINGS AND PRECAUTIONS

##### FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

**Waste disposal:** serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

##### Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

##### Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the devices within 60 min.

1. **Discard devices which show the substrate (well 3) blue colored.**


2. If the sample is added into the well manually, verify that it is perfectly distributed.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the CHORUS TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.  
**The use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<https://www.diesse.it/en/downloads/downloads/strumento:39/>)**
5. Check that the CHORUS TRIO instrument is set up correctly (see CHORUS TRIO Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see CHORUS TRIO Operating Manual).
9. When CHORUS Autosampler (REF 81210) is used, refer to the related Operating Manual
10. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
11. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
12. Do not use the device after the expiry date.
13. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

## 5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

**DD** DEVICES 6 packages each containing 6 devices

Description:

**Position 9:** Space for application of bar code label

**Position 8:** EMPTY

**Position 7:** MICROPLATE WELL

Coated with VDZ

**Position 6:** MICROPLATE WELL

Coated with VDZ

**Position 5:** STOP SOLUTION

Contents: 0.3 M sulfuric acid solution

**Position 4:** CONJUGATE

Contents: HRP-labelled VDZ in buffer containing preservative

**Position 3:** TMB SUBSTRATE

Contents: Tetramethylbenzidine stabilized in citrate buffer

**Position 2:** DILUTION BUFFER

Contents: Buffer solution containing preservative

**Position 1:** EMPTY WELL

in which the sample is transferred.

**Use: equilibrate a package at room temperature**, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

**CALIBRATOR** CALIBRATOR **1 x 0.650 ml**  
Contents: Human serum containing anti-VDZ antibodies, and preservative.  
Liquid, ready for use.

**CONTROL +** POSITIVE CONTROL **1 x 1.100 ml**  
Contents: Human serum containing anti-VDZ antibodies, and preservative.  
Liquid, ready for use.

## MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

## 6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

|                  |                  |
|------------------|------------------|
| DEVICES          | 8 weeks at 2/8°C |
| CALIBRATOR       | 8 weeks at 2/8°C |
| POSITIVE CONTROL | 8 weeks at 2/8°C |

## 7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 5 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be mixed carefully before use.

It is strongly recommended to avoid any harsh treatment of the samples.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

## 8. ASSAY PROCEDURE

The test CHORUS Promonitor ANTI-VEDOLIZUMAB (REF. 86711) can be performed simultaneously **ONLY** with the other CHORUS Promonitor tests.

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense the following volumes in well no. 1 of each device:

|                  |               |
|------------------|---------------|
| SAMPLE           | 130 µl/device |
| CALIBRATOR       | 130 µl/device |
| POSITIVE CONTROL | 130 µl/device |

At each change of batch, use a device for the calibrator.

4. Place the devices in the CHORUS TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

## 9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554  
Fax: 0039 0577 366605  
email: scientificsupport@diesse.it

## 10. INTERPRETATION OF THE RESULTS

The CHORUS TRIO instrument expresses the result in AU/ml, calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is  $\geq 9.7$  AU/ml

NEGATIVE: when the result is  $\leq 8.0$  AU/ml

EQUIVOCAL: when the result is between 8.1 and 9.6 AU/ml

Equivocal results are obtained in case of values which ranges between Limit of detection ( $< 8.0$  AU/ml) and Low Limit of Quantitation (9.7 AU/ml).

## 11. LIMITATIONS

1. Trough antibody levels should be measured. The optimum sample collection time is immediately before the VDZ infusion.
2. This assay measures free anti-VDZ antibodies. It does not measure immunocomplexes formed between VDZ and anti-VDZ antibodies.
3. This assay is not intended for use as an aid in diagnosis, neither as a companion diagnostic test. The clinician should interpret test results in the context of other variables like sampling window, methodology, VDZ levels and other clinical findings like disease activity.

4. The assay performance characteristics have not been established for matrices other than serum.
5. This product should only be used by appropriately trained laboratory professional personnel.

## 12. CALIBRATION RANGE

Calibration range : 8.0-684.0 AU/ml

For samples  $> 684.0$  AU/ml retest the diluted sample in the Negative Control/Sample Diluent (REF 83607 - not supplied with the kit).

## 13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, are negative, with values  $< 8.0$  AU/ml.

## 14. ANALYTICAL SPECIFICITY

3 samples (1 Negative, 1 Low Positive and 1 High Positive) were spiked with the following potentially interfering factors and then tested:

Bilirubin (4.5 mg/dl – 45 mg/dl)

Triglycerides (250 mg/dl – 1500 mg/dl)

Hemoglobin (2.5 mg/ml - 10 mg/ml)

The presence of the above interfering substances in the test serum does not alter the test result (with the exception of hemoglobin at concentration higher than 5 mg/ml).

2 Rheumatoid Factor (RF) positive samples (200 IU/ml) were spiked with Anti-Vedolizumab:

High level – 146.9 AU/ml

Low level – 68.6 AU/ml

Negative –  $< 8.0$  AU/ml

The presence of RF in the test sera does not alter the test result.

## 15. CROSS-REACTIONS

10 samples, not containing anti- vedolizumab antibodies, but containing anti-IFX antibodies (4000.0 AU/ml) (5 samples) and anti-ADL antibodies (4000.0 AU/ml) (5 samples) were tested.

No cross-reactions were found.

## 16. METHOD COMPARISON

In an experimentation 164 samples have been tested with Diesse kit and with a competitor kit (Promonitor antiVDZ ELISA kit).

Data are summarized in the following table:

|        |      | Reference |     |      |
|--------|------|-----------|-----|------|
|        |      | +         | -   | Tot. |
| Diesse | +    | 40        | 0   | 40   |
|        | -    | 0         | 124 | 124  |
|        | Tot. | 40        | 124 | 164  |

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% IC<sub>95%</sub>: 91.2-99.9

Percent Negative Agreement (~Diagnostic Specificity):

100.0% IC<sub>95%</sub>: 97.0-100.0

Positive Predictive Value (PPV): 100.0% IC<sub>95%</sub>: 100.0-100.0

Negative Predictive Value (NPV): 100.0% IC<sub>95%</sub>: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

## 17. PRECISION AND REPEATABILITY

| Sample | Within run   |     | Between run  |      |
|--------|--------------|-----|--------------|------|
|        | Mean (AU/ml) | CV% | Mean (AU/ml) | CV%  |
| 1      | 610.8        | 2.8 | 618.1        | 5.1  |
| 2      | 209.3        | 6.8 | 216.2        | 11.4 |
| 3      | 40.1         | 8.0 | 44.3         | 5.6  |
| 4      | 24.1         | 4.6 | 28.3         | 8.1  |
| 5      | <8.0         | -   | <8.0         | -    |
| 6      | <8.0         | -   | <8.0         | -    |

| Sample | Between lots |     | Between instruments |      |
|--------|--------------|-----|---------------------|------|
|        | Mean (AU/ml) | CV% | Mean (AU/ml)        | CV%  |
| 1      | 576.2        | 7.7 | 642.7               | 5.8  |
| 2      | 233.6        | 5.5 | 239.8               | 13.4 |
| 3      | 45.2         | 7.6 | 43.3                | 10.6 |
| 4      | 27.0         | 7.8 | 26.7                | 11.2 |
| 5      | <8.0         | -   | <8.0                | -    |
| 6      | <8.0         | -   | <8.0                | -    |

## 18. LINEARITY

The linearity of the test was evaluated with 3 samples with known concentration of anti-Vedolizumab antibodies. For each sample 8 serial dilutions were done; the linearity was demonstrated in the range 8.0-675.5 AU/ml.

## 19. BIBLIOGRAPHY

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.
4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical

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### CHORUS Promonitor ANTI-VEDOLIZUMAB

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## 19. BIBLIOGRAFIE

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.
4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

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#### 19. LITERATUR

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.

2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.

3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.

4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

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1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. *J Crohns Colitis*. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." *Inflamm Bowel Dis*. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. *N Engl J Med*. 2013 Aug 22;369(8):711-21.
4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. *Clin Pharmacokinet*. 2017 Nov;56(11):1287-1301.

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## 19. BIBLIOGRAFIA

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.

3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 2013 Aug 22;369(8):711-21.

4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

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## 19. BIBLIOGRAFIE

1. Gisbert JP, Chaparro M. Predictors of Primary Response to Biologic Treatment [Anti-TNF, Vedolizumab, and Ustekinumab] in Patients With Inflammatory Bowel Disease: From Basic Science to Clinical Practice. J Crohns Colitis. 2020 Jun 19;14(5):694-709.
2. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
3. Sandborn WJ, Feagan BG, Rutgeerts P, Hanauer S, Colombel JF, Sands BE, Lukas M, Fedorak RN, Lee S, Bressler B, Fox I, Rosario M, Sankoh S, Xu J, Stephens K, Milch C, Parikh A; GEMINI 2 Study Group. Vedolizumab as induction and maintenance therapy for Crohn's disease. N Engl J Med. 22 Aug 2013;369(8):711-21.
4. Rosario M, Dirks NL, Milch C, Parikh A, Bargfrede M, Wyant T, Fedyk E, Fox I. A Review of the Clinical Pharmacokinetics, Pharmacodynamics, and Immunogenicity of Vedolizumab. Clin Pharmacokinet. 2017 Nov;56(11):1287-1301.

|                                                                                     |                                  |                                                                                                                                                                                                                                                                         |                            |                                                                                                                                                                                                                                         |
|-------------------------------------------------------------------------------------|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Date of manufacture<br>Data di fabbricazione<br>Herstellungsdatum<br>Fecha de fabricación<br>Datum výroby<br>Data di fabbricazione                                                                                                                                      | FR<br>GR<br>PT<br>PL<br>RO | Date de fabrication<br>Ημερομηνία Παραγωγής<br>Data de fabrico<br>Data produkcji<br>Data fabricației                                                                                                                                    |
|    | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Use By<br>Utilizzare entro<br>Verwendbar bis<br>Fecha de caducidad<br>Použijte do<br>Utilizzare entro                                                                                                                                                                   | FR<br>GR<br>PT<br>PL<br>RO | Utiliser jusque<br>Ημερομηνία λήξης<br>Prazo de validade<br>Zużyć do dnia<br>A se folosi până la                                                                                                                                        |
|    | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Do not reuse<br>Non riutilizzare<br>Nicht wieder verwenden<br>No reutilizar<br>Nepoužívejte znovu<br>Nu reutilizați                                                                                                                                                     | FR<br>GR<br>PT<br>PL<br>RO | Ne pas réutiliser<br>Μην κάνετε επαναληπτική χρήση<br>Não reutilizar<br>Nie używać ponownie<br>A nu se refolosi                                                                                                                         |
|    | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Caution, consult accompanying documents<br>Attenzione, vedere le istruzioni per l'uso<br>Achtung, die Gebrauchsanleitung lesen<br>Atención, ver instrucciones de uso<br>Pozor, řiďte se návodem k použití<br>Attenzione, vedere le istruzioni per l'uso                 | FR<br>GR<br>PT<br>PL<br>RO | Attention voir notice d'instructions<br>Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα<br>Atenção, consulte a documentação incluída<br>Uwaga, patrz instrukcja obsługi<br>Atenție, consultați documentele însoțitoare                   |
|    | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Manufacturer<br>Fabbicante<br>Hersteller<br>Fabricante<br>Výrobce<br>Fabbicante                                                                                                                                                                                         | FR<br>GR<br>PT<br>PL<br>RO | Fabricant<br>Κατασκευαστής<br>Fabricante<br>Producuent<br>Producător                                                                                                                                                                    |
|  | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Contains sufficient for <n> tests<br>Contenuto sufficiente per "n" saggi<br>Inhalt reicht für „n“ Tests<br>Contenido suficiente para <n> ensayos<br>Dostatečný obsah pro „n“ testů<br>Contenuto sufficiente per "n" saggi                                               | FR<br>GR<br>PT<br>PL<br>RO | Contenu suffisant pour "n" tests<br>Περιεχόμενο επαρκές για «n» εξετάσεις<br>Conteúdo suficiente para "n" ensaios<br>Wystarczająca ilość materiału do „(liczby)” badań<br>Conținut suficient pt <n> teste                               |
|  | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Temperature limitation<br>Limiti di temperatura<br>Temperaturgrenzwerte<br>Limite de temperatur<br>Limity teploty<br>Limiti di temperatura                                                                                                                              | FR<br>GR<br>PT<br>PL<br>RO | Limites de température<br>Περιορισμοί θερμοκρασίας<br>Limites de temperatura<br>Granice temperatury<br>Limită de temperatură                                                                                                            |
|  | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Consult Instructions for Use<br>Consultare le istruzioni per l'uso<br>Die Gebrauchsanleitung lesen<br>Consulte las instrucciones de uso<br>Řiďte se návodem k použití<br>Consultare le istruzioni per l'uso                                                             | FR<br>GR<br>PT<br>PL<br>RO | Consulter les instructions d'utilisation<br>Συμβουλευτείτε τις οδηγίες χρήσης<br>Consulte as instruções de utilização<br>Należy zapoznać się z instrukcją obsługi<br>Pentru utilizare consultați instrucțiunile                         |
|  | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Catalogue number<br>Numero di catalogo<br>Katalognummer<br>Número de catálogo<br>Číslo katalogu<br>Numero di catalogo                                                                                                                                                   | FR<br>GR<br>PT<br>PL<br>RO | Référence du catalogue<br>Αριθμός καταλόγου<br>Referência de catálogo<br>Numer katalogowy<br>Număr de catalog                                                                                                                           |
|  | EN<br>IT<br>DE<br>ES<br>CS<br>RO | In Vitro Diagnostic Medical Device<br>Dispositivo medico-diagnostico in vitro<br>Medizinisches In-vitro-Diagnostikum<br>Producto sanitario para diagnóstico in vitro<br>Zdravotnický nástroj k diagnostice in vitro<br>Dispozitiv medical pentru diagnosticare in vitro | FR<br>GR<br>PT<br>PL<br>RO | Dispositif médical de diagnostic in vitro<br>In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν<br>Dispositivo médico para diagnóstico in vitro<br>Wyrób medyczny do diagnostyki in vitro<br>Dispozitiv medical pentru diagnosticare in vitro |
|  | EN<br>IT<br>DE<br>ES<br>CS<br>RO | Batch code<br>Codice del lotto<br>Chargennummer<br>Código de lote<br>Kód šarže<br>Codice del lotto                                                                                                                                                                      | FR<br>GR<br>PT<br>PL<br>RO | Code du lot<br>Αριθμός Παρτίδας<br>Código do lote<br>Kod partii<br>Lot                                                                                                                                                                  |
|  | EN<br>IT<br>DE<br>ES<br>CS<br>RO | CE marking of conformity<br>Marcatura CE di conformità<br>CE-Konformität Skenneichnung<br>Marcado CE de conformidad<br>Označení shody CE<br>Marcajul de conformitate CE                                                                                                 | FR<br>GR<br>PT<br>PL<br>RO | Marquage de conformité CE<br>Σημανση συμμορφωσης CE<br>Marcação CE de conformidade<br>Oznakowanie zgodności CE<br>Marcajul de conformitate CE                                                                                           |