# CHORUS VITAMIN EXTRACTION

**REF 86500**

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## **INSTRUCTIONS FOR USE**

### **CHORUS VITAMIN EXTRACTION**

#### **Kit for extraction of samples to be tested with the CHORUS VIT B12 and CHORUS FOLATE kits**

#### **For *in vitro* diagnostic use only**

##### **1. INTENDED USE**

Kit for the extraction of samples to be tested with the CHORUS VIT B12 (REF. 86400) and CHORUS FOLATE (REF. 86402) tests.

It must only be used by professional laboratory personnel.

##### **2. INTRODUCTION**

The kit allows the correct treatment for the extraction of the samples needed to perform the CHORUS VIT B12 and CHORUS FOLATE tests.

##### **3. PRINCIPLE OF THE METHOD**

The kit consists of three reagents to be used for the extraction of the samples on which to dose vitamin B12 and/or folate.

The extraction process allows these vitamins to free themselves from the carrier proteins to which they are normally bound for their absorption and functionality so that they can be correctly detected.

After performing the sample extraction procedure, it is ready to be dosed.

For analysis of the extracted sample, refer to the instructions for use of the CHORUS VIT B12 (REF. 86400) and/or CHORUS FOLATE (REF. 86402) test.

##### **4. PRECAUTIONS**

#### **FOR *IN VITRO* DIAGNOSTIC USE ONLY.**

Wait for the serum samples to reach room temperature (18-30°C) before performing the extraction procedure.

The used, deteriorated or expired product must be treated as potentially infected residue, therefore, disposed of in accordance with the provisions of current laws.

#### **Personnel safety warnings**

- Any material of human origin must be considered potentially infected. All samples must be handled in accordance with the safety rules normally adopted in the laboratory.

- Use single-use gloves and eye protection when handling samples.
- Wash your hands thoroughly after use.
- Refer to the Safety Data Sheet ([www.diesse.it](http://www.diesse.it)) for the safety characteristics of the reagents contained in the kit.
- Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% sodium hypochlorite, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried. All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste. Do not autoclave materials containing sodium hypochlorite.

#### **Analytical warnings**

- Collect the sample taking care not to contaminate the product.
- Do not use the reagents after the expiry date.
- Take all the necessary precautions to avoid contamination of the product during its use.

##### **5. KIT COMPOSITION**

The package is sufficient for 72 extractions.

#### **RELEASING AGENT**

RELEASING BUFFER

**1 x 6,0 ml**

Content: Solution containing sodium hydroxide (<1.5%) and potassium cyanide (<0.03%).

Liquid, to be mixed with the stabilizing buffer.

#### **STABILIZING AGENT**

STABILIZING BUFFER

**1 x 0,220 ml**

Content: Solution containing tris-2 carboxyethyl phosphine (TCEP) (<15%).

Liquid, to be mixed with the releasing buffer.

#### **NEUTRALIZING AGENT**

NEUTRALIZING BUFFER

**1 x 5,0 ml**

Content: Monobasic sodium phosphate salt buffer (<8.5%) able to reduce the pH of the extracted sample, with red dye (0.1%).

Liquid, ready to use.

#### **OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED:**

- Vortex
- Test tubes in which to do sample extraction
- Test tube shaker
- Centrifuge
- Micropipettes capable of accurately withdrawing volumes of 1-1000 µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials

## 6. STORAGE METHOD AND REAGENT STABILITY

The reagents must be stored at 2-8°C.

The expiry date is indicated on the outer label of the package.

The reagents have limited stability after opening:

|                    |         |
|--------------------|---------|
| RELEASING AGENT    | 8 weeks |
| STABILIZING AGENT  | 8 weeks |
| NEUTRALIZING AGENT | 8 weeks |

## 7. SAMPLE TYPE AND STORAGE

The sample is serum obtained from blood collected by normal venous blood sampling and handled with the precautions required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

Fresh serum can be kept for 5 days at 2-8°C; for longer storage periods, freeze at temperatures < -20°C.

The sample can be thawed for a maximum of 3 times. Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

The extracted sample can be used within 8 hours from extraction if kept at 4°C.

It is strongly recommended to avoid any intensive sample processing.

Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

## 8. PROCEDURE

1. Dispense 100 µl of the sample to be extracted into a test tube (not provided).
2. Mix the stabilizing buffer (Stabilizing Agent) and the releasing buffer (Releasing Agent) to prepare a 1/40 extracting solution (Stabilizing/Releasing).  
For example, for the extraction of 6 samples, mix 10 µl of stabilizing buffer and 390 µl of releasing buffer.  
The 1/40 solution is stable at room temperature for 8 hours after preparation.
3. Dispense the previously prepared extraction reagent into the 50 µl test tube.
4. Using the vortex, mix 3 times for 1-2 seconds **immediately after adding the extraction reagent.**
5. Incubate for 15 minutes at room temperature.
6. Dispense the neutralizing buffer into the 50 µl test tube.
7. Using the vortex, mix 3 times for 1-2 seconds **immediately after adding the neutralizing buffer.**
8. Use the extracted sample as indicated in the IFU of the kit to be used.

Should you have any problems, contact Scientific Support:

Phone: 0039 0577 319554  
Fax: 0039 0577 366605  
email: [scientificsupport@diesse.it](mailto:scientificsupport@diesse.it)

## 9. INCIDENT REPORTING

If a serious incident has occurred in relation to this device in the European Union market territory, please report it to the manufacturer and the competent authority of your Member State without delay.

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diesse.it

## FOLATE & VITAMIN B12 SAMPLE PREPARATION

### PURPOSE:

The following illustrated steps were created to help train technicians and simplify the process of sample preparation

### SAMPLES:

Prepare all patients, calibrators and controls before use in the assay

### SAMPLE STABILITY:

Fresh serum can be kept up to 5 days at 2-8°C; for longer periods, please freeze at -20°C or lower. Extracted sample can be used within 8 hours since the extraction, when maintained at 4°C.

### EXTRACTION STEPS

#### STEP 1

- Add 100µl of sample (patient, calibrator or control) into a tube
- Hint: touch side of tube with pipette while dispensing at an angle near the bottom of the tube
- Hint: prepare all samples

#### STEP 2

- Preparation of Extraction Reagent: mix the Stabilizing Agent and the Releasing Agent to prepare a 1/40 solution (Stabilizing/Releasing).  
Ex. For 6 samples, mix 10µl of Stabilizing Agent and 390µl of Releasing Agent
- Add 50µl of Extraction Reagent into sample tubes
- Hint: touch side of tube with pipette while dispensing at an angle slightly above the sample (do not touch the sample to avoid contamination)
- Mix on vortex 1-2 seconds using multiple touch technique (3 x 1-2 seconds), immediately after the addition of the Extraction Reagent
- Hint: mix before moving onto the next tube

#### STEP 3

- Incubate 15 minutes at room temperature

#### STEP 4

- Add 50µl of Neutralizing Agent into sample tubes
- Hint: touch side of tube with pipette while dispensing at an angle slightly above the sample (do not touch the sample to avoid contamination)
- Mix on vortex 1-2 seconds using multiple touch technique (3 x 1-2 seconds), immediately after adding the Neutralizing Agent
- Hint: mix before moving onto the next tube

#### STEP 5

- Use prepared samples accordingly to IFU procedure

diesse.it

|                                                                                     |    |                                                  |    |                                                  |
|-------------------------------------------------------------------------------------|----|--------------------------------------------------|----|--------------------------------------------------|
|    | IT | Data di fabbricazione                            | ES | Fecha de fabricación                             |
|                                                                                     | EN | Date of manufacture                              | FR | Date de fabrication                              |
|                                                                                     | CS | Datum výroby                                     | PT | Data de fabrico                                  |
|                                                                                     | DE | Herstellungsdatum                                | RO | Data fabricației                                 |
|                                                                                     | EL | Ημερομηνία Παραγωγής                             |    |                                                  |
|    | IT | Utilizzare entro                                 | ES | Fecha de caducidad                               |
|                                                                                     | EN | Use By                                           | FR | Utiliser jusque                                  |
|                                                                                     | CS | Použitelné do                                    | PT | Prazo de validade                                |
|                                                                                     | DE | Verwendbar bis                                   | RO | A se folosi până la                              |
|                                                                                     | EL | Ημερομηνία λήξης                                 |    |                                                  |
|    | IT | Nu reutilizați                                   | ES | No reutilizar                                    |
|                                                                                     | EN | Do not reuse                                     | FR | Ne pas réutiliser                                |
|                                                                                     | CS | Nepoužívejte opakovaně                           | PT | Não reutilizar                                   |
|                                                                                     | DE | Nicht wieder verwenden                           | RO | A nu se refolosi                                 |
|                                                                                     | EL | Μην κάνετε επαναληπτική χρήση                    |    |                                                  |
|    | IT | Attenzione, vedere le istruzioni per l'uso       | ES | Atención, ver instrucciones de uso               |
|                                                                                     | EN | Caution, consult accompanying documents          | FR | Attention voir notice d'instructions             |
|                                                                                     | CS | Pozor, čtěte přiložené dokumenty                 | PT | Atenção, consulte a documentação incluída        |
|                                                                                     | DE | Achtung, die Gebrauchsanleitung lesen            | RO | Atenție, consultați documentele însoțitoare      |
|                                                                                     | EL | Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα   |    |                                                  |
|    | IT | Fabbricante                                      | ES | Fabricante                                       |
|                                                                                     | EN | Manufacturer                                     | FR | Fabricant                                        |
|                                                                                     | CS | Výrobce                                          | PT | Fabricante                                       |
|                                                                                     | DE | Hersteller                                       | RO | Producător                                       |
|                                                                                     | EL | Κατασκευαστής                                    |    |                                                  |
|    | IT | Contenuto sufficiente per "n" saggi              | ES | Contenido suficiente para <n> ensayos            |
|                                                                                     | EN | Contains sufficient for <n> tests                | FR | Contenu suffisant pour "n" tests                 |
|                                                                                     | CS | Obsah stačí na <n> testů                         | PT | Conteúdo suficiente para "n" ensaios             |
|                                                                                     | DE | Inhalt reicht für „n“ Tests                      | RO | Conținut suficient pt <n> teste                  |
|                                                                                     | EL | Περιεχόμενο επαρκές για «n» εξετάσεις            |    |                                                  |
|   | IT | Limiti di temperatura                            | ES | Límite de temperatura                            |
|                                                                                     | EN | Temperature limitation                           | FR | Limites de température                           |
|                                                                                     | CS | Teplotní omezení                                 | PT | Limites de temperatura                           |
|                                                                                     | DE | Temperaturgrenzwerte                             | RO | Limită de temperatură                            |
|                                                                                     | EL | Περιορισμοί θερμοκρασίας                         |    |                                                  |
|  | IT | Consultare le istruzioni per l'uso               | ES | Consulte las instrucciones de uso                |
|                                                                                     | EN | Consult Instructions for Use                     | FR | Consulter les instructions d'utilisation         |
|                                                                                     | CS | Čtěte návod k použití                            | PT | Consulte as instruções de utilização             |
|                                                                                     | DE | Die Gebrauchsanleitung lesen                     | RO | Pentru utilizare consultați instrucțiunile       |
|                                                                                     | EL | Συμβουλευτείτε τις οδηγίες χρήσης                |    |                                                  |
|  | IT | Rischio biologico                                | ES | Riesgo biológico                                 |
|                                                                                     | EN | Biological risks                                 | FR | Risques biologiques                              |
|                                                                                     | CS | Biologická rizika                                | PT | Risco biológico                                  |
|                                                                                     | DE | Biologisches Risiko                              | RO | Risc biologic                                    |
|                                                                                     | EL | Βιολογικοί κίνδυνοι                              |    |                                                  |
| <b>REF</b>                                                                          | IT | Numero di catalogo                               | ES | Número de catálogo                               |
|                                                                                     | EN | Catalogue number                                 | FR | Référence du catalogue                           |
|                                                                                     | CS | Katalogové číslo                                 | PT | Referência de catálogo                           |
|                                                                                     | DE | Katalognummer                                    | RO | Număr de catalog                                 |
|                                                                                     | EL | Αριθμός καταλόγου                                |    |                                                  |
| <b>IVD</b>                                                                          | IT | Dispositiv medical pentru diagnosticare in vitro | ES | Producto sanitario para diagnóstico in vitro     |
|                                                                                     | EN | In Vitro Diagnostic Medical Device               | FR | Dispositif médical de diagnostic in vitro        |
|                                                                                     | CS | Lékařské vybavení pro diagnostiku in vitro       | PT | Dispositivo médico para diagnóstico in vitro     |
|                                                                                     | DE | Medizinisches In-vitro-Diagnostikum              | RO | Dispositiv medical pentru diagnosticare in vitro |
|                                                                                     | EL | In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν     |    |                                                  |
| <b>LOT</b>                                                                          | IT | Codice del lotto                                 | ES | Código de lote                                   |
|                                                                                     | EN | Batch code                                       | FR | Code du lot                                      |
|                                                                                     | CS | Kód šarže                                        | PT | Código do lote                                   |
|                                                                                     | DE | Chargennummer                                    | RO | Lot                                              |
|                                                                                     | EL | Αριθμός Παρτίδας                                 |    |                                                  |
| <b>CE</b>                                                                           | IT | Marcajul de conformitate CE                      | ES | Marcado CE de conformidad                        |
|                                                                                     | EN | CE marking of conformity                         | FR | Marquage de conformité CE                        |
|                                                                                     | CS | Oznakowanie zgodności CE                         | PT | Marcação CE de conformidade                      |
|                                                                                     | DE | CE-Konformität Skeneichnung                      | RO | Marcajul de conformitate CE                      |
|                                                                                     | EL | Σημάνση συμμορφωση CE                            |    |                                                  |