W.F.-DAT

REF 26020

DIESSE
DIESSE

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INSTRUCTIONS FOR USE

W.F.-DAT

For the determination of anti-Proteus specific agglutinating antibodies

For in vitro Diagnostic use Only

1. INTENDED USE

Method for the determination of specific agglutinating antibodies against Proteus in human serum.

2. INTRODUCTION

Stained suspensions of Proteus are used for the diagnosis of rickettsiae infection (Weil-Felix test).

The test is based on the fact that the somatic antigens of some strains of Proteus are identical to those of some species of Rickettsie. The sera from subjects infected by some forms of Rickettsie agglutinate suspensions of Proteus strains.

3. PRINCIPLE OF THE METHOD

The serum sample is mixed with a bacterial suspension at the optimal concentration for the semi-automated test on the AUTO-DAT instrument and / or in slide.

A visible agglutination will appear if specific antibodies against bacterial antigens at significant titer are present. Positive samples must then be tested in the tube test, using one drop of the same suspension used for the slide.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

All material of biological origin must be considered potentially infectious. All reagents and samples must be handled in accordance with the safety precautions adopted in the laboratory.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

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1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to

clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to room temperature (18-30°C) before use.

1. After use, return reagents to 2-8°C.
2. Do not use the reagents after the expiry date.
3. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
4. Do not modify the Test Procedure.

5. KIT COMPOSITION AND REAGENT PREPARATION

Each bacterial suspension is sufficient for 60 tests.

BACTERIAL SUSPENSION BACTERIAL SUSPENSION 1 x 1.5 mL

Of the following suspensions:

- Proteus OX2
- Proteus OX19
- Proteus OXK

Contents: Suspension of stained, killed bacteria at an optimized concentration to perform the test. Contains sodium azide 0.09% as preservative.

CONTROL + POSITIVE CONTROL 1 x 0.5 mL

Contents: Polyvalent positive control serum of animal origin containing 0.09% sodium azide as preservative. Contains antibodies that agglutinate all bacterial suspensions of the kit.

BRING THE REAGENTS TO ROOM TEMPERATURE BEFORE USE.

OTHER MATERIALS:

- Slides (3)
- Mixing sticks (100)

MATERIALS REQUIRED BUT NOT PROVIDED:

- Pipettes with disposable needles
- Physiological saline
- 96-well microplates with flat bottom, Auto-DAT Microplate **REF** 26005
- AUTO-DAT Instrument (**REF** 26000)

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the kit label.

The stability of the reagents does not change after opening the bottle, provided that the user is careful to keep the product safe from possible microbial contamination.

7. SPECIMEN COLLECTION AND STORAGE

Fresh serum may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C; it can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Use immediately after thawing.

8. ASSAY PROCEDURE

Bring the reagent to room temperature and mix the vial of bacterial suspension carefully but not vigorously.

INSTRUMENTAL APPLICATION WITH AUTO-DAT

Pipette 25 µL of sample and 25 µL of reagent into a microplate well (REF 26005). Insert the microplate into the instrument and press START. The instrument, at the end of the reaction (which lasts about 15 minutes) will print the result in titer.

SCREENING TEST

If the AUTO-DAT test cannot be performed, reagents may be used for a screening test.

1. Distribute 25 µL of serum in a circle on the slide.
2. Add 25 µL of suspension and mix.
3. Rotate the slide.
4. Read the results after exactly 1 minute.

Positive samples will give agglutination while negative samples will remain homogeneous.

TUBE TITRATION

1. Dilute the sample 1:20 with physiological saline (0.1 mL serum + 1.9 mL saline).
2. Take 1 mL from the first tube and make serial dilutions with 1 mL of saline, as shown in the scheme.
3. Add one drop of suspension to each tube, mix and incubate overnight at 37°C (see the scheme shown below).
4. The blank is performed by mixing the suspension with saline.

Tube	1	2	3	4	5	6	7
Saline solution (ml)	1.9	1	1	1	1	1	1
Serum (ml)	0.1						
Mix and transfer (ml)		1	1	1	1	1	1
							
Dilution/Titer	1/20	1/40	1/80	1/160	1/320	1/640	1/1280

9. TEST VALIDATION ON AUTO-DAT

Using 25 µL of **Positive Control** serum, perform the test on AUTO-DAT and compare the result with the titers reported on the certificate of analysis.
+/- one titer dilution deviations are accepted.

If the result is different from the expected, contact the Scientific Support:

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

Titers 1:80 are considered suspect; higher titers are diagnostically significant:

Disease	Agent	Weil-Felix Reaction		
		OX19	OX2	OXK
Exanthematic typhoid fever	R.prowazekii	>1:80	1:80	-
Tifo murino	R.mooseri	>1:80	1:80	-
Rocky Mountain spotted fever	R.rickettsii	1:80	>1:80	1:40
Boutonneuse fever	R.conori	1:80	>1:80	-
Tsutsunga-mushi fever	R.orientalis	-	-	>1:80

11. LIMITATIONS

Not all rickettsie give a cross-reaction with Proteus. R. Burnetii, the etiological agent of Q-Fever, does not react.

12. METHOD COMPARISON – DIAGNOSTIC SENSITIVITY AND SPECIFICITY

The three reagents were compared with other commercial products, using sera from a hospital routine practice (89), 16 of which had a request for the Weil Felix test.

Only a sample was positive, in agreement with the other method; all the others were negative in both methods.

Diagnostic specificity was therefore 100%. The sensitivity was also 100%, but this data is only of relative value due to the low number of positive samples found during the experimentation.

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