ESR CONTROL

REF 10430

REF 10434

DIESSE

	Capitolo Section Κεφάλαιο Kapitel Capítulo Capítulo Chapitre Rozdział Capitol Fejezet Poglavlje
Modifiche introdotte nella revisione corrente Changes introduced in the current revision Αλλαγές που εισήχθησαν στην τρέχουσα αναθεώρηση Änderungen in der aktuellen Revision eingeführt Cambios introducidos en la revisión actual Alterações introduzidas na revisão atual Modifications introduites dans la révision actuelle Zmiany wprowadzone w aktualnej rewizji Modificări introduse în actuala revizuire A jelenlegi felülvizsgálatban bevezetett változások Promjene uvedene u trenutačnoj reviziji	3

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10. BIBLIOGRAFIA

- Clinical and Laboratories Standard Institute, H26-A2, Validation, Verification and Quality assurance of automated Haematology analyser. Approved standard – second edition.

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REF 10430 4 x 9 ml (2 x Level I + 2 x Level II)

REF 10434 2 x 9 ml (1 x Level I + 1 x Level II)

**For monitoring the precision of instruments
used for the determination of Erythrocyte
Sedimentation Rate (ESR)**

For *in vitro* diagnostic use only

1. INTENDED PURPOSE

Dual-level control material to be used for monitoring the precision of DIESSE instruments used for the determination of Erythrocyte Sedimentation Rate (ESR).

It must only be used by professional laboratory personnel.

2. INTRODUCTION AND PRINCIPLE OF THE METHOD

It is good laboratory practice to use control materials to check the precision of the instruments and methods used.

ESR Control is a control material to be used in the Mini-VES, Ves-Matic and VES-MATIC Original instruments for monitoring their precision

3. REAGENTS

Contents: Suspension of stabilized human red blood cells in a preservative medium.

Preparation: ready for use

Stability: Store the product at 2-10°C.

1. The unopened product is stable up to the expiry date and is stable for 95 days at 2-30°C after opening.
Close the vial carefully after use and store in an upright position.
2. Once prepared in the tube, the product is stable for 7 days at 18-30°C.

4. PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and found negative with FDA-approved methods, for the presence of HbsAg and for the presence of anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee of the absence of infective agents, all materials of human origin must be considered potentially infectious.

All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: material of biological origin once used must be treated as infectious residuals and eliminated according to law.

5. WARNINGS

1. Check that the batch number on the vial corresponds to the lot number reported on the Certificate of Analysis (CoA) enclosed with the present Instructions for Use.
2. Avoid exposure to direct sunlight. After use, clean the screw part of the cap and the vial, and close the cap tightly; store in an upright position.
3. **DO NOT FREEZE AND DO NOT EXPOSE TO TEMPERATURES OVER 30°C.**
4. Good function of the product is guaranteed only if it is stored and used as described in the present Instructions for Use.
5. **It may occur that in some laboratories the results obtained do not fall within the declared range. In these cases it is advisable to perform several tests following the laboratory's own method and thus establish a personal range of acceptability. The laboratory can thus consider this range valid until there is a change of lot, or a change in the execution method for the test.**
6. The supernatant in the unmixed vial may appear clear red; this is normal and does not indicate that the product is unsuitable for use.
7. A loss of colour of the product may be a sign of degradation. **Do not use the product if a loss of quality is suspected.**


6. PROCEDURE

IMPORTANT: Incomplete resuspension can influence the sample preparation and the remaining product left in the vial.

1. Remove the vial from the refrigerator and take to room temperature (20-30 minutes).
2. To mix:
 - a. Prior to first time use carefully resuspend the product contained in the unopened vial using vortex (up to 60 seconds).
 - b. Hold the vial in vertical position and roll it between the palms of the hands for 15-20 seconds.

- c. Keep mixing the vial, turning it completely and rapidly upside down, for 20 times at least.

- d. Check the bottom of the vial to verify that the product is completely resuspended and repeat steps 2b and 2c if necessary.
Allow bubbles to disperse and mix again by inversion 8-10 times immediately before the sampling.
 - e. Divide the sample in aliquots immediately after mixing.
3. Transfer 1 mL (± 0.2) of suspension into a VES-TEC, VACU-TEC or VACU-CODE cuvette without emptying it of the anticoagulant (Mini-VES and Ves-Matic instruments).
 4. Transfer 1 mL (± 0.2) of suspension into a VES-TEC or VACU-TEC cuvette (as indicated in the user manual) without emptying it of the anticoagulant (VES-MATIC Original instrument).
 5. Always apply the barcode label relating to the ESR Control used on the VES-TEC or VACU-TEC cuvette: VES-MATIC Original Level I/Level II. For proper application of the label, refer to the user manual of the VES-MATIC Original instrument.
 6. Insert the cuvette or tube in the instrument like any other sample (see the operating manual of the instrument).
 7. Compare the found result with the value reported on the CoA.

7. EXPECTED RESULTS

The ranges of acceptability for the different instruments are reported on the enclosed CoA.

If the results of the control are not within the range reported on the CoA but are reproducible, the relative working conditions may be the cause of the values not falling within the given range. In this case it is advisable to perform several runs following the laboratory's usual method and thus determine the relative acceptable range.

If the results of the control are out of the range reported in the CoA and are not reproducible, it is advisable to repeat the control using a new batch.

It must be remembered that the ESR Control is for use above all to check the precision of the system (see the paragraph on WARNINGS) and not as a reference standard.

8. PERFORMANCE

The minimum and maximum values reported in the table on the CoA are derived from replicate analysis on automated methods.

9. LIMITATIONS

The ESR Control is not intended for use as a reference standard. The ESR Control cannot be used to control any hematological parameters other than the ESR.

10. REFERENCES

1. Clinical and Laboratories Standard Institute, H26-A2, Validation, Verification and Quality assurance of automated Haematology analyzer. Approved standard – second edition.

11. INCIDENT REPORTING

If a serious incident has occurred in relation to this device in the European Union market territory, please report it to the manufacturer and the competent authority of its Member State without delay.

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